

TN-1383

LC-MS/MS Analysis of Benzodiazepines in Urine with Room Temperature Hydrolysis and One-Step Clean Up Using Beta™-Gone

Andre Sutka¹, Janet Jones², Elías Villalobos², Crystalyn Consolacion³, Stephanie J. Marin³

¹United State Drug Testing Laboratory, ²Kura Biotec®, ³Phenomenex, Inc.



Introduction

Benzodiazepines are a class of depressants that are used for sedation, to relieve anxiety, and reduce muscle spasms and seizures. Commonly prescribed benzodiazepines include Valium® (Diazepam) Xanax® (alprazolam), Halcion® (triazolam), Ativan® (lorazepam) and Klonopin® (clonazepam). These controlled substances can be habit forming. Misuse and dependence are a concern with these drugs. Toxicology labs quantitate benzodiazepines for drug monitoring, drug facilitated crime and impairment investigations, and cause of death determinations. Many labs use LC-MS/MS for analysis of benzodiazepines, due to its sensitivity and specificity. Modern LC-MS/MS instrumentation and chromatography using superficially porous “core-shell” LC columns can generate fast methods for high-throughput drug analysis workflows.

Urine is one of the most commonly used specimens for drug testing methods. Benzodiazepines are extensively metabolized via methylation and hydroxylation. They are also heavily conjugated prior to excretion in the urine. Therefore, most urine methods hydrolyze the urine prior to extraction and analysis by LC-MS/MS to convert glucuronide conjugates back to the parent drug or metabolite. Sample clean up to remove residual hydrolysis enzyme and other matrix components improves quantitation by LC-MS/MS. SPE and other extraction methods require multiple steps and significant method development time. They also require more time in a drug analysis workflow. Beta-Gone is a simple, single step extraction that cleans up hydrolyzed urine specimens. It requires no method development and can be easily incorporated into a urine toxicology drug testing workflow.

Sample Preparation

Eight calibrators at 20, 40, 100, 200, 500, 1000, 2000 and 2500 ng/mL, and three QC samples at 50, 125, and 800 ng/mL prepared in drug free urine were analyzed. Internal standard concentrations were 100 ng/mL for all analytes. All calibration curves used a quadratic fit, except lorazepam which used a linear calibration fit. All calibration curves were weighted 1/x. MRM parameters are listed in **Table 1**. Oxazepam glucuronide was used to measure hydrolysis efficiency.

For each sample, 50 µL of internal standard solution were added to 100 µL of urine. Next, 200 µL Kura® B-One® enzyme were added, and the samples were incubated for 30 min at room temperature. Samples were filtered through Beta-Gone 1 mL cartridges (part number [8B-S139-TAK](#)) into a vial prior to analysis by LC-MS/MS.

LC Conditions

Column: Kinetex™ Phenyl-Hexyl
Dimensions: 50 x 4.6 mm, 2.6 µm
Part No.: [00B-4495-E0](#)
Mobile Phase: A= 10mM Ammonium Formate and 0.05% Formic Acid
 B= Methanol and 0.05% Formic Acid

Gradient:	Time (min)	% B
	0.0	10
	4.5	100
	5.5	100
	6.1	10
	7.5	10

Flow Rate: 0.6 mL/min
Injection Volume: 3 µL
Temperature: 40 °C
LC System: Agilent® 1200 Series LC
Detection: LC-MS/MS
Detector: Sciex® 5500 Triple Quad™

MS/MS Conditions

Ion Source: Turbo Spray
Polarity: ESI Positive Mode
Source Temperature: 650 °C
GS1: 60
GS2: 50
CUR: 40
IS: 2500 V



Table 1. MRM Conditions Retention Times, Recovery and Matrix Effect for Benzodiazepines


Analyte or Internal Standard	Retention Time (min)	% Recovery	Matrix Effect	Q1	Q3	Dwell	DP	CE	CXP
7-aminoclonazepam-d4	3.69	-	-	290.10	254.00	15.00	65.00	29.00	13.00
				290.10	226.00	15.00	65.00	26.00	13.00
7-aminonitrazepam	3.48	107%	44%	252.10	121.00	15.00	40.00	23.00	35.00
				252.10	146.00	15.00	40.00	25.00	53.00
7-aminoclonazepam	3.7	103%	41%	286.10	250.00	15.00	80.00	19.00	13.00
				286.10	222.00	15.00	80.00	27.00	13.00
7-aminoflunitrazepam	3.94	94%	29%	284.10	135.10	15.00	51.00	23.00	39.00
				284.00	227.20	15.00	51.00	23.00	49.00
alpha-OH Alprazolam-d5	4.74	-	-	330.10	302.20	15.00	51.00	30.00	31.00
				330.10	284.00	15.00	51.00	36.00	31.00
alpha-OH Midazolam	4.8	96%	30%	342.10	203.10	15.00	60.00	26.00	10.00
				342.10	168.10	15.00	60.00	43.00	11.00
alpha-OH Triazolam	4.62	106%	49%	359.10	331.10	15.00	90.00	27.00	12.00
				359.10	239.00	15.00	90.00	50.00	12.00
alpha-OH Alprazolam	4.75	103%	34%	325.10	297.10	15.00	60.00	30.00	15.00
				325.10	299.10	15.00	60.00	25.00	15.00
Oxazepam-d5	4.74	-	-	292.20	246.20	15.00	30.00	20.00	35.00
				292.20	274.00	15.00	30.00	19.00	35.00
Oxazepam	4.75	100%	27%	287.10	241.10	15.00	76.00	31.00	18.00
				287.10	269.20	15.00	76.00	14.00	22.00
Triazolam-d4	4.75	-	-	347.20	312.20	15.00	90.00	30.00	50.00
				347.20	243.00	15.00	90.00	40.00	50.00
Lorazepam	4.68	109%	34%	321.10	275.10	15.00	80.00	30.00	10.00
				321.10	229.10	15.00	80.00	30.00	10.00
Temazepam-d5	4.91	-	-	306.00	260.10	15.00	70.00	15.00	35.00
				306.00	177.00	15.00	70.00	40.00	35.00
2-OH Ethylflurazepam	4.74	100%	29%	333.10	211.20	15.00	51.00	25.00	53.00
				333.10	109.20	15.00	56.00	29.00	41.00
Temazepam	4.92	102%	34%	301.10	255.10	15.00	70.00	15.00	8.00
				301.10	177.10	15.00	70.00	40.00	10.00
Nordiazepam-d5	4.95	-	-	276.10	140.00	15.00	100.00	30.00	15.00
				276.10	213.10	15.00	100.00	30.00	15.00
Nordiazepam	4.96	106%	22%	271.10	140.10	15.00	100.00	30.00	15.00
				271.10	165.10	15.00	100.00	30.00	15.00

Results and Discussion

Analyte retention time, matrix effect and recovery are reported in **Table 1**. Matrix effect was between 29% and 49%. Recovery ranged from 96% to 109% with %CVs of $\pm 15\%$ for analytes with a corresponding isotopically labeled internal standard, and within $\pm 20\%$ for analytes without a matching internal standard ($n = 12$). **Figure 1** shows the separation of 11 benzodiazepines using a the Phenomenex Kinetex™ Phenyl-Hexyl LC column. Hydrolysis efficiency was 100%.

Conclusions

The fast 30-minute room temperature hydrolysis using Kura B-One® and a simple one step extraction using Beta-Gone™ produced very clean samples with high recovery, low matrix effect and excellent hydrolysis efficiency. The fast 5-minute LC separation and LC-MS/MS analysis using a SCIEX® 5500 MS/MS system provided accurate and reliable quantitation of benzodiazepines in urine.

Figure 1. Separation of benzodiazepines on a Kinetex Phenyl-Hexyl column

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale. Kinetex, Beta-Gone, and BE-HAPPY are trademarks of Phenomenex. B-One is a trademark of Kura Biotech. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademark, and Triple Quad is a trademark of AB SCIEX Pte. Ltd. XANAX and Halcion are registered marks of Pharmacia & Upjohn Company LLC, which is now part of Pfizer; Valium and Klonopin are registered marks of Hoffmann-La Roche; and Ativan is a registered mark of Valeant Pharmaceuticals International, Inc. Phenomenex is in no way affiliated with Agilent Technologies, Inc. or the pharmaceutical companies for Xanax, Halcion, Valium, Klonopin, and Ativan. Comparative separations may not be representative of all applications. *FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.* © 2025 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists

