


TN-1384

Separation of Irinotecan Hydrochloride and its Chiral Impurity per USP Monograph under Allowable Adjustments

Rajesh Jangam¹, Rajesh Babu Dandamudi, PhD¹, and Sean Orłowicz²

¹ India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India

² Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Introduction

Irinotecan hydrochloride is a chemotherapeutic agent classified as a DNA topoisomerase I inhibitor. Irinotecan is frequently administered alongside other drugs for the treatment of metastatic colon and rectal cancer, and occasionally for small cell lung cancer. The USP monograph specifies the use of an L40 column (4.6-mm x 25-cm; 10-μm), which is a Cellulose tris (3,5 dimethylphenylcarbamate) polymer column. In this technical note we demonstrate the utilization of Lux™- Cellulose-1, (5 μm 150 x 4.6 mm) as an alternative column for the separation of Irinotecan Hydrochloride and its chiral impurity per the allowable adjustments outlined in USP General Chapter <621>. The system Suitability per USP monograph is resolution NLT 2.5 between Irinotecan related compound D and Irinotecan, a percent relative standard deviation (%RSD) of NMT 5.0%, for standard solution and the Irinotecan related compound D peak should be visible for the sensitivity solution. All the system suitability requirements have been met by this column.

LC Conditions

Column: Lux™ Cellulose-1, 5 μm (Part No: [00F-4459-E0](#))
Dimensions: 150 x 4.6 mm
Mobile Phase-A: Hexane: Dehydrated Alcohol: Diethylamine in ratio of 250:250:1 ml, v/v/v.
Diluent: Dehydrated Alcohol: Diethylamine in ratio of 250:1
Isocratic Elution:

Run Time (min)	%	Mobile phase-A
25		100

Flow Rate: 1.5 mL/min
Injection Volume: 12 μL
Temperature: 25°C
LC System: Waters® Arc HPLC with PDA
Detection: UV @ 370 nm

Standard and System suitability solutions

Standard stock solution: 10 mg/ml of USP Irinotecan Related compound-D RS in Diluent.

System suitability solution: 0.1mg/mL Each Irinotecan Hydrochloride and Irinotecan Related compound-D in Diluent

Identification Solution: 1mg/mL of USP-Irinotecan Hydrochloride RS in Diluent.

Standard solution: 1.5 μg/mL of USP Irinotecan Related compound-D in Diluent, from the Standard stock solution.

Sensitivity solution: 0.5 μg/mL of USP Irinotecan Related compound-D reference standard in Diluent, from the Standard stock solution.

Sample solution: 1mg/mL Irinotecan Hydrochloride in Diluent.

Standard and System suitability solutions

System suitability solution and standard solution were run on the Lux Cellulose-1, 5 μm 150 x 4.6 mm column.(Figure-1 on Page no:2)

Figures 1, 3 and 5 show system suitability solution chromatograms and Peak tables 1,3,5 on page no:2 show system suitability criteria Retention time, RRT, Resolution and S/N ratio values from sensitivity solution.

Figures 4 & 5, show standard solution chromatogram and Figure-5 table shows relative standard deviation (% RSD) of the six replicate injection of standard solution.

All solutions were prepared as indicated in the USP monograph for the Irinotecan Hydrochloride

USP Irinotecan HCl (Catalog No.1347609), USP Irinotecan Related compound-D (Catalog No.1347653) were purchased from USP.

Allowable column Adjustments: L/dp Ratio- 25 % to 50%


Column	Length (mm)	ID (mm)	dp (μm)	L/dp	Allowable Range (-25% to +50%)
Original	250	4.6	10	25	18.75 to 37.5
Alternative	150	4.6	5	30	Allowed

Conclusions

This USP monograph outlines the utilization of an L-40 column, measuring 10μm - 250 x 4.6mm. The 10μm particle columns are older generation columns with lower efficiencies compared to the 5μm. Utilizing the USP Allowable Adjustments Calculator, we have used Lux- Cellulose-1 column (5 μm, 150 x 4.6 mm) to evaluate its applicability. Method adjustments have also been considered in accordance with the Allowable Adjustments due to the change in the LC column.


The performance of this column was assessed based on the Relative Retention Time (RRT) of relative compound-D, resolution, signal-to-noise ratio, and the relative standard deviation (RSD) of the standard solution, all of which fell within acceptable limits. This indicates that the Lux-Cellulose-1 column (5 μm, 150 x 4.6 mm) can be used as an alternative to the 10 μm 250 x 4.6 mm column referenced in the monograph.

Figure 1. Irinotecan Hydrochloride System suitability solution




Peak No.	Analyte	Retention Time	Area	% Area	RT Ratio	Resolution	USP Tailing Factor
1	Irinotecan Related Compound-D	4.805	1603142	50.03	0.67		1.2
2	Irinotecan	7.121	1595875	49.97	1.00	2.86	1.2

Figure 2. Irinotecan Hydrochloride USP- Identification solution




Peak No.	Analyte	Retention Time	Area	% Area
1	Irinotecan	7.064	1537204	100.00

Figure 3. Irinotecan Hydrochloride USP- Sensitivity solution(0.5µg/ml)




Peak No.	Analyte	Retention Time	Area	S/N Ratio
1	Irinotecan Related Compound-D	4.826	9941	10

Figure 4. Irinotecan Hydrochloride USP- Standard solution (1.5µg/ml)



Peak No.	Analyte	Retention Time	Area
1	Irinotecan Related Compound-D	4.817	26688

Figure 5. Irinotecan Hydrochloride USP- Standard solution (1.5µg/mL), replicate injections (N=6, overlaid) chromatogram.



Injections (N=6)	Retention Time	Area
Mean	4.824	26608.66
SD	0.005	94.8
%RSD	0.1	0.35

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale. Lux and BE-HAPPY are trademarks of Phenomenex. Waters and Waters Acquity Arc are registered trademarks of Waters Corporation. Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Waters Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2025 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists

