


TN-1395

Aggregate Analysis of Semaglutide Using a Biozen™ dSEC-1 SEC Column

Sujatha Chilakala, Ph.D.¹ and Lucia Geis-Asteggianti, Ph.D.²

¹ Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

² Phenomenex, Inc., Via M. Serenari 15/D, Castel Maggiore, BO 40013 ITALY



Introduction


Semaglutide is an agonist of the glucagon-like peptide-1 (GLP-1) receptor, utilized in the management of type 2 diabetes and obesity. Semaglutide functions by increasing insulin secretion, lowering glucagon levels, and delaying gastric emptying. These actions result in better blood sugar regulation and facilitate weight loss. It emulates the effects of the endogenous hormone GLP-1, which is responsible for appetite control and glucose metabolism¹. As a synthetic peptide therapeutic, Semaglutide is prone to aggregation, which can compromise its safety, efficacy, and shelf-life.

Size Exclusion Chromatography (SEC) is a key analytical technique employed to monitor aggregation by separating molecular species based on size. However, SEC of peptides and small proteins can be challenging due to their marked susceptibility to non-specific interactions. Hydrophobic and/or ionic interactions may compromise chromatographic performance, leading to issues such as peak tailing or analyte adsorption. Due to Semaglutide's hydrophobic nature, SEC analysis often requires organic solvents and acidic conditions to ensure accurate resolution.

In this technical note, we demonstrate the successful performance of Biozen dSEC-1 SEC columns for aggregate analysis of Semaglutide. The hydrophilic nature of the stationary phase imparts inertness to the media, reducing the need for organic solvents in the mobile phase. Thus, improving chromatographic results and assisting in robust and reproducible method development.

Sample Preparation

Semaglutide: OZEMPIC® 8 mg in 3 mL contains the following 30 amino acids, glucagon-like peptide-1 analogue drug.



*Aib = α-aminoisobutyric acid

Image modified from <https://pdb101.rcsb.org/global-health/diabetes-mellitus/drugs/incretins/drug/semaglutide/semaglutide>

Sample aggregation was further induced by exposure to ambient light for 3 days. All samples were injected before and after light exposure, injecting the commercial solution directly at volumes of 7.4 µL to achieve a load of 20 µg of the Semaglutide on column. All analyses were done in triplicate.

LC Conditions

Column: Biozen 1.6 µm dSEC-1, 90 Å
Dimensions: 150 x 4.6 mm
Part No.: [00F-4801-E0](#)
Mobile Phase: (60:40, v/v) 1X PBS*:Acetonitrile
Isocratic: 12 minutes
Flow Rate: 0.40 mL/min
Injection Injection volume selected to introduce 20 µg on
Volume: column
Temperature: 40 °C
LC System: Waters ACQUITY® UPLC H-Class
Detection: UV @ 280 nm

*1X PBS (Phosphate Buffered Saline) contains 137 mM NaCl, 2.7 mM KCl, 10 mM Na2HPO4, and 1.8 mM KH2PO4.

Results and Discussion

In this study, we evaluated key chromatographic performance parameters for both the monomer and aggregates components in a commercial Semaglutide injection solution, using a Biozen dSEC-1 SEC column. Analyses were conducted on both the original product and a light-stressed sample exposed for three days. Effective separation of the monomer and aggregate peaks was achieved. A summary of the main chromatographic results is presented in **Table 1**.

Table 1. Summary of observed retention times, peak areas, % of the species present, resolution between dimer and monomer for Semaglutide before and after photoinduced aggregation (n=3).

Semaglutide		Biozen dSEC-1 (150 x 4.6 mm, 1.6 µm)		
		Aggregate 1	Aggregate 2	Monomer
Before exposure to light	Retention Time (min)	-	2.71	3.03
	Area (%)	-	0.18	99.8
	Resolution	-	-	3.16
After exposure to light	Retention Time (min)	2.61	2.71	3.03
	Area (%)	0.07	0.88	99.0
	Resolution	-	1.21	3.36

Using a Biozen dSEC-1 SEC column (150 × 4.6 mm, 1.6 µm), three peaks corresponding to the monomer and 2 aggregates were successfully separated in a 12-minute run. Resolution of 3.16 and 3.36 was achieved between the monomer and closest eluting aggregate for samples before and after photolytic stress, respectively (**Figure 1**). For both samples, retention time and peak area reproducibility are observed with RSD of no more than 0.28 % and 4.02 %, respectively. Good sensitivity is clearly achieved, being able to detect and quantify the aggregate in final product with S/N of 26, 26126 for aggregate 2 and monomer, respectively. The S/N ratios are even higher for the stressed sample, reaching 14, 99 and 19212 for aggregate 1, aggregate 2 and monomer, respectively.




Figure 1. Chromatogram of Semaglutide analyzed using a Biozen dSEC-1 SEC column (150 x 4.6 mm, 1.6 μ m). Traces of representative injections were made for samples before (blue) and after (red) exposure to light.

Conclusion

We demonstrated that Biozen dSEC-1 SEC columns are suitable size exclusion columns for Semaglutide aggregation studies, showing good resolution between aggregate and monomer and excellent retention time and peak area reproducibility. This was achieved without the need of exposing the sample to mobile phases with high organic content or added acidic modifiers.

References

1 Sorli, C., Harashima, S. I., Tsoukas, G. M., Unger, J., Karsbøl, J. D., Hansen, T., & Bain, S. C. (2017). Efficacy and safety of once-weekly semaglutide monotherapy versus placebo in patients with type 2 diabetes (SUSTAIN 1): a double-blind, randomised, placebo-controlled, parallel-group, multinational, multicentre phase 3a trial. *The lancet Diabetes & endocrinology*, 5(4), 251-260. [https://doi.org/10.1016/S2213-8587\(17\)30013-X](https://doi.org/10.1016/S2213-8587(17)30013-X)



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Biozen is a trademark of Phenomenex. ACQUITY is a registered trademark of Waters Corp. OZEMPIC is a registered trademark of Novo Nordisk A/S.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2025 Phenomenex, Inc. All rights reserved

