
TN-1413

Apixaban Assay, Organic Impurities and Limit of Apixaban Related Compounds -H, -G, and -F per USP monograph

 Trivikram Reddy Gundala, PhD¹, Rajesh Babu Dandamudi, PhD¹, and Sean Orlowicz²
¹India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India

²Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

Apixaban is an anticoagulant medication used to treat and prevent blood clots and stroke in people with nonvalvular atrial fibrillation through directly inhibiting factor Xa.

The USP monograph is available for Apixaban Assay, Organic Impurities, and Related Compounds. The column mentioned in the monograph for the determination of Assay and Limit of Apixaban Related compound-H, Apixaban Related compound-G and Apixaban Related compound-F is L1, 250 x 4.6 mm, 5 μm packing. The column mentioned in the monograph for determination of Organic Impurities is L7 column, 250 x 4.6 mm, 5 μm packing.

In this technical note, we demonstrate the usefulness of the fully porous Luna C18(2), the thermally modified Luna Omega C18, and the superficially porous Kinetex C18 columns as alternate columns in determination of Assay and Limit of Apixaban Related compound-H, Apixaban Related compound-G and Apixaban Related compound-F per Apixaban USP monograph. For the separation of Organic Impurities in Apixaban, the fully porous Luna C8(2) and the superficially porous Kinetex C8 were evaluated.

System suitability for Apixaban assay is a tailing factor no more than (NMT) 2.0 and a %RSD of NMT 0.73 % from the Standard Solution. For Organic Impurities, system suitability is a resolution no less than (NLT) 2.0 between Apixaban Related Compound-D and Related Compound-E from the System Suitability Solution, a tailing factor NMT 2.0 and a %RSD NMT 5.0 % from the Standard Solution, and a Signal-to-Noise ratio (S/N) NLT 10 from the Sensitivity Solution. The system suitability for the Limits of the Related Compounds is a resolution NLT 1.5 between Related Compound-H and Related Compound-G, and between Related Compound-G and Related Compound-F, and a %RSD NMT 5.0 % for each impurity from the Standard Solution.

All the solutions were prepared as indicated in the monograph of Apixaban. The following Reference Standards were purchased locally from Clearysynth Labs Limited.

- Apixaban (CAS Number: 503612-47-3)
- Apixaban Related Compound-D (CAS Number: 1351611-14-7)
- Apixaban Related Compound-E (CAS Number: 1074365-84-6)
- Apixaban Related Compound-F (CAS Number: 536759-91-8)
- Apixaban Related Compound-G (CAS Number: 881386-12-5)
- Apixaban Related Compound-H (CAS Number: 1421823-20-2)

LC Conditions: Assay

Columns: Luna 5 μm C18(2) (00G-4252-E0)
Luna Omega 5 μm C18 (00G-4785-E0)
Kinetex 5 μm C18 (00G-4601-E0)

Dimensions: 250 mm x 4.6 mm

Mobile Phase: Buffer / Acetonitrile 48:52 v/v
(Buffer: 0.77g Ammonium Acetate in Water)

Diluent: Water / Acetonitrile 50:50 v/v

Gradient: Isocratic

Flow Rate: 1.0 mL/min

Injection Volume: 10 μL

Temperature: 45 °C

LC System: Waters Arc HPLC with PDA

Detection: UV @ 280 nm

LC Conditions: Organic Impurities

Columns: Luna 5 μm C8(2) (00G-4249-E0)
Kinetex 5 μm C8 (00G-4608-E0)

Dimensions: 250 mm x 4.6 mm

Mobile Phase: A: 0.1 % Trifluoroacetic Acid in Water
B: Acetonitrile

Diluent: Water / Acetonitrile (60:40, v/v)

Gradient:	Time(min)	%B
	0	30
	10	30
	30	70
	35	70
	40	30
	50	30

Flow Rate: 1.0 mL/min

Injection Volume: 10 μL

Temperature: 25 °C

Autosampler Temperature: 5 °C

LC System: Waters Arc HPLC with PDA

Detection: UV @ 280 nm

LC Conditions: Limit of Related Compounds

Columns: Luna 5 μm C18(2) (00G-4252-E0)
Luna Omega 5 μm C18 (00G-4785-E0)
Kinetex 5 μm C18 (00G-4601-E0)

Dimensions: 250 mm x 4.6 mm

Mobile Phase: MP-A: Buffer: Acetonitrile 80:20
MP-B: Acetonitrile
(Buffer: 0.77g Ammonium Acetate in Water, adjust pH 6.6 with Ammonia Water or Acetic Acid)

Diluent: Acetonitrile

Gradient:	Time(min)	%B
	0	20
	5	20
	25	80
	27	20
	40	20

Flow Rate: 1.0 mL/min


Injection Volume: 10 μL

Temperature: 30 °C

LC System: Waters Arc HPLC with PDA

Detection: UV @ 280 nm

Figure 1. Structure of Apixaban



Results and Discussion

The USP method for assay of Apixaban was successfully verified on the Luna C18(2), the Luna Omega C18, and Kinetex C18 columns (Figure 2). All the system suitability parameters for Assay were well within the acceptance criteria in all the columns in terms of % RSD and Tailing. Apixaban peak was observed at similar RT ~3.5 minutes both the Luna C18(2) column and the Luna Omega C18 column. Whereas on Kinetex C18 column, the peak was eluted at RT ~2.8 minutes. %RSD for Apixaban peak was <0.05 % and provided excellent reproducibility on all columns. In addition, the tailing factor for the Apixaban peak was observed slightly higher on the Kinetex C18 column than the Luna C18(2) and Luna Omega C18 columns.


The USP method for Organic Impurities in Apixaban was successfully verified on the Luna C8(2) and the Kinetex C8 columns. Figure 3 shows the effective separation of Apixaban Related Compound-D and Apixaban Related Compound-E from Apixaban on the Luna C8(2) and Kinetex C8 columns. Figure 4 corresponds to chromatograms of the Sensitivity Solution showing the excellent S/N ratio on both the columns. Figure 5 shows the Apixaban Standard Solution overlay of 6 injections on the Luna C8(2) column and Kinetex C8 column. All the system suitability parameters for Organic Impurities were well within the acceptance criteria on both the columns in terms of %RSD, tailing factor, and resolution. The S/N ratio for Apixaban peak from the Sensitivity Solution was observed 5-8 times of the acceptance criteria in both columns. The reproducibility and % RSD for Apixaban is well within the acceptance

criteria in both the columns. The resolution between the Apixaban RC-D and Apixaban RC-E was observed more in Kinetex C8 column than Luna C8(2).

The USP method for Limit of Apixaban Related Compound-H, Related Compound-G, and Related Compound-F was verified on the Luna C18(2), Luna Omega C18, and Kinetex C18 columns. All the system suitability parameters for Limit of Apixaban Related Compound-H, Related Compound-G, and Related Compound-F were well within the acceptance criteria on all the columns in terms of %RSD and resolution. The reproducibility and %RSD for all the impurities was well within the acceptance criteria on all the columns. The retention times and the resolution observed between the Apixaban Related Compound-G and Apixaban Related Compound-F were relatively similar on the Luna C18(2) column and the Luna Omega C18 as shown in Figure 6. Additionally, the reproducibility of the method was confirmed through six replicate injections with the % RSD <0.5 % for the Luna C18(2) column and Luna Omega C18 column. The Kinetex C18 column showed lower resolution compared to the Luna C18(2) and Luna Omega C18 columns.


Figure 2. Standard Solution – Assay

Luna 5 μm C18(2) Column




Peak No.	Analyte	Retention Time	Area	USP Tailing	USP Plate Count	Area % RSD (N=6)
1	Apixaban	3.743	2125863	1.02	12712	0.05

Kinetex 5 μm Polar C18 Column



Peak No.	Analyte	Retention Time	Area	USP Tailing	USP Plate Count	Area % RSD (N=6)
1	Apixaban	2.853	2125052	1.12	9966	0.04


Luna Omega 5 μm C18 Column



Peak No.	Analyte	Retention Time	Area	USP Tailing	USP Plate Count	Area % RSD (N=6)
1	Apixaban	3.651	2127201	1.03	13884	0.04

Figure 3. System Suitability Solution – Organic Impurities


Luna 5 μm C8(2) Column



App ID 28749

Peak No.	Analyte	Retention Time	Area	% Area	RT Ratio	USP Resolution
1	Apixaban	14.150	15292548	99.71		
2	Apixaban Related Compound-D	20.545	23597	0.15	1.452	19.01
3	Apixaban Related Compound-E	21.537	20322	0.13	1.522	4.45

Kinetex 5 μm C8 Column




App ID 28752

Peak No.	Analyte	Retention Time (min.)	Area	% Area	RT Ratio	USP Resolution
1	Apixaban	8.834	15393949	99.72		
2	Apixaban Related Compound-D	16.409	23332	0.15	1.858	32.69
3	Apixaban Related Compound-E	17.668	20478	0.13	2.000	6.06

Figure 4. Sensitivity Solution – Organic Impurities


Luna Omega 5 μm C8(2) Column



App ID 28750

Peak No.	Analyte	Retention Time	Area	Height	USP S/N
1	Apixaban	14.147	7233	464	67.18

Kinetex 5 μm C8 Column




App ID 28753

Peak No.	Analyte	Retention Time	Area	Height	USP S/N
1	Apixaban	8.822	7951	874	88.28

Figure 5. Standard Solution – Organic Impurities


Luna Omega 5 μm C8(2) Column



App ID 28751

Peak No.	Analyte	Retention Time	Area	USP Tailing	USP Plate Count	Area % RSD (N=6)
1	Apixaban	14.119	150939	1.04	16430	0.15

Kinetex 5 μm C8 Column




App ID 28754

Peak No.	Analyte	Retention Time	Area	USP Tailing	USP Plate Count	Area % RSD (N=6)
1	Apixaban	8.831	153121	0.99	21685	0.07


Figure 6. Standard Solution – Limit of Related Compounds

Luna 5 μm C18(2) Column



App ID 28755


Luna 5 μm C18(2) Column – Overlay Chromatogram



App ID 28756


Peak No.	Analyte	Retention Time	Area	USP Resolution	Area % RSD (N=6)
1	Apixaban Related Compound-H	21.176	32948		0.25
2	Apixaban Related Compound-G	21.802	31286	3.41	0.33
3	Apixaban Related Compound-F	22.295	28612	2.68	0.22

Luna Omega 5 μm C18 Column



App ID 28757


Luna Omega 5 μm C18 Column – Overlay Chromatogram



App ID 28758


Peak No.	Analyte	Retention Time	Area	USP Resolution	Area % RSD (N=6)
1	Apixaban Related Compound-H	20.421	33452		0.53
2	Apixaban Related Compound-G	20.998	31837	3.45	0.37
3	Apixaban Related Compound-F	21.463	29210	2.77	0.30

Kinetex 5 μm C18 Column



App ID 28759

Kinetex 5 μm C18 Column – Overlay Chromatogram



App ID 28760

Peak No.	Analyte	Retention Time	Area	USP Resolution	Area % RSD (N=6)
1	Apixaban Related Compound-H	18.017	39277		0.37
2	Apixaban Related Compound-G	18.622	36870	3.94	0.48
3	Apixaban Related Compound-F	18.916	34425	1.90	0.48

Conclusion

The determination of Assay and Limit of Apixaban Related compound-H, Apixaban Related compound-G and Apixaban Related compound-F were verified on a fully porous Luna C18(2), a thermally modified Luna Omega C18, and a superficially porous Kinetex C18 column as per Apixaban USP monograph. All three columns showed excellent reproducibility for both methods. All the system suitability parameters were well within the acceptance criteria for the methods on all three columns. Based on the experimental findings, all three columns can serve as an alternative for the analysis of Assay and Limit of Apixaban Related compound-H, Apixaban Related compound-G and Apixaban Related compound-F in accordance with the USP monograph.

Additionally, the effective separation of organic impurities in Apixaban as per USP monograph verified using a fully porous Luna C8(2) and a superficially porous Kinetex C8 column. S/N for Apixaban peak from the Sensitivity Solution was observed 5-8 times of the acceptance criteria in both columns and observed a 2 times greater peak height in the Kinetex C8 column. Based on the experimental findings, both the columns can serve as an alternative for the analysis of Apixaban organic impurities in accordance with the USP monograph.



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions
Corporate Office USA**

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale. Phenomenex, the P icon, and the Phenomenex product and service marks mentioned herein are trademarks or registered trademarks of Phenomenex, Inc. in the United States and/or other countries. All other trademarks are the property of their respective owners. One or more authors are affiliated with Phenomenex, Inc. FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2026 Phenomenex, Inc. All rights reserved.

