

TN-1426

Robust and Reproducible LC-MS/MS Amino Acids Analysis with an Amide-Based HILIC Column

 Emanuele De Vita, PhD¹, Dr. Raffaella Salatiello¹, Luigi Margarucci, PhD²; Helen Harrington, PhD³

1. Istituto Diagnostico Varelli S.r.l., Via Cornelia dei Gracchi, 65 - 80126 Napoli, Italy

2. Phenomenex s.r.l. Via Marino Serenari 15D, 40013, Castel Maggiore Bologna, Italy

3. Phenomenex Ltd, Glasshouse, Alderley Park, Macclesfield, SK10 4AT, UK



Introduction

The accurate and reproducible analysis of amino acids is essential across a wide range of applications including biopharmaceutical development, food and nutritional sciences, and metabolomics. Due to their high polarity and zwitterionic nature, amino acids present well-known challenges for reversed phase liquid chromatography, often requiring derivatization or specialized detection strategies. As analytical demands continue to increase for sensitivity, robustness, and workflow efficiency, alternative chromatographic approaches have gained significant attention.

Amide-based HPLC columns operated under hydrophilic interaction liquid chromatography (HILIC) conditions provide an effective solution for the direct separation of underivatized amino acids. These stationary phases offer strong retention of polar analytes through a combination of hydrophilic partitioning and hydrogen-bonding interactions, enabling improved selectivity and resolution across structurally similar amino acids. When paired with high-organic mobile phases, amide columns also deliver enhanced compatibility with mass spectrometric detection, supporting sensitive and reproducible quantitative analysis. This application note describes the applicability of an amide-based stationary phase for the separation of 22 amino acids (Table 1) highlighting the key method considerations such as mobile phase composition, buffer selection and pH control.

UHPLC-MS Conditions

Column: Kinetex Amide

Dimensions: 100 x 2.1mm

Part No.: [00D-4812-AN](#)
Mobile Phase: A: H₂O +0,1 % FA
 B: ACN +0,1 % FA

Gradient:	Time (min)	% B
	0	95
	0,5	95
	6	30
	8	30
	8,1	95
	10	95

Flow Rate: 0.40mL/min

Injection Volume: 2µL

Temperature: 40°C


Detection: Shimadzu 8060 NX (MRM)

MS Conditions

	Ion Source:	ESI (+)
	Heating Gas Flow (L/min)	13
	Nebulizing Gas Flow (L/min)	2
	Drying Gas Flow (L/min)	10
	Capillary spray voltage (kV)	3,0
	DL temperature (°C)	250
	Heat Block temperature (°C)	400

Table 1. 22 Amino Acids dissolved in H₂O:ACN (20:80) +0.1% Formic Acid

Analyte	Q1 [M+H] ⁺	Q3	Dwell Time (ms)	CE
Alanine	90,1	44,0	21	16
		70,1		27
Arginine	175,1	60,1	5	19
		116,0		14
		74,0		19
Aspartic Acid	134,0	88,0	5	15
		116,0		11
		102,0		15
Glutamic Acid	148,0	84,0	9	25
		76,1		2
Glycine	76,1	30,0	9	19
		110,0		19
Histidine	156,1	83,0	9	34
		86,0		15
Isoleucine	132,1	56,0	5	51
		69,0		23
		86,0		15
Leucine	132,1	44,0	5	31
		43,0		35
		130,0		13
Lysine	147,1	84,0	5	30
		133,0		13
Methionine	150,0	104,0	3	15
		56,0		22
		120,0		19
Phenylalanine	166,0	103,0	5	35
		77,0		51
		70,0		23
Proline	116,1	43,0	5	41
		60,0		15
Serine	106,0	42,0	5	30
		103,2		25
Threonine	120,0	165,0	9	12
		136,0		17
		91,0		38
Tyrosine	182,1	72,0	3	16
		57,0		38
		55,0		27
Valine	118,1	152,1	5	19
		74,0		20
		188,0		13
Cystine	241,1	146,0	3	23
		118,0		36
		132,0		30
Tryptophan	205,1	115,0	3	40
		77,0		31
Serotonin	177,1	76,0	5	19
		130,0		21
Cysteine	122,0	74,0	9	24
		74,0		21
		87,0		15
Glutamine	147,1	84,0	5	15
		116,0		13
Asparagine	133,0	74,0	5	21
		87,0		15
		116,0		13

Figure 1. Elution profile of the selected Amino Acids under HILIC conditions.


The Kinetex Amide HPLC column shows effective retention and separation of amino acids under HILIC modality. All tested amino acids display clean and symmetrical peaks, indicating good efficiency and stable retention behavior (**Figure 1**). The amino acid mixture was prepared by diluting appropriate amounts of each stock solution in a water–acetonitrile (20:80, v/v) solution containing 0.1% formic acid. All solutions were freshly prepared prior to analysis and thoroughly mixed to ensure homogeneity.

As illustrated in **Figure 2**, depicting representative elution profiles for neutral, aromatic, acidic, and basic amino acid, the Kinetex Amide column delivers consistent selectivity and performance across diverse amino acid classes.

Despite their different chemical properties, all analytes follow the characteristic HILIC elution pattern and exhibit good peak shape, confirming the Kinetex Amide as a well-performing column for non-derivatized amino acid analysis.

Conclusion

The column's performance was evaluated in terms of retention, resolution, and reproducibility, confirming its suitability for both routine and advanced amino acid analyses. The results demonstrate that amide-based HILIC columns can streamline amino acid workflows while delivering robust, reliable, and high-quality chromatographic performance.

Figure 2. Chromatographic profile of Proline (a), Valine (b), Cysteine (c), Phenylalanine (d), Methionine (e), and Threonine (f)

Ordering Information

2.6 μ m Kinetex Amide Columns (mm)				Cartridges [‡]
Phases	50 x 2.1	100x 2.1	150x 2.1	3/pk
Amide Polyol	00B-4812-AN	00D-4812-AN	00F-4812-AN	AJ0-9800
for ID: 2.1mm				

[‡] SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 All other countries/regions Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

**BE-HAPPY™
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale. Phenomenex, the P icon, and the Phenomenex product and service marks mentioned herein are trademarks or registered trademarks of Phenomenex, Inc. in the United States and/or other countries. All other trademarks are the property of their respective owners. One or more authors are affiliated with Phenomenex, Inc. FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.
© 2026 Phenomenex, Inc. All rights reserved.

