

# Analysis of Underivatized Anionic and Cationic Pesticides in Reversed Phase and HILIC Modes Using a Single Mixed-mode HPLC Column


Zara Jalali, Ramkumar Dhandapani, PhD, Richard Jack, PhD, and Bryan Tackett, PhD Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA,



## **Overview**

Analysis of polar pesticides presents multiple challenges including adequate retention, separation of critical pairs, and reproducibility, to name a few. Traditionally, multiple columns are utilized for analyzing different class of polar pesticides, which results in dedicating two instruments, two different HPLC column selectivities, and reduced lab throughput. In this application note we show an efficient approach to a unique HPLC selectivity that provides optimal separation of various anionic and cationic polar pesticide classes using a single HPLC column. The Luna Polar Pesticides HPLC column was effective in retaining both anionic and cationic pesticides, showing that the same column can be used for both positive and negative polarity modes for analysis. The versatile selectivity of this column enables Reversed Phase and HILIC retention of anionic and cationic pesticides as evident from Figure 1 and Figure 2, respectively.

Figure 1. Underivatized Anionic Pesticide Analysis on a Luna Polar Pesticides HPLC Column.



Peak No.	Analyte	Retention Time (min)	Q1 (m/z)	Q3 (m/z)	Peak No.	Analyte	Retention Time (min)	Q1 (m/z)	Q3 (m/z)
1	АМРА	0.97	110 110	79 63	9	N-Acetyl Glufosinate-D <sub>3</sub>	1.9	225	62
2	AMPA- <sup>13</sup> C, <sup>15</sup> N	0.97	112	79	10	Phosphonic Acid	2.44	81 81	63 79
3	Maleic Hydrazide-D <sub>2</sub>	1.16	113	42	11	Ethephon-D <sub>4</sub>	2.55	147	79
4	Maleic Hydrazide	1.17	111 111	82 55	12	Ethephon	2.56	143 145	107 107
5	Glufosinate	1.3	180 180	63 85	13	Chlorate	3.07	83 85	67 69
6	МРРА	1.7	151 151	63 107	14	Fosetyl Al-D <sub>15</sub>	3.41	114	82
7	Glyphosate	1.71	168 168	63 79	15	Fosetyl Al	3.43	109 109	81 63
8	N-Acetyl-Glufosinate	1.9	222 222	136 134	16	Perchlorate	3.78	99 101	83 85

## **LC Conditions**

Column: Luna 3 μm Polar Pesticides

**Dimensions:** 100 x 2.1 mm **Part No.:** 00D-4798-AN

Guard: SecurityGuard™ ULTRA for Polar

Pesticides: AJO-8789

Mobile Phase: A: 0.3 % Formic Acid in Water

B: 0.3 % Formic Acid in Acetonitrile

Gradient: Time (min) %B
0 2
0.5 2
6 20
7 90
9 90
9.1 2


Flow Rate: 0.3 mL/minInjection Volume:  $1 \mu L$ Temperature:  $40 \, ^{\circ} C$ Detection: MRM
Detector:  $SCIEX^{\circ} 7500$ 

Analyte Concentrations: 50 ng/mL

# **MRM Conditions**

Polarity: Negative
Gas Temperature: 450 °C
GS1: 40 psi
GS2: 30 psi
CUR: 40 psi
IS: -4500 V

Figure 2. Underivatized Cationic Pesticide Analysis on a Luna™ Polar Pesticides HPLC Column.



Peak No.	Analyte	Retention Time (min)	Q1 (m/z)	Q3 (m/z)
1	Ethyl Thiourea	0.77	103	44
1	Ethyl Milourea	0.77	103	86
2	Mepiquat	3.2	189	144
2	iviepiquat		189	102
3	Dranamaaarh	3.28	114	98
3	Propamocarb		114	70
4	Overnatrine	4.09	265	136
4	Oxamatrine		265	205
5	Davaguat	5.64	185	170
5	Paraquat		185	169

# **LC Conditions**

Column: Luna 3 µm Polar Pesticides

**Dimensions:** 100 x 2.1 mm **Part No.:** 00D-4798-AN

 $\textbf{Guard:} \ \ \mathsf{SecurityGuard^{\mathsf{TM}}} \ \ \mathsf{ULTRA} \ \mathsf{for} \ \mathsf{Polar}$ 

Pesticides: AJO-8789

Mobile Phase: A: 50 mM Ammonium Formate

in Water, adjust pH to 3 with

Formic Acid B: Acetonitrile

Gradient:	Time (min)	%В
	0	97
	0.5	97
	4	70
	5	40
	6	40
	6.1	97

Flow Rate: 0.5 mL/min

Detector: SCIEX® 7500

Analyte Concentrations: 40 ng/mL

# **MRM Conditions**

Polarity: Positive
Gas Temperature: 450 °C
GS1: 40 psi
GS2: 40 psi
CUR: 50 psi
IS: 4500 V

# Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

t: +61 (0)2-9428-6444 auinfo@phenomenex.com

# Austria

t: +43 (0)1-319-1301 anfrage@phenomenex.com

## Belaium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

**Canada** t: +1 (800) 543-3681 info@phenomenex.com

t: +86 400-606-8099 cninfo@phenomenex.com

## Czech Republic

t: +420 272 017 077 cz-info@phenomenex.com

## Denmark

t: +45 4824 8048 nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

t: +49 (0)6021-58830-0 anfrage@phenomenex.com

## **Hong Kong**

t: +852 6012 8162 hkinfo@phenomenex.com

## India

t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

## Indonesia

t: +62 21 5019 9707 indoinfo@phenomenex.com

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

## Japan

t: +81 (0) 120-149-262 jpinfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226 tecnicomx@phenomenex.com

# The Netherlands

t: +31 (0)30-2418700 nlinfo@phenomenex.com

# **New Zealand**

t: +64 (0)9-4780951 nzinfo@phenomenex.com

**Norway** t: +47 810 02 005 nordicinfo@phenomenex.com

## Poland

t: +48 22 104 21 72 pl-info@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

Singapore t: +65 6559 4364 sginfo@phenomenex.com

**Slovakia** t: +420 272 017 077 sk-info@phenomenex.com

## Spain

t: +34 91-413-8613 espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

## Switzerland

t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

## Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

## Thailand

t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

# **United Kingdom**

t: +44 (0)1625-501367 ukinfo@phenomenex.com

t: +1 (310) 212-0555 www.phenomenex.com/chat

## All other countries/regions Corporate Office USA

t: +1 (310) 212-0555 www.phenomenex.com/chat

# www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

# Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <a href="www.phenomenex.com/phx-terms-and-conditions-of-sale">www.phenomenex.com/phx-terms-and-conditions-of-sale</a>. Trademarks

Luna, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. SCIEX is a registered trademark of AB SCIEX Pte. Ltd. Disclaimer

Comparative separations may not be representative of all applications SecurityGuard is patented by Phenomenex, U.S. Patent No. 6.162,362

CAUTION: this patent only applies to the analytical-sized award cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges,

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2023 Phenomenex, Inc. All rights reserved.