

Ph. Eur. Monograph 1425: Fosfomycin Trometamol Assay and Related Substances on Luna™ Omega 3 µm SUGAR Column

Zachary Woodward¹, Zeshan Aqeel¹, Bryan Tackett, PhD¹, and Heiko Behr, PhD²

¹Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

²Phenomenex Ltd. Deutschland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany

Overview

Fosfomycin Trometamol is an antibiotic used for bladder infections. In this application note we show the separation of Fosfomycin Trometamol from its related substances following Ph. Eur. Monograph 1425. We used a Luna Omega 3 µm SUGAR column in two different dimensions (250 x 4.6 mm and 150 x 4.6 mm) and compared them to the ZORBAX 5 µm NH₂, 250 x 4.6 mm column originally used in the monograph.


System suitability per Ph. Eur. Monograph for Fosfomycin Trometamol is the same for both Assay and Related Substances: Minimum resolution of 1.5 between the peaks due to Impurity A and Fosfomycin, and minimum peak-to-valley ratio of 1.5 between the peak due to Impurity C and the peak due to Impurity B for reference solution (c). The peak-to-valley ratio is defined as H_p/H_v , where H_p = height above the baseline of the peak due to Impurity C and H_v = height above the baseline of the lowest point of the curve separating this peak from the peak due to Impurity B.

The ZORBAX 5 µm NH₂ column described in the Ph. Eur. monograph for Fosfomycin Trometamol presents challenges when evaluating the peak-to-valley criteria between Impurity B and Impurity C. Specifically, the noisy baseline from unresolved sample compromises the consistent evaluation of the peak-to-valley parameter. The resulting measurement is then subjective from one System Suitability Standard to the next.

The Luna Omega 3 µm SUGAR columns underwent only about 20 minutes of mobile phase equilibration prior to triplicate injections of the blank and displayed stable baselines across these blank injections (data not shown). Complete baseline resolution between Impurity B and Impurity C was achieved using both the prescribed dimensions (250 x 4.6 mm) and the shorter dimensions (150 x 4.6 mm) of the Luna Omega 3 µm SUGAR columns. The longer 250 mm column also afforded improved resolution between Fosfomycin and Impurity A, while preserving the expected run time. While using a smaller 3 µm particle size would be expected to generate higher pressures, the use of a shorter 150 mm column, which maintains the L/dp ratio, alleviates these concerns. The resolution between Fosfomycin and Impurity A meets suitability requirements on the shorter Luna Omega 3 µm SUGAR column but may not afford as much flexibility when accounting for prolonged use of the column during routine sample analysis.

All reference solutions were prepared as indicated in Ph. Eur. monograph 1425 for Fosfomycin Trometamol. The following certified reference standards (CRS) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: Allee Kastner CS 30026 F - 67081 Strasbourg (France):

- F0400000, Fosfomycin Trometamol CRS



LC-UV Conditions

Columns: Luna Omega 3 µm SUGAR, 250 x 4.6 mm ([00G-4775-E0](#))

Luna Omega 3 µm SUGAR, 150 x 4.6 mm ([00F-4775-E0](#))

ZORBAX® 5 µm NH₂, 250 x 4.6 mm

Mobile Phase: [Mobile Phase \(Table 1\)](#)

Flow Rate: 1.0 mL/min (Isocratic)

Injection: 5 µL

Temperature: Ambient with column in column heater, off

Detector: Refractive Index Detector at 35 °C

System: Agilent® 1260

Table 1. Preparation of Test and Reference Solutions

Solution	Composition
Mobile Phase	10.89 g/L of Potassium Dihydrogen Phosphate in HPLC Water.
Test Solution	Same as Reference Solution (a)
Reference Solution (a)	Dissolve 0.600 g of Fosfomycin Trometamol CRS in the Mobile Phase and dilute to 5.0 mL with the Mobile Phase .
Reference Solution (b)	Dilute 1.0 mL of Reference Solution (a) to 100.0 mL with the Mobile Phase . Dilute 3.0 mL of this solution to 10.0 mL with the Mobile Phase .
Reference Solution (c)	Wet 0.3 g of Fosfomycin Trometamol CRS with 60 µL of Water and heat in an oven (Shimadzu GC-MS) at 60 °C for 24 hours. Dissolve the residue in the Mobile Phase and dilute to 20.0 mL with the Mobile Phase (Solution A) . Dissolve 0.6 g of Fosfomycin Trometamol CRS in Solution A and dilute to 5.0 mL with the same solution (<i>in situ</i> degradation to obtain impurities A, B, C, and D).

Figure 1. Fosfomycin Trometamol





Figure 2. System Suitability Test for Assay and Related Substances Using Reference Solution (c)

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sales.

Trademarks

Luna and BE-HAPPY are trademarks of Phenomenex. ZORBAX and Agilent are registered trademarks of Agilent Technologies Inc.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Agilent Technologies Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

