

Separation of Tacrolimus and its Organic Impurities per USP Monograph

Rajesh Jangam¹, Rajesh Babu Dandamudi, PhD¹, and Bryan Tackett, PhD²

¹India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India

²Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview


Tacrolimus is an inhibitor of calcineurin, a phosphatase that regulates production of interleukin-2. In this application note, we report the separation of Tacrolimus and its related organic impurities per the USP monograph using a Luna 3 µm C18 column, a Luna 3 µm C18(2) column, and a Luna Omega 3 µm C18 column.

System suitability per USP Monograph for Tacrolimus Assay is a resolution no less than (NLT) 3.0 between Ascomycin and Tacrolimus and a percent relative standard deviation (%RSD) no more than (NMT) 1.0 % for 5 replicates for the sum of the responses of Tacrolimus, Tacrolimus open ring, and Tacrolimus 19-epimer. System suitability per USP Monograph for Tacrolimus Organic Impurities, Procedure 2, is resolution NLT 3.0 between Ascomycin and Tacrolimus and a percent relative standard deviation (%RSD) no more than (NMT) 10.0 % for 6 replicates for the sum of the responses of Tacrolimus and Tacrolimus 19-epimer.

All system suitability requirements for Tacrolimus Assay and Organic Impurities were met by all columns. However, the newer generation Luna Omega C18 column, which delivers an L1 column with increased inertness and reproducibility, provided the greatest resolution and lowest % RSD making it the ideal column selection for this monograph.

All solutions were prepared as indicated in the USP Monograph for Tacrolimus. USP Tacrolimus RS (Catalog No. 1642802), USP Tacrolimus Related Compound A RS (Catalog No. 1642813), USP Tacrolimus System Suitability Test Mix RS (Catalog No. 1642857), USP Tacrolimus 8-epimer RS (Catalog No. 1642824), and USP Tacrolimus 8-propyl analog (Catalog No. 1642835) were purchased from USP.

Figure 1. Tacrolimus Structure.



LC-UV Conditions

Columns: Luna™ 3 µm C18 ([00F-4114-E0](#))
Luna 3 µm C18(2) ([00F-4251-E0](#))
Luna Omega 3 µm C18 ([00F-4784-E0](#))

Dimensions: 150 x 4.6 mm

Mobile Phase: A: **Solution A** / **Solution B** (4:1, v/v)
B: **Solution A** / **Solution B** (1:4, v/v)

Gradient:	Time (min)	%B
	0	28
	30	28
	53	85
	54	28
	60	28

Flow Rate: 1.5 mL/min

Injection Volume: 20 µL

Temperature: 60 °C

Detector: UV @ 220 nm

System: Waters® ACQUITY Arc® HPLC


Table 1. Preparation of Test and Reference Solutions

Solution	Composition
Solution A	6 mM Phosphoric Acid
Solution B	Acetonitrile / <i>tert</i> -Butyl Methyl Ether (81:19, v/v)
Diluent	Acetonitrile / Water (7:3, v/v).
System Suitability Solution (Assay and Organic Impurities)	3 mg/mL of USP Tacrolimus System Suitability Mixture RS in Diluent . Allow the solution to stand for 3 h at ambient temperature before use. Protect from light by using low-actinic glassware.
Standard Solution (Assay)	3 mg/mL of USP Tacrolimus RS in Diluent . Allow the solution to stand for 3 h at ambient temperature before use. Protect from light by using low-actinic glassware.
Standard Solution (Organic Impurities, Procedure 2)	30 µg/mL of USP Tacrolimus RS in Diluent . Allow the solution to stand for 3 h at ambient temperature before use. Protect from light by using low-actinic glassware.



Figure 2. System Suitability Solution – Assay and Organic Impurities


Luna™ 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	RT Ratio	Area	Resolution
1	Tacrolimus 19-epimer	16.82	0.63	1573240	-
2	Ascomycin	23.26	0.87	846856	3.01
3	Tacrolimus	26.59	1.00	17646968	
4	Tacrolimus 8-epimer	34.27	1.29	382551	-
5	Tacrolimus 8-propyl analog	35.78	1.35	45686	-

N=5 Injections


Luna 3 µm C18(2) Column



Peak No.	Analyte	Retention Time (min)	RT Ratio	Area	Resolution
1	Tacrolimus 19-epimer	16.61	0.63	1560723	-
2	Ascomycin	22.94	0.88	838941	3.23
3	Tacrolimus	26.20	1.00	17172567	
4	Tacrolimus 8-epimer	33.93	1.30	376096	-
5	Tacrolimus 8-propyl analog	35.43	1.35	48686	-

N=5 Injections

Luna Omega 3 µm C18 Column




Peak No.	Analyte	Retention Time (min)	RT Ratio	Area	Resolution
1	Tacrolimus 19-epimer	20.503	0.63	428618	-
2	Ascomycin	28.327	0.88	760690	3.89
3	Tacrolimus	32.370	1.00	16292286	
4	Tacrolimus 8-epimer	37.382	1.15	349065	-
5	Tacrolimus 8-propyl analog	38.406	1.19	41012	-

N=5 Injections



Figure 3. Standard Solution – Assay


Luna™ 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Total Area	Total Area %RSD
1	Tacrolimus open ring	13.93	467785	16347732	0.30
2	Tacrolimus 19-epimer	16.98	1308637		
3	Tacrolimus	26.79	14571310		

N=5 Injections


Luna 3 µm C18(2) Column



Peak No.	Analyte	Retention Time (min)	Area	Total Area	Total Area %RSD
1	Tacrolimus open ring	13.50	497129	16615330	0.78
2	Tacrolimus 19-epimer	16.43	1328070		
3	Tacrolimus	25.74	14790131		

N=5 Injections

Luna Omega 3 µm C18 Column




Peak No.	Analyte	Retention Time (min)	Area	Total Area	Total Area %RSD
1	Tacrolimus open ring	15.62	521842	17348029	0.13
2	Tacrolimus 19-epimer	18.93	1408543		
3	Tacrolimus	29.98	15417644		

N=5 Injections




Figure 4. Standard Solution – Assay

Luna™ 3 µm C18 Column




Peak No.	Analyte	Retention Time (min)	Area	Total Area	Total Area %RSD
1	Tacrolimus 19-epimer	17.80	14774	170120	0.91
2	Tacrolimus	28.51	155346		
N=6 Injections					

Luna 3 µm C18(2) Column



Peak No.	Analyte	Retention Time (min)	Area	Total Area	Total Area %RSD
1	Tacrolimus 19-epimer	17.00	14145	166737	1.27
2	Tacrolimus	27.09	152592		
N=6 Injections					

Luna Omega 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Total Area	Total Area %RSD
1	Tacrolimus 19-epimer	19.35	14204	161680	1.60
2	Tacrolimus	31.02	147476		
N=6 Injections					



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY Arc are registered trademarks of Waters Technology Corporation.

Disclaimer

Comparative separation may not be representative of all applications. Phenomenex is in no way affiliated with Waters Technology Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists