

# Limit of Valsartan Related Compound A per USP Monograph for Amlodipine, Valsartan, and Hydrochlorothiazide Tablets

Alyssa Nodland and Bryan Tackett, PhD  
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

## Overview

Amlodipine, Valsartan, and Hydrochlorothiazide tablets contain a combination of drugs used to treat high blood pressure. Amlodipine is a calcium channel blocker (CCB) that affects calcium movement into the cells of the heart and blood vessels, Valsartan is an angiotensin II receptor blocker (ARB) that relaxes blood vessels, and Hydrochlorothiazide is a thiazide diuretic that reduces the amount of water in the body by increasing the flow of urine and lowering blood pressure. In this application note, we report the limit of Valsartan Related Compound A using a Lux 5 µm Cellulose-1 column compared with a CHIRALCEL 5 µm OD-H column according to the USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets, which specifies a column containing L40 (Cellulose tris-(3,5-dimethylphenylcarbamate)) packing with dimensions of 250 x 4.6 mm.

System suitability per USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets is resolution no less than (NLT) 2.0 between Valsartan and Valsartan related compound A and a percent relative standard deviation (%RSD) of no more than (NMT) 5 % for Valsartan related compound A. All system suitability requirements for the Limit of Valsartan Related Compound A were met by both the Lux Cellulose-1 and CHIRALCEL OD-H.

The Lux Cellulose-1 column is a suitable substitute for the CHIRALCEL OD-H column for the USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets Limit of Valsartan Related Compound A method.

All solutions were prepared as indicated in the USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets. USP Valsartan RS (Catalog No. 1708762) and USP Valsartan Related Compound A RS (Catalog No. 1708773) were purchased from USP.

## LC-UV Conditions

**Columns:** Lux™ 5 µm Cellulose-1 ([00G-4459-E0](#))  
CHIRALCEL® 5 µm OD-H

**Dimensions:** 250 x 4.6 mm

**Mobile Phase:** N-Hexane / 2-Propanol / Trifluoroacetic Acid  
(850:150:1, v/v/v)

**Flow Rate:** 0.8 mL/min (Isocratic)

**Injection Volume:** 20 µL

**Temperature:** 30 °C


**LC System:** Agilent® 1260

**Detection:** UV @ 230 nm

**Table 1.** Preparation of Test Solutions


Solution	Composition
Standard Solution	0.001 mg/mL of USP Valsartan Related Compound A RS in Mobile Phase.
System Suitability Solution	0.04 mg/mL each of USP Valsartan Related Compound A and USP Valsartan RS in Mobile Phase.

**Figure 1.** Drug Structures.



**Figure 2.** System Suitability Solution.


CHIRALCEL® 5 µm OD-H® Column



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Valsartan Related Compound A	8.40	3225.7	3.51
2	Valsartan	13.11	3570.6	

N=5 Injections

Lux™ 5 µm Cellulose-1 Column




Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Valsartan Related Compound A	9.81	5349.74	3.60
2	Valsartan	13.94	5290.46	

N=5 Injections

**Figure 3.** Standard Solution.


CHIRALCEL 5 µm OD-H Column



Peak No.	Analyte	Retention Time (min)	Area	%RSD
1	Valsartan	13.13	101.50	3.55

N=6 Injections

Lux 5 µm Cellulose-1 Column



Peak No.	Analyte	Retention Time (min)	Area	%RSD
1	Valsartan	15.26	153.15	1.19

N=6 Injections



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 3952 5747  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 51 02 180  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
**GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Lux and BE-HAPPY are trademarks of Phenomenex. CHIRALCEL and OD-H are registered trademarks of DAICEL Corporation. Agilent is a registered trademark of Agilent Technologies, Inc.

**Disclaimer**

Comparative separation may not be representative of all applications.

Phenomenex is in no way affiliated with DAICEL Corporation or Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat) to get in touch with one of our Technical Specialists