

Meeting System Suitability for USP Amlodipine, Valsartan, and Hydrochlorothiazide Tablets Assay and Organic Impurities

Rajesh Jangam¹, Rajesh Babu Dandamudi, PhD¹, and Bryan Tackett, PhD²

¹India Phenologix Lab, Phenomenex India, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Holbi, Jala Taluka, Bengaluru 562149, India

²Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview


Amlodipine, Valsartan, and Hydrochlorothiazide tablets contain a combination of drugs used to treat high blood pressure. Amlodipine is a calcium channel blocker (CCB) that affects calcium movement into the cells of the heart and blood vessels, Valsartan is an angiotensin II receptor blocker (ARB) that relaxes blood vessels, and Hydrochlorothiazide is a thiazide diuretic that reduces the amount of water in the body by increasing the flow of urine and lowering blood pressure. In this application note, we report the separation of Amlodipine, Valsartan, and Hydrochlorothiazide and their related organic impurities per the USP monograph using a Luna Omega 3 μ m C18 column as an alternative L1 column to the ACE 3 μ m C18 column originally used to elucidate the monograph.

System suitability per USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets Assay is a symmetry factor no more than (NMT) 2.0 for Amlodipine, Valsartan, and Hydrochlorothiazide and a percent relative standard deviation (%RSD) of NMT 2 % for Amlodipine, Valsartan, and Hydrochlorothiazide. System suitability per USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets organic impurities is resolution no less than (NLT) 2.0 between any adjacent peaks of Benzothiadiazine related compound A, Hydrochlorothiazide, Amlodipine related compound A, Amlodipine, Valsartan related compound B, and Valsartan, a %RSD NMT 5.0 % for Amlodipine related compound A, Benzothiadiazine related compound A, Amlodipine, Valsartan, and Hydrochlorothiazide, and a signal-to-noise (S/N) ratio NLT 10 for Amlodipine, Valsartan, and Hydrochlorothiazide. All system suitability requirements for both Assay and Organic Impurities were met by both the Luna Omega C18 column and the ACE C18 column.

The Luna Omega C18 column gave comparable results and is therefore a suitable substitute for the ACE C18 column for the USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide Tablets Assay and Organic Impurities methods.

All solutions were prepared as indicated in the USP monograph for Amlodipine, Valsartan, and Hydrochlorothiazide tablets. USP Amlodipine Besylate RS (Catalog No. 1029501), UPS Amlodipine Related Compound A RS (Catalog No. 1029512), USP Benzothiadiazine Related Compound A RS (Catalog No. 1057507), USP Hydrochlorothiazide RS (Catalog No. 1314009), USP Valsartan RS (Catalog No. 1708762) and USP Valsartan Related Compound B RS (Catalog No. 1708784) were purchased from USP.

Figure 1. Drug Structures.



LC-UV Conditions

Columns: Luna™ Omega 3 μ m C18 ([00F-4784-E0](#))
ACE® 3 μ m C18

Dimensions: 150 x 4.6 mm

Mobile Phase: A: Acetonitrile / Water / Phosphoric Acid (50:950:1, v/v/v)
B: Acetonitrile / Water / Phosphoric Acid (950:50:1, v/v/v)

Gradient	Time (min)	%B
	0.00	5
	3.00	50
	6.00	60
	10.0	95
	10.1	5
	15.0	5

Flow Rate: 1.5 mL/min

Injection Volume: 10 μ L

Temperature: 40 °C

LC System: Waters® ACQUITY® H-Class

Detection: UV @ 225 nm


Table 1. Preparation of Test Solutions. Use Amber Glassware for All Solutions Containing Drug Substances.

Solution	Composition
Diluent	Acetonitrile / Water (1000:1, v/v)
Standard Solution (Assay)	0.14 mg/mL of USP Amlodipine Besylate RS, 0.064 mg/mL of USP Valsartan RS, and 0.025 mg/mL of USP Hydrochlorothiazide RS in Diluent.
Standard Solution (Organic Impurities)	0.0005 mg/mL of USP Amlodipine Related Compound A RS, 0.0001 mg/mL of USP Benzothiadiazine Related Compound A RS, 0.0003 mg/mL of USP Amlodipine Besylate RS, 0.00015 mg/mL of USP Valsartan RS, and 0.00005 mg/mL of USP Hydrochlorothiazide RS in Diluent.
System Suitability Solution (Organic Impurities)	0.02 mg/mL each of USP Benzothiadiazine Related Compound A RS and USP Valsartan Related Compound B RS, 0.005 mg/mL of USP Amlodipine Related Compound A RS, 0.14 mg/mL of USP Amlodipine Besylate RS, 0.064 mg/mL of USP Valsartan RS, and 0.025 mg/mL of USP Hydrochlorothiazide RS in Diluent.
Sensitivity Solution (Organic Impurities)	0.14 μ g/mL of USP Amlodipine Besylate RS, 0.064 μ g/mL of USP Valsartan RS, and 0.025 μ g/mL of USP Hydrochlorothiazide RS in Diluent.



Figure 2. Standard Solution – Assay.


ACE® 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	%RSD	Symmetry Factor
1	Hydrochlorothiazide	2.49	1198889	1.65	0.9
2	Amlodipine	4.14	1378275	1.63	2.0
3	Valsartan	5.73	1430802	1.72	1.0

N=5 Injections

Luna™ Omega 3 µm C18 Column




Peak No.	Analyte	Retention Time (min)	Area	%RSD	Symmetry Factor
1	Hydrochlorothiazide	2.68	1220034	0.05	0.9
2	Amlodipine	3.97	1404294	0.02	1.8
3	Valsartan	5.95	1455076	0.05	1.0

N=5 Injections


Figure 3. System Suitability Solution – Organic Impurities.

ACE 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Benzothiadiazine related compound A	2.23	1003193	4.2
2	Hydrochlorothiazide	2.47	1148293	
3	Amlodipine related compound A	3.81	50605	6.8
4	Amlodipine	4.11	1317665	
5	Valsartan related compound B	5.05	382561	10.7
6	Valsartan	5.68	1372433	


Luna Omega 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Benzothiadiazine related compound A	2.43	1014133	4.4
2	Hydrochlorothiazide	2.67	1155816	
3	Amlodipine related compound A	3.70	49951	5.4
4	Amlodipine	3.94	1325691	
5	Valsartan related compound B	5.24	385825	10.8
6	Valsartan	5.90	1384201	

Figure 4. Standard Solution – Organic Impurities.


ACE 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	%RSD
1	Benzothiadiazine related compound A	2.23	6857	1.36
2	Hydrochlorothiazide	2.47	4096	2.49
3	Amlodipine related compound A	3.81	6031	1.35
4	Amlodipine	4.15	3208	0.92
5	Valsartan	5.68	3832	2.99

N=6 Injections


Luna Omega 3 µm C18 Column




Peak No.	Analyte	Retention Time (min)	Area	%RSD
1	Benzothiadiazine related compound A	2.43	6695	1.83
2	Hydrochlorothiazide	2.67	4212	1.94
3	Amlodipine related compound A	3.70	6107	0.96
4	Amlodipine	3.98	2937	3.32
5	Valsartan	5.90	3639	3.65

N=6 Injections



Figure 5. Sensitivity Solution – Organic Impurities.

Peak No.	Analyte	Retention Time (min)	Area	S/N
1	Hydrochlorothiazide	2.47	1381	10
2	Amlodipine	4.15	1447	19
3	Valsartan	5.68	1813	12



Peak No.	Analyte	Retention Time (min)	Area	S/N
1	Hydrochlorothiazide	2.67	1648	13
2	Amlodipine	3.98	1469	18
3	Valsartan	5.90	1580	11



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and BE-HAPPY are trademarks of Phenomenex. ACE is a registered trademark of VWR International. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation.

Disclaimer

Comparative separation may not be representative of all applications.

Phenomenex is in no way affiliated with VWR International or Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists