

Application Handbook



**Chromatography Solutions
for API Manufacturers**


Quality filtration products from Phenomenex offer a convenient way to clean up your samples - for improved analytical results with a simple method. Phenex syringe filters are designed for efficient and rapid filtration of almost any solution prior to analysis, and are optimized for superior flow rates and high throughput. Phenex offers a wide variety of membranes ideal for any application.

- High quality filtration product
- Increased column lifetime
- Less system downtime

The housing attaches to any standard Luer lock syringe, so the sample can easily be pushed through the membrane with minimal pressure. The result is a particulate-free eluent that is ready for use with HPLC, GC, or other analytical techniques.

A Sample Preparation Solution


Filtration is one of the most popular sample preparation methods due to the ability to remove particulates before injection onto the column, the ability to use it in parallel with another sample preparation technique, and requires no method development.



Avoid Unnecessary Problems

Using low-quality, low-cost filters can lead to sample contamination by introducing filter membrane and housing extractables directly into your sample. This will negatively affect chromatography and can cause ghost peaks for future runs.

All Phenex filters are made with the highest grade materials available to remove any unwanted interferents.



Phenex Helps Reduce

- | | |
|--------------------|-------------------|
| Extraneous peaks | Sample co-elution |
| False quantitation | Instrument damage |

Ensure Quality

- A Identification**
Membrane type and pore size are clearly marked on individual syringe filters
- B Luer Lock Inlet Tip**
Secures connections to prevent "blow off"
- C Sample Distribution Rings**
Creates even sample distribution for high sample flow rates
- D Medical Grade Polymer Housing**
Offers the most inert syringe filter and helps eliminate unwanted secondary interactions with the filter housing
- E Ultrasonically Welded**
Ensures robust housing - filter integrity

To learn more about sample preparation options, go to

www.phenomenex.com/samplepreparation

Which Filter Membrane Is Right for Me?



Phenex syringe filters are offered in a variety of chemically compatible membranes that are ideal for any application. Proper membrane and size selection are the keys to choosing the best product to maintain the integrity of your sample components as well as to protect your system from particulate contamination.



Select your filter in three EASY steps:

1. What is your sample volume?

≤ 2 mL Sample Volume	2 to 10 mL Sample Volume	10 to 100 mL Sample Volume
4mm Diameter	15mm Diameter	25 - 28mm Diameter

2. What is your LC column particle size?

> 3 μm	< 3 μm	OR
0.45 μm	0.20 μm	Viscous samples such as serum, plasma, or other biological matrices. Solutions with high particulate load (e.g., some environmental or food and beverage applications).
		Glass Fiber Filter with 0.45 μm membrane

3. What type of sample are you working with?

	Aqueous			Solvents	
	Solvents Mixtures	Tissue Culture Media, Buffers	Protein Analysis/ Biological Samples	Non-Aqueous	Aqueous Mixtures
	↓	↓	↓	↓	↓
	Hydrophobic/ Strong Acids			Hydrophilic	
	↓			↓	
	RC (Regenerated Cellulose)	CA (Cellulose Acetate)	PES (Polyethersulfone)	PES (Polyethersulfone)	RC (Regenerated Cellulose)

RC (Regenerated Cellulose)

For Aqueous and Mixed Organic Solutions

A broad range of aqueous and mixed-organic solutions

Fast-flow and ultra-low protein and non-specific binding characteristics

Broadly recommended as an excellent general purpose/high-performance sample filter for most applications

PTFE, Teflon® (Polytetrafluoroethylene)

For 100% Organic Solutions

Well-suited for the clarification of non-aqueous samples

Hydrophobic membrane, excellent for filtration of organic-based, highly acidic or basic samples and solvents

A hydrophobic membrane, that can be made hydrophilic by wetting with alcohol and then flushing with deionized water



Additional Syringe Filter Membranes

Membrane Types	Recommended Uses
PES (Polyethersulfone)	Polyethersulfone membranes exhibit very fast-flow and ultra-low protein binding characteristics. Phenex-PES membranes are typically broadly recommended for filtering critical biological samples, tissue culture media, additives and buffers.
NY (Nylon)	Nylon has inherent hydrophilic characteristics and works well for filtration of many aqueous and mixed-organic samples. In combination with a glass pre-filter (Phenex-GF/NY), this membrane is excellent for the filtration of particle-laden samples, such as foods and beverages, environmental, biofuels, and dissolution samples.
CA (Cellulose Acetate)	Cellulose Acetate (CA) membranes exhibit ultra-low protein binding and are broadly used in the filtration of biological samples. In combination with a glass pre-filter (Phenex-GF/CA), this membrane is excellent for filtration of tissue culture media, general biological sample filtration and clarification.
GF (Glass Fiber)	Glass Fiber (GF) filters are made of inert borosilicate glass and have a nominal 1.2 µm pore size. They are commonly used with highly viscous samples or samples containing high concentrations of particulate matter (e.g., food analysis, biological samples, soil samples, fermentation broth samples, removal of yeasts, molds, etc.).
PVDF (Polyvinylidene Fluoride)	Hydrophilic PVDF membrane provides high flow rates and throughput, low extractables, and broad chemical compatibility. This membrane binds less protein than nylon or PTFE membranes.

Luna Omega is a fully porous silica re-engineered for the 21st century. The thermally modified particles exhibit higher efficiency and greater inertness, providing highly efficient, high resolution separations with excellent peak shape.

The manufacturing process was specifically designed to control batch to batch reproducibility, providing a stable platform for method development and implementation of methods in QC environments across the lifespan of pharmaceutical products.


Thermal Modified Pore Structure



Luna Omega C18

Excellent first choice C18 column demonstrating high levels of inertness

- [Gemifloxacin Mesylate and organic impurities](#)
- [Levothyroxine sodium and related substances](#)




Luna Omega Polar C18

Novel C18 phase providing increased retention for both polar and non-polar analytes

- [Paliperidone and organic impurities](#)


Phenomenex offers a wide range of columns, both high efficiency core-shell and fully porous materials. Our Kinetex range of core-shell selectivities are available in a range of particle sizes to match with the instrumentation you have, and have been proven to separate drug substances, related compounds and impurities as the resources below demonstrate.



The high efficiency delivered by Kinetex core-shell particles allow for high resolution separations, enabling the user to separate complex mixtures of API and impurities.


Kinetex F5 – a pentafluorophenyl stationary phase bonded to our high efficiency core-shell material. The phase provides excellent selectivity for halogenated, conjugated, isomeric or highly polar compounds

- [Ezetimibe and organic impurities](#)




Kinetex Biphenyl – a 100% aqueous stable biphenyl phase bonded to core-shell particles. This stationary phase provides reversed phase selectivity for hydrophobic, aromatic and polar compounds

- [Diazoxide and organic impurities](#)



Kinetex Phenyl-Hexyl – this is a more hydrophobic phenyl phase due to the presence of the C6 linker. It provides aromatic interactions and hydrophobic interactions, providing an alternative selectivity to C8 or C18 phases.


- [Miconazole Nitrate and organic impurities](#)



Core-shell vs. Fully Porous


High Efficiency, High Density Particle

This van Deemter plot indicates that Kinetex core-shell particles exhibit higher plate count (lower HETP) than fully porous materials, allowing for the separation of complex mixtures.




Kinetex C8 – a C8 (octyl) stationary phase, ideally suited for very hydrophobic compounds which might otherwise elute late on a C18 column.

- [Furosemide and related substances](#)




Kinetex C18 – All purpose C18 with excellent methylene selectivity

- [Lamivudine and related substances](#)
- [Acetylsalicylic acid, impurities and related substances](#)




Kinetex XB-C18 – core-shell particles bonded with C18 chains with additional protective butyl side-chains

- [Related substances of Clarithromycin on a Kinetex 3.5 μm XB-C18 column](#)
- [Albuterol organic impurities](#)
- [Atorvastatin related substances](#)



Kinetex EVO – core-shell technology combined with organo-silica grafting, bonded with C18 to provide a robust stationary phase resistant to extremes of pH

- [Ibuprofen related substances on fully porous and core-shell organo-silica HPLC columns](#)



The USP chapter 467 outlines the considerations and testing required by manufacturers of API and drug products. Testing is conducted using static headspace extraction followed by testing using a G43 or G16 stationary phase GC column.

- G43 – 6% cyanopropylphenyl-94% dimethylpolysiloxane – ZB-624*PLUS*
- G16 – Polyethylene glycol compound (average molecular weight about 15,000). A high molecular weight compound of polyethylene glycol and a diepoxide linker – ZB-Wax*PLUS*

We have developed high performing columns with industry leading levels of inertness for this application:

- [ZB-624*PLUS*™](#)
- [ZB-WAX*PLUS*™](#)




These Plus GC solutions provide improved performance thanks to increases in:

- Stability
- Reproducibility
- Inertness
- Low Bleed

Classification of Residual Solvents by Risk Assessment

Residual Solvent Class	Risk Assessment	Appropriate Analytical Procedure
Class 1	<ul style="list-style-type: none"> • Known/strongly suspected human carcinogens • Environmentally hazardous • Solvents to be avoided 	If Class 1 solvents are used or produced and are not removed by a process, then these solvents should be identified (Procedure A & B) and quantified (Procedure C)
Class 2	<ul style="list-style-type: none"> • Nongenotoxic animal carcinogens • Possible irreversible toxicity • Suspected reversible toxicity • Solvents to be limited 	If Class 2 solvents are present at greater than the Concentration Limit specified by <467>, Option 1 or 2 limits, then these solvents should be identified (Procedure A & B) and quantified (Procedure C)
Class 3	<ul style="list-style-type: none"> • Solvents with low toxic potential • No health-based exposure limit [Class 3 residual solvents have PDEs of 50 mg or more per day] 	If Class 3 solvent limit in an individual monograph is greater than 50 mg per day, then that solvents should be identified (Procedure A & B) and quantified (Procedure C)

USP Class 1 Standard Solution on Zebtron™ ZB-624PLUS™ GC Column



Same conditions for all separations:

Column: Zebtron ZB-624PLUS
Dimensions: 30 meter x 0.32 mm x 1.80 µm
Part No.: [ZHM-G040-31](#)
Injection: Split 5:1 @ 140 °C, 1 µL
Recommended Liner: Zebtron PLUS Straight Z-Liner™
Liner Part No.: [AG2-0A03-05](#) (for Agilent® & Thermo Scientific® systems)

Carrier Gas: Helium @ 2.2 mL/min (constant flow)
Oven Program: 40 °C for 20 min, to 240 °C at 10 °C/min
Detector: FID @ 250 °C
Sample

1. 1,1-Dichloroethene
2. 1,1,1-Trichloroethane
3. Carbon tetrachloride
4. Benzene
5. 1,2-Dichloroethane

Nitrosamines are known to be carcinogenic and have the potential to be intermediates in organic synthesis. Due to their potent genotoxicity, nitrosamines have been a serious cause for concern. It has now become a requirement to accurately quantitate this group of compounds in pharmaceuticals during drug development and manufacturing. The formation of Nitrosamine Drug Substance Related Impurities (NDSRI) is also of concern during the API manufacturing process, and their presence presents additional separation challenges. We have developed resources to aid you in the development of LC-MS quantitation.

Nitrosamines, Meeting Regulatory Requirements

Application Handbook

Chromatography Solutions for API Manufacturers

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 All other countries/regions

Corporate Office USA
t: +1 (310) 212-0555
www.phenomenex.com/chat



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions. Kinetex, Luna, Phenex, Zebtron ZB-624PLUS, and ZB-WAXPLUS are trademarks of Phenomenex. SCIEX is a registered trademark of SCIEX Ptd. Ltd. Agilent is a registered trademark of Agilent Technologies, Inc. Thermo Scientific is a registered trademark of Thermo Fisher Scientific Inc. Comparative separations may not be representative of all applications. **FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.** © 2024 Phenomenex, Inc. All rights reserved.