



# Kinetex Core-Shell LC Columns

When Innovation Meets Quality



- Award Winning Particle Technology
- Reproducible, Consistent Results
- High Efficiency, Resolution, and Sensitivity



[www.phenomenex.com/Kinetex](http://www.phenomenex.com/Kinetex)

 **phenomenex™**

# Table of Contents



**Kinetex™**  
Core-Shell Technology

## Industry Applications

## LC Accessories

# Multi-Award Winning Technology

**Kinetex Core-Shell** is a multi-award winning technology, first introduced in 2009 by Phenomenex R&D experts, that has proven to deliver dramatic improvements in efficiency over conventional fully porous media providing advantages such as increased resolution, higher productivity, reduced solvent consumption, and decreased overall costs.

Whether you are running HPLC or UHPLC methods, the Kinetex core-shell family delivers dramatically improved performance over your current column.

Phenomenex designs and manufactures its own silica and organo-silica core-shell particles to ensure quality control throughout the entire manufacturing process. The combination of a consistent, solid high density core along with proprietary column packing technologies ensures optimum bed structure and high column performance.

- Performance gains on ANY LC system
- System-to-system and lab-to-lab method transferability
- Improve the productivity of dated, established legacy methods

## Multi-Award-Winning Technology




# Strong History of Continuous Innovation



**Kinetex™**  
Core-Shell Technology

Kinetex has over 15 years of proven innovation, consistently delivering reliability and top performance across phases with robust retention mechanisms, meeting analytical needs across multiple industries.



Find your core-shell solution at  
[www.phenomenex.com/Kinetex](http://www.phenomenex.com/Kinetex)



# Strong History of Continuous Innovation



**KINETEX<sup>®</sup>**  
 Core-Shell Technology

- 2 Particle Sizes
- 3 Phases

2009

Chromatographers in the following industries can benefit from the advantages of Core-Shell technology:

-  **Agriculture**
-  **Forensics**
-  **Clinical**
-  **Life Science**
-  **Environmental**
-  **Pharmaceutical**
-  **Food and Beverage**
-  **Chemical/Industrial**
-  **Consumer Products**

- 5 Particle Sizes
- 11 Phases
- Available in Micro, Analytical, Semi Prep and Prep Columns



**Kinetex<sup>™</sup>**  
 Core-Shell Technology

NOW

| Fully Porous                                                                        |    | Kinetex Core-Shell                                                                    |  | Average Efficiency Gain with Kinetex* |
|-------------------------------------------------------------------------------------|----|---------------------------------------------------------------------------------------|--|---------------------------------------|
|  | vs |    |  | <b>90 % Higher</b>                    |
|  | vs |    |  | <b>85 % Higher</b>                    |
| Fully Porous                                                                        |    | Kinetex Core-Shell                                                                    |  | Average Efficiency Gain with Kinetex* |
|  | vs |  |  | <b>20 % Higher</b>                    |
|  | vs |  |  | <b>50 % Higher</b>                    |


\* May not be representative of all applications

# Core-Shell Advantage



By using sol-gel processing techniques that incorporate nano-structuring technology, a durable and homogeneous porous shell is grown on a solid silica core to create a Kinetex Core-Shell particle. This particle morphology reduces broadening contributions resulting in extremely high efficiencies when compared to columns featuring fully porous particles.

## Kinetex 2.6 $\mu\text{m}$ Core-Shell Particle




- Higher throughput without sacrificing resolution
- Easy method transfer across LC system platforms
- Reduce solvent consumption with faster analysis
- Achieve lower levels of detection and quantitation



Kinetex Core-Shell



Fully Porous





# Impact of Selectivity on Resolution




**Kinetex™**  
Core-Shell Technology

Selectivity ( $\alpha$ ) has the greatest impact on chromatographic resolution ( $R$ ) when compared to other parameters. Often, the simplest and most effective way to improve your chromatographic results is to change your column's phase or solid support. Phenomenex offers a wide breadth of stationary phase chemistries across multiple solid supports for simplified method development and optimization.

## The Impact of Selectivity on Resolution





Find your Ideal  
Kinetex Selectivity  
on pages 14-15



# Improved Performance

Kinetex 5  $\mu\text{m}$  columns have efficiencies and peak capacities on par with traditionally fully porous 3  $\mu\text{m}$  columns while still only requiring backpressures seen with traditional 5  $\mu\text{m}$  columns in both isocratic and gradient methods.



## Advantages of Low Backpressures with Kinetex 5 $\mu\text{m}$ Columns


- Longer column lifetime
- Higher throughput
- Increased system compatibility and method transferability

## Instantly Improve 5 $\mu\text{m}$ and 3 $\mu\text{m}$ Methods



Immediately improve resolution, productivity, and sensitivity of your current 3  $\mu\text{m}$  and 5  $\mu\text{m}$  HPLC methods with **Kinetex 5  $\mu\text{m}$**  core-shell technology. This core-shell particle was specifically developed for use on standard or older model HPLC systems that may have low pressure limitations.

## Resolution of Core-shell vs. Fully Porous Higher Resolution with no Pressure Increase

### Waters® XBridge® 5 $\mu\text{m}$ C18



### Waters XBridge 3.5 $\mu\text{m}$ C18




Waters and XBridge are registered trademarks of Waters Corporation. Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

# Ultra-High Performance on Both U/HPLC Systems

Dramatically improve the productivity and performance of your existing methods with the use of shorter Kinetex columns, while decreasing your solvent usage! On a low volume HPLC or UHPLC system **Kinetex 2.6µm** columns will provide up to 3x the efficiency of 5µm and potentially double the efficiency of 3µm fully porous media.

## Performance with Kinetex 2.6µm

### HPLC System




#### Conditions for all columns

Columns: Kinetex 2.6µm C18  
Part No.: [00D-4462-E0](#)  
Dimensions: 100 x 4.6 mm  
Mobile A: Water with 0.1% TFA  
Phase: B: Acetonitrile with 0.1% TFA  
Gradient: Time (min) % B  
0 10  
20 70  
Flow Rate: 1.2 mL/min  
Temperature: Ambient  
Detection: UV @ 210 nm  
Sample: Mupirocin degradants

\*Agilent 1100 was optimized with the Core-Shell Performance Enhancement Kit [AQ0-8892](#).  
Comparative separations may not be representative of all applications.

### UHPLC System




Learn more about Kinetex LC solutions at  
[www.phenomenex.com/Kinetex](http://www.phenomenex.com/Kinetex)




# Ultra-High Efficiency with Core-Shell Particles

The band broadening (wide peaks) and lengthy retention times of traditional fully porous products can limit your results. Kinetex Core-Shell ultra-high efficient performance achieves shorter run times, higher levels of sensitivity, and overall better HPLC or UHPLC results.

## Core-Shell vs. Fully Porous Efficiency Levels (plates/m)



## Core-Shell Performance Gains



# Scalable Particle Platform


Easily develop and transfer your method between systems



**Kinetex™**  
Core-Shell Technology

With Kinetex 5  $\mu\text{m}$ , 3.5  $\mu\text{m}$ , 2.6  $\mu\text{m}$ , 1.7  $\mu\text{m}$ , and 1.3  $\mu\text{m}$  Core-Shell Technology, you are no longer restricted from developing high-performance LC methods. These five scalable Kinetex particle sizes allow to develop and transfer your method effortlessly from system to system.

## Complete Scalable Solution from UHPLC to HPLC to PREP LC



Incredible UHPLC efficiency and performance gains

20 % higher efficiency than fully porous 1.7  $\mu\text{m}$  columns


Achieve sub-2  $\mu\text{m}$  performance on HPLC and UHPLC systems

Instantly improve your pharmacopoeia (Ph. Eur. & USP) monographs that require 3.5  $\mu\text{m}$  particle size

3  $\mu\text{m}$  or better efficiencies at 5  $\mu\text{m}$  pressures for HPLC and PREP LC methods

## Scalability: UHPLC/HPLC/PREP LC

### Gingerols




App ID 21433, 21434, 21435, 21436

# Scalable Particle Platform



**Kinetex™**  
Core-Shell Technology

## Kinetex 5µm C18 on Shimadzu® LC-20A



Columns: Kinetex 5µm C18  
Dimension: 250 x 4.6 mm  
Part No.: [00G-4601-E0](#)

Conditions are the same except as noted:

Mobile Phase: Water/Acetonitrile/  
Phosphoric Acid  
(600:400:2)


Flow Rate: 1 mL/min

Temperature: Ambient

Detection: UV @ 237 nm


Sample: 1. Impurity A  
2. Impurity B  
3. Impurity C  
4. Acetylsalicylic acid  
5. Impurity D

## Kinetex 2.6µm C18 on Agilent® 1100



Columns: Kinetex 2.6µm C18  
Dimension: 150 x 4.6 mm  
Part No.: [00F-4462-E0](#)

## Kinetex 1.7µm C18 on Agilent 1290



Columns: Kinetex 1.7µm C18  
Dimension: 50 x 3.0 mm  
Part No.: [00B-4475-Y0](#)  
Mobile Phase: 680:320:2

## Particle Size Recommended Based on System

|         | 5µm | 3.5µm | 2.6µm | 1.7µm | 1.3µm |
|---------|-----|-------|-------|-------|-------|
| UHPLC   |     |       |       |       |       |
| HPLC    |     |       |       |       |       |
| PREP LC |     |       |       |       |       |


To find Kinetex applications go to  
<https://www.phenomenex.com/applications>



# Adaptability and Method Transfers


UHPLC methods developed with fully porous sub-2  $\mu\text{m}$  columns often generate higher backpressure that only certain systems can run. With the **Kinetex 2.6  $\mu\text{m}$**  particle performance you are no longer restricted by system limitations for your HPLC or UHPLC method development.

## Kinetex 2.6 $\mu\text{m}$ , 4.6 mm ID on Agilent® 1100




Column: Kinetex 2.6  $\mu\text{m}$  C18  
Dimensions: 50 x 4.6 mm  
Part No.: [00B-4462-E0](#)  
Mobile Phase: Acetonitrile / Water (50:50)  
Flow Rate: 2.35 mL/min\*  
Temperature: Ambient  
Detection: UV @ 254 nm  
Sample: 1. Uracil  
2. Acetophenone  
3. Toluene  
4. Naphthalene

## Kinetex 2.6 $\mu\text{m}$ , 3.0 mm ID on Shimadzu® Prominence UFLC®




Column: Kinetex 2.6  $\mu\text{m}$  C18  
Dimensions: 50 x 3.0 mm  
Part No.: [00B-4462-Y0](#)  
Mobile Phase: Acetonitrile / Water (50:50)  
Flow Rate: 1.0 mL/min\*  
Temperature: Ambient  
Detection: UV @ 254 nm  
Sample: 1. Uracil  
2. Acetophenone  
3. Toluene  
4. Naphthalene

## Kinetex 2.6 $\mu\text{m}$ , 2.1 mm ID on Agilent 1200



Column: Kinetex 2.6  $\mu\text{m}$  C18  
Dimensions: 50 x 2.1 mm  
Part No.: [00B-4462-AN](#)  
Mobile Phase: Acetonitrile / Water (50:50)  
Flow Rate: 0.49 mL/min\*  
Temperature: Ambient  
Detection: UV @ 254 nm  
Sample: 1. Uracil  
2. Acetophenone  
3. Toluene  
4. Naphthalene

## Kinetex 2.6 $\mu\text{m}$ , 2.1 mm ID on Waters® ACQUITY® UPLC®



Column: Kinetex 2.6  $\mu\text{m}$  C18  
Dimensions: 50 x 2.1 mm  
Part No.: [00B-4462-AN](#)  
Mobile Phase: Acetonitrile / Water (50:50)  
Flow Rate: 0.49 mL/min\*  
Temperature: Ambient  
Detection: UV @ 254 nm  
Sample: 1. Uracil  
2. Acetophenone  
3. Toluene  
4. Naphthalene


\*Please note that the flow rates were scaled to maintain the same linear velocity.

# Recommended Selectivities

Quick guide to easily select your phase


| Acids                                  | Bases                                                | Neutrals              | Aromatics                      |
|----------------------------------------|------------------------------------------------------|-----------------------|--------------------------------|
| C18<br>F5<br>Phenyl-Hexyl              | EVO C18<br>XB-C18<br>Biphenyl<br>Polar C18<br>PS-C18 | C18<br>C8<br>Biphenyl | Biphenyl<br>Phenyl-Hexyl<br>F5 |
| Acids, Bases, and Neutrals             | Highly Polar Compounds                               | High pH               | Isomers                        |
| Polar C18<br>Biphenyl<br>EVO C18<br>F5 | Polar C18<br>F5<br>Biphenyl<br>HILIC                 | EVO C18               | F5                             |

## Separation Mechanisms




### Hydrophobicity

The ability of a phase to hydrophobically interact with carbon groups




### Steric Interaction

The ability of a phase to separate compounds based on structural differences




### Hydrogen Bond Donating Capacity

The ability of a phase to hydrogen-bond with proton accepting groups



### Hydrogen Bond Accepting Capacity

The ability of a phase to hydrogen-bond with proton donating groups



### Cation Selectivity at pH 2.8

The ability of a phase to interact with cation groups at acidic pH

### Cation Selectivity at pH 7.0

The ability of a phase to interact with cation groups at neutral pH

# Select the Right Phase for Your Analysis

Combining the high efficiency of Kinetex Core-Shell Technology with an excellent range of surface chemistries gives you the best opportunity for increased resolution.

## Kinetex Phase Selection




Unbonded silica phase for HILIC conditions to provide selectivity for polar compounds



100 % aqueous stable which allows for excellent reversed phase retention and enhanced polar and aromatic selectivity



All-purpose hydrophobic retention and methylene selectivity



Minimized hydrophobic retention for highly hydrophobic compounds



High pH stability from 1-12 to deliver robust methods and improved peak shape for bases



Highly reproducible pentafluorophenylpropyl phase, exceptional for halogenated, conjugated, isomeric, or highly polar compounds



Greater retention and separation of aromatic compounds



Superior peak shape and enhanced separation of basic compounds under neutral and acidic conditions



A multi-modal C18 column with a unique positive surface modification that displays improved peak shape for basic compounds



Polymerically bonded C18 phase specifically developed for the separation of EU and EPA priority PAHs



Combined C18 and polar modified surface that provide polar and non-polar retention alongside 100% aqueous stability

## Simple Selection of the Suitable Column


|                       | Particle Sizes (µm) | Pore Size (Å) | Effective Surface Area (m <sup>2</sup> /g) | Effective Carbon Load (%) | pH Range   | Pressure Stability (bar/PSI) |
|-----------------------|---------------------|---------------|--------------------------------------------|---------------------------|------------|------------------------------|
| <b>Kinetex Phases</b> |                     |               |                                            |                           |            |                              |
| Kinetex Polar C18     | 2.6                 | 100           | 200                                        | 9                         | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex PS C18        | 2.6                 | 100           | 200                                        | 9                         | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex C18           | 1.3, 1.7, 2.6, 5    | 100           | 200                                        | 12                        | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex EVO C18       | 1.7, 2.6, 5         | 100           | 200                                        | 11                        | 1.5 – 12   | 1,034/15,000                 |
| Kinetex XB-C18        | 1.7, 2.6, 3.5, 5    | 100           | 200                                        | 10                        | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex C8            | 1.7, 2.6, 5         | 100           | 200                                        | 8                         | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex Biphenyl      | 1.7, 2.6, 5         | 100           | 200                                        | 11                        | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex Phenyl-Hexyl  | 1.7, 2.6, 5         | 100           | 200                                        | 11                        | 1.5 – 8.5* | 1,034/15,000                 |
| Kinetex F5            | 1.7, 2.6, 5         | 100           | 200                                        | 9                         | 1.5 – 8.5* | 1,034/15,000                 |

† Shipping conditions may vary slightly in terms of organic to aqueous ratio, depending on column dimensions.

\* pH stability under gradient conditions. pH stability is 1.5-10 under isocratic conditions.


\*\*Pressure limits are stable for all Kinetex columns 4.6 mm ID and under. For 10 mm ID Kinetex columns pressure > 413 bar/6000 psi may compromise column longevity. For 21.1, 30, and 50 mm ID Kinetex columns pressure > 241 bar/3500 psi may compromise column longevity.

# Industry Applications



**Kinetex™**  
Core-Shell Technology

- Pharmaceutical
- Clinical/Toxicology
- PFAS Analysis
- Environmental
- Food & Beverages
- Cannabis




# Pharmaceutical Applications




**Kinetex™**  
Core-Shell Technology


## Ibuprofen Tablet USP Dissolution: A Rapid HPLC Alternative to the Traditional UV Method




## Separation of Benzalkonium Chloride Homologs (C12, C14, and C16)



## NSAID Active Ingredient Nimesulide from a Topical Gel Formulation Using a Kinetex 2.6 $\mu$ m Polar C18 Column



## Doxylamine Impurity Profile Using a Kinetex EVO C18 Core-Shell LC Column




### LC-UV Conditions

Columns: Kinetex 2.6  $\mu$ m EVO C18  
 Dimensions: 50 x 2.1 mm  
 Part No.: [00B-4725-AN](#)  
 Mobile Phase: A: 20 mM Ammonium bicarbonate (pH 10)  
 B: Acetonitrile  
 Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 5   |
| 5          | 09  |

  
 Flow Rate: 0.7 mL/min  
 Temperature: Ambient  
 Detector: UV @ 254 nm  
 Sample: Impurity profile of Doxylamine  
 Concentration: 25 mg/mL API in DMSO

## Cold Medicine API Screen Using a Kinetex EVO C18 Core-Shell LC Column



### LC-UV Conditions


Columns: Kinetex 2.6  $\mu$ m EVO C18  
 Dimensions: 100 x 4.6 mm  
 Part No.: [00D-4725-E0](#)  
 Mobile Phase: A: 0.1 % TFA  
 B: 0.1 % TFA in Acetonitrile  
 Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 2   |
| 6          | 30  |
| 11         | 100 |

  
 Flow Rate: 1.25 mL/min  
 Temperature: Ambient  
 Detector: UV @ 254 nm  
 Sample: 

|                     |                          |
|---------------------|--------------------------|
| 1. Maleic acid      | 10. Brompheniramine      |
| 2. Fumaric acid     | 11. Acetylsalicylic acid |
| 3. Phenylephrine    | 12. 4-Nitrophenol        |
| 4. Acetaminophen    | 13. Impurity             |
| 5. Pheniramine      | 14. Dextromethorphan     |
| 6. Doxylamine       | 15. Diphenhydramine      |
| 7. Pyrilamine       | 16. Clemastine           |
| 8. Chlorpheniramine | 17. Ibuprofen            |
| 9. Guaifenesin      |                          |

## Forced Degradation of Cefaclor Analysis Using Kinetex EVO C18 Core-Shell LC Column



### LC-UV Conditions

Columns: Kinetex 2.6  $\mu$ m EVO C18  
 Dimensions: 100 x 4.6 mm  
 Part No.: [00D-4725-E0](#)  
 Mobile Phase: A: 0.1 % H3PO4 in Water  
 B: 0.1 % H3PO4 in Acetonitrile  
 Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 5   |
| 6          | 40  |


  
 Flow Rate: 1.75 mL/min  
 Temperature: Ambient  
 Detector: UV @ 254 nm  
 Sample: 1. Cefaclor

# Pharmaceutical Applications



**Kinetex™**  
Core-Shell Technology

## Pheniramine and Halogenated Derivatives Separated with a Kinetex 2.6 µm EVO C18 Core-Shell Column




### LC-UV Conditions

Columns: Kinetex 2.6 µm EVO C18  
Dimensions: 100 x 4.6 mm  
Part No.: [00D-4725-E0](#)  
Mobile Phase: A: 10 mM Ammonium bicarbonate (pH 10.2)  
B: Acetonitrile  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 50  |
| 4          | 60  |

  
Flow Rate: 1.25 mL/min  
Temperature: Ambient  
Detector: UV @ 210 nm  
Sample: 1. Pheniramine  
2. Chlorpheniramine  
3. Brompheniramine


## Analysis of Common Statin Drugs Using a Kinetex F5 Core-Shell LC Column



### LC-UV Conditions

Columns: Kinetex 2.6 µm F5  
Dimensions: 50 x 2.1 mm  
Part No.: [00B-4723-AN](#)  
Mobile Phase: A: 0.1 % TFA in Water  
B: Acetonitrile  
Gradient: 20 to 90 % B in 1 min. Hold for 4 min  
Flow Rate: 0.5 mL/min  
Temperature: Ambient  
Detector: UV @ 248 nm  
Sample: 1. Lovastatin  
2. Atorvastatin  
3. Pravastatin  
4. Simvastatin

## 100% Aqueous Stability of Kinetex EVO C18 Core-Shell LC Columns



### LC-UV Conditions


Columns: Kinetex 5 µm EVO C18  
Dimensions: 50 x 4.6 mm  
Part No.: [00B-4633-E0](#)  
Mobile Phase: 0 mM Sodium phosphate pH 2.5  
Flow Rate: 1 mL/min  
Temperature: 30 °C  
Detector: UV @ 270 nm  
Sample: 1. Norepinephrine  
2. Epinephrine  
3. Dopamine

# Pharmaceutical Applications



**Kinetex™**  
Core-Shell Technology

## Ph. Eur. Monograph 2217: Lamivudine Related Substances with Ph. Eur. Method Modernization



App ID 27631

### LC-UV Conditions

Columns: Kinetex 2.6  $\mu$ m C18  
Dimensions: 150 x 3.0 mm  
Part No.: [00F-4462-Y0](#) – Method 4  
Mobile Phase: Mobile Phase Table 1  
Flow Rate: 0.82 mL/min – Method 4  
Injection Volume: 3  $\mu$ L – Method 4  
Temperature: 35 °C  
Detector: UV @ 277 nm  
Sample: Waters® ACQUITY Arc® HPLC – Method 4



**Click Each!**



### Tech Note

European Pharmacopoeia Paracetamol Monograph Draft Method: Achieving Improved Sensitivity, Resolution, and Separation for Paracetamol and All 14 Related Impurities using Kinetex 5  $\mu$ m C18 Core-Shell Columns



### Tech Note

LC-MS/MS Quantitative Analysis of NDMA in Ranitidine Active Pharmaceutical Ingredient (API) and Drug Product using the SCIEX® 4500 QTRAP™



### Tech Note

Ph. Eur. Monograph 2217: Lamivudine Related Substances with Ph. Eur. Method Modernization



### Tech Note

Assay of Doxepin Hydrochloride According to USP Monograph Using Three Different HPLC Columns



### Tech Note

USP Butamben Assay and Organic Impurities by LC-UV using the Kinetex 5  $\mu$ m XB-C18 Core-Shell HPLC Column



### Tech Note

Meeting and Surpassing System Suitability for USP Sildenafil Citrate Assay and Organic Impurities



### Tech Note

Meeting and Surpassing System Suitability for USP Fluconazole and Related Impurities



### Tech Note

USP Assay and Organic Impurities (LC-UV) for Chloroquine Phosphate



### Tech Note


USP Assay (LC-UV) for Lopinavir and Ritonavir Tablets

# Clinical/Toxicology Applications



**Kinetex™**  
Core-Shell Technology

## Total Ion Chromatogram (TIC) in Negative Ion Mode of Steroid Isomers Fully Separated in an 8-minute Method Using a Kinetex C18 Column



### LC-MS/MS Conditions

Columns: Kinetex 2.6  $\mu$ m C18  
Dimensions: 50 x 3.0 mm  
Part No.: [00B-4462-Y0](#)  
Mobile Phase: A: 0.5 mM Ammonium Fluoride in Water  
B: Methanol  
Gradient: 


| Time (min) | % B |
|------------|-----|
| 0          | 40  |
| 2          | 50  |
| 4.5        | 75  |
| 5          | 95  |
| 6          | 95  |
| 6.5        | 40  |
| 8          | 40  |

  
Flow Rate: 0.8 mL/min  
Injection Volume: 5  $\mu$ L  
Temperature: 30 °C  
LC System: Agilent® 1290 Infinity  
Detection: MS/MS  
Detector: SCIEX® 7500 Triple Quad™

### MS/MS Conditions

Ion Source: ESI  
Polarity: Negative  
Source Temperature: 700 °C  
GS1: 60 psi  
GS2: 60 psi  
CUR: 40 psi  
CAD: 10  
IS: -3000V  
EP: -10 V

## Total Ion Chromatogram (TIC) in Positive Ion Mode of Steroid Isomers Fully Separated in an 8-minute Method Using a Kinetex C18 Column



### LC-MS/MS Conditions

Columns: Kinetex 2.6  $\mu$ m C18  
Dimensions: 50 x 3.0 mm  
Part No.: [00B-4462-Y0](#)  
Mobile Phase: A: 0.5 mM Ammonium Fluoride in Water  
B: Methanol  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 40  |
| 2          | 50  |
| 4.5        | 75  |
| 5          | 95  |
| 6          | 95  |
| 6.5        | 40  |
| 8          | 40  |

  
Flow Rate: 0.8 mL/min  
Injection Volume: 5  $\mu$ L  
Temperature: 30 °C  
LC System: Agilent® 1290 Infinity  
Detection: MS/MS  
Detector: SCIEX® 7500 Triple Quad™

### MS/MS Conditions

Ion Source: ESI  
Polarity: Positive  
Source Temperature: 700 °C  
GS1: 60 psi  
GS2: 60 psi  
CUR: 40 psi  
CAD: 10  
IS: 3000V  
EP: 10 V

To view the full analyte table and MS transitions, view the full technical note




# Clinical/Toxicology Applications



**Kinetex™**  
Core-Shell Technology

## Separation of 39 Pain Panel Analytes Using Protein Precipitation



### LC-MS/MS Conditions

Columns: Kinetex 2.6 µm Biphenyl  
Dimensions: 50 x 3.0 mm  
Part No.: [00B-4622-Y0](#)  
Mobile Phase: A: 0.1 % Formic Acid in Water  
B: 0.1 % Formic Acid in Methanol  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 15  |
| 3.5        | 95  |
| 5          | 95  |
| 5.01       | 15  |
| 7          | 15  |

  
Flow Rate: 0.5 mL/min  
Injection Volume: Ambient  
Temperature: 25 °C  
LC System: Agilent® 1200 Series  
Detection: MS/MS  
Detector: SCIEX® 4500


| Peak No.                         | Analyte Name     | RT (min) | Q1    | Q3    | Protein Precipitation |            | Phree PLR |            |
|----------------------------------|------------------|----------|-------|-------|-----------------------|------------|-----------|------------|
|                                  |                  |          |       |       | % Rec.                | % CV (N=4) | % Rec.    | % CV (N=4) |
| 1                                | Alprazolam       | 4.8      | 309.1 | 281.1 | 93                    | 4.3        | 107       | 5.6        |
| 2                                | Amphetamine      | 2.3      | 136.1 | 91.1  | 88                    | 9.6        | 81        | 8.8        |
| 3                                | Benzoylcegonine  | 3.3      | 290.1 | 168.1 | 87                    | 3.9        | 91        | 13.8       |
| 4                                | Codeine          | 2.6      | 300.2 | 152.1 | 84                    | 7.5        | 91        | 15.2       |
| 5                                | Diazepam         | 4.9      | 285   | 193.2 | 85                    | 9.9        | 91        | 12         |
| 6                                | MDMA             | 2.9      | 194.1 | 105.1 | 75                    | 0.9        | 95        | 12.4       |
| 7                                | Methamphetamine  | 2.6      | 150.1 | 91    | 78                    | 13.5       | 85        | 15         |
| 8                                | Norbuprenorphine | 3.6      | 414.3 | 83.2  | 95                    | 10.3       | 93        | 13.4       |
| 9                                | Oxazepam         | 4.4      | 287   | 241   | 93                    | 4.1        | 100       | 5.1        |
| 10                               | Oxymorphone      | 2        | 302.1 | 227   | 72                    | 11.7       | 73        | 4.8        |
| 11                               | PCP              | 4        | 244.3 | 91    | 80                    | 7.4        | 95        | 10.3       |
| 12                               | Propoxyphene     | 4        | 340.3 | 266.3 | 93                    | 12.5       | 100       | 10.6       |
| 13                               | Sufentanil       | 4.1      | 387.2 | 238.1 | 92                    | 14.3       | 90        | 3.2        |
| 14                               | 6MAM             | 2.57     | 328.1 | 165.1 | 71                    | 2.9        | 80        | 13.5       |
| 15                               | Buprenorphine    | 3.9      | 468.3 | 55.2  | 93                    | 12.5       | 97        | 10.6       |
| 16                               | Carisoprodol     | 3.9      | 261.1 | 176.2 | 103                   | 11.9       | 98        | 11         |
| 17                               | Clonazepam       | 4.4      | 316.1 | 270.1 | 94                    | 2          | 96        | 6.5        |
| 18                               | EDDP             | 4.2      | 278.2 | 234.2 | 84                    | 4.1        | 88        | 4.8        |
| 19                               | Fentanyl         | 3.9      | 337.3 | 105.1 | 103                   | 9.3        | 105       | 11.6       |
| 20                               | Flunitrazepam    | 4.7      | 314.1 | 268.2 | 86                    | 1.9        | 97        | 3.3        |
| 21                               | Flurazepam       | 4        | 388.2 | 315.2 | 112                   | 13.3       | 104       | 3.8        |
| 22                               | Hydrocodone      | 2.8      | 300.2 | 199   | 85                    | 8.6        | 94        | 14.2       |
| 23                               | Hydromorphone    | 2.1      | 286.1 | 185.1 | 78                    | 7.2        | 82        | 12.6       |
| 24                               | Lorazepam        | 4.3      | 321   | 275   | 100                   | 4.1        | 109       | 12.3       |
| 25                               | MDA              | 2.7      | 180.1 | 133   | 78                    | 7          | 93        | 8.8        |
| 26                               | MDEA             | 3        | 208.2 | 163   | 81                    | 9.7        | 91        | 15.3       |
| 27                               | Meperidine       | 3.4      | 248.2 | 220.2 | 98                    | 4.9        | 95        | 7.2        |
| 28                               | Methadone        | 4.4      | 310   | 265   | 99                    | 14.9       | 99        | 7.1        |
| 29                               | Midazolam        | 4.1      | 326.1 | 291.1 | 107                   | 11.1       | 93        | 8.6        |
| 30                               | Morphine         | 1.9      | 286.1 | 152.1 | 92                    | 15.7       | 94        | 5.4        |
| 31                               | Naloxone         | 2.56     | 328.2 | 212   | 76                    | 15.1       | 104       | 8.5        |
| 32                               | Naltrexone       | 2.8      | 342.2 | 267.1 | 94                    | 14.1       | 92        | 10.5       |
| 33                               | Nordiazepam      | 4.64     | 271   | 140   | 98                    | 5.3        | 98        | 8.9        |
| 34                               | Norfentanyl      | 3.2      | 233.2 | 84.1  | 99                    | 6.5        | 94        | 11.2       |
| 35                               | Normeperidine    | 3.4      | 234.1 | 160.1 | 93                    | 8.5        | 86        | 13.1       |
| 36                               | Norpropoxyphene  | 4.1      | 308.2 | 100.1 | 99                    | 15         | 97        | 4.4        |
| 37                               | Oxycodone        | 2.8      | 316.1 | 241.2 | 95                    | 2.4        | 110       | 7.6        |
| 38                               | Temazepam        | 4.7      | 301.1 | 255.1 | 90                    | 4.6        | 89        | 5.8        |
| 39                               | Tramadol         | 3.2      | 264.1 | 58.1  | 87                    | 10.2       | 89        | 13.9       |
| % Recovery range for 39 analytes |                  |          |       |       | 71-112%               |            | 73-110%   |            |

# Clinical/Toxicology Applications



**Kinetex™**  
Core-Shell Technology

## Representative Chromatogram of 32 Pain Panel Analytes Extracted from Oral Fluid Matrix Utilizing a Kinetex 2.6 $\mu$ m Biphenyl Column



### LC-MS/MS Conditions

Column: Kinetex 2.6  $\mu$ m Biphenyl  
Dimensions: 50 x 4.6 mm  
Part No.: [00B-4622-E0](#)  
Mobile Phase: A: 10 mM Ammonium Formate  
B: Methanol  
Gradient: 


| Time (min) | % B |
|------------|-----|
| 0          | 15  |
| 1          | 70  |
| 3.5        | 95  |
| 5.5        | 85  |
| 5.51       | 15  |
| 7          | 15  |

  
Flow Rate: 0.6 mL/min  
Injection Volume: 5  $\mu$ L  
Temperature: Ambient  
LC System: Agilent® 1260 Infinity  
Detection: MS/MS  
Detector: SCIEX® 4500 Triple Quad™

### MS/MS Conditions

Ion Source: ESI  
Polarity: Positive  
Source Temperature: 650 °C  
GS1: 70  
GS2: 70  
CUR: 25  
IS: 5000

## Separation of THC and Metabolites using $\beta$ -Gone™ Sample Preparation and the Kinetex 2.6 $\mu$ m C18 Column



### LC-MS/MS Conditions

Column: Kinetex 2.6  $\mu$ m C18  
Dimensions: 50 x 201 mm  
Part No.: [00B-4462-AN](#)  
Mobile Phase: A: 0.1 % Formic acid in Water  
B: 0.1 % Formic acid in Methanol  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0.25       | 68  |
| 5.25       | 70  |
| 7.75       | 80  |
| 7.95       | 100 |
| 8.95       | 100 |
| 9.15       | 68  |
| 11.0       | 68  |

  
Flow Rate: 3 mL/min  
Injection Volume: 10  $\mu$ L  
Temperature: 40 °C  
LC System: Agilent® 1260 Infinity  
Detection: MS/MS  
Detector: SCIEX API 4000 QTRAP®  
Analytes: 1. D9 THC-OH  
2. D9 THC-COOH  
3. D9 THC  
4. D8 THC

### MS/MS Conditions

Polarity: Positive or Negative  
Source Temperature: 600 °C  
GS1: 50  
GS2: 50  
CUR: 10  
IS: +5500 or -4500

# Clinical/Toxicology Applications



**Kinetex™**  
Core-Shell Technology



Click Each!

Additional applications available on  
<https://www.phenomenex.com/applications>



## Tech Note

25-OH Vitamin D2 and D3 from Serum



## Tech Note

A Fast Approach of a Supported Liquid Extraction (SLE) Method to Determine 25-OH Vitamin D2/D3 in Human Serum Using LC-MS/MS ([TN-0128](#))



## Tech Note

Analysis of Fat-Soluble Vitamins using the Kinetex Core-Shell Biphenyl LC Column



## Tech Note

Analyzing Testosterone in Human Serum by UHPLC using High Efficiency Kinetex 1.7µm C18 Core-Shell Columns



## Tech Note

Comparison of Two Clean-up Techniques and a Quick Analysis of 39 Pain Management Drugs from Serum ([AN-1019](#))



## Tech Note

Evaluation of Newer HPLC/UHPLC Technologies for the Analysis of Steroids ([TN-1097](#))



## Tech Note

Increased Sensitivity of THC and Metabolites using B-Gone and Kinetex 2.6 µm C18 Column by LC-MS/MS



## Tech Note

LC-MS/MS Analysis of a Comprehensive Steroid Panel from Serum for Using Supported Liquid Extraction



## Tech Note

LC-MS/MS based Quantification of 47 Therapeutic Drug Monitoring Compounds in Serum: A Simple Sample Preparation Strategy for Efficient Analysis ([TN-0159](#))



## Tech Note

Quantitation of Clinical Research Steroid Analytes from Serum Utilizing Solid Phase Extraction with LC-MS/MS



## Tech Note

Quantitative Analysis of Illicit Drugs of Abuse in Whole Blood by LC-MS/MS



## Tech Note

Steroid Panel by LC-MS/MS




## Tech Note

Simplifying Urine Drug Testing by Combining IMCSzyme® RT and β-Gone™ Plus β-Glucuronidase Removal




## TDCA and PFOS Separation on Kinetex C18



### LC-MS/MS Conditions

Column: Kinetex 1.7  $\mu$ m C18  
Dimensions: 50 x 2.1 mm  
Part No.: 00B 4475 AN (Kinetex)  
Mobile Phase: A: Acetonitrile  
B: 2 mM Ammonium Acetate in Water / Acetonitrile (95:5, v/v)  
Gradient: Time (min) % B Flow Rate ( $\mu$ L/min)  
0 98 350  
0.2 98 350  
4 70 400  
7 45 400  
9 25 400  
10 5 400  
10.4 98 400  
11.8 98 400  
12 98 350  
15 98 350  
Injection Volume: 5  $\mu$ L  
Temperature: 40 °C  
Instrument: Agilent® 1260 Quaternary  
Detection: MS/MS  
Detector: SCIEX® Triple Quad™ 4500


## Perfluorinated Compounds in Aqueous Matrices



### LC-MS/MS Conditions

Column: Kinetex 5  $\mu$ m EVO C18 100Å  
Dimensions: 100 x 2.1 mm  
Part No.: 00D-4633-AN  
SecurityGuard™: AJO-9298  
Mobile Phase: A: 0.4 % v/v Ammonium hydroxide in Water  
B: Methanol  
Gradient: Time (min) % B  
0 90  
3.1 20  
4.5 20  
6.1 90  
11 90  
14 90  
Flow Rate: 0.3 mL/min  
Injection: 5  $\mu$ L  
Temperature: Ambient  
Instrument: Thermo TSQ Quantum® Ultra QQQ  
HPLC System: Thermo Accela™ 1250  
Sample: Analyte Retention Time (Min)  
1. PFBA 5.69  
2. PFPeA 6  
3. PFBS 6.3  
4. PFHxA 6.83  
5. PFHpA 7.39  
6. PFHxS 7.5  
7. 6:2 FTS 7.68  
8. PFOA 7.7  
9. PFOS 7.7  
10. PFNA 7.92  
11. PFDS 7.9  
12. 8:2 FTS 8.08  
13. PFDA 8.09  
14. N-MeFOSAA 8.15  
15. PFDA 8.2  
16. N-EtFOSAA 8.23  
17. PFUnDA 8.25  
18. PFDoA 8.4  
19. PFTeDA 8.56  
20. PFTeDA 8.74


## Method Validation for the California Expanded List of PFAS Compounds



### HPLC Conditions (for both)

Column: Kinetex EVO C18  
100 x 2.1mm (00D-4633-AN)  
Delay Column: Kinetex EVO C18  
50 x 2.1mm (00B-4633-AN)  
Mobile Phase: A: 20 mM Acetic Acid  
B: 25 mM Ammonium Hydroxide in Methanol  
Gradient: Time (min) % B  
0 5  
1.2 45  
3.6 65  
11 90  
13 90  
13.01 5  
17 5  
Column Temp: 350°C  
Injection Volume: 10 µL  
HPLC System: Agilent® 1100

## Matrix Studies Chromatographic Separation of Spiked PFAS Compounds in a Wastewater Sample






## Polycyclic Aromatic Hydrocarbons (PAHs) in Water by LC/UV


### LC-UV Conditions

Column: Kinetex 2.6  $\mu$ m C18  
Dimensions: 100 x 4.6 mm  
Part No.: [00D-4462-E0](#)  
SecurityGuard™: [AJ0-8768](#)  
Mobile Phase: A: Water  
B: Acetonitrile  
Gradient: Time (min) % B  
0 70  
10 100  
Flow Rate: 1.5 mL/min  
Temperature: 30 °C  
Detection: UV @ 254 nm (ambient)


Sample: 1. Naphthalene  
2. Acenaphthylene  
3. Fluorene  
4. Acenaphthene  
5. Phenanthrene  
6. Anthracene  
7. Floranthene  
8. Pyrene  
9. Chrysene  
10. Benz[a] anthracene  
11. Benzo[b] fluoranthene  
12. Benzo[k] fluoranthene  
13. Benzo[a] pyrene  
14. Dibenzo[a,h] anthracene  
15. Indeno[1,2,3-cd] pyrene  
16. Benzo [g,h,i] Perylene



## Separation of 16 PAHs in EPA Method 610 Using a Core-Shell Kinetex 3.5 $\mu$ m PAH 100 X 4.6 mm Column




## Improved Analysis of the EPA Method 610 PAH Mixture a Core-Shell Kinetex 3.5 $\mu$ m PAH 100 x 4.6 mm Column Under Adjusted Gradient Conditions






## Multi-Residue Pesticides




|                        |                          |                   |                         |                   |                      |                         |                              |                          |
|------------------------|--------------------------|-------------------|-------------------------|-------------------|----------------------|-------------------------|------------------------------|--------------------------|
| 1. Methamidophos       | 21. Aldicarb             | 41. Acetamiprid   | 61. Carbetamide         | 81. Oxadixyl      | 101. Clofentezine    | 121. Diniconazole       | 141. Clethodim Z             | 161. Azoxystrobin        |
| 2. Methomyl            | 22. Promecarb            | 42. Mexacarbate   | 62. Oxamyl              | 82. Metalaxyl     | 102. Fenpropimorph   | 122. Benalaxyl          | 142. Etoxazole               | 162. Trifloxystrobin     |
| 3. Fenuron             | 23. Aminocarb            | 43. Dioxacarb     | 63. Dicrotophos         | 83. Penconazole   | 103. Buprofezin      | 123. Dimoxystrobin      | 143. Benzoximate             | 163. Benfuracarb         |
| 4. Cyromazine          | 24. Propoxur             | 44. Mepanipirim   | 64. 3-Hydroxycarbofuran | 84. Vamidothion   | 104. Fenazaquin      | 124. Diclobutrazol      | 144. Flufenacet              | 164. Spirodiclofen       |
| 5. Acephate            | 25. Ethirimol            | 45. Monocrotophos | 65. Pirimicarb          | 85. Myclobutanil  | 105. Quinoxifen      | 125. Pencycuron         | 145. Pyridaben               | 165. Carfentrazone-ethyl |
| 6. Fuberidazole        | 26. Acibenzolar-S-methyl | 46. Bendiocarb    | 66. Prometryn           | 86. Chloroxuron   | 106. Tebuconazole    | 126. Epoxiconazole      | 146. Picoxystrobin           | 166. Mandipropamid       |
| 7. Propamocarb         | 27. Chlortoluron         | 47. Mevinphos E   | 67. Terbutryn           | 87. Thiamethoxam  | 107. Diflubenzuron   | 127. Fenarimol          | 147. Propargite              | 167. Fenpyroximate       |
| 8. Tricyclazole        | 28. Omethoate            | 48. Mevinphos Z   | 68. Forchlorfenuron     | 88. Paclobutrazol | 108. Fenamidone      | 128. Tebufenpyrad       | 148. Tetraconazole           | 168. Fluoxastrobin       |
| 9. Butocarboxim        | 29. Simetryn             | 49. Terbutmeton   | 69. Linuron             | 89. Triadimefon   | 109. Hexaconazole    | 129. Ipconazole         | 149. Spirotetramat           | 169. Temephos            |
| 10. Carbenazim         | 30. Metribuzin           | 50. Ethiofencarb  | 70. Clothianidin        | 90. Amitraz       | 110. Kresoxim-methyl | 130. Zoxamide           | 150. Fluquinconazole         | 170. Flufenoxuron        |
| 11. Isoprocarb         | 31. Monolinuron          | 51. Methiocarb    | 71. Thiacloprid         | 91. Triadimenol   | 111. Nuairimol       | 131. Fenbuconazole      | 151. Prochloraz              | 171. Butafenacil         |
| 12. Cymoxanil          | 32. Pymetrozine          | 52. Prometon      | 72. Imidacloprid        | 92. Imazalil      | 112. Flusilazole     | 132. Propiconazole      | 152. Bromuconazole cis       | 172. Novaluron           |
| 13. Cycluron           | 33. Pyracarbolid         | 53. Secbumeton    | 73. Thiobencarb         | 93. Spiroxamine   | 113. Bupirimate      | 133. Boscalid           | 153. Bromuconazole trans     | 173. Hydramethylnon      |
| 14. Pyrimethanil       | 34. Thiofanox            | 54. Ametryn       | 74. Metobromuron        | 94. Mefenacet     | 114. Desmedipham     | 134. Thiophanate-methyl | 154. Furathiocarb            | 174. Metaflumizone       |
| 15. Carbaryl           | 35. Thidiazuron          | 55. Tebuthiuron   | 75. Fludioxinil         | 95. Bifenazate    | 115. Triticonazole   | 135. Triflumizole       | 155. Pyraclostrobin          | 175. Chlorfluazuron      |
| 16. Thiabendazole      | 36. Formetanate          | 56. Dimethoate    | 76. Diethofencarb       | 96. Phenmedipham  | 116. Metconazole     | 136. Hexythiazox        | 156. Dimethomorph (isomer 2) | 176. Spinosyn A          |
| 17. Dinotefuran        | 37. Methabenzthiazuron   | 57. Diuron        | 77. Mepronil            | 97. Fenhexamid    | 117. Iprovalicarb    | 137. Tebufenozide       | 157. Dimethomorph (isomer 1) | 177. Spinosyn D          |
| 18. Aldicarb-sulfoxide | 38. Carbofuran           | 58. Fluometuron   | 78. Nitenpyram          | 98. Flutriafol    | 118. Pyriproxyfen    | 138. Piperonyl butoxide | 158. Rotenone                | 178. Spinetoram A        |
| 19. Isoprotruron       | 39. Aldicarb-sulfone     | 59. Siduron       | 79. Methoprotirne       | 99. Furalaxyl     | 119. Flutolanil      | 139. Triflumuron        | 159. Ethiprole               | 179. Spinetoram B        |
| 20. Fenobucarb         | 40. Butoxycarboxim       | 60. Carboxin      | 80. Neburon             | 100. Fenoxycarb   | 120. Cyazofamid      | 140. Clethodim E        | 160. Alanycarb               |                          |

## Multi-Toxin Screening in Cereal Based Foods by LC-MS/MS

### Negative Polarity



### Positive Polarity



TIC of all analytes with negative and positive fast polarity switching.

#### LC-MS/MS Conditions


Column: Kinetex 2.6µm XB-C18  
Dimensions: 50 x 2.1 mm  
Part No.: [00B-4496-AN](#)  
Mobile Phase: A: Water with 5mM Ammonium acetate and 0.5 % Acetic acid  
B: Methanol with 5mM Ammonium acetate and 0.5 % Acetic acid  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 2   |
| 2          | 2   |
| 5          | 80  |
| 5.2        | 98  |
| 8          | 98  |

  
Flow Rate: 450 µL/min  
Temperature: Ambient  
Detection: Tandem Mass Spectrometer (MS/MS) 550°C  
Instrument: SCIEX® API 5500™  
Sample: 1. 15-Acetyldeoxynivalenol (15-AcDON)  
2. 3-Acetyldeoxynivalenol (3-AcDON)  
3. Aflatoxin B1 (AFB1)  
4. Aflatoxin B2 (AFB2)  
5. Aflatoxin G1 (AFG1)  
6. Aflatoxin G2 (AFG2)  
7. Deoxynivalenol (DON)  
8. Diacetoxyscirpenol (DAS)  
9. Fumonisin B1 (FB1)  
10. Fumonisin B2 (FB2)  
11. Fusarenon X (FUS X)  
12. HT-2 toxin (HT2)  
13. Monoacetoxyscirpenol (MAS)  
14. Nivalenol (NIV)  
15. Ochratoxin A (OTA)  
16. T-2 toxin (T2)  
17. Zearalenone (ZON)

## Neonicotinoid Pesticide Residues in Honey by QuEChERS and LC-MS/MS

roQ™ QuEChERS AOAC kits were used for extraction and clean up.



Note: 400 ng/g extracted neonicotinoids spiked in honey.

#### LC-MS/MS Conditions

Column: Kinetex 2.6µm Biphenyl  
Dimensions: 50 x 2.1 mm  
Part No.: [00B-4622-AN](#)  
SecurityGuard™ ULTRA Cartridges: [AJ0-8782](#)  
SecurityGuard ULTRA Cartridge [AJ0-9000](#)  
Holder:  
Mobile Phase: A: 0.1 % Formic acid in Water  
B: 0.1 % Formic acid in Acetonitrile  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0.1        | 5   |
| 10         | 70  |
| 12         | 70  |
| 12.1       | 5   |
| 14.1       | 5   |

  
Flow Rate: 0.5 mL/min  
Injection Volume: 25 µL  
Temperature: 40 °C  
Detection: Tandem Mass Spec (MS/MS)  
Instrument: SCIEX API 4500 QTRAP®  
Sample: 1. Nitenpyram  
2. Thiamethoxam  
3. Clothianidin  
4. Imidacloprid  
5. Acetamiprid  
6. Azoxystrobin  
7. Coumaphos  
8. Chlorpyrifos



## Antibiotics in Meat by LC-MS/MS


### LC-MS/MS Conditions

Column: Kinetex 2.6µm C18  
Dimensions: 50 x 2.1 mm  
Part No.: [00B-4462-AN](#)  
Mobile A: 0.1 % Formic Acid in Water  
Phase: B: 0.1 % Formic Acid in Methanol  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 2   |
| 0.30       | 80  |
| 7.27       | 99  |
| 7.37       | 99  |
| 8.27       | 2   |

  
Flow Rate: 0.5mL/min  
Temperature: 40 °C  
Inj. Volume: 10µL  
Detection: Mass Spectrometer (MS) (300 °C)  
Sample: 

|                            |                        |
|----------------------------|------------------------|
| 1. Sulfanilamide           | 17. Spiramycin         |
| 2. Amoxicillin             | 18. Chlorotetracycline |
| 3. Lincomycin              | 19. Sulfadoxine        |
| 4. Sulfadiazine            | 20. Clindamycin        |
| 5. Sulfathiazole           | 21. Tilmicosin         |
| 6. Ampicillin              | 22. Chloramphenicol    |
| 7. Thiamphenicol           | 23. Sulfadimethoxine   |
| 8. Sulfamerazine           | 24. Sulfamethoxazole   |
| 9. Tetracycline            | 25. Erythromycin       |
| 10. Ciprofloxacin          | 26. Tulosin            |
| 11. Enrofloxacin           | 27. Josamycin          |
| 12. Danofloxacin           | 28. Penicillin G       |
| 13. Sulfamethazine         | 29. Cloxacillin        |
| 14. Sarafloxacin           | 30. Flunixin           |
| 15. Sulfamethoxypyridazine |                        |
| 16. Florfenicol            |                        |



**30 Antibiotics in  
under 8 minutes**

App ID 18555

## Nitrates and Nitrites by Ion Pair LC-UV


### LC-UV Conditions

Columns: Kinetex 2.6µm C18  
Dimensions: 100 x 3.0 mm  
Part No.: [00D-4462-YO](#)  
Mobile Phase: 5mM Tetrabutylammonium Hydrogen Sulfate + 5mM Ammonium Bicarbonate (pH adjusted to 6.5 with Sodium Hydroxide)  
Gradient: 

| Time (min) | % B |
|------------|-----|
| 0          | 90  |
| 71         | 71  |
| 4          | 49  |
| 6          | 6   |

  
Flow Rate: 0.6mL/min  
Injection: 20 µL (in DI water)  
Temperature: Ambient  
Detection: UV/Vis @ 218 nm (Bandwidth 4 nm)  
HPLC System: Agilent® 1100 HPLC with Quaternary Pump  
Analytes: 


|            |
|------------|
| 1. Nitrite |
| 2. Nitrate |



App ID 23940



## Polyphenols in Wine by LC-UV



Click Each!



### Tech Note

Multi-Class Antibiotics From Ground Meat



### Tech Note

Fast and Robust Analysis of Organic Acids from Wine using HPLC-UV



### Tech Note


Extraction and Separation of Phenolic Antioxidants from Edible Oils using a Kinetex 2.6  $\mu$ m Polar C18 Column

To find more Kinetex Applications for Food, go to  
[www.phenomenex.com/food](http://www.phenomenex.com/food)






## Alpha Beta Gamma Delta Tocopherols and Tocotrienols on Kinetex 2.6 µm F5 150x4.6mm



## Catechins and Phenolic acids on a Kinetex 2.6 µm F5 100 x 2.1



# Cannabis Applications



**Kinetex™**  
Core-Shell Technology


## 18 Cannabinoids for Potency Testing by LC-UV

### LC-UV Conditions

Column: Kinetex 2.6  $\mu$ m C18  
Dimensions: 150 x 4.6 mm  
Part No.: [00F-4462-E0](#)  
Mobile Phase: A: 20 mM Ammonium Formate, pH 2.9 with Formic Acid  
B: Acetonitrile  
Isocratic: Isocratic 24:76 (A/B)  
Flow Rate: 1.5 mL/min  
Injection: 2  $\mu$ L  
Volume:  
Back Pressure: ~260 bar  
Temperature: 40 °C  
Flow Rate: 1.5 mL/min  
Detection: UV @ 228 nm

Analytes: 1. CBDVA Cannabidivarinic acid  
2. CBDV Cannabidivarin  
3. CBDA Cannabidiolic acid  
4. CBGA Cannabigerolic acid  
5. CBG Cannabigerol  
6. CBD Cannabidiol  
7. THCV Tetrahydrocannabivarin  
8. THCVA Tetrahydrocannabivarinic acid  
9. CBNA Cannabinolic acid

10. CBN Cannabinol  
11. EXO-THC Exo-tetrahydrocannabinol  
12. D9-THC  $\Delta^9$ -Tetrahydrocannabinol  
13. D8-THC  $\Delta^8$ -Tetrahydrocannabinol  
14. THCA-A Tetrahydrocannabinolic acid A  
15. CBCA Cannabichromenic acid  
16. CBL Cannabicycol  
17. CBC Cannabichromene  
18. CBLA Cannabicyclic acid




## 10 Cannabinoids by LC-MS/MS

### LC-MS/MS Conditions


Column: Kinetex 2.6  $\mu$ m Biphenyl  
Dimension: 150 x 4.6 mm  
Part No.: [00F-4622-E0](#)  
Mobile Phase: A: Water + 5 mM Ammonium acetate + 0.1% Formic acid  
B: Methanol/Water (98:2) + 5 mM Ammonium acetate

| Gradient | Time (min) | % B |
|----------|------------|-----|
|          | 3.5        | 100 |
|          | 4          | 100 |
|          | 4.1        | 90  |
|          | 5          | 90  |

Flow Rate: 1 mL/min  
Injection: 2  $\mu$ L  
Temperature: Ambient  
Detector: SCIEX QTRAP® 6500+  
Detection: ESI Positive - MS/MS  
Sample: 1. CBDV  
2. CBDA  
3. CBGA  
4. CBG  
5. CBD  
6. THCV  
7. CBN



## Pesticides on Kinetex 2.6 $\mu$ m Biphenyl 50 x 4.6 mm



### LC-MS/MS Method Parameters

Column: Kinetex 2.6  $\mu$ m Biphenyl 50 x 4.6 mm  
 Part No.: [00B-4662-E0](#) (Kinetex)  
[00B-4760-E0](#) (Luna Omega)  
 Mobile Phase: A: 5mM Ammonium formate + 0.1 %  
 Formic acid in Water  
 B: 5mM Ammonium formate + 0.1 %  
 Formic acid in 98:2

| Gradient: | Time (min) | % B |
|-----------|------------|-----|
|           | 0          | 10  |
|           | 1          | 10  |
|           | 4.3        | 80  |
|           | 8.7        | 95  |
|           | 10.5       | 95  |
|           | 10.6       | 10  |
|           | 14         | 10  |

Injection: 10  $\mu$ L  
 Flow Rate: 0.4 mL/min  
 Temperature: 40°C  
 Detection: MS/MS (ESI+)  
 Instrument: SCIEX® API 4000™  
 Detector: SCIEX 3500 Triple Quad

Note: The 16 minute LC gradient chromatographically separates all pesticide residues of interest as well as observed endogenous cannabis flower interferences, and the same gradient and run-time can be used on 150mm length columns for greater capacity and lifetime.

### Analytes:

|                   |                         |                   |                     |
|-------------------|-------------------------|-------------------|---------------------|
| 1. Acephate       | 11. Carbaryl            | 21. Dimethoate    | 31. Hexythiazox     |
| 2. Acequinocyl    | 12. Carbofuran          | 22. Ethoprophos   | 32. Imazalil        |
| 3. Acetamiprid    | 13. Chlorantraniliprole | 23. Etofenprox    | 33. Imidacloprid    |
| 4. Aldicarb       | 14. Chlorpyrifos        | 24. Etoxazole     | 34. Kresoxim-methyl |
| 5. Avermectin B1a | 15. Clofentezine        | 25. Fenoxycarb    | 35. Malathion A     |
| 6. Avermectin B1b | 16. Cyfluthrin          | 26. Fenpyroximate | 36. Metalaxyl       |
| 7. Azoxystrobin   | 17. Cypermethrin        | 27. Fipronil      | 37. Methiocarb      |
| 8. Bifenazate     | 18. Daminozide          | 28. Fipronil NH4  | 38. Methomyl        |
| 9. Bifenthrin     | 19. Diazinon            | 29. Flonicamid    |                     |
| 10. Boscalid      | 20. Dichlorvos          | 30. Fludioxinil   |                     |



Click Each!



#### Guide

Cannabis Testing Guide



#### Tech Note

Comprehensive Cannabis Analysis from One Extract Using One Column and One Solvent System



#### Tech Note

7 Primary Terpenes in Cannabis by LC-MS/MS



#### Tech Note

Determination of Pesticide Residues in Cannabis by LC-MS/MS

For more applications go to  
[www.phenomenex.com/cannabis](http://www.phenomenex.com/cannabis)



# Optimize Your Analysis and Column Lifetime



## The Phenomenex Security Portfolio

1. SecurityCAP™: Minimize solvent contamination and increase lab safety
2. SecurityLINK™: Get better HPLC/UHPLC results with fingertight, zero dead-volume connections
3. SecurityGuard™ and SecurityGuard ULTRA: Increase LC column lifetime


# SecurityCAP™

Protect Any U/HPLC Solvent from Contamination  
and Staff from Hazardous Solvents Vapors.



**Kinetex™**  
Core-Shell Technology

The SecurityCAP mobile phase safety filters prevent hazardous solvent vapors and gases from leaving the solvent reservoir negatively impacting the health of laboratory workers and visitors. The integrated filter membrane captures dust and other contaminants preventing solvent contamination or bacterial growth, which could negatively impact both your chromatography and HPLC/UHPLC system.




# SecurityLINK™

Simplify Your System and Column Connections with SecurityLINK Fingertight Connections




**Kinetex™**  
Core-Shell Technology


The SecurityLINK fingertight fitting connector simplifies your system and column connections using Torque Limiting Technology that prevents column damaging overtightening.




## Tubing Options



PEEKsil™




Stainless Steel



PEEK-Lined Stainless Steel

## What is Torque Limiting Technology?

Once the perfect connection has been made through finger tightening, the SecurityLINK fitting offers a haptic “click” to confirm that optimum torque has been reached. This ensures a consistent connection each and every time and prevents over or under tightening that may cause column or performance issues.



Order now at  
[www.phenomenex.com/SecurityLINK](http://www.phenomenex.com/SecurityLINK)




# SecurityGuard™ and SecurityGuard ULTRA

Protect your LC column. Protect your research outcomes.

SecurityGuard is a cost convenient universal column protection system that effectively protects your U/HPLC columns from the damaging effects of chemical contaminants without altering your chromatographic results.

- Protects chromatographic results
- Extends UHPLC, HPLC and PREP column lifetime
- Does not alter chromatography
- Simple to use
- Saves time and money
- Compatible with most U/HPLC columns

## With SecurityGuard ULTRA



(24000 times magnification)


## SecurityGuard ULTRA Cartridge System

The SecurityGuard ULTRA cartridge system protects ultra-high performance columns, like Kinetex, from damaging contaminants and microparticulates.

High Pressure  
Rated Format

- Extend Kinetex column lifetime
- Simple to use
- Pressure rated to 20000 psi (1378 bar)
- Fits virtually all manufacturers' columns 2.1 to 4.6 mm ID

## Without SecurityGuard ULTRA




(24000 times magnification)

# Shop Online in 6 Quick Steps



**Kinetex™**  
Core-Shell Technology



**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.  
[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

# Kinetex™ Core-Shell LC Columns

## Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
info@phenomenex.com

## China

t: +86 400-606-8099  
cninfo@phenomenex.com

## Czech Republic

t: +420 272 017 077  
cz-info@phenomenex.com

## Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

## Hong Kong

t: +852 6012 8162  
hkinfo@phenomenex.com

## India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

## Indonesia

t: +62 21 3952 5747  
indoinfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

## Japan

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

## Poland

t: +48 22 51 02 180  
pl-info@phenomenex.com

## Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

## Singapore

t: 800-852-3944  
sginfo@phenomenex.com

## Slovakia

t: +420 272 017 077  
sk-info@phenomenex.com

## Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

## Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

## Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

## Thailand

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
info@phenomenex.com

## 🌐 All other countries/regions Corporate Office USA

t: +1 (310) 212-0555  
www.phenomenex.com/chat



## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

### Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Kinetex, SecurityGuard, SecurityCAP, SecurityLINK,  $\beta$ -Gone, roQ, Luna, and BE-HAPPY are trademarks of Phenomenex. SCIEX is a registered trademark and API 4000 is a trademark of AB SCIEX Pte. Ltd. Agilent and ZORBAX are registered trademarks of Agilent Technologies, Inc. ACQUITY and Waters are registered trademarks of Waters Technologies Corporation. AIT is a registered trademark of AIT France. Canary-Safe is a trademark of Analytical Sales and Services. S.C.A.T. is a registered trademark of S.C.A.T.-Europe GmbH. VICI Jour is a registered trademark of Valco Instruments Co. L.P. VapLock is a trademark of Western Fluids, LLC. DURAN is a registered trademark of DWK Life Sciences GmbH. Timestrip is a registered trademark of Timestrip UK Ltd. Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with authors.

### Disclaimer

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.