

# HPLC Troubleshooting Mini Guide

# Baseline Issues



## Content

1. Overview
2. Baseline Issues
3. Key Problem Areas and Preventive Maintenance
4. Protect Your HPLC Column



# Overview

# 1

## Locating and Correcting the Problem

A systematic approach to identifying the problem is the best path to troubleshooting your HPLC system. This guide is organized by four major categories of symptoms to help you quickly identify the source of the problem(s) you are encountering:

- pressure abnormalities
- leaks
- peak problems
- baseline issues



When you have corrected the problem, record the incident in the system recordbook to help with future problems.

## Prevention

Many LC problems can be prevented with routine preventive maintenance such as replacing pump seals regularly. Consistent preventive maintenance practices will enhance lab productivity, avoid system critical damage, equipment downtime and costly repairs.


## Where to Get Additional Help

1. Chat with Phenomenex technical experts. Phenomenex has experienced technical consultants who can assist you with any chromatography issue in real time. To chat now go to [www.phenomenex.com/chat](http://www.phenomenex.com/chat).
2. The operator's and service manuals for the instrument should be consulted. These contain exploded diagrams, troubleshooting procedures for specific models, and part numbers to help you order replacement parts.
3. Other people in the lab may have had experience solving a problem which is giving you trouble; they can be a helpful resource.
4. The manufacturer of your instrument can help you. Most LC manufacturers offer free technical support to their customers.
5. Phenomenex offers seminars on HPLC/UHPLC. Join PhenoAcademy for specialty troubleshooting webinars, [www.phenomenex.com/phenoacademy](http://www.phenomenex.com/phenoacademy)
6. Other resources:
  - J.W. Dolan and L.R. Snyder, **Troubleshooting LC Systems**, Humana Press, NJ (1989).
  - L.R. Snyder and J.J. Kirkland, **Introduction to Modern Liquid Chromatography**, 2nd ed., Wiley, NY (1979).
  - D.J. Runser, **Maintaining and Troubleshooting HPLC Systems - A User's Guide**, Wiley, NY (1981).
  - J.W. Dolan, "LC Troubleshooting", LC/GC Magazine. This is a monthly column.

# Baseline Issues


# 2

## Baseline drift


POSSIBLE CAUSE	SOLUTION
1. Column temperature fluctuation. (Even small changes cause cyclic baseline rise and fall. Most often affects refractive index and conductivity detectors, or UV detectors at high sensitivity or in direct photometric mode.)	1. Control column and mobile phase temperature, use heat exchanger before detector  
2. Nonhomogeneous mobile phase. (Drift usually to higher absorbance, rather than cyclic pattern from temperature fluctuation)	2. a. Prepare fresh mobile phase with HPLC grade reagents b. Gently swirl or mix any premixed mobile phases that have been idle for longer than a day c. Ensure mobile phase is being properly degassed; check degasser.
3. Contaminant or air buildup in detector cell	3. Flush cell with methanol or other strong solvent. If necessary, clean cell with 1N HNO <sub>3</sub> (never with HCl). Add a back pressure regulator after the detector to prevent out gassing in the detector flow cell.
4. Plugged outlet line after detector. (High pressure cracks cell window, producing noisy baseline)	4. Unplug or replace line. Refer to detector manual to replace flow cell.
5. Mobile phase mixing problem or change in flow rate	5. Correct composition / flow rate. To avoid, routinely monitor composition and flow rate
6. Slow column equilibration, especially when changing mobile phase	6. a. Flush with intermediate strength solvent, run 10-20 column volumes of new mobile phase before analysis b. Ion-Pairing & HILIC mobile phases may take up to 60 column volumes for equilibration
7. Mobile phase contaminated, deteriorated, or prepared from low quality materials	7. Check make-up of mobile phase. Use highest grade chemicals and HPLC solvents. Prepare fresh mobile phase.
8. Strongly retained materials in sample (high k') can elute as very broad peaks and appear to be a rising baseline. (Gradient analyses can aggravate problem)	8. Use guard column SecurityGuard™ or SecurityGuard ULTRA is recommended. If necessary, flush column with strong solvent between injections or periodically during analysis
9. Mobile phase recycled but detector not adjusted	9. Reset baseline. Use new mobile phase when dynamic range of detector is exceeded
10. Detector (UV) not set at absorbance maximum but at slope of curve	10. Change wavelength to UV absorbance maximum

# Baseline Issues (continued)

## Baseline noise (regular)

POSSIBLE CAUSE	SOLUTION
1. Air in mobile phase, detector cell, or pump	1. Degas mobile phase. Flush system to remove air from detector cell or pump
2. Leak	2. Check system for loose fittings. Check pump for leaks, salt build-up, unusual noises. Change pump seals if necessary
	
3. Incomplete mobile phase mixing	3. Mix mobile phase by hand or use less viscous solvent
4. Temperature effect (column at high temperature, detector unheated)	4. Reduce differential or add heat exchanger
5. Other electronic equipment on same line	5. Isolate LC, detector or recorder to determine if source of problem is external. Correct as necessary
6. Pump pulsations	6. Incorporate pulse dampener into system

## Baseline noise (irregular)

POSSIBLE CAUSE	SOLUTION
1. Leak	1. Check for loose fittings. Check pump for leaks, salt build-up, unusual noises. Change seals if necessary. Check for detector cell leak
	
2. Mobile phase contaminated, deteriorated, or prepared from low quality materials	2. Check make-up of mobile phase
3. Mobile phase solvents immiscible	3. Select and use only miscible solvents
4. Detector/recorder electronics	4. Isolate detector and recorder electronically. Refer to instruction manual to correct problem
5. Air trapped in system	5. Flush system with strong solvent
6. Air bubbles in detector	6. Purge detector. Install backpressure device after detector
7. Detector cell contaminated (even small amounts of contaminants can cause noise)	7. Clean cell by flushing with 1N HNO <sub>3</sub> (never with HCl)
8. Weak detector lamp	8. Replace lamp
9. Column leaking silica or packing material	9. Replace column
10. Mobile phase mixer inadequate or malfunctioning	10. Repair or replace the mixer or mix off-line if isocratic

# Key Problem Areas and Preventive Maintenance



The chart below lists the most common problems that occur with each LC module. In the right-hand column are listed preventive maintenance practices that can reduce the failure rate. The numbers in parentheses are suggested intervals between maintenance. The operator's and service manuals for your LC may have additional suggestions for preventive maintenance of your model of LC.

## Reservoir

POSSIBLE CAUSE	SOLUTION
1. Blocked inlet frit	1. a. Replace (3-6 mo.) b. Filter mobile phase, 0.5µm filter
2. Gas bubbles	2. Degas mobile phase

## Pump

POSSIBLE CAUSE	SOLUTION
1. Air bubbles	1. Degas mobile phase
2. Pump seal failure	2. Replace (3 mo.)
3. Check valve failure	3. Filter mobile phase, use inlet-line frit. Keep spare.

## Injector

POSSIBLE CAUSE	SOLUTION
1. Rotor seal wear	1. a. Don't overtighten b. Filter samples

## Column

POSSIBLE CAUSE	SOLUTION
1. Blocked frit	1. a. Filter mobile phase b. Filter samples c. Use in-line filter and/or guard column
2. Void at head of column	2. a. Avoid mobile phase pH at or near column max pH limit b. Install column at low flow rates to avoid pressure shock c. Use guard column

# Protect Your HPLC Column

# 4


## Contaminants Can Cause the Following:

- High Backpressure
- Split Peaks
- Broad Peaks
- Baseline Noise
- Baseline Drift
- Loss of Resolution
- Irreversible Column Damage
- System Damage

## Protect Your HPLC Column. Protect Your Results.

The SecurityGuard™ and SecurityGuard ULTRA cartridge systems effectively protect analytical columns from the damaging effect of contaminants that could impact results and data quality. Either cartridge system is designed to trap contaminants without altering your chromatography.

How It Works:



SecurityGuard and SecurityGuard ULTRA standard can adjust to fit any manufacturer's female/inverted endfitting.



Additional information can be found at  
[www.phenomenex.com/securityguard](http://www.phenomenex.com/securityguard)

# Phenex™ Syringe Filters

## For Sample and Solvent Filtration Prior to Chromatography

- Less system downtime
- More consistent, reproducible results
- Increased column lifetime

### Features and Benefits

- Low protein adsorption
- Broad chemical compatibility
- Minimized extractables
- Excellent flow rate
- High total throughput
- Low hold-up volume
- Certified quality
- 100 % integrity tested



#### MEMBRANE TYPES

RC (Regenerated Cellulose)

PTFE, Teflon® (Polytetrafluoroethylene)

PES (Polyethersulfone)

PVDF (Polyvinylidene Fluoride)

NY (Nylon)

CA (Cellulose Acetate)

GF (Glass Fiber)



To order and find out more about this product  
visit [www.phenomenex.com/phenex](http://www.phenomenex.com/phenex)

# HPLC Troubleshooting Mini Guide

# Baseline Issues

## BE-HAPPY™ GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.  
[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

### Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Czech Republic

t: +420 272 017 077  
cz-info@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### Hong Kong

t: +852 6012 8162  
hkinfo@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Indonesia

t: +62 21 3952 5747  
indoinfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Japan

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

### Poland

t: +48 22 51 02 180  
pl-info@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: 800-852-3944  
sginfo@phenomenex.com

### Slovakia

t: +420 272 017 077  
sk-info@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

### Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

### Thailand

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions**  
**Corporate Office USA**  
t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)



[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

No part of this booklet may be copied without prior written permission from Phenomenex, Inc. USA. While every attempt has been made to ensure the accuracy of the information contained in this guide, Phenomenex assumes no responsibility for its use. We welcome any additions or corrections for incorporation into future editions. Phenex, SecurityGuard and BE-HAPPY are trademarks of Phenomenex. Teflon is a registered trademark of E.I. du Pont de Nemours and Co. SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362. CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges. © 2024 Phenomenex, Inc. All rights reserved.