

# Ph. Eur. Monograph 2733: Escitalopram Oxalate Related Substances on Luna® 5 µm C18(2) and Luna Omega 5 µm C18

Janja Vasilić<sup>1</sup>, Amra Perva Uzunalić<sup>1</sup>, Dirk Hansen<sup>2</sup> and Heiko Behr<sup>2</sup>

<sup>1</sup>University of Maribor, Faculty of Chemistry and Chemical Engineering Group for Separation Analysis, Smetanova 17, 2000 Maribor, Slovenia

<sup>2</sup>Phenomenex Ltd. Deutschland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany

## Overview

Escitalopram Oxalate is mainly used for the treatment of major depressive disorder or generalized anxiety disorder. It is administered orally.

In this application note we are showing the effective separation of Escitalopram Oxalate from its impurities and related substances according to Ph. Eur. Monograph 2733. The fully porous Luna 5 µm C18(2) and Luna Omega 5 µm C18 columns met the system suitability requirement for the peak-to-valley ratio between Escitalopram and impurity D ( $\geq 5.0$ ). Both columns provided a baseline separation between Escitalopram and impurity D. Finally, we are showing the consistency in analytical performance of Luna Omega 5 µm C18 across columns from multiple media batches.

All reference solutions were prepared as indicated in Ph.Eur. monograph 2733 for Escitalopram Oxalate. The following certified reference standards (CRS) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) – Council of Europe; Postal address: Allee Kastner CS 30026 F - 67081 Strasbourg (France):

- Y0001796 Escitalopram Oxalate CRS
- Y0001797 Escitalopram for system suitability CRS

## Experimental

**Table 1:** Preparation of buffer solution, mobile phase solutions, test and reference solutions

Solution	Composition
Buffer solution	Dissolve 3.4 g $\text{KH}_2\text{PO}_4$ in 900 mL water for chromatography. Adjust with phosphoric acid to pH 3.0. Dilute to 1000 mL with water for chromatography.
Mobile phase A	Acetonitrile:Buffer solution 10:90 (V/V)
Mobile phase B	Buffer solution:Acetonitrile 35:65 (V/V)
Test solution	Dissolve 25 mg of the substance to be examined (for this study: Escitalopram oxalate CRS) in mobile phase A and dilute to 50.0 mL with mobile phase A.
Reference solution (a)	Dilute 5 mg of Escitalopram for system suitability CRS (containing impurity D) in mobile phase A and dilute to 10.0 mL with mobile phase A.
Reference solution (b)	Dilute 1.0 mL of the test solution to 100.0 mL with mobile phase A. Dilute 1.0 mL of this solution to 10.0 mL with mobile phase A.

## LC-UV Conditions

**Columns:** Luna 5 µm C18(2) ([00G-4252-E0](#))

Luna Omega 5 µm C18 ([00G-4785-E0](#))

**Dimension:** 250 x 4.6 mm

**Mobile Phase:** A = Acetonitrile:Buffer solution 10:90 (V/V)

B = Buffer solution:Acetonitrile 35:65 (V/V)

Gradient:	Time (min)	%B	Flow Rate (mL/min)
	0	5	1.0
	2	5	1.0
	37	35	1.0
	47	100	1.0
	47.1	100	2.0
	62	100	2.0
	62.1	5	1.0
	70	5	1.0

**Injection:** 20 µL


**Temperature:** 45 °C

**Detector:** UV @ 237 nm and 254 nm

**System:** Agilent® 1260 Infinity II




## Chromatographic Results

**Figure 1:** Reference solution (a) on Luna® 5 µm C18(2)**Table 2:** Results reference solution (a) on Luna 5 µm C18(2) at 237 nm and 254 nm


No.	Escitalopram Rt (min) at 237 nm	Peak area Escitalopram at 237 nm	Resolution factor (Escitalopram/Impurity D) at 237 nm	Hp/Hv at 237 nm	Impurity D Rt (min) at 254 nm	Impurity D RRt at 254 nm	Peak area impurity D at 254 nm
1	37.917	21090773	1.063	-*	37.331	0.985	25176
2	37.912	21120286	1.057	-*	37.327	0.985	25314
3	37.915	21135675	1.059	-*	37.331	0.985	25119
4	37.925	21132248	1.058	-*	37.341	0.985	24968
5	37.918	21116763	1.060	-*	37.332	0.985	25127
6	37.918	21068527	1.063	-*	37.332	0.985	25143
Average	37.917	21110712	1.06	-*	37.332	0.985	25141
% RSD	0.011	0.12	0.24	-*	0.012	0.003	0.44

\*No Hp/Hv ratio could be calculated because of baseline separation between Escitalopram and impurity D.


**Figure 2:** Reference solution (b) on Luna 5 µm C18(2)**Table 3:** Results reference solution (b) on Luna 5 µm C18(2) at 237 nm and 254 nm

No.	Rt (min) at 237 nm	Area of Escitalopram peak at 237 nm	Tailing factor (T) at 237 nm	Rt (min) at 254 nm	Area of Escitalopram peak at 254nm	Tailing factor (T) at 254 nm
1	38.677	20348	1.064	38.681	4778	1.000
2	38.667	20219	1.053	38.669	4912	1.090
3	38.653	20285	1.033	38.654	4951	1.030
4	38.643	20425	1.040	38.642	4701	0.973
5	38.628	20253	1.041	38.629	5020	1.113
6	38.611	20409	1.069	38.601	4863	1.082
Average	38.646	20323	1.050	38.646	4871	1.048
% RSD	0.063	0.42	1.39	0.075	2.4	5.3



**Figure 3:** Test solution on Luna® 5 µm C18(2)**Table 4:** Results test solution on Luna 5 µm C18(2) at 237 nm and 254 nm


No.	Escitalopram Rt (min) at 237 nm	Peak area Escitalopram at 237 nm	Unknown Rt (min) at 254 nm	Unknown RRt at 254 nm	Peak area unknown at 254 nm
1	37.928	20558317	34.957	0.922	25174
2	37.935	20551239	34.959	0.922	25671
3	37.925	20552239	34.950	0.922	25841
4	37.924	20550636	34.950	0.922	25997
5	37.907	20549039	34.937	0.922	26402
6	37.906	20546436	34.929	0.921	26634
Average	37.921	20551318	34.947	0.922	25953
% RSD	0.031	0.019	0.034	0.008	2.015

**Figure 4:** Reference solution (a) on Luna Omega 5 µm C18**Table 5:** Results Reference solution (a) on Luna Omega 5 µm C18 at 237 nm and 254 nm


No.	Escitalopram Rt (min) at 237 nm	Peak area Escitalopram at 237 nm	Resolution factor (Escitalopram/Impurity D) at 237 nm	Hp/Hv at 237 nm	Impurity D Rt (min) at 254 nm	Impurity D RRt at 254 nm	Peak area impurity D at 254 nm
1	39.103	20949728	0.955	-*	38.548	0.986	25058
2	39.092	20940230	0.954	-*	38.538	0.986	25220
3	39.076	20963085	0.952	-*	38.522	0.986	25222
4	39.100	20977182	0.953	-*	38.546	0.986	25077
5	39.102	20959523	0.953	-*	38.548	0.986	25039
6	39.133	20953808	0.955	-*	38.579	0.986	24995
Average	39.101	20957259	0.954	-*	38.547	0.986	25102
% RSD	0.048	0.060	0.128	-*	0.048	0.001	0.382

\*No Hp/Hv ratio could be calculated because of baseline separation between Escitalopram and impurity D.




**Figure 5:** Reference solution (b) on Luna® Omega 5 µm C18**Table 6:** Results reference solution (b) on Luna Omega 5 µm C18 at 237 nm and 254 nm

No.	Rt (min) at 237 nm	Area of Escitalopram peak at 237 nm	Tailing factor (T) at 237 nm	Rt (min) at 254 nm	Area of Escitalopram peak at 254nm	Tailing factor (T) at 254 nm
1	40.029	20360	1.025	40.025	4886	1.075
2	40.028	20331	1.024	40.031	4899	1.057
3	40.029	20346	1.009	40.030	4798	1.040
4	40.023	20298	1.023	40.026	4916	1.037
5	40.019	20455	1.024	40.023	4918	0.935
6	40.010	20261	1.038	40.005	4797	1.046
Average	40.023	20342	1.024	40.023	4869	1.032
% RSD	0.019	0.32	0.89	0.023	1.2	4.8

**Figure 6:** Test solution on Luna Omega 5 µm C18**Table 7:** Results Test solution on Luna Omega 5 µm C18 at 237 nm and 254 nm

No.	Escitalopram Rt (min) at 237 nm	Peak area Escitalopram at 237 nm	Unknown Rt (min) at 254 nm	Unknown RRT at 254 nm	Peak area unknown at 254 nm
1	39.434	20770725	36.539	0.927	33943
2	39.419	20745765	36.523	0.927	33187
3	39.391	20739482	36.496	0.927	31757
4	39.388	20748513	36.491	0.926	32421
5	39.376	20718286	36.483	0.927	32994
6	39.399	20681838	36.499	0.926	32724
Average	39.401	20734101	36.505	0.927	32838
% RSD	0.054	0.15	0.058	0.008	2.2



**Figure 7:** Overlap Reference solution (a) on 3 Luna® Omega 5 µm C18 columns from 3 different batches at 237 nm**Table 8:** Results reference solution (a) on 3 Luna Omega 5 µm C18 columns from 3 different batches at 237 nm

No.	Batch 1			Batch 2			Batch 3		
	Escitalopram Rt (min)	Impurity D Rt (min)	Resolution	Escitalopram Rt (min)	Impurity D Rt (min)	Resolution	Escitalopram Rt (min)	Impurity D Rt (min)	Resolution
1	39.314	38.763	0.956	39.366	38.814	0.953	39.103	38.548	0.955
2	39.309	38.759	0.956	39.310	37.758	0.953	39.092	38.538	0.954
3	39.309	38.757	0.959	39.321	38.769	0.952	39.076	38.522	0.952
4	39.304	38.753	0.958	39.320	38.767	0.956	39.100	38.546	0.953
5	39.298	38.747	0.957	39.305	38.755	0.953	39.102	38.548	0.953
6	39.267	38.716	0.958	39.308	38.757	0.953	39.133	38.579	0.955
Average	39.300	38.749	0.956	39.322	38.770	0.953	39.101	38.547	0.954
%RSD	0.044	0.044	0.11	0.057	0.057	0.13	0.048	0.048	0.13

## Conclusion

The Luna 5 µm C18(2) (**Figure 1** and **Table 2**) and Luna Omega 5 µm C18 (**Figure 4** and **Table 5**) columns tested in this study fulfilled the system suitability requirements of Ph. Eur. Monograph 2733. In addition, the high batch-to-batch reproducibility of the Luna Omega 5 µm C18 column has been verified by the comparison of separation results for 3 columns from 3 different batches (**Figure 7** and **Table 8**).



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/LiveChat](http://www.phenomenex.com/LiveChat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 5010 9707  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 104 21 72  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: +65 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**☎ All other countries/regions****Corporate Office USA**

t: +1 (310) 212-0555  
info@phenomenex.com

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Luna is registered trademark and BE-HAPPY is a trademark of Phenomenex. Agilent is a registered trademark of Agilent Technologies.

**Disclaimer**

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Agilent Technologies.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

