


Separation of Solifenacin Succinate Stereoisomers per USP Monograph using Lux® 5 µm Amylose-1 Column

Lauren Nakasone, Zeshan Aqeel, and Dr. Bryan Tackett
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501-1430, USA



Overview

Solifenacin Succinate is a chiral drug which is commercialized as pure stereoisomer. It is an anticholinergic drug used for the treatment of patients with overactive bladder. The current trend is to commercialize drugs formulated in their pure stereoisomeric form, which means containing only the single stereoisomer responsible for pharmaceutical activity. The (R,S)-stereoisomer of Solifenacin Succinate is the active form and therefore isolating it from the other stereoisomers is of key importance.

In this application note, we report the stereoisomeric separation between Solifenacin Succinate and its stereoisomers using the Lux Amylose-1 chiral column according to the USP monograph for Solifenacin Succinate. The USP Monograph suggests an ambient running temperature, but we used a lower temperature of 15 °C to help with peak shape, staying within the allowable adjustments per USP <621>.

System suitability per USP Monograph for Solifenacin Succinate is for resolution of not less than (NLT) 1.5 between the Solifenacin (R,S)-stereoisomer and the (S,S)-stereoisomer peaks, and NLT 2.0 between the (R,R)-stereoisomer and Solifenacin Succinate peaks. The separation of the stereoisomers was achieved and peaks corresponding to the (R,S)- and (S,S)-stereoisomers had a resolution of 7.92, and the peaks corresponding to the (R,R)-stereoisomer and Solifenacin Succinate had a resolution of 13 (**Figure 1**). These far surpassed the required resolution to meet system suitability. The higher resolution indicates the powerful chiral recognition ability of the Lux 5 µm Amylose-1 chiral column.

The percent relative standard deviation (%RSD) of the six replicate injections of the Standard Solution must be no more than (NMT) 10.0 % to meet system suitability requirements. All parameters had a %RSD less than 10.0 %, meeting system suitability requirements (**Figure 2**).

All reference solutions were prepared as indicated in the USP Monograph for Solifenacin Succinate. Solifenacin Identification Mixture B RS (Catalog No. 1615322), and Solifenacin Succinate RS (Catalog no. 1615300) were purchased from USP.

LC-UV Conditions

Column: Lux 5 µm Amylose-1 (250 x 4.6 mm)

Part No.: [00G-4732-E0](#)

Mobile Phase: Chromatographic n-Heptane:Absolute Alcohol:Diethylamine (800:200:1, v/v/v)

Pressure (bar): 76

Flow Rate: 0.8 mL/min

Injection: 10 µL

Temperature: 15 °C


Detector: UV @ 220 nm

System: Agilent® 1260 Binary UHPLC

Table 1. Preparation of Solutions


Solution	Composition
System Suitability Solution	Dissolve 2 mg of USP Solifenacin Identification Mixture B RS in 1 mL mobile phase
Standard Solution	Dissolve 1 mg of USP Solifenacin Succinate RS in 1 mL of mobile phase and making a 1:500 dilution of the 1 mg/mL solution to produce a 0.002 mg/mL solution



Figure 1. System Suitability Solution

Peak	Analyte	Resolution
1	(R,S)-Stereoisomer	7.92
2	(S,S)-Stereoisomer	
3	(R,R)-Stereoisomer	13.00
4	Solifenacin Succinate	

App ID 26770

Figure 2. Standard Solution

App ID 26771

	Time	Area	Height	Width
AVG	22.655	21.450	0.330	1.093
STD	0.188	1.912	0.029	0.100
%RSD	0.830	8.913	8.783	9.141

Number of injections = 6



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Lux is a registered trademark and BE-HAPPY is a trademark of Phenomenex.

Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help!
Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists