


Performance of Surface Modified Luna® Omega and Kinetex® PS Columns for the Separation of Tricyclic Antidepressants

Zeshan Aqeel, Dr. Jason Anspach, Dr. Bryan Tackett, and Anania Assefa Amenu
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501-1430, USA



Overview

Tricyclic Antidepressants (TCAs) are a class of medications that are used primarily to treat depression. TCAs generally have a three-ring structure, with an attached secondary or tertiary amine, and they are classified as weak bases ($pK_a \sim 8.5$).

Many C18 columns show strong ionic interactions of basic amines with residual unbonded silanol groups under low pH mobile phase conditions. However, the Luna Omega 1.6 μm PS C18 and the Kinetex 2.6 μm PS C18 columns are designed with a positively charged functional group on the surface. This reduces or eliminates the ionic interactions of aromatic amines and other basic compounds. This translates into a dramatic improvement in peak shape for basic analytes, even when using a low ionic strength 0.1% formic acid buffer system.

In this work, we compared the performance of six columns in the separation of three TCA analytes run under the same HPLC conditions. Compared to the traditionally bonded C18 columns with no positive ligand on the surface, the Luna and Kinetex PS C18 columns showed superior separation producing excellent peak shape and overall resolution (**Figure 1, (b) and (c)**). Both Luna Omega and Kinetex PS columns provided narrow peak widths with minimal peak tailing compared to the results obtained using competitor C18 columns with similar surface modification technology (**Figure 1, (e) and (f)**).

Having the PS chemistry on both high performance fully porous and core-shell silica allows for the optimization of high column efficiency and selectivity.

LC-UV Conditions

Column: Luna Omega 1.6 μm PS C18 ([00B-4752-AN](#))
Luna Omega 1.6 μm C18 ([00B-4742-AN](#))
Kinetex 2.6 μm PS C18 ([00B-4780-AN](#))
Kinetex 2.6 μm C18 ([00B-4462-AN](#))
Waters® CORTECS™ 2.7 μm C18+
Waters ACQUITY™ CSH 1.7 μm C18

Dimensions: 50 x 2.1 mm

Mobile Phase: A = 0.1% Formic Acid in Water
B = 0.1% Formic Acid in Acetonitrile

Pressure (bar): 420 (Luna Omega 1.6 μm PS C18)
480 (Luna Omega 1.6 μm C18)
270 (Kinetex 2.6 μm PS C18)
251 (Kinetex 2.6 μm C18)
207 (Waters Cortecs 2.7 μm C18+)
450 (Waters ACQUITY CSH 1.7 μm C18)

Flow Rate: 0.5 mL/min


Injection: 1 μL

Temperature: 30 °C

Detector: UV @ 254 nm

System: Agilent® 1260 Binary UHPLC

Kinetex PS C18



Luna Omega PS C18







Figure 1. Results and Observations**(a). Luna® Omega 1.6 µm C18**

Peak #	Analyte
1	Desipramine
2	Imipramine
3	Amitriptyline

(b). Luna Omega 1.6 µm PS C18**(c). Kinetex® 2.6 µm PS C18****(d). Kinetex 2.6 µm C18****(e). Waters® CORTECS™ 2.7 µm C18+****(f). Waters ACQUITY™ CSH 1.7 µm C18**

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and Kinetex are registered trademarks and BE-HAPPY is a trademark of Phenomenex.

Agilent is a registered trademark of Agilent Technologies, Inc.

Waters is a registered trademark and CORTECS and ACQUITY are trademarks of Waters Technologies Corporation.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Agilent Technologies, Inc. or Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.

