Pphenomenex

Ph. Eur. Monograph 2232: Losartan Potassium Related Substances on Luna® 5 μm C18(2)

Laura Nakasone¹, Heiko Behr², and Dr. Bryan Tackett¹

¹Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

²Phenomenex Ltd. Deutschland, Zeppelinstr. 5, 63741 Aschaffenburg, Germany



Overview

Losartan potassium belongs to a class of drugs called angiotensin receptor blockers (ARBs). It is an oral medication beneficial for treating hypertension. There is an increase in interest recently in the effectiveness of Losartan potassium in lowering high blood pressure.

In this application note we demonstrate the effective identification and separation of Losartan potassium from its impurities and related substances according to Ph.Eur. Monograph 2232. To meet system suitability requirements, the peak-to-valley ratio must be a minimum of 2 between Impurity A and Impurity G for reference solution (b). The peak-to-valley ratio is defined as H_p/H_v , where H_p = height above the baseline of the peak due to impurity M and H_v = height above the baseline of the lowest point of the curve separating this peak from the peak due to impurity G.

We used the fully porous Luna 5 μ m C18(2) column and compared it to the Kromasil 5 μ m C18 column originally used in the monograph. Both columns met system suitability and all peaks due to impurities J, K, L, M, and G were identified (**Figure 4**), with the Luna C18(2) column providing a larger peak-to-valley ratio, indicative of increased separation.

All reference solutions were prepared as indicated in Ph. Eur. Monograph 2232 for Losartan potassium. Losartan potassium Certified Reference Standard (CRS) (catalog no. Y0001062), Losartan for System Suitability CRS (catalog no. Y0001076), Losartan Impurity D CRS (catalog no. Y0001072), and Triphenylmethanol (Impurity G) (catalog no. Y0001654) were purchased from the European Directorate for the Quality of Medicines & HealthCare (EDQM) Council of Europe.

LC-UV Conditions

Column: Kromasil® 5 μm C18, 250 x 4.6 mm

Luna 5 µm C18(2), 250 x 4.6 mm (00G-4252-E0)

Pressure: 239 bar (Kromasil) 237 bar (Luna)

Mobile Phase: A: Dilute 1.0 mL of phosphoric acid to 1000 mL with

water

B: Acetonitrile

Gradient: Tim	e (min)	%В
	0	25
	5	25
	30	90
	40	90

Flow Rate: 1.3 mL/min Injection: 10 μL Temperature: 35 °C

Detector: UV @ 220 nm


System: Agilent® 1260 Binary UHPLC


Figure 1. Structure of Losartan Potassium

Table 1. Preparation of Test and Reference Solutions

Solution	Composition
Test Solution	Dissolve 30.0 mg of the Losartan Potassium CRS in methanol and dilute to 100.0 mL with the same solvent.
Reference Solution (a)	Dilute 1.0 mL of the Test Solution to 100.0 mL with methanol. Dilute 1.0 mL of this solution to 10.0 mL with methanol.
Reference Solution (b)	Dissolve 6.0 mg of triphenylmethanol (Impurity G) in 100.0 mL of methanol. Dilute 1.0 mL of the solution to 100.0 mL with methanol. Use 1.0 mL of this solution to dissolve the contents of a vial of Losartan for System Suitability CRS (containing impurities J, K, L and M) and sonicate for 5 min.
Reference Solution (c)	Dissolve 3.0 mg of Losartan Impurity D CRS in methanol and dilute to 100.0 mL with the same solvent. Dilute 1.5 mL of this solution to 100.0 mL with methanol.

Figure 2. Chromatograms for Test Solution







Peak No.	Analyte	Retention Time (min)	Area	Height	Symmetry Factor
1	Losartan Potassium	13.528	2664.1	448.3	1.114

Peak No.	Analyte	Retention Time (min)	Area	Height	Symmetry Factor
1	Losartan Potassium	13.578	3346.4	566.4	1.158

Figure 3. Chromatograms for Reference Solution (a)




Peak No.	Analyte	Retention Time (min)	Area	Height	Symmetry Factor	
1	Losartan Potassium	13.554	2.3	3.7e-1	1.067	




Peak No.	Analyte	Retention Time (min)	Area	Height	Symmetry Factor
1	Losartan Potassium	13.564	2.5	4.4e-1	1.038

Figure 4. Chromatograms for Reference Solution (b)




Peak No.	Analyte	Retention Time (min)	RRT to Losartan	Area	Height	Symmetry Factor	Peak-to- Valley Ratio
1	Impurity J	18.90	1.40	3.8	6.3e-1	1.04	-
2	Impurity K	20.48	1.51	7.8	1.3	1.04	-
3	Impurity L	21.88	1.62	3.7	6.0e-1	1.04	-
4	Impurity M	24.51	1.81	3.7	5.7e-1	0.97	37.78
5	Impurity G	25.04	1.85	7	1.1	1.06	-




Peak No.	Analyte	Retention Time (min)	RRT to Losartan	Area	Height	Symmetry Factor	Peak-to- Valley Ratio
1	Impurity J	19.09	1.41	3.9	6.5e-1	1.0	-
2	Impurity K	20.64	1.52	8.0	1.3	1.06	-
3	Impurity L	21.98	1.62	3.8	5.9e-1	0.96	-
4	Impurity M	24.54	1.81	3.7	5.5e-1	1.02	55.00
5	Impurity	25.07	1.85	718.2	2.8	1.04	-

Figure 5. Chromatograms for Reference Solution (c)



Peak No.	Analyte	Retention Time (min)	RRT to Losartan	Area	Height	Symmetry Factor
1	Impurity	11.84	0.88	9.0e-1	1.3e-1	1.30



Peak No.	Analyte	Retention Time (min)	RRT to Losartan	Area	Height	Symmetry Factor
1	Impurity	11.84	0.88	9 Oe-1	1 3e-1	1 30

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belaium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681 info@phenomenex.com

t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077 cz-info@phenomenex.com

Denmark

t: +45 4824 8048 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162 hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Indonesia

t: +62 21 5010 9707 indoinfo@phenomenex.com

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262 jpinfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951 nzinfo@phenomenex.com

Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72 pl-info@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944 sginfo@phenomenex.com

Slovakia t: +420 272 017 077 sk-info@phenomenex.com

t: +34 91-413-8613 espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367 ukinfo@phenomenex.com

t: +1 (310) 212-0555 info@phenomenex.com

All other countries/regions Corporate Office USA

t: +1 (310) 212-0555 info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions

Trademarks

Luna is a registered trademark and BE-HAPPY is a trademark of Phenomenex. Kromasil is a registered trademark of Nouryon, Separation Products. Agilent is a registered trademark of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Nourvon, Separation Products or Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2021 Phenomenex, Inc. All rights reserved.

