

Separation of Zolmitriptan and its Organic Impurities per USP Monograph

Trivikram Reddy Gundala¹, Krishna Chaitanya Routhu¹, Swetha Kotikalaupdi¹, and Bryan Tackett, PhD²

¹India PhenoLogix, Mahadeva Kodigehalli, Hobli, Jala Taluka, Bengaluru 562149, India

²Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Zolmitriptan is a selective serotonin receptor agonist of the 1B and 1D subtypes and its action on 5-HT 1B/1D receptors causes vasoconstriction in intracranial blood vessels. The development of a quick and efficient analysis of Zolmitriptan and its related organic impurities is of interest for generic drug manufacturers. In this application note, we report on the separation of Zolmitriptan and its related organic impurities using three batches of Luna Omega 5 μ m C18 columns according to the USP monograph for Zolmitriptan.

The use of the 150 mm length Luna Omega 5 μ m C18 column is an allowable adjustment per USP <621> since the L/dp ratio ($150/5 = 30,000$) is within the allowable range of -25 to $+50$ % of the L/dp ratio ($125/5 = 25,000$) for the original 125 mm length, 5 μ m column used to elucidate the assay method.

System suitability per USP Monograph for the Zolmitriptan Assay requires resolution no less than (NLT) 5.0 between Zolmitriptan and Zolmitriptan Related Compound E, a percent relative standard deviation (%RSD) of no more than (NMT) 0.73 % for Zolmitriptan, and a symmetry factor NMT 2.0 for Zolmitriptan. All system suitability requirements for Zolmitriptan Assay were met with all three batches of Luna Omega C18 columns (Figure 2).

System suitability per USP Monograph for the Zolmitriptan Organic Impurities requires resolution NLT 5.0 between Zolmitriptan and Zolmitriptan Related Compound E and a symmetry factor NMT 3.0 for Zolmitriptan. These requirements for system suitability for Organic Impurities were met with all three batches of Luna Omega C18 columns (Figure 3).

Across the three batches of Luna Omega C18 columns tested, the retention time reproducibility for Assay was 1.33 % and 1.79 % for Zolmitriptan and Related Compound E, respectively (Table 2).

All solutions were prepared as indicated in the USP Monograph for Zolmitriptan. USP Zolmitriptan RS (Catalog No. 1727009) and USP Zolmitriptan Related Compound E RS (Catalog No. 1727064) were purchased from USP.

LC-UV Conditions

Column: Luna™ Omega 5 μ m C18

Dimensions: 150 X 4.6 mm

Part No.: [00F-4785-EQ](#)

Mobile Phase: **Mobile Phase** (Table 1)

Flow Rate: 1.3 mL/min

Injection Volume: 20 μ L

Temperature: Ambient

Detector: UV @ 225 nm

System: Waters® Arc HPLC

Table 1. Preparation of Solutions

Solution	Composition
Mobile Phase	Acetonitrile / Water / Trifluoroacetic Acid / Triethylamine (135:865:1:0.25, v/v/v/v)
Diluent	Mobile Phase
Standard Solution – Assay	0.025 mg/mL of USP Zolmitriptan RS in Diluent
System Suitability Solution – Assay	0.12 μ g/mL of USP Zolmitriptan Related Compound E RS and 25 μ g/mL of USP Zolmitriptan RS in Diluent
System Suitability Solution – Organic Impurities	0.5 μ g/mL of USP Zolmitriptan Related Compound E RS and 0.1 mg/mL of USP Zolmitriptan RS in Diluent

Figure 1. Zolmitriptan Structure










Figure 2. System Suitability Solution – Assay

Figure 3. System Suitability Solution – Organic Impurities

Peak No.	Analyte	Retention Time (min)	Resolution	Symmetry Factor
1	Zolmitriptan	7.08	11.0	1.48
2	Zolmitriptan Related Compound E	11.03		-



Peak No.	Analyte	Retention Time (min)	Resolution	Symmetry Factor
1	Zolmitriptan	7.01	9.86	1.61
2	Zolmitriptan Related Compound E	10.87		-
N = 6 Injections				



Peak No.	Analyte	Retention Time (min)	Resolution	Symmetry Factor
1	Zolmitriptan	6.90	10.64	1.41
2	Zolmitriptan Related Compound E	10.68		-
N = 6 Injections				



Table 2. Batch-to-Batch Retention Time Reproducibility of Luna™ Omega 5 µm C18 Columns for Assay

Batch No.	Zolmitriptan	Zolmitriptan Related Compound E
1	7.22	11.14
2	7.12	10.91
3	7.03	10.75
Mean	7.12	10.93
STD	0.10	0.20
%RSD	1.33	1.79



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Luna and BE-HAPPY are trademarks of Phenomenex. Waters is a registered trademarks of Waters Technologies Corporation.

Disclaimer

Comparative separations may not be representative of all applications
Phenomenex is in no way affiliated with Waters Technologies Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.

