




Advanced Oligonucleotide Analysis

Increased Recovery and Reproducibility
with **NEW bioZen Oligo**

- BioTi™ Hardware Reduces Sample Loss and Adsorption
- Robustness at High pH and Temperature
- Core-Shell Advantage for High Efficiency



www.phenomenex.com/bioZenOligo

 **phenomenex®**

Inside the bioZen Oligo Biocompatible Hardware Difference




The use of bio-inert hardware not only improves the chromatographic performance and consistency of oligonucleotides, but also provides improvements in sensitivity, enabling both quantitation and characterization.




BioTi LC Hardware Minimizes Adsorptive Interactions and Sample Loss

Biocompatible Column Hardware





Stainless Steel Column Hardware



BioTi versus Traditional Stainless Steel Hardware

Oligos can chelate to trace heavy metals in stainless steel column hardware, leading to poor recovery, inconsistent chromatography and problematic carryover. The bioZen Oligo bio-inert hardware provides greater sensitivity as well as improved recovery, demonstrating this column's optimal utility for oligonucleotide characterization and quantitation.




BioTi ensures method robustness and consistency from injection to injection!

Patented Technology and Advanced Core-Shell Particle Chemistry


The bioZen™ Oligo Columns utilize organo-silica core-shell particles with a highly consistent morphology that minimizes band broadening associated with diffusion and mass transfer, leading to higher efficiency and minimal peak widths, which is critical for the separation of closely eluting impurities associated to synthetic oligonucleotides.

Patented Core-Shell Particle Chemistry




High Efficiency Core-Shell Particle

After meticulous core construction, a uniform porous silica layer is grown around the spherical solid silica core. This unique combination of precise particle architecture and particle size provides dramatic leaps in performance.



Better Performance than Fully Porous Particles

Core-Shell Technology provides extremely high efficiencies for both low and high oligo concentrations. Industry leading column packing technology in combination with high particle consistency and density helps create optimal bed structure which reduces band broadening effects of Eddy Diffusion to produce highly reproducible columns that generate greater performance compared to fully porous particles. This ultra-high efficiency can be leveraged to achieve increased resolution, improved sensitivity, and higher productivity.



bioZen Oligo Applications

The newest addition to the bioZen™ portfolio brings a unique combination of core-shell versatility and high pH ruggedness necessary for oligonucleotide separations. Additionally, bioZen Oligo is packed in a unique bio-inert titanium hardware designed to minimize the sample loss and adsorption issues typically seen with stainless steel hardware.

Explore the following bioZen Oligo Resources:



Column Hardware Impact on Oligonucleotides

[The Effect of Column Hardware on the Analysis of Synthetic Oligonucleotides by LC-MS](#)



Chromatographic Comparison using bioZen Oligo

[Comparison of Fully Porous and Superficially Porous Particles for the Analysis of Synthetic Oligonucleotides](#)

[Comparison of Polymer and Silica Based Particles for the Analysis of Synthetic Oligonucleotides](#)



Method Optimization for Oligonucleotides

[Effect of Temperature on Single Stranded Oligonucleotide Analysis](#)

[Effect of Perfluorinated Alcohols in Impurity Analysis for Synthetic Oligonucleotides](#)

[Effect of Alkylamine Based Ion-pairing Reagents in Analysis of siRNA](#)

[Impurity Analysis of a 22mer DNA Phosphorothioate](#)

Interested in learning more about these applications?

Chat with our Oligo Experts



www.Phenomenex.com/OligoChat

Introducing bioZen™ Oligo


A New Solution for Oligonucleotide Analysis and Characterization

With a single innovative product line, bioZen separation products provide enhanced characterization over an incredibly wide range of techniques.

Oligonucleotide

NEW

bioZen Oligo
1.7 μ m and 2.6 μ m
Organo-silica core-shell bonded particle with a pH 1-12 stable C18 stationary phase that delivers speed, robust methods, and improved peak shape for Oligo analysis.



9 Particle Chemistries

Intact/Subunit	Size Exclusion (SEC)	Peptide	Glycan
<p>bioZen Intact XB-C8 3.6μm Large pore core-shell particle for fast intact and subunit biologic entry. C8 provides highly useful moderate hydrophobic selectivity.</p>	<p>bioZen SEC-2 1.8μm Extremely inert, high density fully porous particle with high efficiency and low molecular weight (LMW) separation range of 1k–450 kDa.</p>	<p>bioZen Peptide PS-C18 1.6μm and 3μm Excellent retention by combined positively charged surface ligand and C18 ligand.</p>	<p>bioZen Glycan 2.6μm Provides optimal combination of high efficiency and selectivity for released glycans.</p>
<p>bioZen WidePore C4 2.6μm Core-shell particle with butyl stationary phase and optimal wide pore size distribution for better resolution of large biologics, including monoclonal antibodies and subunit analysis.</p>	<p>bioZen SEC-3 1.8μm Extremely inert, high density fully porous particle with high efficiency and high molecular weight (HMW) separation range of 10k–700 kDa.</p>	<p>bioZen Peptide XB-C18 1.7μm and 2.6μm Overall retention of both acidic and basic peptides through C18 stationary phase with di-isobutyl side chains.</p>	<p>bioZen WCX 6μm Monosized particles grafted with linear polycarboxylate chains to envelop and separate proteins from acidic/basic variants</p>
Ion-Exchange			

Ordering Information

bioZen™ Products - Powered by Biocompatible Hardware

bioZen Columns (mm)								Biocompatible Guard Cartridges		
	50 x 2.1	100 x 2.1	150 x 2.1	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	for 2.1 mm	for 4.6 mm	Holder
bioZen 2.6 µm Glycan	00B-4773-AN	00D-4773-AN	00F-4773-AN	—	—	—	—	AJ0-9800	—	AJ0-9000
bioZen 1.6 µm Peptide PS-C18	00B-4770-AN	00D-4770-AN	00F-4770-AN	—	—	—	—	AJ0-9803	—	AJ0-9000
bioZen 3 µm Peptide PS-C18	00B-4771-AN	—	00F-4771-AN	00B-4771-E0	—	00F-4771-E0	—	AJ0-7605	AJ0-7606	KJ0-4282
bioZen 1.7 µm Peptide XB-C18	00B-4774-AN	00D-4774-AN	00F-4774-AN	—	—	—	—	AJ0-9806	—	AJ0-9000
bioZen 2.6 µm Peptide XB-C18	00B-4768-AN	00D-4768-AN	00F-4768-AN	00B-4768-E0	—	00F-4768-E0	—	AJ0-9806	AJ0-9808	AJ0-9000
bioZen 2.6 µm WidePore C4	00B-4786-AN	00D-4786-AN	00F-4786-AN	00B-4786-E0	00D-4786-E0	00F-4786-E0	00G-4786-E0	AJ0-9816	AJ0-9818	AJ0-9000
bioZen 3.6 µm Intact XB-C8	00B-4766-AN	00D-4766-AN	00F-4766-AN	00B-4766-E0	—	00F-4766-E0	—	AJ0-9812	AJ0-9814	AJ0-9000
50 x 2.1 100 x 2.1 150 x 2.1 250 x 2.1 50 x 4.6 100 x 4.6 150 x 4.6 250 x 4.6 300 x 4.6								for 4.6 mm	Holder	
bioZen 1.8 µm SEC-2	00B-4769-AN	—	00F-4769-AN	—	—	—	00F-4769-E0	—	00H-4769-E0	AJ0-9850
bioZen 1.8 µm SEC-3	00B-4772-AN	—	00F-4772-AN	—	—	00D-4772-E0	00F-4772-E0	—	00H-4772-E0	AJ0-9851
50 x 2.1 100 x 2.1 150 x 2.1 250 x 2.1 50 x 4.6 100 x 4.6 150 x 4.6 250 x 4.6 300 x 4.6								for 4.6 mm	Holder	
bioZen 6 µm WCX	00B-4777-AN	00D-4777-AN	00F-4777-AN	00G-4777-AN	00B-4777-E0	00D-4777-E0	00F-4777-E0	00G-4777-E0	—	AJ0-9400
50 x 2.1 100 x 2.1 150 x 2.1 250 x 2.1 50 x 4.6 100 x 4.6 150 x 4.6 250 x 4.6								for 2.1 mm	for 4.6 mm	Holder
bioZen 2.6 µm Oligo	00B-4790-AN	00D-4790-AN	00F-4790-AN	00B-4790-E0	00D-4790-E0	00F-4790-E0	—	AJ0-9820	AJ0-9822	AJ0-9000
bioZen 1.7 µm Oligo	00B-4791-AN	00D-4791-AN	00F-4791-AN	—	—	—	—	AJ0-9820	AJ0-9822	AJ0-9000


Sample Preparation

bioZen Solid Phase Extraction	Format	Sorbent Mass	Part Number	Unit
bioZen N-Glycan Clean-Up	Microelution 96-Well Plate	5 mg/well	8M-S009-NGA	1/box



bioZen MagBeads Streptavidin Coated

Formats	Part No.	Concentration	Bead Size
25 mg (~50 samples)	K50-9531		
50 mg (~100 samples)	K50-9532		
500 mg (~1000 samples)	K50-9533	20 mg/mL	1.0 µm



BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Have Questions or
Need a Quote?

Chat with our Oligo Experts

chat now

Quotes, Methods, Tips... We're here to help

www.Phenomenex.com/OligoChat



bioZen Oligo



Advanced Oligonucleotide Analysis

Australia
t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria
t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium
t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada
t: +1 (800) 543-3681
info@phenomenex.com

China
t: +86 400-606-8099
cniinfo@phenomenex.com

Czech Republic
t: +420 272 017 077
cz-info@phenomenex.com

Denmark
t: +45 4824 8048
nordicinfo@phenomenex.com

Finland
t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France
t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany
t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong
t: +852 6012 8162
hkinfo@phenomenex.com

India
t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland
t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy
t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg
t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico
t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands
t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand
t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway
t: +47 810 02 005
nordicinfo@phenomenex.com

Poland
t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal
t: +351 221 450 488
ptinfo@phenomenex.com

Singapore
t: +65 800-852-3944
sginfo@phenomenex.com

Slovakia
t: +420 272 017 077
sk-info@phenomenex.com

Spain
t: +34 91-413-8613
espinfo@phenomenex.com

Sweden
t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland
t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan
t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand
t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom
t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA
t: +1 (310) 212-0555
info@phenomenex.com

All other countries/regions
Corporate Office USA
t: +1 (310) 212-0555
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen, BioTi, and BE-HAPPY are trademarks of Phenomenex.

bioZen Oligo is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2021 Phenomenex, Inc. All rights reserved.