#### Allowable Adjustments to

# United States Pharmacopeia (USP) Methods

#### LC - Isocratic Elution

Component	United States Pharmacopeia (USP)		
Stationary Phase	No change of the identity of the substituent permitted.  No replacement of C18 by C8  Other physico chemical characteristics must be similar  Change from totally porous (TPP) to superficially porous particles (SPP) is allowed		
Particle Size and Column Length	Particle size (dp) and/or column length (L) can be adjusted, if the L/dp ratio remains constant or in the range between -25 % and +50 %. When switching to Core-Shell particles, other combinations of L and dp can be used, if the number of theoretical plates (N) is between -25 % and +50 % of the original column.		
Column Internal Diameter	Can be adjusted as wanted.  If a smaller internal diameter (ID) or particle size is used, extra column band broadening may need to be minimized by factors such as instrument connection, detector cell volume, sampling rate and injection volume		
Flow Rate	$^{\pm}$ 50 % (after the adjustment due to changes in column ID and particle size)  • When particle size or column internal diameter is changed (e.g. from 4.6 mm ID and 5 μm at 1.0 mL/min to 3.0 mm ID and 3 μm), the flow rate is adjusted using the following equation: $F_2 = F_1 \times \frac{dc_2^2 \times dp_1}{dc_1^2 \times dp_2} = 1.0 \times \frac{3.0^2 \times 5}{4.6^2 \times 3} = 0.71 \ mL/min$		
Column Temperature	$\pm$ 10 °C (where the operating temperature is specified).		
Composition of the Mobile Phase	$\pm$ 30 % (relative) for minor components, but no component is altered by more than 10 % absolute. A minor component comprises less or equal than (100/n) %, n being the total number of components of the mobile phase.		
Mobile Phase pH	± 0.2		
Concentration of Salts in the Buffer	± 10 %		
Detector Wavelength	No adjustment permitted.		
Injection Volume	When the column dimensions are changed, the following equation may be used for adjusting the injection volume: $V_{Inj2} = V_{Inj1} \times \frac{L_2 \times dc_2^2}{L_1 \times dc_1^2}$ Even in the absence of any column dimension change, it may be varied if the system suitability criteria remain within their established acceptability limits.		

Source: European Pharmacopeia 11.0, Chapter 2.2.46. Chromatographic separation techniques, p. 86-96 and USP-NF 2022, General Chapter <621> Chromatography.



#### Allowable Adjustments to

# United States Pharmacopeia (USP) Methods

#### **LC – Gradient Elution**

Component	United States Pharmacopeia (USP)		
Stationary Phase	No change of the identity of the substituent permitted.		
Particle Size and Column Length	Particle size (dp) and/or column length (L) can be adjusted, if the L/dp ratio remains constant or in the range between -25 % and +50 %. When switching to Core-Shell particles, other combinations of L and dp can be used, if the ratio ( $t_R/W_p$ ) <sup>2</sup> is between -25 % and +50 % of the original column, for each peak used to test the system suitability.		
Column Internal Diameter	Can be adjusted as wanted.		
Flow Rate	Is adjusted if column internal diameter and particle size is changed.  • When changing the column from 4.6 mm ID with 5 µm particle size at 2.0 mL/min to 2.1 mm ID with 3 µm particle size, the flow rate is adjusted using the following equation: $F_2 = F_1 \times \frac{dc_2^2 \times dp_1}{dc_1^2 \times dp_2} = 2.0 \times \frac{2.1^2 \times 5}{4.6^2 \times 3} = 0.7 \ mL/min$		
Gradient Volume	When changing the column dimension, each gradient segment volume is adjusted by means of the gradient time using the following equation: $t_{G2} = t_{G1} \times \frac{F_1}{F_2} \times \frac{L_2 \times dc_2^2}{L_1 \times dc_1^2} = 3 \ min \times \frac{2.0}{0.7} \times \frac{100 \times 2.1^2}{150 \times 4.6^2} = 3 \ min \times 0.4 = 1.2 \ min$		
Column Temperature	$\pm$ 5 °C (where the operating temperature is specified).		
Composition of the Mobile Phase + Gradient	Adjustments of the composition of the mobile phase and the gradient are acceptable, if:  The system suitability criteria are fulfilled  The principal peak(s) elute(s) within ± 15 % of the indicated retention time(s). This requirement does not apply when the column dimensions are changed  The first peaks are sufficiently retained and the last peaks are eluted		
Mobile Phase pH	± 0.2		
Concentration of Salts in the Buffer	± 10 %		
Dwell Volume	Gradient time points (t in min) can be adapted to compensate differences in dwell volume between the system used for method development ( $D_0$ in mL) and that actually used (D in mL). The adapted time points ( $t_c$ in min) at the current flow rate (F in mL/min) can be calculated using the following equation: $t_C = t - \frac{(D-D_0)}{F} = 1 \ min - \frac{(1.0 \ mL - 0.5 \ mL)}{1 \ mL/min} = 0.5 \ min$		
Detector Wavelength	No adjustment permitted.		
Injection Volume	When the column dimensions are changed, the following equation may be used for adjusting the injection volume: $V_{Inj2} = V_{Inj1} \times \frac{L_2 \times dc_2^2}{L_1 \times dc_1^2}$ • Even in the absence of any column dimension change, injection volume may be varied if the system suitability criteria remain within their established acceptability limits • When the injection volume is decreased, special attention is given to (limit of) detection and repeatability of the peak response(s) • An increase is permitted, if the linearity and resolution of the peak(s) to be determined remain satisfactory		
	, at any case to permitted, it are integrity and recordation of the pearloy to be determined formant satisfactory		

Source: European Pharmacopeia 11.0, Chapter 2.2.46. Chromatographic separation techniques, p. 86 - 96 and USP-NF 2022, General Chapter <621> Chromatography.



	Column Classification	Recommended Phenomenex Column	Particle Shape
L1	Octadecyl silane chemically bonded to porous or non-porous silica or ceramic microparticles, 1.5 to 10 µm in diameter, or a monolithic rod.	Kinetex "C18 Kinetex EVO C18 Kinetex PS C18 Kinetex PS C18 Kinetex YB C18 Kinetex YB-C18 Luna" C18(2) Luna Omega C18 Luna Omega PS C18 Luna Omega PS C18 Luna Omega Polar C18 Gemini" NX-C18 Gemini "NX-C18 Gemini" Tysion-RP Synergi Fusion-RP Biozen" Peptide PS-C18 Onyx" C18 Jupiter" C18 Clarity" Oligo-MS Clarity Oligo-MS	Core-Shell Core-Shell Core-Shell Core-Shell Core-Shell Spherical Core-Shell Core-Shell Core-Shell Core-Shell
L2	Octadecyl silane chemically bonded to silica gel of a controlled surface porosity that has been bonded to a solid spherical core, $30\ to\ 50\ \mu m$ in diameter.		
L3	Porous silica particles, 1.5 to $10\mu m$ in diameter, or a monolithic silica rod.	Kinetex HILIC Luna Silica(2) Onyx Silica	Core-Shell Spherical Monolith
L4	Silica gel of controlled surface porosity bonded to a solid spherical core, 30 to $50\mu m$ in diameter.		
L5 L6	Alumina of controlled surface porosity bonded to a solid spherical core, 30 to 50 µm in diameter.  Strong cation-exchange packing: sulfonated fluorocarbon polymer coated on a solid spherical core, 30 to 50 µm in diameter.		
L7	Octyl silane chemically bonded to totally or superficially porous silica particles, 1.5 to $10\mu m$ in diameter, or a monolithic silica rod.	Kinetex C8 Luna C8(2) Onyx C8 Biozen Intact XB-C8	Core-Shell Spherical Monolith Core-Shell
L8	An essentially monomolecular layer of aminopropyl-silane chemically bonded to totally porous silica gel support, 1.5 to $10\mu m$ in diameter, or a monolithic silica rod.	Luna NH₂ Luna Omega SUGAR	Spherical Spherical
L9	Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to $10\mu m$ in diameter.	Luna SCX	Spherical
L10	Nitrile groups chemically bonded to porous silica particles or superficially porous particles, 1.5 to 10 $\mu m$ in diameter, or a monolithic silica rod.	Luna CN	Spherical
L11	Phenyl groups chemically bonded to porous silica particles or superficially porous particles, 1.5 to $10\mu m$ in diameter, or a monolithic silica rod.	Kinetex Biphenyl Kinetex Phenyl-Hexyl Synergi Polar-RP Luna Phenyl-Hexyl Gemini C6-Phenyl Prodigy PH-3	Core-Shell Core-Shell Spherical Spherical Spherical Spherical
L12	Strong anion-exchange packing made by chemically bonding a quaternary amine to a solid silica spherical core, 30 to $50\mu m$ in diameter.	-	·
L13	Trimethylsilane chemically bonded to porous silica particles, 3 to 10 µm in diameter.	Develosil® TMS-UG (C1) 130 Å	Spherical
L14	Silica gel having a chemically bonded, strongly basic quaternary ammonium anion-exchange coating, 5 to 10 µm in diameter.	PhenoSphere™ SAX	Spherical
L15	Hexyl silane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter.	PhenoSphere C6	Spherical
L16 L17	Dimethyl silane chemically bonded to porous silica particles, 5 to 10 µm in diameter.  Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 6 to 12 µm in diameter.	Rezex™ RHM-Monosaccharide Rezex ROA-Organic Acid	Spherical Spherical
L18 L19	Amino and cyano groups chemically bonded to porous silica particles, 3 to 10 µm in diameter.  Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the	Rezex RCM-Monosaccharide	Spherical
L13	calcium form, 5 to 15 µm in diameter.	Rezex RCU-Sugar Alcohols	Spherical
L20	Dihydroxypropane groups chemically bonded to porous silica or hybrid particles, 1.5 to $10\mu m$ in diameter, or a monolithic silica rod.	Luna HILIC BioSep™-SEC-S Biozen dSEC-2	Spherical Spherical Spherical
L21	A rigid, spherical styrene-divinylbenzene copolymer, 3 to $30\mu\text{m}$ in diameter.	PolymerX™ RP-1 Phenogel™ 100 Å	Spherical Spherical
L22 L23	A cation-exchange resin made of porous polystyrene gel with sulfonic acid groups, 5 to 15 µm in diameter. An anion-exchange resin made of porous polymethacrylate or polyacrylate gel with quaternary ammonium groups, 7-12 µm in size.	Rezex ROA-Organic Acid Shodex® IEC QA-825	Spherical Spherical
L24 L25	Polyvinylalcohol chemically bonded to porous silica particles, 5 µm in diameter.  Packing having the capacity to separate compounds with a MW range from 100 to 5000 daltons (as determined by polyethylene oxide), applied to neutral, anionic, and cationic water-soluble polymers.  A polymethacrylate resin base, crosslinked with poly-hydroxylated ether (surface contained some residual carboxyl functional groups) was found suitable.	PolySep <sup>™</sup> -GFC-P2000 Shodex OHpak SB-802.5HQ	Spherical Spherical

1105-	N. L. and C. and	B	D. I'. I. 81
	Column Classification	Recommended Phenomenex Column	
L26 L27	Butyl silane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter.	Jupiter 300 C4 Biozen WidePore C4 Sepra Silica	Spherical Core-Shell
L27	Porous silica particles, 30 to 50 µm in diameter.  A multifunctional support, which consists of a high purity, 100 Å, spherical silica substrate that has been	Sepra Silica	Irregular
LZO	bonded with anionic exchanger, amine functionality in addition to a conventional reversed phase C8 functionality.		
L29	Gamma alumina, reversed phase, low carbon percentage by weight, alumina-based polybutadiene spherical particles, 5 µm diameter with a pore volume of 80 Å.		
L30	Ethyl silane chemically bonded to a totally porous silica particle, 3 to 10 µm in diameter.		
L31	A hydroxide-selective, strong anion-exchange resin-quaternary amine bonded on latex particles attached to a core of $8.5\mu m$ macroporous particles having a pore size of 2000Å and consisting of ethylvinylbenzene cross-linked with $55\%$ divinyl benzene.		
L32	A chiral ligand-exchange resin packing-L-proline copper complex covalently bonded to irregularly shaped silica particles, 5 to $10\mu m$ in diameter.		
L33	Packing having the capacity to separate dextrans by molecular size over a range of 4,000 to 500,000 daltons. It is spherical, silica-based and processed to provide pH stability.	Yarra SEC-2000 / SEC-3000 BioSep-SEC-S2000 / SEC-S3000	Spherical Spherical
L34	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the lead form, 7 to $9\mu m$ in diameter.	Rezex RPM-Monosaccharide	Spherical
L35	A zirconium-stabilized spherical silica packing with a hydrophilic (diol-type) molecular monolayer bonded phase having a pore size of 150Å.	(BioSep-SEC-S2000 or Yarra SEC-2000 may be used)	Spherical Spherical
L36	3,5-dinitrobenzoyl derivative of L-phenylglycine covalently bonded to 5 μm aminopropyl silica.		
L37	Polymethacrylate gel packing having the capacity to separate proteins by molecular size over a range of 2,000 to 40,000 daltons.	PolySep-GFC-P3000 Shodex OHpak SB-803HQ	Spherical Spherical
L38	Methacrylate-based size-exclusion packing for water-soluble samples.	PolySep-GFC-P series Shodex OHpak SB-800HQ	Spherical Spherical
L39	Hydrophilic polyhydroxymethacrylate gel of totally porous spherical resin.	PolySep-GFC-P series Shodex OHpak SB-800HQ series Shodex RSpak DM-614	Spherical Spherical Spherical
L40	Cellulose tris-3,5-dimethylphenylcarbamate coated porous silica particles, 3 µm to 20 µm in diameter.	Lux Cellulose-1	Spherical
L41	Immobilized $\alpha_1$ -acid glycoprotein on spherical silica particles, 5 $\mu$ m in diameter.		
L42 L43	Octylsilane and octadecylsilane groups chemically bonded to porous silica particles, 5 µm in diameter.  Pentafluorophenyl groups chemically bonded to silica particles or superficially porous particles, by a propyl spacer, 1.5 to 10 µm in diameter.	Kinetex F5 Luna PFP(2)	Core-Shell Spherical
L44	A multifunctional support, which consists of a high purity, 60 Å, spherical silica substrate that has been bonded with a cationic exchanger, sulfonic acid functionality in addition to a conventional reversed phase C8 functionality.	Lulia I I I (2)	орненса
L45	Beta cyclodextrin, R, S-hydroxypropyl ether derivative, bonded to porous silica particles, 3 to 10 µm in diameter	Shiseido Chiral CD-Ph	Spherical
L46	Polystyrene/divinylbenzene substrate agglomerated with quaternary amine functionalized latex beads, about 9 to 11 $\mu$ m in diameter.		
L47	High capacity anion-exchange microporous substrate, fully functionalized with a trimethylamine group, 8 µm in diameter.		
L48	Sulfonated, cross-linked polystyrene with an outer layer of submicron, porous, anion-exchange microbeads, 5 to 15 µm in diameter.		
L49	A reversed phase packing made by coating a thin layer of polybutadiene on to spherical porous zirconia particles, 3 to 10 μm in diameter.		
L50	Multifunction resin with reversed phase retention and strong anion-exchange functionalities.  The resin consists of ethylvinylbenzene, 55 % cross-linked with divinylbenzene copolymer,  3 to 15 µm in diameter, and a surface area of not less than 350 m²/g. Substrate is coated with quaternary ammonium functionalized latex particles consisting of styrene cross-linked with divinylbenzene.		
L51 L52	Amylose tris-3,5-dimethylphenylcarbamate-coated, porous, spherical, silica particles, 3 to 10 µm in diameter.  A strong cation-exchange resin made of porous silica with sulfopropyl or sulfoethyl groups, 1 to 10 µm in diameter.	Lux Amylose-1	Spherical
L53	Weak cation-exchange resin consisting of ethylvinylbenzene, 55 % cross-linked with divinylbenzene copolymer, 3 to 15 µm diameter. Substrate is surface grafted with carboxylic acid and/or phosphoric acid functionalized monomers. Capacity not less than 500 µEq/column.		
L54	A size exclusion medium made of covalent bonding of dextran to highly cross-linked porous agarose beads, 5 to 15 µm in diameter.		
L55	A strong cation-exchange resin made of porous silica coated with polybutadiene-maleic acid copolymer, about 5 µm in diameter.		
L56 L57	Propyl silane chemically bonded to totally or superficially porous silica particles, 3 to 10 µm in diameter.  A chiral-recognition protein, ovomucoid, chemically bonded to silica particles, about 5 µm in diameter, with a pore size of 120 Å.	Ultron ES-0VM	Spherical
L58	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the sodium form, about 6 to 30 $\mu m$ in diameter.	Rezex RNM-Carbohydrate	Spherical
L59	Size-exclusion separations of proteins (separation by molecular weight) over the range of 5 to 7000 kDa. Spherical (1.5 to $10\mu m$ ), silica or hybrid packing with a hydrophilic coating.	Yarra SEC-2000 BioSep-SEC-S2000 Yarra SEC-3000 BioSep-SEC-S3000	Spherical Spherical Spherical Spherical
L60	Spherical, porous silica gel, 10 µm or less in diameter, surface has been covalently modified with alkyl amide groups and endcapped.		
L61	Hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of 13 µm microporous particles, pore size less than 10Å, and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene with a latex coating composed of 85 nm diameter microbeads bonded with alkanol quaternary ammonium ions (6%).		

JSP (	olumn Classification	Recommended Phenomenex Column	n Particle Shap
.62	C30 silane bonded phase on a fully porous spherical silica or superficially porous particles, 3 to 15 µm in diameter.	Develosil RP-Aqueous Develosil RP-Aqueous-AR	Spherical Spherical
.63	Glycopeptide teicoplanin linked through multiple covalent bonds to a 100Å spherical silica.	Bovoloon in Aquoddo An	орпопош
.64	Strongly basic anion-exchange resin consisting of $8\%$ crosslinked styrene divinylibenzene copolymer with a quaternary ammonium group in the chloride form, $45$ to $180\mu m$ in diameter.		
.65	Strongly acidic cation-exchange resin consisting of 2 $\%$ sulfonated crosslinked styrene divinylbenzene copolymer with a sulfonic acid group in the hydrogen form, 63 to 250 $\mu m$ in diameter.		
66	A crown ether coated on a 5 µm particle size silica gel substrate. The active site is (S)-18-crown-6-ether.	Acchinels ODD 50	Cabariaal
67 68	Porous vinyl alcohol copolymer with a C18 alkyl group attached to the hydroxyl group of the polymer, 2 to 10 µm in diameter.  Spherical, porous silica, 10 µm or less in diameter, the surface of which has been covalently modified with alkyl	Asahipak ODP-50	Spherical
69	Sprietrian, portuge since, it of the result and the surface of which has been covariently into metal and armide groups and not endcapped.  Ethylvinylbenzene/divinylbenzene substrate agglomerated with quaternary amine functionalized 130 nm latex		
70	beads, about 6.5 µm in diameter.  Cellulose tris (phenyl carbamate) coated on 5 µm silica.		
71	A rigid, spherical polymethacrylate 4 to 6 µm in diameter.	Shodex RSpak DE-413	Spherical
70	(C) should be size and 0.5 distance line was linked a social state banded to silice	Shodex RSpak DE-613	Spherical
2	(S)-phenylglycine and 3,5-dinitroaniline urea linkage covalently bonded to silica.  A rigid, spherical polydivinylbenzene particle 5 to 10 µm in diameter.		
74	A strong anion-exchange resin consisting of a highly cross-linked core of 7 µm macroporous particles having a 100 Å average pore size and consisting of ethylvinylbenzene cross-linked with 55 % divinylbenzene and an anion-exchange layer grafted to the surface, which is functionalized with alkyl quaternary ammonium ions.		
5	A chiral-recognition protein, bovine serum albumin (BSA), chemically bonded to silica particles, about 7 $\mu$ m in diameter, with a pore size of 300 Å.		
76	Silica-based weak cation-exchange material, 5 μm in diameter. Substrate is surface polymerized polybutadiene-maleic acid to provide carboxylic acid functionalities. Capacity not less than 29 μEq/column.		
77	Weak cation-exchange resin consisting of ethylvinylbenzene, $55\%$ cross-linked with divinylbenzene copolymer, 6 to $9\mu m$ diameter. Substrate is surface grafted with carboxylic acid functionalized groups. Capacity not less than $500~\mu Eq/column$ ( $4mm \times 25cm$ ).		
8	A silane ligand that consists of both reversed phase (an alkyl chain longer than C8) and anion-exchange (primary, secondary, tertiary, or quaternary amino groups) functional groups chemically bonded to porous or non-porous or ceramic micro-particles, 1.0 to 50 µm in diameter or a monolithic rod.		
9	A chiral-recognition protein, human serum albumin (HSA), chemically bonded to silica particles, about 5 µm in diameter.		
0	Cellulose tris(4-methylbenzoate)-coated, porous, spherical, silica particles, 5 to 20 µm in diameter.	Lux Cellulose-3	Spherical
81	A hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of 9 µm porous particles having a pore size of 2000 Å units and consisting of ethylvinylbenzene cross-linked with 55 % divinylbenzene with a latex coating composed of 70 nm diameter microbeads (6 % crosslinked) bonded with alkanol quaternary ammonium ions.		
32	Polyamine chemically bonded to cross-linked polyvinyl alcohol polymer, 4 - 5 µm in diameter	Asahipak NH <sub>2</sub> P-50	Spherical
33	A hydroxide-selective, strong anion-exchange resin-quaternary amine bonded on latex particles attached to a core of 10.5 µm microporous particles having a pore size of 10Å and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene.		
34	Weak cation-exchange resin consisting of ethylvinylbenzene, 55 % cross-linked with divinylbenzene copolymer, 5 µm diameter. Substrate is surface grafted with carboxylic acid functionalized groups. Capacity not less than 8400 µEq column (5 mm x 25 cm).		
35	A silane ligand that consists of both reversed phase (an alkyl chain longer than C8) and weak cation-exchange (carboxyl groups) functional groups chemically bonded to porous or non-porous particles, 1.0 to 50 µm in diameter.		
36	Fused core particle with a highly polar ligand possessing multiple hydroxyl groups tethered to the silica gel outer layer.		
37	Dodecyl silane chemically bonded to porous or superficially porous silica particles, 1.5 to $10\mu m$ in diameter.	Synergi Max-RP	Spherical
8	Glycopeptide vancomycin linked through multiple covalent bonds to 100 Å spherical silica.		
89	Packing having the capacity to separate compounds with a molecular weight range from 100 - 3000 dalton (as determined by polyethylene oxide), applied to neutral and anionic water-soluble polymers. A polymethacrylate resin base, cross-linked with polyhydroxylate ether (surface contains some residual cationic functional groups).	Shodex OHpak SB-802.5 HQ	Spherical
90	Amylose tris-[(S)-alpha-methylbenzylcarbamate] coated on porous, spherical silica particles, 3 to 10 µm in diameter.		
91	Strong anion-exchange resin consisting of monodisperse porous polystyrene/divinylbenzene beads coupled with quaternary amine. Bead size is 3 to $10\mu m$ .		
32	A strong anion-exchange resin consisting of a highly cross-linked core of 5-9 µm macroporous particles having a 100 Å average pore size and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene and an anion-exchange layer grafted to the surface, which is functionalized with alkanol quaternary ammonium ions.		
93	Cellulose tris (3,5-dimethylphenylcarbamate) reversed phase chiral stationary phase coated on 3 or $5\mu m$ silica gel particles.	Lux Cellulose-1	Spherical
14	A strong anion-exchange resin consisting of highly cross-linked 15 $\mu m$ microporous particles functionalized with very low cross-linked latex (0.5 %) to provide alkanol quaternary ammonium ion-exchange sites.		
95	Highly polar alkyl ligand comprising five hydroxyl groups that are chemically bonded to totally porous or superficially porous silica, or a monolithic silica rod.		
96	Alkyl chain, reversed phase bonded to totally or superficially porous silica designed to retain hydrophilic and other polar compounds when using highly aqueous mobile phases, including 100 % aqueous, 1.5 $\mu$ m to 10 $\mu$ m in diameter.	Kinetex Polar C18 Luna Omega Polar C18 Synergi Hydro-RP Synergi Fusion-RP	Core-Shell Core-Shell Spherical Spherical Core-Shell Spherical Spherical

USP (	column Classification	Recommended Phenomenex Column	Particle Shape
L97	Weak cation-exchange resin consisting of a highly cross-linked core of 5.5 µm porous particles having a pore size of 2000Å and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene. Substrate is surface grafted with carboxylic acid functionalized groups. Capacity not less than 2400 µEq/column (4 mm x 25 cm).	Tiecommended i nenomenea condimi	Tarticle Shape
L98	Weak cation-exchange resin consisting of a highly cross-linked core of $8\mu m$ microporous particles having an average pore size of $10\text{Å}$ and consisting of ethylvinylbenzene cross-linked with $55\text{\%}$ divinylbenzene. Substrate is surface grafted with carboxylic acid functionalized groups. Capacity not less than $46\mu \text{Eq/column}$ ( $4\text{mm} \times 5\text{cm}$ ).		
L99	Amylose tris-(3,5- dimethylphenylcarbamate), immobilized on porous, spherical, silica particles, 3 to 5 µm in diameter	Lux i-Amylose-1	Spherical
L100	A 55% cross-linked, microporous, hydrophobic resin core ( $9\mu m$ microporous particles having a pore size of $10\text{Å}$ ) that consists of a bilayer of anion and cation-exchange latex. The first layer is fully sulfonated ( $140n m$ ) and the second layer is fully aminated ( $76n m$ ).		
L101	Cholesteryl groups chemically bonded to porous or non-porous silica or ceramic micro-particles, 1.5 to $10\mu m$ in diameter, or a monolithic rod.		
L102	(Naproxen, (S,S)Whelk-O 1) 1-(3,5- dinitrobenzamido)-1,2,3,4- tetrahydrophenanthrene covalently bonded to porous spherical silica particles, 5 to 10 µm in diameter.		
L103	A hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of 7.5 µm porous particles having a pore size of 2000Å and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene electrostatically bonded with hyperbranched alkanol quaternary ammonium ions.		
L104	Triazole groups chemically bonded to porous silica particles, 1.5 to 10 $\mu m$ in diameter.		
L105	A strong anion-exchange resin consisting of a highly cross-linked 9 µm supermacroporous (2000 Å) particles functionalized with very low cross-linked latex (0.2 %) to provide alkyl quaternary ammonium ion sites.		
L106	Weak cation-exchange resin consisting of ethylvinylbenzene, $55\%$ cross-linked with divinylbenzene copolymer, $5-8\mu m$ diameter, macroporous particles having an average pore size of $100\text{Å}$ units. Substrate is surface grafted with carboxylic acid and phosphonic acid functional groups. Capacity not less than $2800\mu Eq/column$ (4 mm x $25cm$ ).		
L107	Cellulose tris(4-methylbenzoate)-coated porous spherical particles, 3 to $5\mu m$ in diameter, for use with reversed phase mobile phases.	Lux Cellulose-3	Spherical
L108	A chiral-recognition protein, cellobiohydrolase (CBH), chemically bonded to silica particles, about $5\mu m$ in diameter.		
L109	Spherical particles of porous graphitic carbon, 1.5 to 30 µm in diameter.		
L110	A strong anion-exchange resin consisting of a highly cross-linked 13 µm microporous (less than 10 Å) particles coated with very low cross-linked latex (0.5%) to provide alkanol quaternary ammonium ion-exchange sites.		
L111	Polyamine chemically bonded to porous spherical silica particles, 5 µm in diameter.		
L112	A hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of $8.5\mu m$ porous particles having a pore size of $2000\text{Å}$ units and consisting of ethylvinylbenzene cross-linked with $55\%$ divinylbenzene with a latex coating composed of $65\text{nm}$ diameter microbeads ( $5\%$ cross-linked) bonded with alkanol quaternary ammonium ions.		
L113	A hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of 7.5 µm porous particles having a pore size of 2000Å and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene with a latex coating composed of 65 nm diameter microbeads (5% crosslinked) bonded with alkanol quaternary ammonium ions.		
L114	Sulfobetaine graft-polymerized to totally or superficially porous silica, 1.5 to 10 µm in diameter, or a monolithic rod. Packing having densely bonded zwitterionic groups with 1:1 charge balance.		
L115	Ethylvinylbenzene/divinylbenzene substrate (55% cross-linked) agglomerated with quaternary amine functionalized 275 nm latex microbeads (6% cross-linked), about 8.5 µm in diameter.		
L116	Sulfonated ethylvinylbenzene/divinylbenzene substrate agglomerated with hydrophilic quaternary amine functionalized glycidyl-derivative methacrylate microbeads, approximately 2 to 50 µm in diameter.		
L117	A crown ether coated on a $5\mu m$ particle size silica gel substrate. The active site is (R)-18-crown-6-ether.		
L118	Aqueous polymerized C18 groups on silica particles, 1.2 to $5\mu m$ in diameter.	Kinetex PAH	Core-Shell
L119	Cellulose tris-(3,5-dichlorophenylcarbamate), immobilized on porous, spherical, silica particles, 3 to $5\mu m$ in diameter.	Lux i-Cellulose-5	Spherical
L120	A hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of $13\mu m$ microporous particles having a pore size of less than $10\text{\r{A}}$ units and consisting of ethylvinylbenzene cross-linked with $55\%$ divinylbenzene with a latex coating composed of $65\text{nm}$ diameter microbeads ( $8\%$ cross-linked) bonded with alkanol quaternary ammonium ions. Capacity not less than $10\mu \text{Cq/column}$ ( $4\text{mm}$ x $5\text{cm}$ ).		
L121	A hydroxide-selective, strong anion-exchange resin consisting of a highly cross-linked core of 11 µm porous particles having a pore size of less than 10 Å units and consisting of ethylvinylbenzene cross-linked with 55 % divinylbenzene electrostatically bonded with hyperbranched alkanol quaternary ammonium ions.		
L122	Sulfobetaine graft-polymerized to totally or superficially porous hydrophilic polymer particles, 1.0 to 10 µm in diameter, or a monolithic rod. Packing having densely bonded zwitterionic groups with 1:1 charge balance.		
L123	$\label{eq:control} \textbf{Cellulose tris} (3\text{-chloro-4-methylphenylcarbamate}) \ coated \ porous \ silica \ particles, 3 \ to \ 20 \ \mu m \ in \ diameter.$	Lux Cellulose-2	Spherical
L124	Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the silver form, average 9 µm in diameter.	Rezex RSO-Oligosaccharide	Spherical
L125	, , ,	Shodex IC YS-50	Spherical
L126	Amylose tris (3-chlorophenylcarbamate), immobilized on porous, spherical, silica particles, 1 $\mu$ m to 20 $\mu$ m in diameter.		
L127	A crown ether chemically bonded to a $5\mu m$ particle size silica gel substrate. The active site is (S)- pseudo-18-crown-6-ether.	Sumichiral 0A-8000	Spherical
L128	Porous particles of polystyrene divinyl benzene with linear molecular weight operating range from 200 to 2,000,000 g/mol (polystyrene equivalent), 5 µm in diameter.		
L130	Cellulose tris(3,5-dimethylphenylcarbamate) coated on 10 Âm silica gel particles.	Lux Cellulose-1	Spherical

**Australia** 

t: +61 (0)2-9428-6444 auinfo@phenomenex.com

**Austria** t: +43 (0)1-319-1301 anfrage@phenomenex.com

**Belgium** t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

**Canada** t: +1 (800) 543-3681 info@phenomenex.com

China

t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic t: +420 272 017 077

cz-info@phenomenex.com

**Denmark** t: +45 4824 8048 nordicinfo@phenomenex.com

**Finland** t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

**Germany** t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong t: +852 6012 8162

hkinfo@phenomenex.com

**India** t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Indonesia t: +62 21 3952 5747 indoinfo@phenomenex.com

t: +353 (0)1 247 5405 eireinfo@phenomenex.com

**Italy** t: +39 051 6327511 italiainfo@phenomenex.com

**Japan** t: +81 (0) 120-149-262 jpinfo@phenomenex.com

**Luxembourg** t: +31 (0)30-2418700

nlinfo@phenomenex.com

t: 01-800-844-5226 tecnicomx@phenomenex.com

**The Netherlands** t: +31 (0)30-2418700 nlinfo@phenomenex.com

**New Zealand** 

t: +64 (0)9-4780951 nzinfo@phenomenex.com

**Norway** t: +47 810 02 005 nordicinfo@phenomenex.com

**Poland** t: +48 22 51 02 180 pl-info@phenomenex.com

**Portugal** t: +351 221 450 488 ptinfo@phenomenex.com

**Singapore** t: 800-852-3944 sginfo@phenomenex.com

Slovakia

t: +420 272 017 077 sk-info@phenomenex.com

**Spain** t: +34 91-413-8613 espinfo@phenomenex.com

**Sweden** t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

**Switzerland** t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

**Thailand** t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

**United Kingdom** t: +44 (0)1625-501367

ukinfo@phenomenex.com

**USA** t: +1 (310) 212-0555 www.phenomenex.com/chat

All other countries/regions
 Corporate Office USA

www.phenomenex.com/chat



#### www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

Trademarks

Biozen, Gemini, Kinetex, Luna, Jupiter, Clarity, and Lux are registered trademarks and Aeris, Yarra, Synergi, Onyx, Prodigy, PhenoSphere, Rezex, BioSep, PolymerX, PolySep, Sepra, and Phenogel are trademarks of Phenomenex. Develosil is a registered trademark of Nomura Chemical Co. Shodex and Asahipak are registered trademarks of Showa Denko K.K. Develosil is a registered trademark of Nomura Chemical Co. Shodex and Asahipak are registered trademarks of Showa Denko K.K. Ultron is a registered trademark of Shinwa Chemical Industries. Chirex is a trademark of Chirex, Inc., licensed to Phenomenex. Shiseido is a registered trademark of Shiseido Company, Ltd.

Gemini and Kinetex EVO are patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.