TN-0154 Robust Analysis of Nitrosamine in a Losartan Drug Substance using Strata[™] Activated Carbon Extraction

Zandra Baja, Ramkumar Dhandapani, PhD, Philip Leong, and Bryan Tackett, PhD Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

The presence of nitrosamine as an impurity has been detected in drug substances and drug products, such as Losartan and other "-sartan" drugs. Losartan is an angiotensin II receptor blocker (ARB) and is used in the treatment of high blood pressure. N-Nitrosodimethylamine (NDMA) and N-Nitrosodiethylamine (NDEA) are Class 2A carcinogens that have been found to be present as impurities in several different generic "-sartan" drug substances and drug products. NDMA is highly toxic and is a known carcinogen in lab animals and, along with NDEA, has been classified as a probable human carcinogen. Accurate identification and quantification of these impurities is essential to assure quality and safety of the drug substance and drug product.

Here, a method that employs activated carbon for the extraction of polar analytes such as N-Nitrosodimethylamine is used. While dealing with polar molecules like nitrosamines, traditional sample preparation using Solid Phase Extraction (SPE) might not be efficient. Since the analyte is polar, a specialized extraction technique is needed. Activated carbon contains porous carbon with defined pore volume and high surface area that is an appropriate choice for polar analytes like nitrosamines.

In this technical note, we are providing the optimized extraction protocol for the extraction of nitrosamines from a Losartan drug substance using Strata Activated Carbon followed by GC-MS analysis. Traditional methods for analysis of nitrosamines require GC-MS/MS with a headspace capability to meet the low-level detection limits. With Strata Activated Carbon sample preparation and pseudo large volume injections, GC-MS was sufficient to achieve low-level detection of nitrosamines.

Sample Preparation

Condition:	Strata Activated Carbon, 400 mg cartridge (Part No. <u>CS0-9210</u>) with 2 washes of 3 mL Methylene Chloride, 3 washes of 3 mL Methanol
Equilibrate:	Cartridges with 6 washes of 3 mL Water
Load:	100 mL sample onto cartridges
Dry:	Cartridges for 10 min
Elute:	With 1 mL Methylene Chloride and repeat elution 3 more times (3 mL total volume), add 10 μL Ethylene Glycol
Water Removal:	Pass eluent through Sodium Sulfate tubes, 1 g/6 mL (Part No. <u>88-5124-JCH</u>) that is prewetted with Methylene Chloride and wash with 3 mL Methylene Chloride
Evaporate:	Solvent to approximately 1 mL in a water bath at 20-25 °C under a gentle stream of Nitrogen
Adjust:	Final volume to 1 mL with Methylene Chloride



GC-MS Conditions

т

Column:	Zebron™ ZB-WA	(PLUS [™]		
Dimension:	30 meter x 0.32 r	nm x 1.00 μm	1	
Part No.:	7HM-G013-22			
Injection:	Pulse Splitless (0.	5 min @ 140	Kpa) @ 200 °C,	8 μL
Recommended Liner:	Zebron PLUS Line	er for Shimadz	u® 2010, 3.4 m	m ID Single Taper
Liner Part No.:	AG2-4B10-05			
C	Helium @ 1.8 mL/min (Constant Flow)			
Carrier Gas:	Hellum @ 1.8 mL	/min (Consta	nt Flow)	
	Ramp(°C/min)		Time(min)	
			-	
		Temp (°C)	Time(min)	
	Ramp(°C/min)	Temp (°C) 45	Time(min) 3.0	
	Ramp(°C/min) - 10	Temp (°C) 45 130	Time(min) 3.0 3.0	
	Ramp(°C/min) - 10 15 40	Temp (°C) 45 130 190	Time(min) 3.0 3.0 0.0	
Oven Program:	Ramp(°C/min) 10 15 40 SIM	Temp (°C) 45 130 190	Time(min) 3.0 3.0 0.0	

Table 1. GC-MS SIM Parameters

Peak No.	Analyte Name	m/z
1	N-Nitrosodimethylamine-d6 (NDMA-d6)	80
2	N-Nitrosodimethylamine (NDMA)	74
3	N-Nitrosomethylethylamine (NMEA)	88
4	N-Nitrosodiethylamine (NDEA)	102
5	N-Nitrosodipropylamine (NDPA)	70
6	N-Nitrosodibutylamine (NDBA)	84
7	N-Nitrosopiperidine (NPIP)	114
8	N-Nitrosopyrrolidine (NPYR)	100

Page 1 of 7

Ramkumar Dhandapani, PhD

Ramkumar has earned a PhD in Analytical Chemistry and has over 17 years troubleshooting and hands-on experience in chromatography. He loves to write poems, watch, and read Shakespeare's plays.

Results and Discussion

Figure 1 is a representation of the Strata™ Activated Carbon extraction cartridges. The porous carbon is optimized for reproducible extraction of polar analytes. Although the 400 mg format was employed in this application, there is also a 2 g/6 mL format available for larger sample sizes.

All nitrosamines tested in this study were extracted and reconstituted with Dichloromethane (DCM), so a DCM blank was run to demonstrate there were no nitrosamines in the interface (**Figure 2**). As presented in **Figures 3**, **4**, and **5** the individual standards, including the internal standard (IS) N-Nitrosodimethylamine-d6 (NDMA-d6), were identified after extraction using the Strata Activated Carbon extraction cartridges at 0.5, 1, and 5 μ g/mL respectively. This shows accuracy to low level concentrations of nitrosamines. The analytical method proved to be linear for the extraction and analysis at this concentration range.

Losartan and other members of the "-sartan" family of drugs are known to contain nitrosamine impurities or produce nitrosamines as they are broken down. **Figure 6** is a representative chromatogram of a Losartan sample after extraction using the Strata Activated Carbon extraction cartridge. For analysis, the Zebron[™] ZB-WAX_{PLUS[™]} was utilized. The high polar selectivity, extensive cross-linkage in the stationary phase, and the highly deactivated fused silica in this column provide MS compatibility and an inert environment for active analytes. This results in low baseline noise and symmetric peaks for even the most challenging analytes like nitrosamines.

To establish accuracy of analysis, Losartan sample was spiked with nitrosamines at 0.5 μ g/mL (Figure 7). Repeatability of the analysis was demonstrated using multiple injections of the sample spike. Finally, to clearly identify nitrosamines, and to show low detection in real samples, the Losartan sample was overlapped with nitrosamine standard at 5.0 μ g/mL concentration (Figure 8).




Figure 1. Strata Activated Carbon Extraction Cartridges

Table 2. % Recovery for Extracted Losartan Sample


Analyte Name	Conc (µg/mL)	% Recovery
NDMA-d6	0.9	-
NDMA	0.525	70
NMEA	0.516	103
NDEA	0.529	106
NDPA	0.498	7
NDBA	0.548	83
NPIP	0.514	68
NPYR	0.536	77

Figure 2. DCM Blank




Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists

Figure 3. Nitrosamine Standard at 0.5 µg/mL after Strata™ Activated Carbon Extraction




Peak No.	Analyte	RT (min)	% RSD
1	NDMA-d6 (IS)	11.10	-
2	NDMA	11.11	0.3
3	NMEA	11.86	0.0
4	NDEA	12.38	0.5
5	NDPA	15.29	2.7
6	NDBA	18.03	6.4
7	NPIP	18.53	3.5
8	NPYR	18.85	3.5
N = 3 Injec	ctions		

Figure 4. Nitrosamine Standard at 1.0 $\mu\text{g}/\text{mL}$ after Strata Activated Carbon Extraction



Peak No.	Analyte	RT (min)	% RSD
1	NDMA-d6 (IS)	11.10	-
2	NDMA	11.11	0.4
3	NMEA	11.86	0.3
4	NDEA	12.38	0.7
5	NDPA	15.29	3.2
6	NDBA	18.03	6.4
7	NPIP	18.53	4.3
8	NPYR	18.85	4.0
N = 3 Injec	ctions		

Figure 5. Nitrosamine Standard at 5.0 µg/mL after Strata Activated Carbon Extraction



Peak No.	Analyte	RT (min)	% RSD
1	NDMA-d6 (IS)	11.11	-
2	NDMA	11.12	0.1
3	NMEA	11.86	0.1
4	NDEA	12.38	0.2
5	NDPA	15.29	0.3
6	NDBA	18.04	2.6
7	NPIP	18.53	0.9
8	NPYR	18.85	1.2
N = 3 Injec	ctions		

?

Figure 6. Losartan Sample after Strata™ Activated Carbon Extraction







Figure 7. Chromatogram of 0.5 µg/mL Nitrosamine Spiked Losartan Sample after Strata Activated Carbon Extraction




Peak No.	Analyte	RT (min)	% RSD
1	NDMA-d6 (IS)	11.10	-
2	NDMA	11.12	0.8
3	NMEA	11.86	0.6
4	NDEA	12.38	0.9
5	NDPA	15.29	4.2
6	NDBA	18.03	9.4
7	NPIP	18.52	6.3
8	NPYR	18.85	6.1
N = 6 Injec	ctions		





Peak No.	Analyte	RT (min)
1	NDMA-d6 (IS)	11.11
2	NDMA	11.12
3	NMEA	11.86
4	NDEA	12.38
5	NDPA	15.29
6	NDBA	18.04
7	NPIP	18.53
8	NPYR	18.85



Conclusions

Strata[™] Activated Carbon extraction provided reproducible extraction of several nitrosamines in a Losartan sample matrix. The results demonstrated the range of detection of nitrosamines after extraction using the Strata Activated Carbon extraction cartridge and the ability to identify the presence or absence of nitrosamine in a Losartan sample. The combination of sample prep extraction using Strata Activated Carbon and pseudo large volume injection with pressure pulse helped achieve lower detection limits with simple GC-MS detection.

Strata Tubes and Cartridges Ordering Information

Strata Activated Carbon				
Part	Sorbent Mass	Volume	Format	Units/pk
<u>CS0-9209</u>	2 g	6 mL	Tube	30
<u>CS0-9210</u>	400 mg	Pass through	Cartridge	50

Zebron[™] ZB-WAX_{PLUS}[™] GC Columns Ordering Information

			<u> </u>
ID(mm)	df(µm)	Temp. Limits °C	Part No.
10-Meter			
0.10	0.10	20 to 250/260	7CB-G013-02
15-Meter			
0.25	0.25	20 to 250/260	<u>7EG-G013-11</u>
0.53	1.00	20 to 230/240	7EG-G013-22
20-Meter			
0.18	0.18	20 to 250/260	7FD-G013-08
30-Meter			
0.25	0.25	20 to 250/260	<u>7HG-G013-11</u>
0.25	0.50	20 to 250/260	<u>7HG-G013-17</u>
0.32	0.25	20 to 250/260	7HM-G013-11
0.32	0.50	20 to 250/260	7HM-G013-17
0.32	1.00	20 to 230/240	7HM-G013-22
0.53	0.25	20 to 250/260	<u>7HK-G013-11</u>
0.53	3.00	20 to 230/240	<u>7HK-G013-22</u>
60-Meter			
0.25	0.15	20 to 250/260	<u>7KG-G013-05</u>
0.25	0.25	20 to 250/260	<u>7KG-G013-11</u>
0.25	0.50	20 to 250/260	<u>7KG-G013-17</u>
0.32	0.25	20 to 250/260	<u>7KM-G013-11</u>
0.32	0.50	20 to 250/260	<u>7KM-G013-17</u>
0.53	1.00	20 to 230/240	<u>7KK-G013-22</u>

Note: If you need a 5 in. cage, contact Technical support via Phenomenex.com/chat or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

?



Description	Application	Inlet Style	Dimensions	Deactivation	Part No.	Unit
Zebron™ PLUS Liners fo	or 17A, 2014 and 2025	Models				
Single Taper Z-Liner™	Pesticides	S/SL	3.4 x 95	PLUS Inert	<u>AG2-3B13-01</u> <u>AG2-3B13-05</u> <u>AG2-3B13-25</u>	ea 5/pk 25/pk
Straight Z-Liner	Dirty samples, Volatiles, High initial oven temperatures	S/SL	3.4 x 95	PLUS Inert	AG2-3B03-01 AG2-3B03-05 AG2-3B03-25	ea 5/pk 25/pk
Zebron PLUS Liners for	2010 Models					
Single Taper	Dirty samples, Volatiles, High initial oven temperatures	S/SL	3.4 x 95	PLUS Inert	<u>AG2-4B10-01</u> <u>AG2-4B10-05</u> <u>AG2-4B10-25</u>	ea 5/pk 25/pk
Single Taper Z-Liner	Pesticides	S/SL	3.4 x 95	PLUS Inert	<u>AG2-4B13-01</u> <u>AG2-4B13-05</u> <u>AG2-4B13-25</u>	ea 5/pk 25/pk
Straight	Volatiles	S/SL	3.4 x 95	PLUS Inert	<u>AG2-4B00-01</u> <u>AG2-4B00-05</u> <u>AG2-4B00-25</u>	ea 5/pk 25/pk
Straight Z-Liner	Dirty samples, Volatiles, High initial oven temperatures	S/SL	3.4 x 95	PLUS Inert	AG2-4B03-01 AG2-4B03-05 AG2-4B03-25	ea 5/pk 25/pk

Sodium Sulfate Tubes Ordering Information

Strata™			
Format	Sorbent Mass	Part Number	Unit
Tube	1 g	<u>8B-S124-JCH</u>	6 mL (30/box)
Giga™ Tube	5 g	<u>8B-S124-LEG</u>	20 mL (20/box)



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on <u>www.phenomenex.com</u>, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal <u>www.phenomenex.com/Chat</u>.

India

Indonesia

Ireland

Japan t: +81 (0) 120-149-262

Mexico

t: +62 21 5019 9707

t: +353 (0)1 247 5405

Italy t: +39 051 6327511

Luxembourg t: +31 (0)30-2418700

t: 01-800-844-5226

The Netherlands

t: +31 (0)30-2418700

New Zealand

Norway t: +47 810 02 005

t: +48 22 104 21 72

Poland

t: +64 (0)9-4780951

t: +91 (0)40-3012 2400

indiainfo@phenomenex.com

indoinfo@phenomenex.com

eireinfo@phenomenex.com

italiainfo@phenomenex.com

jpinfo@phenomenex.com

nlinfo@phenomenex.com

nlinfo@phenomenex.com

nzinfo@phenomenex.com

nordicinfo@phenomenex.com

pl-info@phenomenex.com

tecnicomx@phenomenex.com

Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belgium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada t: +1 (800) 543-3681 info@phenomenex.com

China t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic t: +420 272 017 077 cz-info@phenomenex.com

Denmark t: +45 4824 8048 nordicinfo@phenomenex.com

Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com

Hong Kong t: +852 6012 8162 hkinfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right. www.phenomenex.com/behappy

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

> Singapore t: +65 6559 4364 sginfo@phenomenex.com

Slovakia t: +420 272 017 077 sk-info@phenomenex.com

Spain t: +34 91-413-8613 espinfo@phenomenex.com

Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand t: +66 (0) 2 566 0287 thaiinfo@phenomenex.com

United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com

USA t: +1 (310) 212-0555 www.phenomenex.com/chat

 All other countries/regions Corporate Office USA
t: +1 (310) 212-0555
www.phenomenex.com/chat

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at <u>www.phenomenex.com/phx-terms-and-conditions-of-sale</u>. Trademarks

Strata, Zebron, Z-Liner, Giga, and BE-HAPPY are trademarks of Phenomenex. Shimadzu is a registered trademark of Shimadzu Corporation Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Shimadzu Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures. © 2022 Phenomenex, Inc. All rights reserved.

Page 7 of 7

