

# APPLICATIONS



## Choosing the Right UHPLC Column for Peptide Mapping

Zeshan Aqeel and Jeff Layne  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### Overview


Analysts are increasingly relying upon LC-MS to perform their peptide and protein digest analyses, and TFA is not often used in conjunction with LC-MS due to its strong signal-suppression properties. Thus, when performing peptide analysis by LC-MS using typical mobile phases such as water/acetonitrile/formic acid, secondary interactions between basic amino acids on the peptide and the underlying silica can lead to broad peaks and peak tailing, and can interfere with optimal resolution between close-eluting peaks. To obtain optimal performance for peptide analysis using either LC-UV or LC-MS, the ideal media should be highly efficient (narrow peaks) and highly inert (minimal peak tailing when using MS-compatible modifiers like formic acid). In addition, a more retentive media is often favorable, particularly for peptide mapping, to maximize retention of small, polar peptide fragments that may not have a high affinity for a typical C18 phase, particularly in the absence of TFA (trifluoroacetic acid).

Luna<sup>®</sup> Omega HPLC and UHPLC products, ranging in size from 1.6  $\mu$ m for UHPLC work up to 5  $\mu$ m for conventional HPLC as well as preparative HPLC, are the ideal choice for the analysis of small peptide and tryptic mapping. The high surface area and dense C18 bonding provide Luna Omega with excellent retention for small, more polar peptides and fragments. The thermally modified silica gives Luna Omega an efficiency advantage over similarly-sized fully porous silica products, which can translate into narrower peaks, improved sensitivity, and improved resolution in many cases. This thermal treatment process also renders the silica highly inert for minimal secondary interactions, even when using LC-MS mobile phase modifiers such as 0.1% formic acid. Lastly, the unique Luna Omega PS C18 chemistry contains a proprietary, positively-charged functional group on the surface that makes the surface even less likely to exhibit peak tailing when analyzing basic peptides, and can also provide a selectivity quite distinct from a standard C18 phase, which may be useful for separating a target peptide from a matrix interference that co-elutes on a traditional C18.

The three figures inside contain representative chromatograms of a tryptic digest of BSA obtained using three different UHPLC columns (Waters<sup>®</sup> ACQUITY<sup>®</sup> BEH C18 1.7  $\mu$ m, **Figure 1**; Luna Omega 1.6  $\mu$ m C18, **Figure 2**; and Luna Omega 1.6  $\mu$ m PS C18, **Figure 3**). We do recognize that the use of a 50 x 2.1 mm format for analyzing peptide mapping is not ideal, and one would expect a more conventional format to be 100 x 2.1 mm or 150 x 2.1 mm to obtain maximal information. However, our goal in this comparison is to focus on the comparative differences in chromatographic performance between the three columns, as this information can be expected to scale directly to other column dimensions as well.

As you can see, the “conventional” C18 phases (**Figures 1 and 2**) provide very similar performance for this map under these running conditions. Overall retention times are similar, and the profile of the eluting peaks is very similar as well. One could argue that certain fragments are better resolved on one column versus the other, but overall one could conclude that an analyst should be relatively confident that both columns would be equally successful for this given separation. In contrast, with the unique


### PS C18



#### Luna Omega PS C18

Unique, 100% aqueous stable mixed-mode phase that provides both polar and non-polar retention. The surface contains a positive charged ligand which aids in the retention of acidic compounds through ionic interactions, while the C18 ligand promotes general reversed phase retention. The positively charged surface also improves basic compound peaks shape through ionic repulsion.

### C18



#### Luna Omega C18

Rugged and highly efficient C18 with strong focus on hydrophobic retention of non-polar and polar compounds.


Luna Omega PS C18 phase (**Figure 3**), we obtain a completely distinct elution profile for the peptide fragments, particularly in the region highlighted in the yellow box. You can see for that region in particular, the Luna Omega PS C18 appears to provide significantly better separation between the peaks, a shift in selectivity that could prove useful in separating a target peptide fragment from isobaric interferences.

Thus, when developing methods for the tryptic mapping using either LC-UV or LC-MS, Luna Omega C18 and Luna Omega PS C18 should be evaluated as part of any thorough method development process, with the Luna Omega C18 column providing a more conventional elution profile and the Luna Omega PS C18 providing a unique and distinct selectivity solution.




# APPLICATIONS

**Figure 1**  
Waters® ACQUITY® BEH 1.7 µm C18; 420 Bar



## UHPLC Conditions

**Figure 2**  
Luna® Omega 1.6 µm C18; 413 Bar



Conditions same for all columns

Column: As specified

Dimensions: 50 x 2.1 mm

Mobile Phase: A: 0.1% Formic acid in Water

B: 0.1% Formic acid in Acetonitrile

Gradient	Time (min)	% B
	0	3
	20	35

Flow Rate: 0.5 mL/min


Temperature: 40 °C

Injection: 0.5 mL

Detection: MS/MS (SCIEX API 4000™)

Sample: BSA Tryptic Digest

**Figure 3**  
Luna Omega 1.6 µm PS C18; 301 Bar



\*Unique and distinct elution profile of peptides in this region\*

Phenomenex is not affiliated with Waters Technologies Corporation.  
Comparative separations may not be representative of all applications.

# APPLICATIONS

## Luna<sup>®</sup> Omega Ordering Information

1.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
<b>Polar C18</b>	00B-4748-A0	00D-4748-A0	00F-4748-A0
<b>C18</b>	00B-4742-A0	00D-4742-A0	00F-4742-A0

1.6 µm Minibore Columns (mm)					SecurityGuard <sup>™</sup> ULTRA Cartridges <sup>‡</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>Polar C18</b>	00A-4748-AN	00B-4748-AN	00D-4748-AN	00F-4748-AN	AJ0-9505
<b>PS C18</b>	00A-4752-AN	00B-4752-AN	00D-4752-AN	00F-4752-AN	AJ0-9508
<b>C18</b>	00A-4742-AN	00B-4742-AN	00D-4742-AN	00F-4742-AN	AJ0-9502

for 2.1 mm ID

3 µm Minibore Columns (mm)					SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0*
<b>Polar C18</b>	00A-4760-AN	00B-4760-AN	00D-4760-AN	00F-4760-AN	AJ0-7600
<b>PS C18</b>	00A-4758-AN	00B-4758-AN	00D-4758-AN	00F-4758-AN	AJ0-7605

for ID: 2.0 - 3.0 mm

3 µm MidBore <sup>™</sup> Columns (mm)				SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
<b>Polar C18</b>	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJ0-7600
<b>PS C18</b>	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJ0-7605

for ID: 2.0 - 3.0 mm

3 µm Analytical Columns (mm)					SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
<b>Polar C18</b>	00B-4760-E0	00D-4760-E0	00F-4760-E0	00G-4760-E0	AJ0-7601
<b>PS C18</b>	00B-4758-E0	00D-4758-E0	00F-4758-E0	00G-4758-E0	AJ0-7606

for ID: 3.2-8.0 mm

5 µm Minibore Columns (mm)					SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	4 x 2.0*
<b>Polar C18</b>	00A-4754-AN	00B-4754-AN	00D-4754-AN	00F-4754-AN	AJ0-7600
<b>PS C18</b>	00A-4753-AN	00B-4753-AN	00D-4753-AN	00F-4753-AN	AJ0-7605

for ID: 2.0 - 3.0 mm

5 µm MidBore <sup>™</sup> Columns (mm)				SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
<b>Polar C18</b>	00B-4754-Y0	00D-4754-Y0	00F-4754-Y0	AJ0-7600
<b>PS C18</b>	00B-4753-Y0	00D-4753-Y0	00F-4753-Y0	AJ0-7605

for ID: 2.0 - 3.0 mm

5 µm Analytical Columns (mm)					SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
<b>Polar C18</b>	00B-4754-E0	00D-4754-E0	00F-4754-E0	00G-4754-E0	AJ0-7601
<b>PS C18</b>	00B-4753-E0	00D-4753-E0	00F-4753-E0	00G-4753-E0	AJ0-7606

for ID: 3.2-8.0 mm

5 µm Axia <sup>™</sup> Packed Preparative Columns (mm)			SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	150 x 21.2	250 x 21.2	15 x 21.2**
<b>Polar C18</b>	00F-4754-P0-AX	00G-4754-P0-AX	AJ0-7603
<b>PS C18</b>	00F-4753-P0-AX	00G-4753-P0-AX	AJ0-7608

for ID: 18-29 mm

5 µm Axia <sup>™</sup> Packed Preparative Columns (mm)				SecurityGuard <sup>™</sup> Cartridges (mm)
Phases	150 x 30	250 x 30	250 x 50	15 x 30.0♦
<b>Polar C18</b>	00F-4754-U0-AX	00G-4754-U0-AX	00G-4754-V0-AX	AJ0-7604
<b>PS C18</b>	00F-4753-U0-AX	00G-4753-U0-AX	00G-4753-V0-AX	AJ0-7609

for ID: 30-49 mm

‡ SecurityGuard ULTRA Cartridges require holder, Part No.: AJ0-9000

\* SecurityGuard Analytical Cartridges require holder, Part No.: KJ0-4282

\*\* PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8223

♦ PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8277



# APPLICATIONS

## Australia

t: +61 (0)2-9428-6444  
f: +61 (0)2-9428-6445  
auinfo@phenomenex.com

## Austria

t: +43 (0)1-319-1301  
f: +43 (0)1-319-1300  
anfrage@phenomenex.com

## Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: +1 (800) 543-3681  
f: +1 (310) 328-7768  
info@phenomenex.com

## China

t: +86 400-606-8099  
f: +86 (0)22 2532-1033  
phen@agela.com

## Denmark

t: +45 4824 8048  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: +358 (0)9 4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: +33 (0)1 30 09 21 10  
f: +33 (0)1 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: +49 (0)6021-58830-0  
f: +49 (0)6021-58830-11  
anfrage@phenomenex.com

## India

t: +91 (0)40-3012 2400  
f: +91 (0)40-3012 2411  
indiainfo@phenomenex.com

## Ireland

t: +353 (0)1 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: +39 051 6327511  
f: +39 051 6327555  
italiainfo@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: 01-800-844-5226  
f: 001-310-328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## New Zealand

t: +64 (0)9-4780951  
f: +64 (0)9-4780952  
nzinfo@phenomenex.com

## Norway

t: +47 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: +1 (800) 541-HPLC  
f: +1 (310) 328-7768  
info@phenomenex.com

## Spain

t: +34 91-413-8613  
f: +34 91-413-2290  
espinfo@phenomenex.com

## Sweden

t: +46 (0)8 611 6950  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## United Kingdom

t: +44 (0)1625-501367  
f: +44 (0)1625-501796  
ukinfo@phenomenex.com

## USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com

## All other countries Corporate Office USA

t: +1 (310) 212-0555  
f: +1 (310) 328-7768  
info@phenomenex.com



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Luna is a registered trademark and SecurityGuard, Axia, and MidBore are trademarks of Phenomenex. ACQUITY and Waters are registered trademarks of Waters Technologies Corporation. API 4000 is a trademark of AB SCIEX Pte. Ltd. AB SCIEX is being used under license.

### Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Waters Technologies Corporation.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6, 162, 362. CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP or ULTRA holders, or to any cartridges.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 647, 383.

© 2017 Phenomenex, Inc. All rights reserved.