

APPLICATIONS

Determination of Average DAR of an Antibody Drug Conjugate (ADC) Mimic using a bioZen™ 1.8 µm SEC-2 Column

Chris Nortcliffe Ph.D.¹ and Helen Whitby Ph.D.² ¹SCIEX UK, Warrington ²Phenomenex, Ltd., Queens Avenue, Hurdsfield Ind. Est., Macclesfield, Cheshire SK10 2BN UK

Overview

Antibody drug conjugates (ADC) are a rapidly expanding field in the biopharma industry using the specificity of antibodies to deliver cytotoxic payloads. Their appeal lies in the ability to target and kill tumour cells without compromising healthy ones. A stable link between the antibody and the drug is crucial and drugs linked to antibodies through cysteine residues provide a very reliable and reproducible strategy for achieving an even Drug Antibody Ratio (DAR), something which is essential to ensure good immunogenicity and low patient toxicity. Cysteine linked ADCs provide a specific challenge as the disulphide bonds used to hold the protein together are broken making the complexes unable to withstand the conditions of typical reversed phase LC-MS. The protein in this form is held together only through noncovalent interactions and dissociates under reversed phase conditions. Other analytical approaches to these compounds include subunit analysis of ADCs or HIC, however each of these also have their own analytical shortcomings. In this application note we demonstrate that if kept under native conditions it is possible to get Average DAR of a cysteine conjugate. This application note serves to highlight how native MS can be used to get intact mass for non-covalently linked large molecules, including cysteine linked ADCs.

We analyzed a commercially available cysteine linked ADC mimic obtained from Sigma and were able to identify DAR 0 through 8 using the bioZen SEC-2 in tandem with the SCIEX[®] X500B.

Using XIC (**Figure 1**) each drug antibody proteoform was used to calculate DAR. The average DAR was calculated to be 3.4 which corresponded to the reported Sigma value of $4.0 \pm$ 0.8. We see under 'native' conditions (**Figure 2 and 3**) the ADC remains intact allowing distinct DAR species to be observed and upon reconstruction DAR and glycoforms can be labelled.

LC Conditions

Column: bioZen 1.8 μm SEC-2 Dimension: 150 x 4.6 mm Part No.: 00F-4769-E0 Recommended Guard: SecurityGuard™ ULTRA Guard Cartridge Part No.: AJ0-9850 Guard Holder Part No.: AJ0-9000 Mobile Phase: 100 mM Ammonium Acetate Flow Rate: 200 μL/min Temperature: 25 °C Detector: QTOF (SCIEX X500B) Sample: Sigma ADC Mimic (MSQC8), 100 μg



Figure 2







Hav Visit



APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on <u>www.phenomenex.com</u>, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal <u>www.phenomenex.com/LiveChat</u>.

Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch) beinfo@phenomenex.com

Canada t: +1 (800) 543-3681 info@phenomenex.com

China t: +86 400-606-8099 cninfo@phenomenex.com

Denmark t: +45 4824 8048 nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

Germany t: +49 (0)6021-58830-0 anfrage@phenomenex.com

www.phenomenex.com

India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Ireland t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com

Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com Singapore t: +65 800-852-3944 sginfo@phenomenex.com

Spain t: +34 91-413-8613 espinfo@phenomenex.com

Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com

USA t: +1 (310) 212-0555 info@phenomenex.com

All other countries Corporate Office USA t: +1 (310) 212-0555 info@phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

bioZen and SecurityGuard are trademarks of Phenomenex. SCIEX is a registered trademark of AB Sciex Pte. Ltd. AB SCIEX™ is being used under license FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2019 Phenomenex, Inc. All rights reserved.