


# APPLICATIONS

## Limit of Free Salicylic Acid in Aspirin Tablets Using the Kinetex<sup>®</sup> 5 µm C18 Under USP Allowable Adjustments

Zeshan Aqeel, J Preston, Ryan Splitstone, and Phil Koerner  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



### Overview

In this application is the robust separation of Aspirin (Acetylsalicylic Acid) and related impurity Salicylic Acid under the USP conditions. The application demonstrates the potential method improvements that can be achieved per the allowable adjustments outlined in USP General Chapter <621> relative to the original column and conditions referenced in the monograph.

Monograph: Acetylsalicylic Acid  
pK<sub>a</sub>: 3.5  
LogP: 1.2

### USP Monograph: Aspirin Tablet, Limit of Free Salicylic Acid

**Standard Stock Solution** 0.015 mg/mL of USP Salicylic Acid RS in *Diluent*

**Standard Solution** 0.015 mg/mL of USP Salicylic Acid RS and 0.5 mg/mL of USP Aspirin RS in *Diluent*

**Diluent** Mixture of Acetonitrile and Formic Acid (99:1)

### Column

**Size** Method 1: 300 x 3.9 mm, Method 2: 150 x 4.6 mm

**Stationary Phase** Method 1: µBondapak<sup>®</sup> 10 µm C18, Method 2: Kinetex 5 µm C18

**Temperature** 30 °C

**Mobile Phase** 2 g/L of Sodium 1-heptanesulfonate in a mixture of Acetonitrile and water (15:85).  
Adjust with glacial Acetic Acid to a pH of 3.4.

**Isocratic** Premixed

**Flow Rate** 2.0 mL/min

**Detector** UV @ 280 nm

**Injection Volume** 10 µL of System Suitability solution and Standard solution

### System Suitability – System Suitability solution

Sample: System Suitability solution:

- Resolution (Rs): NLT 2.0 between Salicylic Acid and Aspirin for System Suitability solution
- Relative Standard Deviation: NMT 4.0 % for Standard solution (5 replicate injections)

\* Retention times are provided for information only and are not mandatory. The relative retention times for Salicylic Acid and Aspirin are about 0.7 and 1.0, respectively.


# APPLICATIONS

## Method 1

Original Method: System Suitability solution for Aspirin

**Column:**  $\mu$ Bondapak 10  $\mu$ m C18 Fully Porous  
**Dimensions:** 300 x 3.9 mm  
**Flow Rate:** 2.0 mL/min  
**Sample:** 1. Salicylic Acid  
2. Aspirin

**Average Rs:** 6.08 min **Pressure:** 276 bar  
**Rs RSD:** 3.62%




App ID: 25862

## Method 2

Alternative Method: System Suitability solution for Aspirin

**Column:** Kinetex<sup>®</sup> 5  $\mu$ m C18 Core-Shell  
**Dimensions:** 150 x 4.6 mm  
**Part No.:** 00F-4601-E0  
**Flow Rate:** 2.0 mL/min  
**Sample:** 1. Salicylic Acid  
2. Aspirin

**Average Rs:** 6.32 **Pressure:** 235  
**Rs RSD:** 1.16 %



App ID: 25861

# APPLICATIONS

## Adjustments for Meeting System Suitability

Method Parameter	Allowed Adjustments (isocratic elution)	Method 1	Method 2
Mobile Phase pH	± 0.2 units	As specified	As specified
Concentration of Salts in Buffer	± 10 %	As specified	As specified
Composition of the Mobile Phase	± 30 % Relative; cannot exceed ± 10 % Absolute adjustment; cannot be reduced to zero	As specified	As specified
Wavelength of Detector	No deviations permitted	280 nm (as specified)	As specified
Injection Volume	Can be adjusted as much as needed; must be consistent with linearity, precision, and detection requirements	10 µL (as specified)	10 µL (Allowed)
Column Temperature	± 10 °C	30 °C (Allowed)	30 °C (Allowed)
Stationary Phase	No change of the identity of the substituent permitted (e.g. no replacement of C18 by C8)	L1 (as specified)	As specified
Column Length	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	300 mm (as specified)	150 mm (Allowed)
Column Internal Diameter	Can be adjusted so long as linear velocity is maintained	4.0 (as specified)	4.6 mm (Allowed)
Particle Size	Column length (L) to particle size diameter (dp) ratio can be adjusted between -25 % and +50 %*	10 µm (as specified)	5 µm (Allowed)
Flow Rate	± 50 % (at given ID)	2.0 mL/min (as specified)	2.0 mL/min**

\*Alternatively (as for the application of particle size adjustment to superficially porous particles), other L/dp combinations can be used provided that the number of theoretical plates (N) is within -25% to +50%.

\*\*To maintain equivalent linear velocity for 4.6 mm ID column, flow rate should be adjusted to 2.64 mL/min; the flow rate of 2.0 mL/min is within the allowed adjustment of ± 50 %.

## Allowable Column Adjustments: L/dp Ratio -25 % to 50 %

Column	Length (mm)	ID (mm)	dp (µm)	L/dp	Allowable Range
Original	300	4.0	10	30,000	22,500 – 45,000
Alternative	150	4.6	5	30,000	ALLOWED

## Method Summary and Comparison

	Method 1	Method 2
Column	µBondpak® 10 µm C18	Kinetex® 5 µm C18
Aspirin Average Rt (min)	8.5	5.5
Salicylic Acid and Aspirin Average Rs	6.08	10.05
System Suitability solution Rs RSD (n=5)	3.24 %	0.89 %
Backpressure (Bar)	276	235

# APPLICATIONS

## Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/ChatNow](http://www.phenomenex.com/ChatNow).

### Australia

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

### Austria

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

### Belgium

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

### Canada

t: +1 (800) 543-3681  
info@phenomenex.com

### China

t: +86 400-606-8099  
cninfo@phenomenex.com

### Denmark

t: +45 4824 8048  
nordicinfo@phenomenex.com

### Finland

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

### France

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

### Germany

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

### India

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

### Ireland

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

### Italy

t: +39 051 6327511  
italiainfo@phenomenex.com

### Luxembourg

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### Mexico

t: 01-800-844-5226  
tecnicomx@phenomenex.com

### The Netherlands

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

### New Zealand

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

### Norway

t: +47 810 02 005  
nordicinfo@phenomenex.com

### Poland

t: +48 (12) 881 0121  
pl-info@phenomenex.com

### Portugal

t: +351 221 450 488  
ptinfo@phenomenex.com

### Singapore

t: +65 800-852-3944  
sginfo@phenomenex.com

### Spain

t: +34 91-413-8613  
espinfo@phenomenex.com

### Sweden

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

### Switzerland

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

### Taiwan

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

### United Kingdom

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

### USA

t: +1 (310) 212-0555  
info@phenomenex.com

● All other countries/regions  
Corporate Office USA  
t: +1 (310) 212-0555  
info@phenomenex.com

## [www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™**  
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

### Trademarks

Kinetex are registered trademarks and BE-HAPPY is a trademarks of Phenomenex. µBondpack is a registered trademark of Waters Technologies Corporation.

### Disclaimer

Phenomenex is not affiliated with Waters Corporation. Chromatograms may not be representative of all separations.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.