

APPLICATIONS

18 Cannabinoids for Potency Testing By LC-UV

Scott Krepich, Sean Orlowicz, Zeshan Aqeel, and Jeff Layne
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview


Legalization of medical and recreational marijuana and hemp have lead to an expansion of potency testing needs throughout the globe. Routine testing workflows are now benefiting from resolution of additional exotic cannabinoids to keep up with the diverse range of cultivars and formulations reaching the market.

Here we demonstrate a practical chromatographic separation of 18 cannabinoids with flexibility across HPLC and UHPLC platforms in a simple, isocratic workflow. While 5 primary cannabinoids, Δ^9 -THC, CBD, THCA, CBDA, and CBN are the most common and abundant, new indications are continuously being explored with dozens of other cannabinoids.

LC-UV Conditions

Column: Kinetex[®] 2.6 μ m C18
Dimension: 150 x 4.6mm
Part No.: [00F-4462-E0](#)
Mobile Phase: A: 20 mM Ammonium Formate, pH 2.9 with Formic Acid
B: Acetonitrile
Isocratic: Isocratic 24:76 (A/B)
Flow Rate: 1.5 mL/min
Injection Volume: 2 μ L
Back Pressure: ~260 bar
Temperature: 40 $^{\circ}$ C
Detection: UV @ 228 nm
Analytes: 18 Cannabinoids

Retention Time	Abbrev	Full Name
1.731	CBDVA	Cannabidivarinic acid
1.963	CBDV	Cannabidivarin
2.337	CBDA	Cannabidiolic acid
2.508	CBGA	Cannabigerolic acid
2.704	CBG	Cannabigerol
2.813	CBD	Cannabidiol
2.991	THCV	Tetrahydrocannabivarin
3.321	THCVA	Tetrahydrocannabivarinic acid
3.672	CBNA	Cannabinolic acid
3.968	CBN	Cannabinol
4.753	EXO-THC	Exo-tetrahydrocannabinol
4.907	D9-THC	Δ^9 -Tetrahydrocannabinol
5.084	D8-THC	Δ^8 -Tetrahydrocannabinol
5.670	THCA-A	Tetrahydrocannabinolic acid A
5.819	CBCA	Cannabichromenic acid
6.004	CBL	Cannabicycol
6.173	CBC	Cannabichromene
6.615	CBLA	Cannabicyclolic acid



App ID 25749

APPLICATIONS

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 800-852-3944
sginfo@phenomenex.com

Spain

t: +34 91-413-8613
esinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
info@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
guarantee

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.
www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Kinetex is a registered trademark of Phenomenex.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2020 Phenomenex, Inc. All rights reserved.