

## APPLICATIONS

### bioZen<sup>™</sup> 3 µm NX-C18 Fractionation Column has Improved Loading Capacity in Comparison to Competitor Brand

Dr. Jason Anspach, Roxana Eggleston-Rangel, and Dr. Bryan Tackett Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

### **Overview**

High pH fractionation is a strategy commonly used to increase proteome coverage and its main purpose is to maximize protein identification in complex sample matrices. In proteomics it is often the situation where the amount of sample that is available for analysis is limited. This presents the analyst with a challenge to successfully collect all the information needed to advance their research. A powerful tool that can, and often is, used to increase the effectiveness of the MS efficiency is to miniaturize the LC separation. While the mass of the sample available is often small in proteomics applications, the volume that the sample is present in typically is not. The injection of a large volume on a small ID column directly is a very timeconsuming process. Thus, in this technical note we explore the loading capacity of the bioZen 3 µm NX-C18 fractionation column and how it compares to the competitor column.

A mixture of Angiotensin III Peptide (RVYIHPF) was injected at a range of masses on the columns. The injection mass range chosen was well within the limits of detection for this method (Figure 1). The peak area response to the mass load on the column was monitored to investigate the loadability. A non-linear response of peak area to mass would indicate that the loadability of the trap had been exceeded. None of the masses reached the capacity threshold for the bioZen NX-C18 fractionation column as the peak area was consistently increasing, but the Waters XBridge column did not reach the same amount of angiotensin before reaching the same peak area (Figure 2). This data suggests that the bioZen NX-C18 fractionation column has a 50 % higher loading capacity than the Waters XBridge column.

### LC-MS/MS Conditions

-			
Column:	bioZen 3 μm NX-C18 Fractionation		
	Waters <sup>®</sup> XBridge <sup>®</sup> 3.5 µm BEH C18		
Dimension:	150 x 2.1 mm		
	150 x 2.1 mm		
Part No.:	<u>00F-4793-AN</u> (bioZen)		
Pressure (bar):	220 bar		
Mobile Phase:	A: 10 mM Ammonium Formate pH 10 in Water		
	B: 10 mM Ammonium Formate pH 10 in 90 %		
	Acetonitrile, 10 % Water		
Gradient:	Time (min)	%A	%В
	0	50	50
	7	50	50
Flow Rate:	300 μL/min		
Temperature:	Ambient		
LC System:	EASY-nLC™ 1200		
Detection:	UV @ 220 nm		
Injection Volume:	10 μL		



Figure 1. Chromatograms of loading capacity of (A) bioZen<sup>™</sup> NX-C18 Fractionation and (B) Waters<sup>®</sup> XBridge<sup>®</sup> BEH C18.

**\_ICATIONS** 





Have questions or want more details on implementing this method? We would love to help! Visit **www.phenomenex.com/LiveChat** to get in touch with one of our Technical Specialists.

### APPLICATIONS





**Ophenomenex**<sup>®</sup>

...breaking with tradition™

В





Have questions or want more details on implementing this method? We would love to help! Visit **www.phenomenex.com/LiveChat** to get in touch with one of our Technical Specialists.



# PPLICATIONS

#### Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/LiveChat.

Australia t: +61 (0)2-9428-6444 auinfo@phenomenex.com

Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com

**Belgium** t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch)

beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681 info@phenomenex.com

China t: +86 400-606-8099 cninfo@phenomenex.com

Czech Republic t: +420 272 017 077 cz-info@phenomenex.com

Denmark t: +45 4824 8048 nordicinfo@phenomenex.com

Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com

**Germany** t: +49 (0)6021-58830-0 anfrage@phenomenex.com

#### www.phenomenex.com

India t: +91 (0)40-3012 2400 indiainfo@phenomenex.com

Ireland t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

Mexico t: 01-800-844-5226 tecnicomx@phenomenex.com

The Netherlands t: +31 (0)30-2418700 nlinfo@phenomenex.com

New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com

Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Poland t: +48 22 104 21 72 pl-info@phenomenex.com

**Portugal** t: +351 221 450 488 ptinfo@phenomenex.com

Singapore t: +65 800-852-3944 sginfo@phenomenex.com

**Slovakia** t: +420 272 017 077 sk-info@phenomenex.com

Spain t: +34 91-413-8613 espinfo@phenomenex.com

Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland t: +41 (0)61 692 20 20 swissinfo@phenomenex.com

Taiwan t: +886 (0) 0801-49-1246 twinfo@phenomenex.com

Thailand t: +66 (0) 2 566 0287

thaiinfo@phenomenex.com

United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com

USA t: +1 (310) 212-0555 info@phenomenex.com

All other countries/regions Corporate Office USA t: +1 (310) 212-0555

info@phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com



Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right. www.phenomenex.com/behappy

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions. Trademarks

bioZen and BE-HAPPY are trademarks of Phenomenex. Waters and XBridge are registered trademarks of Waters Technologies Corporation.

EASY-nLC is a trademark of Thermo Fisher Scientific. Disclaimer

**Terms and Conditions** 

Comparative separations may not be representative of all applications. Phenomenex is not affiliated with Waters Corporation.