

TN-1359

Separation of Rivaroxaban and its Organic Impurities per USP Monograph

Alyssa Nodland and Bryan Tackett, PhD

Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

Rivaroxaban is a factor Xa inhibitor, an anticoagulant, that decreases the clotting ability of the blood and helps prevent clots from forming in the blood vessels. In this technical note, we report the separation of Rivaroxaban from its organic impurities per the USP monograph using a Luna™ Omega 3 µm C18 column as an alternative L1 column to a ZORBAX® Eclipse® 3.5 µm XDB-C18 column originally used to elucidate the monograph. The allowable adjustments pertaining to gradient separations are highlighted together with the necessary calculations that are required for the adjusted method conditions to remain compliant with the original monograph.

System suitability per USP Monograph for the Rivaroxaban Assay requires a tailing factor no more than (NMT) 2.0 and a percent relative standard deviation (%RSD) of NMT 0.73 % for six replicate injections. System suitability per USP Monograph for the Rivaroxaban Related Organic Impurities requires resolution no less than (NLT) 8.0 between Rivaroxaban related Compound G and Rivaroxaban, %RSD NMT 5.0 %, and a Signal-to-Noise (S/N) ratio NLT 10

According to USP General Chapter <621>, the configuration of the equipment employed may significantly alter the resolution, retention time, and relative retentions described. Differences in system dwell volume can have an impact on the results obtained for gradient methods. Monographs preferably include an isocratic step before the start of the gradient program so that an adaptation can be made to the gradient time points to take account of differences in dwell volume between the system used for analytical procedure development and that actually used for implementation. The systems used in this study had different dwell volumes (1.5 mL for the Agilent® 1260, and 0.3944 mL for the Waters® ACQUITY® H-Class) than described in the monograph (0.92 mL)*, so the time points stated in the gradient table must be replaced by the adapted time points calculated using the following equation:

$$t_c = t - \frac{(D - D_0)}{F}$$

t_c = adapted time point (min)

t = time point indicated in the monograph (min)

D = Dwell volume (mL)

D_0 = Dwell volume used for development of the method (mL)

F = flow rate (mL/min)

The adapted time points for the ZORBAX Eclipse 3.4 µm XDB-C18 column would be:

Agilent 1260	Waters ACQUITY H-Class
$t_c = 2 - \frac{(1.5 - 0.92)}{1}$	$t_c = 2 - \frac{(0.3944 - 0.92)}{1}$
$t_c = 1.42$	$t_c = 2.53$

Each time point in the gradient would be adjusted by -0.58 minutes for the Agilent 1260 system and +0.53 minutes for the Waters H-Class system and is shown in the gradient table of the LC Conditions on the next page.

Adjustments to column dimensions for gradient methods will be allowed provided that the L/dp ratio remains constant or within the range between -25 % to +50 % of the prescribed L/dp ratio indicated in the monograph. In this monograph, the indicated column length was 150 mm, and the particle size was 3.5 µm; therefore, the Luna Omega 3 µm column used here would be an allowed adjustment. When the particle size is changed, the flow rate requires adjustment because smaller-particle columns will require higher linear velocities to deliver the same performance. The flow rate is adjusted for particle size using the following equation:

$$F_2 = F_1 \times \frac{dc_2^2 \times dp_1}{dc_1^2 \times dp_2}$$

F_1 = flow rate indicated in the monograph (mL/min)


F_2 = adjusted flow rate (mL/min)

dc_1 = internal diameter of the column indicated in the monograph (mm)

dc_2 = internal diameter of the column used (mm)

dp_1 = particle size indicated in the monograph (µm)

dp_2 = particle size of the column used (µm)



The adjusted flow rate for the Luna Omega 3 µm column would be:

$$F_2 = 1 \times \frac{3^2 \times 3.5}{3^2 \times 3.0}$$

$$F_2 = 1.17$$

Adjustments to column dimensions, particle size, or mobile phase volumetric flow rate for gradient methods will impact the slope of the gradient, which can impact selectivity. It is therefore important to adjust the gradient times using the following equation:

$$t_{G2} = t_{G1} \times \left(\frac{F_1}{F_2}\right) \times \left(\frac{L_2 \times dc_1^2}{L_1 \times dc_2^2}\right)$$

t_{G1} = gradient time indicated in the monograph, or adjusted for dwell volume (min)

t_{G2} = adjusted gradient time (min)

F_1 = flow rate indicated in the monograph (mL/min)

F_2 = adjusted flow rate (mL/min)

L_1 = column length indicated in the monograph (mm)

L_2 = new column length (mm)

dc_1 = internal diameter of the column indicated in the monograph (mm)

dc_2 = internal diameter of the column used (mm)

For the second gradient segment, the adjusted gradient time would be:

Agilent 1260	Waters ACQUITY H-Class
$t_{G2} = 1.42 \times \left(\frac{1}{1.17}\right) \times \left(\frac{150 \times 3^2}{150 \times 3^2}\right)$	$t_{G2} = 2.53 \times \left(\frac{1}{1.17}\right) \times \left(\frac{150 \times 3^2}{150 \times 3^2}\right)$
$t_{G2} = 1.22$	$t_{G2} = 2.16$

A gradient adjustment factor can be calculated and used to determine the new gradient segment times using:

$$\text{Gradient adjustment factor} = \left(\frac{t_{G2}}{t_{G1}}\right)$$


Agilent 1260	Waters ACQUITY H-Class
$\text{Gradient adjustment factor} = 0.86$	$\text{Gradient adjustment factor} = 0.85$

The new gradient timetable for the Luna Omega 3 µm C18 column is shown in the gradient table of the LC Conditions on the next page.

All requirements for System Suitability for Rivaroxaban Assay and Organic Impurities were met by all columns.

All solutions were prepared as indicated in the USP Monograph for Rivaroxaban. USP Rivaroxaban RS (Catalog No. 1604530), USP Rivaroxaban Related Compound B RS (Catalog No. 1604552), USP Rivaroxaban Related Compound D RS (Catalog No. 1604563), USP Rivaroxaban Related Compound G RS (Catalog No. 1604596), and USP Rivaroxaban Related Compound J RS (Catalog No. 1604574) were purchased from USP.

Figure 1. Rivaroxaban Structure



LC Conditions

Column: ZORBAX® Eclipse® 3.5 µm XDB-C18, 150 x 3.0 mm
Luna™ Omega 3 µm C18, 150 x 3.0 mm ([00F-4784-Y0](#))

Mobile Phase: A: Methanol / **Solution A** (5:95, v/v)
B: Acetonitrile

Gradient	Time (min)	ZORBAX Eclipse	Luna Omega	%B	
		Agilent® 1260	Waters® ACQUITY® H-Class		
	0	0	0	24	
	2	1.42	1.22	24	
	8	7.42	6.36	77	
	25	24.42	20.93	24	
	37	36.42	31.22	24	
	37.1	36.52	31.30	98	NLT 7 min of column equilibration with the initial mobile phase conditions is recommended between injections.
	45	44	38.30	98	

Flow Rate: 1.0 mL/min (ZORBAX Eclipse)
1.17 mL/min (Luna Omega)

Injection Volume: 5 µL

Temperature: 60 °C

Detector: UV @ 250 nm


System: Agilent 1260
Waters ACQUITY H-Class UHPLC

Table 1. Preparation of Solutions

Solution	Composition
Solution A	Dissolve 1.36 g of Potassium Dihydrogen Phosphate, 1 g of Sodium Hexane Sulfonate, and 200 µL of Phosphoric Acid in Water. Dilute with Water to 1 L.
Solution B	Dissolve 1.36 g of Potassium Dihydrogen Phosphate and 200 µL of Phosphoric Acid in Water. Dilute with Water to 1 L.
Diluent	Acetonitrile / Solution B (40:60, v/v).
Standard Solution (Assay)	0.5 mg/mL of USP Rivaroxaban RS in Diluent .
System Suitability Solution (Organic Impurities)	0.5 mg/mL of USP Rivaroxaban RS and 0.5 µg/mL each of USP Rivaroxaban Related Compound B RS, USP Rivaroxaban Related Compound D RS, USP Rivaroxaban Related Compound G RS, and USP Rivaroxaban Related Compound J RS in Diluent .
Standard Solution (Organic Impurities)	0.5 µg/mL of USP Rivaroxaban RS in Diluent .
Sensitivity Solution (Organic Impurities)	0.25 µg/mL of USP Rivaroxaban RS in Diluent from the Standard Solution .

Figure 2. Standard Solution - Assay


ZORBAX Eclipse 3.5 µm XDB-C18 Column, Agilent 1260



Peak No.	Analyte	Retention Time (min)	Area	Area %RSD	Symmetry Factor
1	Rivaroxaban	17.00	6387.75	0.23	0.76

N=6 Injections


Luna Omega 3 µm C18 Column, Agilent 1260



Peak No.	Analyte	Retention Time (min)	Area	Area %RSD	Symmetry Factor
1	Rivaroxaban	15.61	5553.95	0.63	0.86

N=6 Injections

Luna Omega 3 µm C18 Column, Waters ACQUITY H-Class




Peak No.	Analyte	Retention Time (min)	Area	Area %RSD	Symmetry Factor
1	Rivaroxaban	16.90	1136825	0.60	0.77

N=6 Injections



Figure 3. System Suitability Solution – Organic Impurities


ZORBAX® Eclipse® 3.5 µm XDB-C18 Column, Agilent® 1260



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Related Compound B	5.40	10.2	-
2	Related Compound D	10.04	11.1	-
3	Related Compound G	15.26	4.5	10.23
4	Rivaroxaban	17.08	5559.4	
5	Related Compound J	27.97	6.7	-

N=6 Injections


Luna™ Omega 3 µm C18 Column, Agilent 1260



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Related Compound B	5.42	7.4	-
2	Related Compound D	9.29	8.1	-
3	Related Compound G	14.10	3.6	11.36
4	Rivaroxaban	15.66	4779.4	
5	Related Compound J	24.63	5.5	-

N=6 Injections

Luna Omega 3 µm C18 Column, Waters® ACQUITY® H-Class




Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Related Compound B	6.80	2102	-
2	Related Compound D	10.08	1444	-
3	Related Compound G	14.87	2011	9.22
4	Rivaroxaban	16.24	1824246	
5	Related Compound J	25.05	2373	-

N=5 Injections




Figure 4. Standard Solution – Organic Impurities

ZORBAX® Eclipse® 3.5 µm XDB-C18 Column, Agilent® 1260



Luna™ Omega 3 µm C18 Column, Agilent 1260



Luna Omega 3 µm C18 Column, Waters® ACQUITY® H-Class






Figure 5. Sensitivity Solution – Organic Impurities


ZORBAX Eclipse 3.5 µm XDB-C18 Column, Agilent 1260



Luna Omega 3 µm C18 Column, Agilent 1260



Luna Omega 3 µm C18 Column, Waters ACQUITY H-Class



Luna™ Omega Ordering Information

3 µm MidBore™ Columns (mm)		SecurityGuard™ Cartridges (mm)		
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*/10pk
Polar C18	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJ0-7600
PS C18	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJ0-7605
C18	00B-4784-Y0	00D-4784-Y0	00F-4784-Y0	AJ0-7611
SUGAR	—	—	00F-4775-Y0	AJ0-4496

for ID: 2.0 – 3.0 mm

*SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☎ All other countries/regions Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™ GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Luna, Midbore, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. Waters and ACQUITY are registered trademarks of Waters Technologies Corporation. Agilent, ZORBAX, and Eclipse are registered trademarks of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in no way affiliated with Waters Technologies Corporation or Agilent Technologies, Inc.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists

