

TN-0166

Microelution SPE of Anticonvulsants from Serum for Fast, Sustainable Sample Preparation for Analysis by LC-MS/MS

Shahana Wahab Huq, Stephanie Marin, PhD, and Bryan Tackett, PhD
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Introduction

SPE requires more method development and processing time, but by optimizing the pretreatment, wash, and elution conditions to remove matrix interferences and maximize recovery of target analytes, it provides the cleanest samples for LC-MS/MS analysis. Microelution SPE requires less solvent for the wash and elution steps. Typical wash volumes are 100-200 µL and typical elution solvent volumes are 50-100 µL. The low volume of elution solvent means extracted samples can often be simply diluted prior to analysis, skipping a time-consuming evaporation step that requires specialized equipment. Microelution SPE also uses less sample volume and reagents, reducing solvent usage and waste. These advantages can be used to create faster, more sustainable workflows.

First, a Strata™-X Method Development 96-well plate was used to determine the best sorbent and extraction conditions for a panel of Anticonvulsant drugs and metabolites. The Strata-X-CW sorbent, under acidic load and wash conditions and basic elution conditions, provided the maximum % absolute recovery (TN-0163). Previously, we demonstrated that the Luna Omega 3 µm Polar C18 column was the better option for LC separation of Anticonvulsant analytes prior to LC-MS/MS analysis (TN-1357).

In this technical note, we outline optimal microelution SPE conditions for a panel of 11 Anticonvulsant drug analytes using the Strata-X-CW microelution 96 well plate (Table 1). This was combined with a fast LC method using a Luna Omega 3 µm Polar C18 LC column to resolve all target analytes and determine absolute % recovery and % CV. Percent recovery for the extracted samples was calculated as follows:

$$\% \text{ Recovery} = \left(\frac{\text{Pre-spiked serum analyte}}{\text{Post-spiked serum analyte}} \times 100 \right)$$

Precision was determined as % CV with N=4 replicates.

Sample Preparation

Microelution SPE Conditions

Step	Description
Sample Pretreatment:	10 µL human serum was spiked with Anticonvulsants standard mix and internal standards and then diluted with 200 µL of 25 mM Ammonium Formate, pH ~3.5 adjusted.
Condition:	Strata-X-CW Microelution 96-well plate, 2 mg/well (Part No.: 8M-S035-4GA) with 200 µL of Methanol.
Equilibrate:	200 µL Water.
Load:	200 µL diluted pre-treated sample.
Wash 1:	200 µL of Acidic Buffer 25 mM Ammonium Formate, pH ~3.5 adjusted
Wash 2:	200 µL Methanol / Water (2:8, v/v).
Dry:	1 minute at 20-25 in. Hg.
Elute:	2 aliquots of 50 µL of 5 % Ammonium Hydroxide in Methanol.
Dry Down:	Bypass.
Reconstitution:	Bypass.
Dilute:	Dilute with 200 µL mobile phase A (0.1 % Formic Acid in Water) before injection.

LC Conditions

Column: Luna™ Omega 3 µm Polar C18
Dimensions: 50 x 3.0 mm
Part No.: 00B-4760-Y0
Mobile Phase: A: 2 mM Ammonium Acetate
 B: 2 mM Ammonium Acetate in Methanol

Gradient:	Time (min)	% B
	0	20
	0.5	20
	1.5	40
	2.5	80
	3	95
	3.5	95
	3.51	20
	5	20

Flow Rate: 0.8 mL/min
Injection Volume: 5 µL
Temperature: 40 °C
LC System: Agilent® 1260 Infinity
Detection: MS/MS
Detector: SCIEX® 6500 Triple Quad™

MS/MS Conditions

Ion Source: ESI
Polarity: Dual Polarity
Source Temperature: 450° C
GS1: 55 psi
GS2: 60 psi
CUR: 35 psi
IS: +2500 V or -2500V

Table 1. MS Transitions.

Analyte	Q1 Mass (Da)	Q3 Mass (Da)
Pregabalin	160.1	55
Gabapentin	172	137
Levetiracetam	171	126.1
Lamotrigine	256.1	211.1
Felbamate	178.1	117.1
Lacosamide	251.1	108
Zonisamide	212.9	132
Topiramate	338	78
Oxcarbazepine	253	236
Carbamazepine Epoxide	253.1	180.1
Carbamazepine	237.1	194.1



Results and Discussion

A Strata™-X-CW 2 mg/well 96-well microelution plate was employed for extraction of Anticonvulsants in serum under AB conditions (acidic load and basic elution) based on our initial SPE sorbent screening results (utilizing a Strata-X 30 mg/well 96-well SPE Method Development Plate, [TN-0163](#)). The optimized microelution method aimed toward higher recovery of analytes and direct injection of extracted samples that bypasses dry-down and reconstitution. Although the Ethyl Acetate / Isopropanol / Ammonium Hydroxide (7:2:1, v/v/v) elution solvent resulted in good recovery for Anticonvulsants extracted using the Strata-X-C SPE sorbent from the Strata-X Method Development 96-well plate (**Table 2**), it wasn't ideal for direct injection of the extracted samples from the microelution SPE extraction. The elution conditions for the microelution plate was changed to Methanol treated with 5 % Ammonium Hydroxide and showed improved analyte recovery (by 70-80 %) especially for Zwitterionic compounds. As a protic solvent, Methanol was a stronger solvent for polar analytes which increased solubility and resulted in complete elution of Pregabalin and Gabapentin from the weak cation exchange SPE sorbent. The optimized microelution method is shown in **Figure 2** (with representative chromatogram of extracted samples shown in **Figure 1**), compared to the Method Development plate conditions that were previously established.


The 11 anticonvulsants studied here have a wide range of hydrophobicities and acid-base characteristics. The decrease in organic strength in the second wash from 50 to 20 % Methanol in the microelution method significantly improved recovery of the polar and neutral Topiramate ($\log P = 0.12$, $pK_a = 11.09$) Felbamate ($\log P = -0.7$, $pK_a = 3.82$) by 2- to 3-fold. However, the more polar Levetiracetam with a negative $\log P$ of -0.6 and a pK_a of 2.6 couldn't sustain the 20 % organic wash and had poor recovery (<25 %).

Percent recovery for the optimized microelution method at 10 ng/mL was between 80-124 % (**Table 2**) except the polar compounds Topiramate (55 %) and Levetiracetam (<25 %). Percent recovery at 3 different concentrations (10 ng/mL, 40 ng/mL, and 400 ng/mL) was >80 % for all analytes except Levetiracetam at all three concentrations, Oxcarbazepine at 40 ng/mL, and Topiramate at 10 ng/mL. Percent CVs were within ± 20 % for all analytes (**Table 3**). A 100-fold recovery improvement was observed for Topiramate from the lowest concentration (10 ng/mL) to the highest concentration (400 ng/mL). However, this loss of recovery may be an artifact of poor MS signal at the lower concentrations.

Matrix Effect was 34 % to 68 % at 10 ng/mL. It improved to 80 \pm 5 % at 40 ng/mL, except Topiramate (95 %) and Carbamazepine Epoxide (112 %). At 400 ng/mL, Matrix Effect was 100 \pm 10 % except Carbamazepine Epoxide (251 %). Process Efficiency was 33 % to 59 % at 10 ng/mL for all analytes except Carbamazepine Epoxide (85 %). Process Efficiency was 50 \pm 10 % at 40 ng/mL except Oxcarbazepine (62 %), Topiramate (86 %), and Carbamazepine Epoxide (164 %). At 400 ng/mL, Process Efficiency was 100 \pm 20 % except Oxcarbazepine (74 %) and Carbamazepine Epoxide (382 %).

Calibration curves over a 100-fold dynamic range (5 ng/mL to 500 ng/mL) with 1/x weighting demonstrated excellent linearity with R^2 values ≥ 0.993 for all target compounds including the Zwitterions Pregabalin and Gabapentin (**Figure 3**). The only exception Levetiracetam ($R^2 = 0.976$) with <25 % recovery (**Table 2**).

Figure 1. Analysis of Anticonvulsants Extracted from Serum Using Strata-X-CW Microelution Plate, Under Acidic Load and Basic Elution, on a Luna™ Omega 3 μ m Polar C18 Column. Spiked Concentration of Antipsychotics in Serum is 30 ng/mL.



Peak No.	Analyte	Retention Time (min)	Polarity
1	Pregabalin	0.7	Positive
2	Gabapentin	0.8	Positive
3	Levetiracetam	0.88	Positive
4	Zonisamide	1.53	Positive
5	Lacosamide	1.93	Positive
6	Felbamate	2.08	Positive
7	Lamotrigine	2.28	Positive
8	Carbamazepine Epoxide	2.5	Positive
9	Topiramate	2.5	Negative
10	Oxcarbazepine	2.6	Positive
11	Carbamazepine	2.8	Positive

Table 2. Recovery Comparison of Anticonvulsants Extracted from Serum Using the Strata-X-C SPE of the Strata-X Method Development 30 mg/well 96-Well Plate and the Strata-X-CW Microelution 2 mg/well 96-Well Plate.

Analyte	Strata-X-C SPE Sorbent of the Strata-X Method Development 30 mg/well 96-Well Plate *Spiked conc. of serum = 1 ng/mL		Strata-X-CW Microelution 2 mg/well 96-Well Plate *Spiked conc. of serum = 10 ng/mL	
	Elution Solvent: Ethyl Acetate / Isopropanol / Ammonium Hydroxide (7:2:1, v/v/v)		Elution Solvent: 5 % Ammonium Hydroxide in Methanol	
	% Recovery	% CV (N = 4)	% Recovery	% CV (N = 4)
Pregabalin	53	9.0	95	17.3
Carbamazepine	82	3.0	88	12.5
Felbamate	31	17.6	91	11.2
Gabapentin	50	8.9	87	13.6
Lacosamide	31	21.6	80	13.8
Lamotrigine	76	7.6	98	17.3
Zonisamide	25	17.1	85	14.3
Oxcarbazepine	37	17.4	97	19.4
Carbamazepine Epoxide	96	2.9	124	15.8
Topiramate	26	54.1	54.9	20.1
Levetiracetam	<25	N/A	<25	N/A

Figure 2. SPE Method Comparison of Anticonvulsants from Serum Using the Strata™-X-C SPE of the Strata-X Method Development 30 mg/well 96-Well Plate (Left) and the Strata-X-CW Microelution 2 mg/well 96-Well Plate (Right).


Step	Description	Step	Description
Sample Pretreatment:	500 µL human serum was spiked with Anticonvulsants standard mix at a concentration of 1 ng/mL then diluted with 1 mL of 25 mM Ammonium Formate, pH 4-5 adjusted.	Sample Pretreatment:	10 µL human serum was spiked with Anticonvulsants standard mix and internal standards and then diluted with 200 µL of 25 mM Ammonium Formate, pH ~3.5 adjusted.
Condition:	Strata-X Method Development 96-well plate, 30 mg/well (Part No.: KSO-8209) with 1 mL of Methanol.	Condition:	Strata-X-CW Microelution 96-well plate, 2 mg/well (Part No.: 8M-S035-4GA) with 200 µL of Methanol.
Equilibrate:	1 mL Water.	Equilibrate:	200 µL Water.
Load:	About 1.5 mL of pre-treated sample.	Load:	200 µL diluted pre-treated sample.
Wash 1:	1 mL of 25 mM Ammonium Formate, pH 4-5 adjusted.	Wash 1:	200 µL of Acidic Buffer 25 mM Ammonium Formate, pH ~3.5 adjusted
Wash 2:	1 mL Methanol / Water (1:1, v/v).	Wash 2:	200 µL Methanol / Water (2:8, v/v).
Dry:	5-8 minutes at 20-25 in. Hg.	Dry:	1 minute at 20-25 in. Hg.
Elute:	2 aliquots of 300 µL of Ethyl Acetate / Isopropanol / Ammonium Hydroxide (7:2:1, v/v/v). 15-20 minutes at 40 °C under a gentle stream of Nitrogen.	Elute:	2 aliquots of 50 µL of 5 % Ammonium Hydroxide in Methanol.
Dry Down:	Time \approx 25 min Total Dry Time \approx 35 min	Dry Down:	Bypass. Total Dry Time \approx 1 min
Reconstitution:	500 µL on initial mobile phase spiked with 5 ng/mL Internal Standard mix (Carbamazepine-D ₁₀ , Gabapentin-D ₁₀ , Topiramate-D ₁₂). Time \approx 1 min	Reconstitution:	Bypass.
Total Evaporation and Reconstitution Time:	\approx 36 min/plate	Dilute:	Dilute with 200 µL mobile phase A (2 mM Ammonium Acetate) before injection.
Total Reagent Volume:	6.1 mL/well ; 585.6 mL/plate	Total Evaporation and Reconstitution Time:	\approx 1 min/plate
		Total Solvent Volume:	1.3 mL/well ; 124.8 mL/well

Table 3. Method Performance and Qualification Data Utilizing Optimized Conditions and the Strata-X-CW Microelution Plate for the Extraction of Anticonvulsants from Serum.

Analyte	Concentration (ng/mL)	% Recovery	% CV	% Matrix Effect	PE	Precision	Accuracy	Linear regression R ²	Regression Equation
Pregabalin	10	95	17.3	52	49	12.8	102	0.995	y = 0.154 x + 0.0246
	40	98	3	84	50	3.1	101		
	400	105	5.1	98	100	10.8	101		
Carbamazepine	10	88	12.5	54	48	13.8	91	0.997	y = 0.162 x + 0.0656
	40	87	6.6	83	51	7.5	98		
	400	101	5.4	103	104	5.9	103		
Felbamate	10	91	11.2	56	51	5.4	85	0.996	y = 0.39 x + 0.14
	40	88.2	9.6	83	52	9.8	97		
	400	106.2	8.3	107	116	8.2	105		
Gabapentin	10	87	13.6	58	51	15.5	95	0.996	y = 0.454 x + 0.168
	40	94	4.5	85	47	5.8	99		
	400	98.7	8.5	102	98	8.6	101		
Lacosamide	10	80.2	13.8	57	45	5.9	86	0.995	y = 0.498 x + 0.173
	40	86.2	11.7	83	50	11.7	98		
	400	101.7	7	107	112	6.5	104		
Lamotrigine	10	98	17.3	60	59	14.1	95	0.994	y = 0.09 x + 0.040
	40	90	14.3	85	55	14.7	99		
	400	105	8.5	110	119	7.7	104		
Zonisamide	10	85	14.3	58	49	20	78	0.995	y = 0.236 x + 0.0888
	40	89	6.2	82	50	6.3	101		
	400	106	8.2	104	112	7.6	107		
Oxcarbazepine	10	96.5	19.4	34	33	4	85	0.993	y = 0.0987 x + 0.0404
	40	71	11.8	77	62	12	95		
	400	85	8.5	78	74	10.7	104		
Carbamazepine Epoxide	10	124.1	15.8	68	85	14.7	97	0.994	y = 0.502 x - 0.0569
	40	133	11.8	112	164	11.8	103		
	400	152	8.6	251	382	7.8	105		
Topiramate	10	54.9	51	68	37	78	70	0.994	y = 0.126 x + 0.0245
	40	92	31.7	95	86	36	109		
	400	109	22.2	109	106	20	105		
Levetiracetam	10	<25						0.9766	y = 0.124 x + 0.159
	40	<25							
	400	<25							



Figure 3. Calibration Curves for Four Selected Anticonvulsants Over 100-Fold Dynamic Range (5 ng/mL to 500 ng/mL).



Conclusions

The method using the Strata™-X-CW SPE Microelution plate resulted in a simple, rapid extraction for identification and quantitation of 10 Anticonvulsants from human serum. The method reduces time and solvent usage and can be automated for a high throughput, more sustainable sample preparation workflow for LC-MS/MS analysis.



SPE Ordering Information

Strata™-X Method Development 96-Well Plate

Part No.	Description	Unit
K50-8209	Strata-X, -X-C, -X-CW, and -X-AW 30 mg/well each	ea

Strata-X Microelution 96-Well Plates (ea)

Phase	2 mg
Strata-X-AW	8M-S038-4GA
Strata-X-A	8M-S123-4GA
Strata-X	8M-S100-4GA
Strata-X-C	8M-S029-4GA
Strata-X-CW	8M-S035-4GA

Luna™ Omega Ordering Information

3 µm MidBore™ Columns (mm)		SecurityGuard™ Cartridges (mm)		
Phases	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*/10pk
Polar C18	00B-4760-Y0	00D-4760-Y0	00F-4760-Y0	AJ0-7600
PS C18	00B-4758-Y0	00D-4758-Y0	00F-4758-Y0	AJ0-7605
C18	00B-4784-Y0	00D-4784-Y0	00F-4784-Y0	AJ0-7611
SUGAR	—	—	00F-4775-Y0	AJ0-4496

for ID: 2.0 – 3.0 mm

*SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

☎ All other countries/regions Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™ GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Strata, Luna, MidBore, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX is a registered trademark and Triple Quad is a trademark of AB SCIEX Pte. Ltd.

Disclaimer

Phenomenex is in no way affiliated with Agilent Technologies, Inc.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists