

TN-1361

Meeting System Suitability for USP Abacavir and Lamivudine Tablets Assay

Alyssa Nodland and Bryan Tackett, PhD
Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Introduction

Abacavir and Lamivudine in combination are used with other medications to treat Human Immunodeficiency Virus (HIV) infection. In this technical note, we report the separation of Lamivudine from its related chiral and diastereomeric impurities in Abacavir and Lamivudine tablets per the USP monograph using a Luna™ Omega 3 µm C18 column and a Kinetex™ 2.6 µm EVO C18 as alternative L1 columns to a XTERRA® 3.5 µm MS C18 column originally used to elucidate the monograph. The allowable adjustments pertaining to gradient separations are highlighted together with the necessary calculations that are required for the adjusted method conditions to remain compliant with the original monograph.

System suitability per USP Monograph for the Abacavir and Lamivudine Tablets Assay requires resolution no less than (NLT) 1.0 between Lamivudine-S-Oxide and Lamivudine-R-Oxide and NLT 1.0 between Lamivudine diastereomer and Lamivudine, and a percent relative standard deviation (%RSD) of no more than (NMT) 1.5 % each for Abacavir and Lamivudine with five replicate injections.

According to USP General Chapter <621>, adjustments to column dimensions for gradient methods will be allowed provided that the L/dp ratio remains constant or within the range between -25 % to +50 % of the prescribed L/dp ratio indicated in the monograph. In this monograph, the indicated column length was 150 mm, and the particle size was 3.5 µm; therefore, the Luna Omega 3 µm C18 column and the Kinetex 2.6 µm EVO C18 column used here would be an allowed adjustment. When the particle size is changed, the flow rate requires adjustment because smaller-particle columns will require higher linear velocities to deliver the same performance. The flow rate is adjusted for particle size using the following equation:

$$F_2 = F_1 \times \frac{dc_2^2 \times dp_1}{dc_1^2 \times dp_2}$$

F_1 = flow rate indicated in the monograph (mL/min)
 F_2 = adjusted flow rate (mL/min)
 dc_1 = internal diameter of the column indicated in the monograph (mm)
 dc_2 = internal diameter of the column used (mm)
 dp_1 = particle size indicated in the monograph (µm)
 dp_2 = particle size of the column used (µm)

The adjusted flow rate for the columns would be:

Luna Omega C18	Kinetex EVO C18
$F_2 = 1.5 \times \frac{4.6^2 \times 3.5}{4.6^2 \times 3.0}$	$F_2 = 1.5 \times \frac{4.6^2 \times 3.5}{4.6^2 \times 2.6}$
$F_2 = 1.75$	$F_2 = 2.02$

Adjustments to column dimensions, particle size, or mobile phase volumetric flow rate for gradient methods will impact the slope of the gradient, which can impact selectivity. It is therefore important to adjust the gradient times using the following equation:

$$t_{G2} = t_{G1} \times \left(\frac{F_1}{F_2}\right) \times \left(\frac{L_2 \times dc_2^2}{L_1 \times dc_1^2}\right)$$

t_{G1} = gradient time indicated in the monograph, or adjusted for dwell volume (min)
 t_{G2} = adjusted gradient time (min)
 F_1 = flow rate indicated in the monograph (mL/min)
 F_2 = adjusted flow rate (mL/min)
 L_1 = column length indicated in the monograph (mm)
 L_2 = new column length (mm)
 dc_1 = internal diameter of the column indicated in the monograph (mm)
 dc_2 = internal diameter of the column used (mm)

For the second gradient segment, the adjusted gradient time would be:

Luna Omega C18	Kinetex EVO C18
$t_{G2} = 4.00 \times \left(\frac{1.5}{1.75}\right) \times \left(\frac{150 \times 4.6^2}{150 \times 4.6^2}\right)$	$t_{G2} = 4.00 \times \left(\frac{1.5}{2.02}\right) \times \left(\frac{100 \times 4.6^2}{150 \times 4.6^2}\right)$
$t_{G2} = 3.43$	$t_{G2} = 1.98$

A gradient adjustment factor can be calculated and used to determine the new gradient segment times using:

$$\text{Gradient adjustment factor} = \left(\frac{t_{G2}}{t_{G1}}\right)$$

Luna Omega C18	Kinetex EVO C18
Gradient adjustment factor = 0.858	Gradient adjustment factor = 0.495


The new gradient timetables for the Luna Omega 3 µm C18 column and the Kinetex 2.6 µm EVO C18 column are shown in the gradient table of the LC Conditions on the next page.

All requirements for System Suitability for Abacavir and Lamivudine tablets Assay were met by all columns; however, the Kinetex EVO C18 column provided significantly reduced retention times and decreased resolution versus the original column. The Luna Omega 3 µm C18 column provided improved resolution between Lamivudine and its diastereomer.

All solutions were prepared as indicated in the USP Monograph for Abacavir and Lamivudine tablets. USP Abacavir Sulfate RS (Catalog No. 1000408), USP Lamivudine RS (Catalog No. 1356836), USP Lamivudine Resolution Mixture C RS (Catalog No. 1356869) were purchased from USP. USP Lamivudine Resolution Mixture C RS is a mixture of Lamivudine and the following impurities (other impurities may also be present):

- Uracil
- Lamivudine-Uracil derivative
- Cytosine
- Lamivudine-S-Sulfoxide
- Lamivudine-R-Sulfoxide
- Lamivudine Carboxylic Acid
- Lamivudine Diastereomer
- Salicylic Acid

Figure 1. Abacavir and Lamivudine Structures.



LC Conditions

Column: XTERRA® 3.5 µm MS C18, 150 x 4.6 mm

Luna™ Omega 3 µm C18, 150 x 4.6 mm ([00F-4784-E0](#))Kinetex™ 2.6 µm EVO C18, 100 x 4.6 mm ([00D-4725-E0](#))

Mobile Phase: A: Water / Trifluoroacetic Acid (2000:1, v/v)

B: Acetonitrile / Methanol / Trifluoroacetic Acid (1000:1000:1, v/v/v)

	XTERRA MS C18	Luna Omega C18	Kinetex EVO C18	%B
Gradient:	Time (min)	Adjusted Time (min)	Adjusted Time (min)	
	0	0	0	0
	4	3.43	1.98	0
	12	10.29	5.94	30
	12.1	10.37	5.99	60
	13.1	11.23	6.49	60
	13.2	11.31	6.54	0
	20.2	18.31	13.00	0

NLT 7 min of column equilibration with the initial mobile phase conditions is recommended between injections.

Flow Rate: 1.5 mL/min 1.75 mL/min 2.02 mL/min

Injection Volume: 10 µL

Temperature: 40 °C

Detector: UV @ 270 nm


System: Agilent® 1260

Table 1. Preparation of Solutions

Solution	Composition
Diluent	0.1 N Hydrochloric Acid
Standard Solution	0.35 mg/mL of USP Abacavir Sulfate RS and 0.15 mg/mL of USP Lamivudine RS in Diluent. Sonicate to dissolve prior to final dilution. Pass through a nylon filter (Clarify™ Syringe Filter, Part No.: AF8-7707-12).
System Suitability Solution	Dissolve the contents of one vial of USP Lamivudine Resolution Mixture C RS in 2.5 mL of Diluent. Note: One vial of USP Lamivudine Resolution Mixture C RS contains 0.8 mg of USP Lamivudine Resolution Mixture C RS.

Figure 2. Standard Solution.


XTERRA 3.5 µm MS C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Area %RSD
1	Lamivudine	7.08	2628.84	0.23
2	Abacavir	10.80	2640.18	0.35

N=5 Injections


Luna Omega 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Area %RSD
1	Lamivudine	7.20	2393.20	0.16
2	Abacavir	10.78	2453.76	0.20

N=5 Injections

Kinetex 2.6 µm EVO C18 Column




Peak No.	Analyte	Retention Time (min)	Area	Area %RSD
1	Lamivudine	1.67	2099.52	0.15
2	Abacavir	5.21	2132.2	0.21

N=5 Injections



Figure 3. System Suitability Solution.


XTERRA® 3.5 µm MS C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Lamivudine-S-Oxide	2.38	267.14	4.12
2	Lamivudine-R-Oxide	2.80	154.56	
3	Lamivudine Diastereomer	6.91	124.42	3.59
4	Lamivudine	7.18	2674.94	

N=5 Injections


Luna™ Omega 3 µm C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Lamivudine-S-Oxide	3.06	242.40	4.33
2	Lamivudine-R-Oxide	3.51	141.00	
3	Lamivudine Diastereomer	6.86	113.96	4.30
4	Lamivudine	7.15	2447.48	

N=5 Injections

Kinetex™ 2.6 µm EVO C18 Column



Peak No.	Analyte	Retention Time (min)	Area	Resolution
1	Lamivudine-S-Oxide	0.85	225.36	2.42
2	Lamivudine-R-Oxide	0.95	126.98	
3	Lamivudine Diastereomer	2.30	99.60	2.86
4	Lamivudine	2.66	2149.52	

N=5 Injections



Luna™ Omega Ordering Information

3 µm Analytical Columns (mm)			SecurityGuard™ Cartridges (mm)		
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*/10pk
Polar C18	00B-4760-E0	00D-4760-E0	00F-4760-E0	00G-4760-E0	AJ0-7601
PS C18	00B-4758-E0	00D-4758-E0	00F-4758-E0	00G-4758-E0	AJ0-7606
C18	00B-4784-E0	00D-4784-E0	00F-4784-E0	00G-4784-E0	AJ0-7612
SUGAR	—	00D-4775-E0	00F-4775-E0	00G-4775-E0	AJ0-4495

for ID: 3.2 – 8.0 mm

*SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)

Kinetex™ Ordering Information

2.6 µm Analytical Columns (mm)				SecurityGuard ULTRA Cartridges [†]			
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	00A-4725-E0	00B-4725-E0	—	00D-4725-E0	00F-4725-E0	00G-4725-E0	AJ0-9296
PS C18	00A-4780-E0	00B-4780-E0	—	00D-4780-E0	00F-4780-E0	00G-4780-E0	AJ0-8949
Polar C18	00A-4759-E0	00B-4759-E0	—	00D-4759-E0	00F-4759-E0	—	AJ0-9530
Biphenyl	—	00B-4622-E0	—	00D-4622-E0	00F-4622-E0	—	AJ0-9207
XB-C18	—	00B-4496-E0	00C-4496-E0	00D-4496-E0	00F-4496-E0	—	AJ0-8768
C18	00A-4462-E0	00B-4462-E0	00C-4462-E0	00D-4462-E0	00F-4462-E0	—	AJ0-8768
C8	—	00B-4497-E0	00C-4497-E0	00D-4497-E0	00F-4497-E0	—	AJ0-8770
HILIC	—	00B-4461-E0	00C-4461-E0	00D-4461-E0	00F-4461-E0	—	AJ0-8772
Phenyl-Hexyl	—	00B-4495-E0	00C-4495-E0	00D-4495-E0	00F-4495-E0	—	AJ0-8774
F5	00A-4723-E0	00B-4723-E0	—	00D-4723-E0	00F-4723-E0	—	AJ0-9320

for 4.6 mm ID

*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 3952 5747
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 51 02 180
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: 800-852-3944
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
info@phenomenex.com

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Luna, Kinetex, Clarify, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. XTERRA is a registered trademark of Waters Technologies Corporation. Agilent is a registered trademarks of Agilent Technologies, Inc.

Disclaimer

Comparative separations may not be representative of all applications. Phenomenex is in no way affiliated with Waters Technologies Corporation or Agilent Technologies, Inc.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.



Have questions or want more details on implementing this method? We would love to help! Visit www.phenomenex.com/Chat to get in touch with one of our Technical Specialists

