



## TN-1365

# Identification and Sensitive Quantitation of N-nitroso N-desmethyl Orphenadrine Impurity in Orphenadrine Citrate API

Lakshmanan Deenadayalan<sup>1</sup>, Sashank Pillai<sup>1</sup>, Rahul Baghla<sup>2</sup>, Eshani Nandita, PhD<sup>2</sup>, and Bryan Tackett, PhD<sup>3</sup>

<sup>1</sup>SCIEX Lab, Hitech Defence and Aerospace Park Industrial Area, Mahadeva Kodigehalli, Hobli, Jala Taluka, Bengaluru 562149

<sup>2</sup>AB Sciex LLC, 500 Old Connecticut Path, Framingham, MA 01701, USA

<sup>3</sup>Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501, USA

## Introduction

Nitrosamines are highly potent carcinogens classified into various risk categories (class 1-5). Categorization was performed using the Carcinogenic Potency Categorization Approach (CPCA), where the severity was determined based on the acceptable intake, activating, or deactivating features defined structurally. N-nitroso N-desmethyl Orphenadrine (NNDO) can be formed in Orphenadrine Citrate (OC) because of the presence of the secondary amine, placing it under the class 1 category as per the CPCA approach.

Due to increased concerns posed by the presence of nitrosamines in therapeutics, the EU has set a regulation limit of 18 ng/day. Considering the daily dosage and the regulation limit, N-nitroso N-desmethyl Orphenadrine should be analyzed below the 90 ng/mg limit.

This technical note demonstrates a sensitive method for the identification (**Figure 1**) and sensitive quantitation of N-nitroso N-desmethyl Orphenadrine impurity in Orphenadrine Citrate API using the QTRAP version of the SCIEX 7500 system. A limit of quantitation (LOQ) of 2.5 pg/mL was achieved with baseline separation of N-nitroso N-desmethyl Orphenadrine and Orphenadrine Citrate API (**Figure 2**).

## Sample Preparation

**Standards:** Calibration curve dilutions of N-nitroso N-desmethyl Orphenadrine were prepared across a range of concentrations using a Acetonitrile / Water (50:50, v/v) diluent (2.5, 5, 10, 50, 100, 500, 1000, 5000, and 10000 pg/mL).

**Samples:** A 2 mg of Orphenadrine Citrate API was weighed into a suitable vessel. A 20 mL aliquot of diluent (Acetonitrile / Water (50:50, v/v)) was added and vortexed thoroughly to yield a 0.1 mg/mL concentration. A 100 pg/mL solution of N-nitroso N-desmethyl Orphenadrine was spiked in 0.1 mg/mL of Orphenadrine Citrate API solution to achieve a final concentration of 5 pg/mL.

## LC Conditions

**Column:** Kinetex™ 2.6 μm C8  
**Dimensions:** 100 x 2.1 mm  
**Part No.:** [OOD-4497-AN](#)  
**Mobile Phase:** A: 5 mM Ammonium Formate + 0.1 % Formic Acid in Water  
 B: Acetonitrile

Gradient	Time (min)	% B
	0	40
	2	40
	5	98
	8	98
	8.1	40
	12	40

**Flow Rate:** 0.5 mL/min  
**Injection Volume:** 3 μL  
**Temperature:** 40 °C  
**LC System:** SCIEX® ExionLC™  
**Detection:** MS/MS  
**Detector:** SCIEX QTRAP® 7500

**NOTE:** The LC flow was diverted to waste for the first 2 min to prevent Orphenadrine Citrate API from entering the mass spectrometer and after 5 min during column wash.

## MS/MS Conditions

**Ionization Mode:** ESI  
**Polarity:** Positive  
**Source Temperature:** 300 °C  
**GS1:** 45 psi  
**CUR:** 70 psi  
**CAD:** 8  
**IS:** 1700 V  
**Scan Rate:** 10000 Da/s  
**CE ± CE Spread:** 30 ± 15 V  
**Fixed Fill Time:** 50 ms  
**MS/MS Scan Range:** 50-290 Da

**Table 1.** MRM Transitions.

Analyte	Precursor Ion (m/z)	Fragment Ion (m/z)	CE (V)	CXP (V)	DP (V)
N-nitroso N-desmethyl Orphenadrine-01	285.0	181.09	20	6	35
N-nitroso N-desmethyl Orphenadrine-02	285.0	166.11	45	10	35



**Results and Discussion**

Baseline chromatographic separation was achieved between Orphenadrine Citrate and the N-nitroso N-desmethyl Orphenadrine. The N-nitroso N-desmethyl Orphenadrine retained on the Kinetex™ C8 column much longer with a retention time of 3.8 min, while the Orphenadrine Citrate API eluted at a retention time of 1.9 min (Figure 2).

N-nitroso N-desmethyl Orphenadrine was analyzed across the concentration range of 2.5 to 10000 pg/mL. To evaluate reproducibility, each calibration standard was analyzed in triplicate. Linearity was achieved across concentrations ranging from 2.5 to 10000 pg/mL with a correlation of determination (R<sup>2</sup>) of >0.999 for both quantifier and qualifier ions (Figure 3). An LDR of 3.6 orders of magnitude was achieved. No interference in the diluent blank was observed in the calibration curve samples (Figure 4).

The specification limit (90 pg/mg) was calculated based on the maximum daily dose of 200 mg/day. The N-nitroso N-desmethyl Orphenadrine was analyzed at 5 pg/0.1 mg of API, which is below the calculated specification limit of 90 pg/mg. Recovery was calculated against the neat solution, where the peak area from N-nitroso N-desmethyl Orphenadrine in the control Orphenadrine Citrate API solution was subtracted from the peak area from spiked N-nitroso N-desmethyl Orphenadrine in Orphenadrine Citrate. The average recovery was 90 % with a %CV of 3.8, evaluated in triplicate (Table 2).

Analytical performance was evaluated based on the criteria that the accuracy of the calculated mean should be between 80 % and 120 % at the LOQ and between 85 % and 115 % at the higher concentrations. In addition, the %CV of the calculated mean of the concentration should be <20 % at the LOQ and <15 % at all higher concentrations. The assay accuracy was within ±9 % of the actual concentration and the %CV was <9. Calculated percent accuracy and %CV values were within the acceptance criteria at each concentration level (Figure 5).

In the Orphenadrine Citrate API sample, a peak was observed at the retention time of N-nitroso N-desmethyl Orphenadrine, around 3.8 min (Figure 1). The unknown impurity was identified by comparing its MS/MS spectra with the standard N-nitroso N-desmethyl Orphenadrine standard. Identification was performed using full scan MS/MS experiments with library searching. Data acquisition was performed using the linear ion trap feature of the QTRAP® 7500 system using an MRM to Enhanced Product Ion (EPI) experiment. In this case, the selected MRM transitions for N-nitroso N-desmethyl Orphenadrine were used to create an EPI survey scan in an Information Dependent Acquisition (IDA) experiment setting. Here, the MS/MS spectra from the unknown impurity was compared to standard N-nitroso N-desmethyl Orphenadrine sample spectra. MS/MS spectra matching identified the impurity as N-nitroso N-desmethyl Orphenadrine (Figure 1), but it was below the calculated specification limit of 90 pg/mg.

Figure 1. Detect and Verify Nitrosamine Impurity Using the MRM > EPI Approach on the QTRAP 7500 System.




Figure 2. Good Chromatographic Separation was Achieved Between N-nitroso N-desmethyl Orphenadrine and Orphenadrine Citrate API.






Figure 3. Calibration Curve for Quantitation of N-nitroso N-desmethyl Orphenadrine Quantifier Ion (285 → 181.09) and Qualifier Ion (285 → 166.11). A Weighing Factor of 1/x was Applied for both Calibration Curves.



**Figure 4 .** Representative Extracted Ion Chromatograms of the Diluent (Left) and LOQ, 2.5 pg/mL (Right).



**Table 2 .** Recovery and Precision Calculation.

ID	Peak Area of NNDO at 5 pg/mL Neat Standard	Peak Area of NNDO in OC Control Sample	Peak Area of Spiked 5 pg/mL of NNDO in 0.1 mg of OC API Solution
Sample_01	13405	6548	18164
Sample_02	13349	5760	18439
Sample_03	12807	5824	18999
Mean	13187	6289	18154
SD	269.7	582.4	696
%CV	2.0	9.3	3.8
Peak Area Response at 0.15 ng/mL after Control Correction			11865
Recovery (%)			100.2

**Figure 5 .** Quantitative Performance of N-nitroso N-desmethyl Orphenadrine Quantifier Ion (285 → 181.09) and Qualifier Ion (285 → 166.11).

Row	Component Name	Actual Concentration	Num. Values	Mean	Standard Deviation	Percent CV	Average Accuracy across Replicates
1	NNDO-01	2.50	3 of 3	2.362	0.037	1.57	94.5
2	NNDO-01	5.00	3 of 3	5.407	0.135	2.49	108.
3	NNDO-01	10.00	3 of 3	10.676	0.452	4.23	107.
4	NNDO-01	50.00	3 of 3	47.976	0.969	2.02	96.0
5	NNDO-01	100.00	3 of 3	98.815	2.410	2.44	98.8
6	NNDO-01	500.00	3 of 3	486.642	3.627	0.745	97.3
7	NNDO-01	1000.00	3 of 3	980.718	11.817	1.20	98.1
8	NNDO-01	5000.00	3 of 3	4968.405	51.023	1.03	99.4
9	NNDO-01	10000.00	3 of 3	10107....	138.167	1.37	101.

Row	Component Name	Actual Concentration	Num. Values	Mean	Standard Deviation	Percent CV	Average Accuracy across Replicates
1	NNDO-02	2.50	3 of 3	2.529	0.045	1.79	101.
2	NNDO-02	5.00	3 of 3	4.736	0.421	8.89	94.7
3	NNDO-02	10.00	3 of 3	10.606	0.433	4.09	106.
4	NNDO-02	50.00	3 of 3	49.755	1.413	2.84	99.5
5	NNDO-02	100.00	3 of 3	99.131	2.717	2.74	99.1
6	NNDO-02	500.00	3 of 3	494.118	3.562	0.721	98.8
7	NNDO-02	1000.00	3 of 3	992.681	10.993	1.11	99.3
8	NNDO-02	5000.00	3 of 3	4974.377	86.185	1.73	99.5
9	NNDO-02	10000.00	3 of 3	10182....	107.332	1.05	102.

**Conclusions**

An LOQ of 2.5 pg/mL was achieved for N-nitroso N-desmethyl Orphenadrine in Orphenadrine Citrate API. Linearity was achieved at concentrations ranging from 2.5 pg/mL to 10000 pg/mL with an R<sup>2</sup> >0.999 for both quantifier and qualifier ions covering a linear dynamic range of 3.6 orders of magnitude. An impurity in Orphenadrine Citrate API was identified as N-nitroso N-desmethyl Orphenadrine by comparing the impurity MS/MS spectra with the N-nitroso N-desmethyl Orphenadrine standard MS/MS spectra. Good quantitative performance was demonstrated with accurate and highly reproducible (%CV <9) results on the QTRAP® 7500 system. The method demonstrated the quantitation of N-nitroso N-desmethyl Orphenadrine impurity below the calculated specification limit (90 pg/mg) in the Orphenadrine Citrate API.



## Kinetex™ Ordering Information

Phases	2.6 µm Minibore Columns (mm)			SecurityGuard™ ULTRA Cartridges (mm)*		
	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	<a href="#">00A-4725-AN</a>	<a href="#">00B-4725-AN</a>	—	<a href="#">00D-4725-AN</a>	<a href="#">00F-4725-AN</a>	<a href="#">AJ0-9298</a>
PS C18	<a href="#">00A-4780-AN</a>	<a href="#">00B-4780-AN</a>	—	<a href="#">00D-4780-AN</a>	<a href="#">00F-4780-AN</a>	<a href="#">AJ0-8951</a>
Polar C18	<a href="#">00A-4759-AN</a>	<a href="#">00B-4759-AN</a>	—	<a href="#">00D-4759-AN</a>	<a href="#">00F-4759-AN</a>	<a href="#">AJ0-9532</a>
Biphenyl	<a href="#">00A-4622-AN</a>	<a href="#">00B-4622-AN</a>	—	<a href="#">00D-4622-AN</a>	<a href="#">00F-4622-AN</a>	<a href="#">AJ0-9209</a>
XB-C18	<a href="#">00A-4496-AN</a>	<a href="#">00B-4496-AN</a>	<a href="#">00C-4496-AN</a>	<a href="#">00D-4496-AN</a>	<a href="#">00F-4496-AN</a>	<a href="#">AJ0-8782</a>
C18	<a href="#">00A-4462-AN</a>	<a href="#">00B-4462-AN</a>	<a href="#">00C-4462-AN</a>	<a href="#">00D-4462-AN</a>	<a href="#">00F-4462-AN</a>	<a href="#">AJ0-8782</a>
C8	<a href="#">00A-4497-AN</a>	<a href="#">00B-4497-AN</a>	<a href="#">00C-4497-AN</a>	<a href="#">00D-4497-AN</a>	<a href="#">00F-4497-AN</a>	<a href="#">AJ0-8784</a>
HILIC	<a href="#">00A-4461-AN</a>	<a href="#">00B-4461-AN</a>	<a href="#">00C-4461-AN</a>	<a href="#">00D-4461-AN</a>	<a href="#">00F-4461-AN</a>	<a href="#">AJ0-8786</a>
Phenyl-Hexyl	<a href="#">00A-4495-AN</a>	<a href="#">00B-4495-AN</a>	<a href="#">00C-4495-AN</a>	<a href="#">00D-4495-AN</a>	<a href="#">00F-4495-AN</a>	<a href="#">AJ0-8788</a>
F5	<a href="#">00A-4723-AN</a>	<a href="#">00B-4723-AN</a>	—	<a href="#">00D-4723-AN</a>	<a href="#">00F-4723-AN</a>	<a href="#">AJ0-9322</a>

for 2.1 mm ID

\*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)

**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 3952 5747  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 51 02 180  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™  
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/phx-terms-and-conditions-of-sale](http://www.phenomenex.com/phx-terms-and-conditions-of-sale).

**Trademarks**

Kinetex, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. SCIEX and QTRAP are registered trademarks and ExionLC is a trademark of AB SCIEX Pte. Ltd.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.  
© 2024 Phenomenex, Inc. All rights reserved.

