

TN-1047

Alternative Selectivity of Chiral Stationary Phases Based on Cellulose tris(3-chloro-4-methylphenylcarbamate) and Cellulose tris(3,5-dimethylphenylcarbamate)

Tivadar Farkas, PhD and Bryan Tackett, PhD
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA



Introduction

Polysaccharide-based chiral stationary phases (CSP) are widely used due to their wide chiral recognition ability. Several cellulose and amylose derivatives are extremely effective in the separation of a wide range of compounds of interest in the pharmaceutical industry. This technical note demonstrates the different chiral recognition capabilities of CSPs based on Cellulose tris(3-chloro-4-methylphenylcarbamate) and Cellulose tris(3,5-dimethylphenylcarbamate). Over 180 racemates of pharmaceutical interest were analyzed on these two phases in normal phase (NP), polar-organic (PO) and reversed phase (RP) separation modes. Numerous examples including important classes of drug compounds as well as statistical data demonstrates that Cellulose tris(3-chloro-4-methylphenylcarbamate) offers a good alternative to the commonly used Cellulose tris(3,5-dimethylphenylcarbamate) in the separation of difficult racemic mixtures.

LC Conditions

Column: Lux™ 5 µm Cellulose-1 ([00G-4459-E0](#))
Lux 5 µm Cellulose-2 ([00G-4457-E0](#))
CHIRALCEL® 5 µm OD-H®

Dimensions: 250 x 4.6 mm

Mobile Phase: See **Table 1**

Flow Rate: 1.0 mL/min (Isocratic)

Injection Volume: 5 - 20 µL (depending on analyte response)

Temperature: Ambient

LC System: Agilent® 1100


Detection: UV @ 220 nm

Sample: 500 µg/mL racemate dissolved in mobile phase

Figure 1. Structures of Chiral Selective Phases.

Lux Cellulose-1

Cellulose tris(3,5-dimethylphenylcarbamate).



Lux Cellulose-2

Cellulose tris(3-chloro-4-methylphenylcarbamate).




Table 1. Mobile Phase Compositions.

Mobile Phase	Basic and Neutral Compounds	Acidic and Neutral Compounds
NP	0.1 % Diethylamine in Hexane / Isopropanol	0.1 % Acetic Acid (or Formic Acid) in Hexane / Isopropanol
PO	0.1 % Diethylamine in Methanol / Isopropanol	0.1 % Acetic Acid (or Formic Acid) in Methanol / Isopropanol
	0.1 % Diethylamine in Acetonitrile / Isopropanol	0.1 % Acetic Acid (or Formic Acid) in Acetonitrile / Isopropanol
RP	0.1 % Diethylamine in Acetonitrile / Water	0.1 % Acetic Acid (or Formic Acid) in Acetonitrile / Water
	0.1 % Diethylamine in Methanol / Water	0.1 % Acetic Acid (or Formic Acid) in Methanol / Water



Results and Discussion

Lux™, the Phenomenex brand of polysaccharide-based chiral selective phases, was introduced into the market in 2008. One phase, Lux Cellulose-1 is based on Cellulose tris(3,5-dimethylphenyl carbamate) like other chiral phases on the market (e.g., CHIRALCEL® OD-H®). Another phase, Lux Cellulose-2, is a new member to the family of polysaccharide based chiral selective phases and uses Cellulose tris(3-chloro-4-methylphenyl carbamate) (**Figure 1**) as a chiral selector; this new chemistry delivers a unique selectivity versus other phases.

In this study over 180 diverse compounds of pharmaceutical interest were screened on the Lux line of chiral selective phases as well as other comparative medias to better characterize the selectivity delivered by each Lux phase.

Table 1 summarizes the screening conditions used for each column; different types of mobile phases (NP or PO) as well as additives used (0.1 % Formic Acid or Acetic Acid for acidic analytes or 0.1 % Diethylamine for basic analytes).

Figures 2-5 show several representative examples of the different selectivities provided by Lux Cellulose-1 and Lux Cellulose-2 in chiral separations across normal phase, polar-organic, and reversed phase separation modes. Representative compounds such as various β -blockers, Warfarin, Sulcanazole, Milnacipran, and Clenbuterol demonstrate the complementary behavior of Lux Cellulose-2 to the commonly used cellulose tris(3,5-dimethylphenylcarbamate) based CSPs (CHIRALCEL OD-H and Lux Cellulose-1) in the separation of difficult racemates.

Figure 2 demonstrates the behavior of the two Lux phases in normal phase separations. While Lux Cellulose-1 generally demonstrates slightly better resolution and increased retention versus CHIRALCEL OD-H, there are several cases where Lux Cellulose-2 is a better separation choice when using normal phase. While Bopindolol is equally well separated on the two Lux phases, Oxprenolol enantiomers are better resolved on Lux Cellulose-2. Toliprolol enantiomers are separated on Lux Cellulose-1 with spectacular resolution but at the expense of extensive retention for one of the

enantiomers. Lux Cellulose-2 barely separates racemic Oxprenolol under similar mobile phase conditions, but with minimal optimization (i.e. a reduction of Isopropanol in the mobile phase) a better separation is achieved with Lux Cellulose-2 with significantly shorter analysis time.

Figure 3 shows additional normal phase enantiomeric separations using Lux Cellulose-1, Lux Cellulose-2, and CHIRALCEL OD-H. Such separations further demonstrate the complementary selectivity offered by Lux Cellulose-2 versus the Cellulose tris(3,5-dimethylphenyl carbamate) phases Lux Cellulose-1 and CHIRALCEL OD-H.


Figures 4 and 5 demonstrate the effect of mobile phase composition on chiral resolution. **Figure 4** demonstrates the dramatic changes in selectivity for each phase when traditional normal phase separation is substituted with polar organic separation mode using either Acetonitrile or Methanol as mobile phase. **Figure 5** shows changes in selectivity are observed in reversed phase mode. Different solvents can alter the steric structure of the polysaccharide backbone and the arrangement of binding sites, providing alternative selectivity for separating difficult chiral compounds. Such mobile phase alteration offers a powerful tool in developing and optimizing chiral separations.

Figure 6 compares the success rates of Lux Cellulose-1 and Lux Cellulose-2 in the analysis of over 180 racemates in normal phase or polar organic separation modes. The number of uniquely baseline resolved racemates are given at the bottom of the bar graph. The same selection criteria was applied to partially separated racemates. For example, Lux Cellulose-2 shows good chiral recognition in Acetonitrile mobile phase with 9 baseline separations of racemates that could not be separated on Lux Cellulose-1. This complementary enantioselectivity of Lux Cellulose-2 over Lux Cellulose-1 is most evident in Acetonitrile / Isopropanol mobile phase mixtures and is less pronounced in standard normal phase mixtures (Hexane / Isopropanol) and Methanol mixtures.




Figure 2. Enantioseparations of β -Blockers in Normal Phase.


Toliprolol on Lux™ Cellulose-1
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (80:20, v/v)




Oxprenolol on Lux Cellulose-1
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (90:10, v/v)




Toliprolol on Lux Cellulose-2
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (80:20, v/v)






Oxprenolol on Lux Cellulose-2
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (90:10, v/v)



Toliprolol on Lux Cellulose-2
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (90:10, v/v)



Bopindolol on Lux Cellulose-1
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (90:10, v/v)



Bopindolol on Lux Cellulose-2
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (90:10, v/v)






Figure 3. Enantioseparations in Normal Phase.


Warfarin on Lux™ Cellulose-1
0.1 % Formic Acid in Hexane / 0.1 % Formic Acid in Isopropanol (60:40, v/v)




Sulcanazole on Lux Cellulose-1
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (60:40, v/v)




Warfarin on Lux Cellulose-2
0.1 % Formic Acid in Hexane / 0.1 % Formic Acid in Isopropanol (60:40, v/v)





Sulcanazole on Lux Cellulose-2
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (60:40, v/v)




Warfarin on CHIRALCEL® OD-H®
0.1 % Formic Acid in Hexane / 0.1 % Formic Acid in Isopropanol (60:40, v/v)



Sulcanazole on CHIRALCEL OD-H
0.1 % Diethylamine in Hexane / 0.1 % Diethylamine in Isopropanol (60:40, v/v)



Warfarin



Sulcanazole



Figure 4. Complementary Enantioselectivity in Normal Phase and Polar-Organic.











Figure 4 Cont'd. Complementary Enantioselectivity in Normal Phase and Polar-Organic.

Chlormezanone on Lux™ Cellulose-1

0.1 % Diethylamine in Acetonitrile / 0.1 % Diethylamine in Isopropanol (60:40, v/v)





Figure 5. Enantioseparations in Reversed Phase.


Clenbuterol on Lux Cellulose-2

0.1 % Diethylamine in Methanol / 0.1 % Diethylamine in Water (80:20, v/v)




Chlormezanone on Lux Cellulose-2

0.1 % Diethylamine in Acetonitrile / 0.1 % Diethylamine in Isopropanol (95:5, v/v)




Clenbuterol on Lux Cellulose-2

0.1 % Diethylamine in Acetonitrile / 0.1 % Diethylamine in Water (60:40, v/v)





Chlormezanone on Lux Cellulose-1

0.1 % Diethylamine in Methanol / 0.1 % Diethylamine in Isopropanol (90:10, v/v)




Clenbuterol on Lux Cellulose-2

0.1 % Diethylamine in Acetonitrile / 0.1 % Diethylamine in Water (40:60, v/v)



Chlormezanone




Clenbuterol

Mobile Phase	Rs
Hexane / Isopropanol	0.00
Methanol / Isopropanol	0.00
Acetonitrile	0.00



Figure 6 . Success Rates for over 180 Racemates on Lux™ Cellulose-1 and Cellulose-2.



Lux™ Ordering Information

5 µm Minibore and Analytical Columns (mm)						SecurityGuard™ Cartridges (mm)	
Phases	50 x 2.0	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 2.0*	4 x 3.0*
						/10pk	/10pk
i-Amylose-1	00B-4762-B0	00B-4762-E0	00D-4762-E0	00F-4762-E0	00G-4762-E0	AJ0-8640	AJ0-8641
i-Amylose-3	—	00B-4779-E0	00D-4779-E0	00F-4779-E0	00G-4779-E0	AJ0-8651	AJ0-8650
i-Cellulose-5	—	00B-4756-E0	00D-4756-E0	00F-4756-E0	00G-4756-E0	AJ0-8631	AJ0-8632
Cellulose-1	—	00B-4459-E0	00D-4459-E0	00F-4459-E0	00G-4459-E0	AJ0-8402	AJ0-8403
Cellulose-2	00B-4457-B0	00B-4457-E0	00D-4457-E0	00F-4457-E0	00G-4457-E0	AJ0-8398	AJ0-8366
Cellulose-3	—	00B-4493-E0	00D-4493-E0	00F-4493-E0	00G-4493-E0	AJ0-8621	AJ0-8622
Cellulose-4	—	—	00D-4491-E0	00F-4491-E0	00G-4491-E0	AJ0-8626	AJ0-8627
Amylose-1	00B-4732-B0	—	00D-4732-E0	00F-4732-E0	00G-4732-E0	AJ0-9337	AJ0-9336

*SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)

for ID: 2.0-3.0 mm 3.2-8.0 mm



Need a different column size or sample preparation format?

No problem! We have a majority of our available dimensions up on www.phenomenex.com, but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal www.phenomenex.com/Chat.

Australia

t: +61 (0)2-9428-6444
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
info@phenomenex.com

China

t: +86 400-606-8099
cninfo@phenomenex.com

Czech Republic

t: +420 272 017 077
cz-info@phenomenex.com

Denmark

t: +45 4824 8048
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
anfrage@phenomenex.com

Hong Kong

t: +852 6012 8162
hkinfo@phenomenex.com

India

t: +91 (0)40-3012 2400
indiainfo@phenomenex.com

Indonesia

t: +62 21 5019 9707
indoinfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
italiainfo@phenomenex.com

Japan

t: +81 (0) 120-149-262
jpinfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
nlinfo@phenomenex.com

Mexico

t: 01-800-844-5226
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
nordicinfo@phenomenex.com

Poland

t: +48 22 104 21 72
pl-info@phenomenex.com

Portugal

t: +351 221 450 488
ptinfo@phenomenex.com

Singapore

t: +65 6559 4364
sginfo@phenomenex.com

Slovakia

t: +420 272 017 077
sk-info@phenomenex.com

Spain

t: +34 91-413-8613
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
nordicinfo@phenomenex.com

Switzerland

t: +41 (0)61 692 20 20
swissinfo@phenomenex.com

Taiwan

t: +886 (0) 0801-49-1246
twinfo@phenomenex.com

Thailand

t: +66 (0) 2 566 0287
thaiinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

🌐 **All other countries/regions**
Corporate Office USA

t: +1 (310) 212-0555
www.phenomenex.com/chat

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at international@phenomenex.com

BE-HAPPY™
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

www.phenomenex.com/behappy

Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/phx-terms-and-conditions-of-sale.

Trademarks

Lux, SecurityGuard, and BE-HAPPY are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. CHIRALCEL and OD-H are registered trademarks of DAICEL Corporation.

Disclaimer

Comparative separations may not be representative of all applications.

Phenomenex is in now way affiliated with Agilent Technologies, Inc. or DAICEL Corporation.

SecurityGuard is patented by Phenomenex. U.S. Patent No. 6,162,362.

CAUTION: this patent only applies to the analytical-sized guard cartridge holder, and does not apply to SemiPrep, PREP, or ULTRA holders, or to any cartridges.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2023 Phenomenex, Inc. All rights reserved.

