

Application Note

Sample Clean-Up After N-Linked Glycan Release and Labeling Using bioZen[™] N-Glycan Clean-Up HILIC Solid Phase Extraction (SPE) Microelution 96-Well Plate

Brian Rivera, Jessica Detsch, and Eric Chapa Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Overview

This application note compares two glycan clean-up solutions for a commercially available in-solution digestion and labeling kit for glycan profiling of glycoproteins. bioZen N-Glycan Clean-Up is a HILIC solid phase extraction product in a microelution 96well plate that has excellent retention and recovery of labeled, released n-glycans. In this study, a trastuzumab biosimilar and a heavily sialylated glycoprotein (alpha-1 acid glycoprotein, AGP) were deglycosylated and glycans were labeled using the Waters® GlvcoWorks[™] Rapifluor-MS[™] workflow. To remove excess dye, the samples were extracted using two different HILIC plates: bioZen N-Glycan Clean-Up microelution plate and the Glyco-Works HILIC µElution plate. For the trastuzumab biosimilar, the recovery of key labeled glycans, as determined by peak areas, were then compared. For AGP, a weighted average of charge (Z-value) was determined. The bioZen N-Glycan Clean-Up microelution plate gave comparable results to the GlycoWorks HILIC µElution plate, providing an alternative SPE product in the glycan workflow offered at a cost savings for labs running glycan clean-up and analysis.

Experimental Conditions

Waters GlycoWorks RapiFluor-MS N-glycan profiling kit was obtained from Waters (Waltham, MA, USA).

Deglycosylation and Labeling: Samples were reconstituted to 2 mg/mL in pure water. Deglycosylation and labeling were performed per vendor recommendations. Each sample was loaded in quadruplicate onto both glycan clean-up plates.

Solid Phase Extraction Protocol

96-Well Plate:	Waters GlycoWorks HILIC µElution Plate
	bioZen N-Glycan Clean-Up Microelution Plate
Part No.:	8M-S009-NGA (bioZen)
Condition:	200 µL Water
Equilibrate:	200 µL Water/Acetonitrile (15:85)
Load:	Labeled sample (~400 µL)
Wash:	2x 600 µL Formic acid/Water/Acetonitrile (1:9:90)
Elute:	3x 30 µL 200 mM Ammonium acetate in Acetonitrile/Water (5:95)

Each sample was analyzed using a generic amide LC column on a Waters ACQUITY[®] UPLC iClass equipped with fluorescence detection using standard methods and conditions (data not shown).

Results

The efficiency and peak areas of the released N-glycans from the trastuzumab biosimilar were compared between the bioZen N-Glycan Clean-Up microelution plate and the Waters Glyco-Works HILIC μ Elution Plate. G1 and G1F isomers were summed for simplicity and the data is shown in **Table 1**.

Figure 1.

Average Peak Areas of Labeled Glycans

Labeled Glycan	bioZen N-Glycan Clean-Up Average Peak Area	Waters GlycoWorks HILIC µElution Average Peak Area	% Difference
G0	132898	128611	103%
G0F	1131650	1082913	105%
G1/G1'	59236	56873	104%
G1F/G1F'	35168	332001	106%

Because AGP is composed mainly of multi-antennae, sialylated glycans, Z-values were used for comparison. Z-value takes a weighted average of contributing charge shares for different levels of sialylation and could reveal a bias or difference in recovery of different forms of sialylated glycans.

Figure 2.

Comparison of Z-value determination for AGP

	bioZen N-Glycan Clean-Up Microelution Plate			Waters GlycoWorks µElution Plate		
Peak Group	Percentage of peak group area (%)	Multiplication Factor	Charge Number Share	Percentage of peak group area (%)	Multiplication Factor	Charge Number Share
Asialo	0.3	0	0.0	0.4	0	0.0
Monosialo-	14.4	1	14.4	14.6	1	14.6
Disialo-	1.8	2	3.6	2.0	2	3.9
Trisialo-	57.6	3	172.7	57.8	3	173.3
Tetrasialo-	25.9	4	103.5	24.8	4	99.2
			Z=294.2			Z=291.0

Conclusion

The bioZen N-Glycan Clean-Up microelution plate is an appropriate replacement for the Waters GlycoWorks µElution plates in the Waters GlycoWorks *Rapi*Flour-MS glycan workflow solution.



Application Note

Ordering Information

Sample Preparation

bioZen [™] Solid Phase Extraction	Format	Sorbent Mass	Part Number	Unit
bioZen N-Glycan Clean-Up	Microelution 96-Well Plate	5 mg/well	8M-S009-NGA	1/box

LC Columns

bioZen Columns (mm)	Biocompatible Guard C			
	100 x 2.1	150 x 2.1	for 2.1 mm	Holder
bioZen 2.6 µm Glycan	00D-4773-AN	00F-4773-AN	AJ0-9800	AJ0-9000



Australia

- t: +61 (0)2-9428-6444 auinfo@phenomenex.com
- Austria t: +43 (0)1-319-1301 anfrage@phenomenex.com
- Belaium t: +32 (0)2 503 4015 (French) t: +32 (0)2 511 8666 (Dutch)
- beinfo@phenomenex.com
- **Canada** t: +1 (800) 543-3681 info@phenomenex.com

China t: +86 400-606-8099 cninfo@phenomenex.com

Denmark t: +45 4824 8048 nordicinfo@phenomenex.com

Finland t: +358 (0)9 4789 0063 nordicinfo@phenomenex.com

- France t: +33 (0)1 30 09 21 10 franceinfo@phenomenex.com
- **Germany** t: +49 (0)6021-58830-0

anfrage@phenomenex.com India t: +91 (0)40-3012 2400

indiainfo@phenomenex.com

Ireland t: +353 (0)1 247 5405 eireinfo@phenomenex.com

Italy t: +39 051 6327511 italiainfo@phenomenex.com

Luxembourg t: +31 (0)30-2418700 nlinfo@phenomenex.com

www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Mexico t: 01-800-844-5226

- tecnicomx@phenomenex.com
- The Netherlands t: +31 (0)30-2418700
- nlinfo@phenomenex.com
- New Zealand t: +64 (0)9-4780951 nzinfo@phenomenex.com
- Norway t: +47 810 02 005 nordicinfo@phenomenex.com

Portugal t: +351 221 450 488 ptinfo@phenomenex.com

Singapore t: +65 800-852-3944 sginfo@phenomenex.com

Spain t: +34 91-413-8613 espinfo@phenomenex.com

Sweden t: +46 (0)8 611 6950 nordicinfo@phenomenex.com

Switzerland t: +41 61 692 20 20 swissinfo@phenomenex.com

United Kingdom t: +44 (0)1625-501367 ukinfo@phenomenex.com

USA t: +1 (310) 212-0555 info@phenomenex.com

All other countries Corporate Office USA





If products in this technical note do not perform at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions. Trademarks

bioZen is a trademark of Phenomenex. Waters and ACQUITY are registered trademarks and *Rapi*Fluor-MS and GlycoWorks are trademarks of Waters Technologies Corp.

Disclaimer

Phenomenex is not affiliated with Waters Technologies Corp. Comparative separations may not be representative of all applications.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2018 Phenomenex, Inc. All rights reserved.