

Introducing PFAS Proficiency Testing EPA methods 1633, 533, and 537.1

ISO 17043 accredited Proficiency Testing and Quality Control Standards ready to ship

A New Standard of Confidence to Your PFAS Analysis

Phenova continues to support laboratories conducting PFAS analysis by providing ISO accredited standards in accordance with EPA Methods 533, 537.1 and now 1633.

See complete list of analytes on the next page.

Order Now



Easily Fit Your PT Studies Into Your Lab's Schedule

- Water Pollution
- Water Supply
- Air & Emissions
- Biosolids
- Solids/Hazardous Waste
- Microbiology
- Fuels/UST
- DWMQA



PFAS Proficiency testing studies are available quarterly.
Proficiency Testing studies for lab accreditation and renewals offered regularly.

Book PT Study

Complete List of Analytes for PFAS Proficiency Testing



Components *	WP	WS	HW/Soil
4,8-Dioxa-3H-perfluorononanoic acid (ADONA)	X	X	X
9-chloro hexadecane fluoro-3-oxanone-1-sulfonic acid (9Cl-PF3ONS)	X	X	X
11-CHLOROEICOSAFLUORO-3-OXAUNDECANE-1-SULFONICACID (11Cl-PF30UdS)	X	X	X
4,4,5,5,6,6,6-Heptafluorohexanoic acid (3:3 FTCA)	X		X
2H,2H,3H,3H-Perfluorooctanoic acid (5:3 FTCA)	X		X
2H,2H,3H,3H-Perfluorodecanoic acid (7:3 FTCA)	X		X
1H,1H,2H,2H-Perfluorohexane sulfonic acid (4:2FTS)	X	X	X
1H,1H,2H,2H-Perfluorooctane sulfonic acid (6:2FTS)	X	X	X
1H,1H,2H,2H-Perfluorodecanesulfonic acid (8:2FTS)	X	X	X
Hexafluoropropylene oxide diamer acid (HFPO-DA)	X	X	X
Perfluorodecanesulfonic acid (PFDoS)	X		X
N-Ethyl perfluorooctanesulfonamidoacetic acid (NEtFOSAA)	X	X	X
N-ethylperfluorooctanesulfonamide (NEtFOSA)	X		X
N-ethylperfluorooctanesulfonamido)-ethanol (NEtFOSE)	X		X
Nonafluoro-3,6-dioxaheptanoic acid (NFDHA)	X	X	X
N-Methyl perfluorooctanesulfonamidoacetic acid (NMeFOSAA)	X	X	X
N-methylperfluorooctanesulfonamide (NMeFOSA)	X		X
N-methylperfluorooctanesulfonamido)-ethanol (NMeFOSE)	X		X
Perfluorobutanoic acid (PFBA)	X	X	X
Nonafluorobutane-1-sulfonic acid (PFBS)	X	X	X
Perfluorodecanoic acid (PFDA)	X	X	X
Perfluorodecanoic acid (PFDoA)	X	X	X
Perfluorodecanesulfonic acid (PFDS)	X		X
Perfluoro(2-ethoxyethane) sulfonic acid (PFEESA)	X	X	X
Perfluoroheptanoic acid (PFHpA)	X	X	X
Perfluoroheptanesulfonic acid (PFHpS)	X	X	X
Perfluorohexanoic acid (PFHxA)	X	X	X
Perfluorohexanesulfonic acid (PFHxS)	X	X	X
Perfluoro(3-methoxypropionic) acid (PFMPA)	X	X	X
Perfluoro(4-methoxybutanoic) acid (PFMBA)	X	X	X
Perfluorononanoic acid (PFNA)	X	X	X
Perfluoro-1-Nonanesulfonic acid (PFNS)	X		X
Pentadecafluorooctanoic acid (PFOA)	X	X	X
Perfluorooctane sulfonic acid (PFOS)	X	X	X
Perfluorooctanesulfonamide (PFOSA)	X		X
Perfluoropentanoic acid (PFPeA)	X	X	X
Perfluoropentanesulfonic acid (PFPeS)	X	X	X
Perfluorotetradecanoic acid (PFTeA)	X	X	X
Perfluorotridecanoic acid (PFTrDA)	X	X	X
Perfluoroundecanoic acid (PFUnA)	X	X	X



Chemical Testing Laboratory
Certificate No. 2427.03

Reference Material Producer
Certificate No. 2427.02

Proficiency Testing Provider
Certificate No. 2427.01

Order Now

Have questions? Reach out to Phenova at:
info@phenova.com or call **303-940-0033**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions. Phenova is a trademark of Phenomenex. For research use only. Not for use in clinical diagnostic procedures. © 2023 Phenomenex, Inc. All rights reserved.

