


# Vitamin Testing

**Top applications and reproducible Phenomenex products to achieve your vitamin testing needs**

- Simple cleanup procedures
- Rapid HPLC analysis
- Accurate, reproducible results

# Here to Support You with Your Clinical Research

## Personalized Support



With Phenomenex, your lab has its very own technical consultant who is available to assist with method development, optimization, troubleshooting, and product recommendations. Feel free to contact us with any questions you may have!

**Support@Phenomenex.com**

## Take Advantage of Our Lab



Our customers asked us for more support and we responded with PhenoLogix, a full-service analytical support laboratory. Let our PhenoLogix clinical research gurus save your lab time and money.

- Method development and optimization in our PhenoLogix laboratory
- On-site method development
- On-site hands-on training

Explore Explore PhenoLogix- a dedicated dedicated team to help with any of your questions and lab needs.

**[www.phenomenex.com/PhenoLogix](http://www.phenomenex.com/PhenoLogix)**

## 24 Hour Access to Information

### Technical Database

Access our vast database of technical notes, applications, webinars, and more. Search by industry, document type, and separation mode to easily find what you need.

**[www.phenomenex.com/MyLinks](http://www.phenomenex.com/MyLinks)**

### Tools

Utilize our web tools to find product recommendations, create customized methods, and more.

**[www.phenomenex.com/Tools](http://www.phenomenex.com/Tools)**



For Research Use Only. Not for use in diagnostic procedures.

© 2022 Phenomenex, Inc. All rights reserved.



# VITAMIN TESTING

## TABLE OF CONTENTS

### 04-14

#### APPLICATIONS

Methylmalonic Acid (MMA) From Plasma .....	4
Vitamin B3 from Plasma: Nicotinic Acid and Nicotinamide .....	5
Vitamin B6 from Plasma: Pyridoxal 5'-Phosphate (PLP), 4-Pyridoxic Acid (PA), and Pyridoxal (PL) .....	6-7
Vitamin C from Plasma: Ascorbic Acid.....	8
25-OH-Vitamin D2 and D3 from Serum.....	9
Vitamin A and E from Human Serum.....	10-11
High pH Vitamin B1 and B6 in Whole Blood....	12-13
25-OH Vitamin D2 and D3 from Plasma.....	14

### 15-16

#### SAMPLE PREPARATION

Protein Precipitation .....	15
Simplified Liquid Extraction (SLE) .....	16

### 20-21

#### HPLC COLUMNS

Complete Scalable Solution from UHPLC to HPLC .....	17-19
Setting the Standard for pH Method Development .....	20-21

# Methylmalonic Acid (MMA) From Plasma

## Sample Preparation

Combine 0.5 mL of 1 % aqueous Acetic Acid and 50  $\mu$ L of internal standard with 100  $\mu$ L blank, standard, or sample

## SPE Method


**Cartridge:** Strata™-X-AW 30 mg/1 mL  
**Part No.:** [8B-S038-TAK](#)  
**Condition:** 1 mL of Methanol  
**Equilibrate:** 1 mL of 1 % Acetic Acid in Water  
**Load:** Pretreated sample  
**Wash:** 0.5 mL Methanol/Water (50:50)  
**Dry:** 5-10 min at max vacuum (or positive pressure)  
**Elute:** 2x 0.6 mL 2 % NH<sub>4</sub>OH in Methanol  
**Evaporate:** Evaporate solvent to dryness @ 45-50 °C under a gentle stream of Nitrogen  
**Reconstitute:** 200  $\mu$ L 0.1 % Formic Acid in Water

## LC Conditions

**Analytical Column:** Luna™ Omega 1.6  $\mu$ m PS C18  
**Dimensions:** 50 x 2.1 mm  
**Part No.:** [00B-4752-AN](#)  
**Guard Cartridge:** SecurityGuard™ ULTRA [AF0-8497](#)  
**Mobile Phase:** A: 0.1 % Formic Acid in Water  
B: 0.1 % Formic Acid in Acetonitrile

Gradient:	Time (min)	B (%)
	0.01	2
	2	90
	3	90
	3.01	2
	5	2

**Flow Rate:** 0.4 mL/min  
**Injection Volume:** 5  $\mu$ L  
**Temperature:** 40 °C



Representative chromatogram of an extracted sample. Pooled human plasma was spiked with standards to 1.5  $\mu$ g/mL of succinic acid and 750 nmol/L of methylmalonic acid above the endogenous concentrations and processed by solid phase extraction. Peaks in order of elution: plasma interference (0.81 min), succinic acid (1.00 min), methyl-D3-malonic acid (1.20 min), and methylmalonic acid (1.23 min)

## MS/MS Source Conditions

**Detector:** SCIEX® 4000 QTRAP®  
**Mode:** Negative Ionization Mode  
**Scan Type:** MRM  
**Curtain Gas (CUR):** 10.0 psi  
**Collision Gas (CAD):** Medium  
**IonSpray Voltage (IS):** -4500 V  
**Temperature (TEM):** 600 °C  
**Ion Source Gas 1 (Gas1):** 50 psi  
**Ion Source Gas 2 (Gas2):** 50 psi  
**Interface Heater (ihe):** On



Have different front end systems  
or require different methods?


Contact us to discuss Sample Prep  
or Column options available.

## Vitamin B3 from Plasma


# Nicotinic Acid and Nicotinamide by LC-MS/MS

### Sample Preparation


1. Add 300  $\mu$ L Acetonitrile to a well of an Impact™ Protein Precipitation Plate.
2. Add 100  $\mu$ L plasma/serum to the well. Mix 3 times with a pipette tip (or vortex the whole plate briefly).
3. Filter under vacuum (5mm Hg) for 5 minutes. Ensure that a collection plate is positioned underneath the Impact Protein Precipitation Plate.
4. Cover collection plate with a sealing mat. Sample is now ready to be injected onto the LC-MS/MS. If sample will not be injected immediately, transfer the filtrate to an amber Verex™ autosampler vial (ambient) to protect from light.



App ID 20909



App ID 20906



### LC-MS/MS Conditions

**Column:** Gemini™ 3  $\mu$ m C18  
**Dimensions:** 100 x 4.6 mm  
**Part No.:** [00D-4439-E0](#)  
**Guard:** SecurityGuard™ cartridge C18, 4 x 3.0 mm  
**Part No.:** [AJO-7597](#) + [KJO-4282](#)  
**Mobile Phase:** A: 0.1 % Formic Acid in Water  
 B: Methanol  
**Gradient:**

Time (min)	% B
0	10
2.5	90
2.6	10
4	10

**Flow Rate:** 0.6 mL/min  
**Temperature:** 22 °C  
**Injection Volume:** 2  $\mu$ L  
**Detector:** MS/MS (SCIEX® API 4000™) ESI+  
**Samples:** 1. Nicotinamide  
 2. Nicotinic Acid

### Data Summary

Analyte	LOD		LOQ	
	ng/mL	S/N Ratio	ng/mL	S/N Ratio
Nicotinamide	2	4.3	10	12.4
Nicotinic Acid	2	3.3	10	10.1



## Vitamin B6 from Plasma

# Pyridoxal 5'-Phosphate (PLP), 4-Pyridoxic Acid (PA), and Pyridoxal (PL) by HPLC-Fluorescence


### Sample Preparation

1. Thaw plasma samples and plasma/serum spiked calibrators or pre-manufactured calibration standards and controls at ambient temperature. Protect from light.
2. Pipette 200  $\mu$ L plasma blank, calibration standards, controls, and plasma samples into appropriately labeled 0.6 mL amber microcentrifuge tubes.
  - Briefly vortex the calibrators and controls immediately prior to sampling.
  - Mix the plasma samples by gentle inversion immediately prior to sampling.
  - Protect the tubes from light.
3. Add 30  $\mu$ L of 100 mg/mL semicarbazide/glycine solution into all the tubes containing samples; cap the tubes, vortex for 15 sec.
4. Incubate in the dark at room temperature for 30 min.
5. Uncap the tubes; add 25  $\mu$ L of 20 % meta-phosphoric acid to the controls and samples.
6. Recap the tubes and vortex for 30 sec.
7. Centrifuge for 5 min at 14,000 rpm at room temperature.  
Note: The relative centrifugal force (RCF) = 16,000 g.
8. Transfer 150  $\mu$ L of supernatant to an amber autosampler glass Verex™ vial.
9. Cover the vial with a screw cap and place it in the autosampler at room temperature.

See more applications at  
[www.phenomenex.com/ClinicalGuides](http://www.phenomenex.com/ClinicalGuides)



## HPLC-Fluorescence Conditions




**Column:** Gemini™ 3µm NX-C18  
**Dimension:** 100 x 4.6 mm  
**Part No.:** [00D-4453-E0](#)  
**Guard:** SecurityGuard™ guard cartridge system  
**Part No.:** [AJ0-8368](#) + [KJ0-4282](#)  
**Mobile Phase:** A: 20 mM Sodium phosphate and 1.0 mL Acetic acid in 1 L DI water, pH 6  
 B: Acetonitrile/Methanol (70:30)  
**Gradient:**

Time (min)	% B
0.0	5
5.0	60
5.1	95
6.0	5
7.0	5

**Flow Rate:** 1 mL/min  
**Column Temp:** 35 °C  
**Injection Volume:** 30 µL  
**Detection:** Fluorescence, Ex 360, Em 450  
**Sample:** 1. Pyridoxal 5'-phosphate  
 2. 4-Pyridoxic Acid  
 3. Pyridoxal


Representative calibration curve for Pyridoxal 5'-phosphate. The method was linear across the concentration range of 6.25-200 nmol/L.



## Data Summary

### Statistical data of PLP, PA, and PL in plasma.

Analyte	LOD	LOQ	Intra Assay Precision (N=3)							
			nmol/L	nmol/L	8 nmol/L		80 nmol/L		160 nmol/L	
					Mean Recovery	RSD %	Mean Recovery	RSD %	Mean Recovery	RSD %
PLP	2	4	86.95	2.50	97.82	0.91	95.60	0.88		
PA	2	4	94.96	4.33	95.14	3.30	95.09	2.15		
PL	2	4	91.59	4.50	97.79	2.76	94.77	4.24		



# Vitamin C from Plasma Ascorbic Acid by LC-UV

## Sample Preparation

### Rapid Protein Precipitation:

Add 300 µL cold 5 % Meta-phosphoric Acid (4 °C) to the wells of an Impact™ Protein Precipitation Plate (see pg. 11 for ordering information).


Add 100 µL plasma/serum directly into the 5 % Meta-phosphoric acid.

Mix 5 times by aspirating with same pipette tip.

Centrifuge\* the Impact plate at 500 g (with collection plate underneath) for 5 min at 4 °C. Purified filtrate is collected in the collection plate.


\* A vacuum manifold may be used however 25mm Hg vacuum pressure must be applied for up to 10 minutes or until sample is completely pulled through the Impact plate

## LC-UV Conditions



**Column:** Kinetex™ 5 µm XB-C18  
**Dimension:** 150 x 4.6 mm  
**Part No.:** [00F-4605-E0](#)  
**Guard:** SecurityGuard™ ULTRA guard cartridge system  
**Part No.:** [AJ0-8768](#) + [AJ0-9000](#)  
**Mobile Phase:** A: 0.1 % Formic Acid in Water  
 B: Acetonitrile  
**Gradient:** Time (min) % B  
 0.0 0  
 3.5 0  
 3.6 100  
 5.0 100  
 5.1 0  
 7.0 0  
**Flow Rate:** 0.8 mL/min  
**Column Temp:** 22 °C  
**Injection Volume:** 30 µL  
**Detection:** UV-Vis @ 245 nm  
**Sample:** 1. Vitamin C (Ascorbic Acid)  
 2. Uric acid

Representative calibration curve for Vitamin C. The method was linear across a concentration range of 0 - 100 µg/mL.



### Recovery of Vitamin C from spiked human plasma

Added Vitamin C (µg/mL)	Observed (µg/mL)	Recovery (%)
0	5.0	-
3	7.9	98.8
5	9.0	90.0
10	14.1	94.0
20	24.8	99.2
30	35.1	100.3
40	46.8	104.0
60	65.5	100.8
Mean		98.2
CV		4.8

### Intra- and inter-day imprecision of plasma Vitamin C analysis

	Intra Imprecision (n=12)			Inter Imprecision (n=6)		
	Mean (µg/mL)	SD	% CV	Mean (µg/mL)	SD	% CV
<b>QC 1</b>	1.1	0.0309	2.90	1.1	0.0532	4.82
<b>QC 2</b>	11.2	0.3552	3.16	11.1	0.3857	3.47
<b>QC 3</b>	34.1	0.5910	1.74	34.9	0.6290	1.80

# 25-OH Vitamin D2, 25-OH-Vitamin D3, and 3-Epi-25-OH-Vitamin D3 from Human Serum


by On-line SPE-LC-MS/MS

## Sample Preparation

Precipitate 100 µL of human serum by adding 200 µL 5:2:1 Methanol: Acetonitrile: 2 % ZnSO<sub>4</sub> and 30 µL of 1 µg/mL working internal standard in water, mix and centrifuge. Inject supernatant.

## On-line SPE-LC-MS/MS Conditions

- On-line SPE Extraction Cartridge:** Strata™ C8 20 µm
- Dimension:** 20 x 20 mm
- Part No.:** [00M-S101-B0-CB](#)
- Column:** Kinetex™ 2.6 µm F5
- Dimension:** 100 x 4.6 mm
- Part No.:** [00D-4723-E0](#)
- Mobile Phase:** A: 0.1 % Formic Acid in Water  
B: 0.1 % Formic Acid in Methanol (A/B 15:85)
- Flow Rate:** 750 µL/min
- Needle Wash:** Wash 1: Methanol/Water (50:50)  
Wash 2: 0.1% Formic Acid in Water
- Injection Volume:** 40 µL
- Detection:** MS/MS (SCIEX® API 4000™ QTRAP®), APCI+
- Sample:**
- 25-OH Vitamin D3
  - 3-epi-25 OH Vitamin D3
  - 25-OH Vitamin D2




App ID 22965

## LC Conditions: Thermo Cohesive® System


Step	Start	Sec	Loading Pump				Eluting Pump				Comments				
			Flow	Grad	%A	%B	%C	%D	Tee	Loop		Flow	Grad	%A	%B
1	0.00	30	0.75	Step	30.0	70.0	-	-	-	Out	0.75	Step	20.0	80.0	Extract sample
2	0.50	5	0.33	Step	15.0	85.0	-	-	-	Out	0.33	Step	20.0	80.0	Slow down pumps
3	0.58	30	0.33	Step	15.0	85.0	-	-	T	In	0.33	Step	15.0	85.0	Transfer analytes
4	1.08	60	1.50	Step		100.0	-	-	-	Out	0.75	Ramp	15.0	85.0	Separate analytes, wash extraction column
5	2.08	120	1.50	Step		100.0	-	-	-	Out	0.75	Step	15.0	85.0	Elute analytes, wash extraction column
6	4.08	30	0.75	Step	15.0	85.0	-	-	-	In	0.75	Step	20.0	80.0	Wash columns and valves
7	4.58	60	0.75	Step	30.0	70.0	-	-	-	In	0.75	Step	20.0	80.0	Fill transfer loop, equilibrate HPLC column
8	5.58	90	0.75	Step	30.0	70.0	-	-	-	Out	0.75	Step	20.0	80.0	Equilibrate columns

Focus Mode Configuration

## Representative Chromatograms of LLOQ (2.5 ng/mL) in Human Serum



App ID 22966




App ID 22967


# Vitamin A and E from Human Serum by LC-MS/MS

## Sample Preparation

1. Dilute 200  $\mu\text{L}$  human serum with 100  $\mu\text{L}$  Isopropanol (IPA) and 150  $\mu\text{L}$  of water. Vortex 30 seconds.
2. Load onto the Novum™ MAX 96-well plate (Part No.: [8E-S138-5GA](#)). Apply a short pulse of vacuum for 10-15 seconds.
3. Wait 5 minutes.
4. Elution: Dispense Ethyl Acetate/Acetone (90:10), 2 x 900  $\mu\text{L}$  (2 aliquots) to elute by gravity (~ 5 minutes) and collect the eluent. Apply vacuum at 5mm Hg for 20-30 secs to complete the extraction.
5. Dry down: Evaporate the final extract to complete dryness under slow stream of nitrogen at 45 °C.
6. Reconstitute: Dried residue in 200  $\mu\text{L}$  of initial mobile phase.



## LC-MS/MS Analysis of Vitamin A and E Using Dual Polarity Technique in MS



<b>Column:</b> Kinetex™ 5 $\mu\text{m}$ EVO C18	<b>Instrument:</b> Agilent® 1260
<b>Dimensions:</b> 100 x 2.1mm	<b>MS/MS Instrument:</b> SCIEX® Triple Quad™ 5000
<b>Part No.:</b> <a href="#">00D-4633-AN</a>	<b>Analyte:</b> 1. Vitamin A (retinol) 2. $\gamma$ -Tocopherol (Vitamin E) 3. $\alpha$ -Tocopherol (Vitamin E)
<b>Mobile Phase:</b> A: Water B: Isopropanol/acetonitrile (1:1)	
<b>Gradient:</b>	
<b>Time (min)</b>	<b>B (%)</b>
0.0	65
3.5	95
4.0	95
4.1	65
5.05	65
<b>Flow Rate:</b> 0.6mL/min	
<b>Injection Volume:</b> 5 $\mu\text{L}$	
<b>Temperature:</b> Ambient	




Have different front end systems  
or require different methods?

Contact us to discuss Sample Prep  
or Column options available.


## Linearity Curve of Vitamin A Extracted Samples on Novum™ MAX

(Matrix: Doubly-Charcoal Stripped Serum)




## Linearity Curve of α-Tocopherol Extracted Samples on Novum MAX

(Matrix: Egg White Albumin)



## Linearity Curve of γ-Tocopherol Extracted Samples on Novum MAX

(Matrix: Egg White Albumin)



## Improved Recovery and Cleanliness (92-110 %; CV=2-9 %) from Novum Max

Analyte	% Recovery
Vitamin A	92
α-Tocopherol (Vit E)	105
γ-Tocopherol (Vit E)	110



# High pH Vitamin B1 and B6 in Whole Blood

## Experimental Conditions

### Optimized Sample Extraction Method

Human whole blood samples were frozen immediately at  $-20^{\circ}\text{C}$  after collection. It is important to freeze the sample for at least 24 hours prior to analysis in order to prevent the analyte from decomposition, especially TDP.

1. Pipette  $100\ \mu\text{L}$  of thawed hemolyzed blood into a  $1.8\ \text{mL}$  centrifuge tube
2. Add  $300\ \mu\text{L}$  of working internal standard ( $20\ \text{ng/mL}$  of Pyridoxine- $\text{D}_2$  and  $50\ \text{ng/mL}$  of Thiamine- $^{13}\text{C}_4$  in DI water) and mix for 30 seconds
3. Add  $30\ \mu\text{L}$  of  $70\%$   $\text{HClO}_4$  and mix for 1 minute to precipitate proteins
4. Centrifuge sample at  $14,000\ \text{rpm}$ s for 10 minutes to pellet the protein
5. Transfer  $200\ \mu\text{L}$  of supernatant into an autosampler vial for LC MS/MS analysis

Note: Since the analytes are light sensitive, the extraction steps were performed in amber color centrifuge tube and were protected from light.


### LC/MS/MS Method Parameters

**Column:** Gemini<sup>®</sup>  $5\ \mu\text{m}$  C18  
**Dimensions:**  $50 \times 4.6\ \text{mm}$   
**Part No.:** 00B-4435-E0  
**SecurityGuard Cartridge:** AJ0-7597  
**Mobile Phase:** A:  $10\ \text{mM NH}_4\text{HCO}_3$  in water,  $\text{pH}8.8$   
B: Methanol


Gradient:	Time (min)	B (%)
	0.01	0
	1.5	0
	5	60
	6.5	60
	6.51	0
	9	0

**Flow Rate:**  $600\ \mu\text{L}/\text{min}$   
**Injection Volume:**  $10\ \mu\text{L}$   
**Instrument:** Agilent<sup>®</sup> 1260 LC  
**Detection:** MS/MS (ESI+) (SCIEX API 4500<sup>™</sup>)  
**Sample:** 1. Pyridoxal 5-phosphate (PLP)  
2. Thiamine Diphosphate (TDP)  
3. Pyridoxine  $\text{D}_2$   
4. Thiamine- $^{13}\text{C}_4$

### Structure of TDP and PLP



TDP  
LogP: -5.80




PLP  
LogP: -2.09

**Accuracy and precision for TDP and PLP. Five replicates at two concentrations were analyzed over two runs (n=10 for each concentration).**

	TDP 100 ng/mL	PLP 100 ng/mL	TFP 200 ng/mL	PLP 200 ng/mL
Mean Conc. Found (ng/mL)	96.3	92.8	224	201
STDV	11.3	11.7	18.2	26.8
CV%	11.8	12.6	8.13	13.3
Accuracy (%)	96.3	92.8	112	101


**Representative Chromatogram in Whole Blood at LLOQ (20 ng/mL)**




# 25-OH Vitamin D2 and D3 from Plasma by LC-MS/MS

## Sample Preparation

1. Add 50  $\mu\text{L}$  of precipitating reagent containing internal standard to a 1.5 mL centrifuge tube.
2. Pipette 100  $\mu\text{L}$  serum into the centrifuge tube.
3. Vortex 20-30 sec.
4. Inspect each tube to ensure no unmixed sample remains in the bottom of the tube. A homogenous mixture is critical. If unmixed sample remains at the bottom of the tube, dislodge by inverting and tapping, then re-vortex.
5. Centrifuge 15 min at 13,000 rpm.
6. Transfer supernatant into a sample vial without disturbing the pellet.



## LC-MS/MS Conditions



Analyte	Q1	Q3
25-OH Vitamin D3-D6	389.3	263.3
25-OH Vitamin D2	395.3	209.3
25-OH Vitamin D3	383.2	257.2

<b>Column:</b> Kinetex™ 2.6 $\mu\text{m}$ C18	<b>Gradient:</b> Time (min) %B
<b>Dimension:</b> 50 x 4.6 mm	0.00 8
<b>Part No.:</b> <a href="#">00B-4462-E0</a>	0.08 8
<b>Guard:</b> SecurityGuard™ ULTRA guard cartridge system	3.40 100
<b>Part No.:</b> <a href="#">AJO-8768</a> + <a href="#">AJO-9000</a>	4.80 100
<b>Mobile Phase:</b> A: 0.05 % Formic Acid in Water	6.00 8
B: 5 mM Ammonium Acetate with 0.1 % Formic Acid in Water	<b>Flow Rate:</b> 1 mL/min
	<b>Column Temp:</b> 35 °C
	<b>Injection Volume:</b> 50 $\mu\text{L}$
	<b>Detection:</b> MS/MS (SCIEX® API 4000™)


# Protein Precipitation

## Impact Rapid Protein Precipitation




### Impact Rapid Protein Precipitation

1. Quickly cleanup sample by passing biological samples through the Impact filter
2. Use Impact filter to quickly clean-up biological sample.
3. Increase reproducibility with the leak-free membrane, preventing premature sample breakthrough and incomplete protein precipitation



### How it Works:



### Ordering Information

#### Impact Protein Precipitation Plates

Part No.	Description	Unit
<b>Impact Precipitation Plates</b>		
<a href="#">CE0-7565</a>	Impact Protein Precipitation, Square Well, Filter Plate, 2 mL	2/box
<a href="#">CE0-7566</a>	Impact Protein Precipitation, Square Well, Filter Plate, 2 mL, Long Drip	2/box
<b>Impact Starter Kit for Protein Precipitation</b>		
<a href="#">CE0-8201</a>	Impact Protein Precipitation Plate (2 ea) Collection Plate 2 mL	2 ea
<a href="#">AHO-8199</a>	Sealing Mat, Santoprene™	2 ea

#### Accessories


Part No.	Description	Unit
<b>Collection Plates (deep well, polypropylene)</b>		
<a href="#">AHO-7192</a>	96-Well Collection Plate 350 µL/well	50/pk
<a href="#">AHO-7193</a>	96-Well Collection Plate 1 mL/well	50/pk
<a href="#">AHO-7194</a>	96-Well Collection Plate 2 mL/well	50/pk
<a href="#">AHO-8635</a>	96-Well Collection Plate, 2 mL Square/Round-Conical	50/pk
<a href="#">AHO-8636</a>	96-Well Collection Plate, 2 mL Round/Round, 8 mm	50/pk
<a href="#">AHO-7279</a>	96-Well Collection Plate, 1 mL/well Round, 7 mm	50/pk
<b>Sealing Mats</b>		
<a href="#">AHO-8597</a>	Sealing Mats, Pierceable, 96-Square Well, Silicone	50/pk
<a href="#">AHO-8598</a>	Sealing Mats, Pre-Slit, 96-Square Well, Silicone	50/pk
<a href="#">AHO-8631</a>	Sealing Mats, Pierceable, 96-Round Well 7 mm, Silicone	50/pk
<a href="#">AHO-8632</a>	Sealing Mats, Pre-Slit, 96-Round Well 7 mm, Silicone	50/pk
<a href="#">AHO-8633</a>	Sealing Mats, Pierceable, 96-Round Well 8 mm, Silicone	50/pk
<a href="#">AHO-8634</a>	Sealing Mats, Pre-Slit, 96-Round Well 8 mm, Silicone	50/pk
<a href="#">AHO-7362</a>	Sealing Tape Pad	10/pk
<b>Vacuum Manifolds</b>		
<a href="#">AHO-8950</a>	96-Well Plate Manifold, Universal with Vacuum Gauge	ea

# Novum™ Simplified Liquid Extraction (SLE) A Variety of Formats to Fit Your Sample and Throughput Requirements



## Tubes

Process samples as small as 100 µL or as large as 1 mL using Novum SLE tubes. Ideal for all types of applications including Bioanalytical, Food Safety, and Environmental.




### Ordering Information

Novum Simplified Liquid Extraction (SLE) Tubes				
Novum SLE Tubes	1 cc	3 cc	6 cc	12 cc
Maximum Sample Volume (before dilution)	200 µL	400 µL	1 mL	2 mL
Recommended Elution Volume	2x 600 µL	2x 900 µL	2x 2.5 mL	2x 5 mL
Part No.	<a href="#">8B-S138-FAK</a>	<a href="#">8B-S138-5BJ</a>	<a href="#">8B-S138-JCH</a>	<a href="#">8B-S138-KDG</a>
Unit	100/pk	50/pk	30/pk	20/pk

## 96-Well Plates

Process 96 samples at once in an easily automatable 96-well plate. Perfect for high-throughput applications.




### Ordering Information

Novum Simplified Liquid Extraction (SLE) 96-Well Plates				
Novum SLE Tubes	MINI	MAX	PRO MINI	PRO MAX
Maximum Sample Volume (before dilution)	300 µL	400 µL	300 µL	400 µL
Recommended Elution Volume	1x 1 mL	2x 900 µL	1x 1 mL	2x 900 µL
Part No.	<a href="#">8E-S138-FGA</a>	<a href="#">8E-S138-5GA</a>	<a href="#">8E-S539-FGA</a>	<a href="#">8E-S539-5GA</a>
Unit	1/pk	1/pk	1/pk	1/pk

For accessories that are compatible with Novum Simplified Liquid Extraction (SLE) Products, see pp. 13

# Thrive with Kinetex Core-Shell Technology! Complete Scalable Solution from UHPLC to HPLC



## Column Characteristics

Kinetex Phases	Shipping Solvent <sup>†</sup>	Particle Sizes (μm)	Pore Size (Å)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	pH Stability	Reversed Phase	Normal Phase	HILIC	100% Aqueous Stable
<b>Polar C18</b>	Acetonitrile/Water (50:50)	2.6	100	200	9	1.5-8.5*	●			●
<b>PS C18</b>	Acetonitrile/Water (50:50)	2.6	100	200	9	1.5-8.5*	●			●
<b>C18</b>	Acetonitrile/Water (50:50)	1.3, 1.7, 2.6, 5	100	200	12	1.5-8.5*	●			
<b>EVO C18</b>	Acetonitrile/Water (45:55)	1.7, 2.6, 5	100	200	11	1-12	●			●
<b>XB-C18</b>	Acetonitrile/Water (50:50)	1.7, 2.6, 3.5, 5	100	200	10	1.5-8.5*	●			
<b>C8</b>	Acetonitrile/Water (45:55)	1.7, 2.6, 5	100	200	8	1.5-8.5*	●			
<b>Biphenyl</b>	Acetonitrile/Water w/ 0.1 % Formic Acid (50:50)	1.7, 2.6, 5	100	200	11	1.5-8.5*	●			●
<b>Phenyl-Hexyl</b>	Acetonitrile/Water (45:55)	1.7, 2.6, 5	100	200	11	1.5-8.5*	●			
<b>F5</b>	Acetonitrile/Water (40:60)	1.7, 2.6, 5	100	200	9	1.5-8.5*	●		●	●
<b>HILIC</b>	Acetonitrile/100 mM Ammonium Formate (93:7)	1.7, 2.6, 5	100	200	0	2.0-7.5		●	●	
<b>PAH</b>	Acetonitrile/Water (65:35)	3.5	—	—	12	1.5-8.5*	●			

<sup>†</sup> Shipping conditions may vary slightly in terms of organic to aqueous ratio, depending on column dimensions.  
\* pH stability under gradient conditions. pH stability is 1.5-10 under isocratic conditions.

# Ordering Information (cont'd)

3.5 µm Minibore and Analytical Columns (mm)						SecurityGuard™ ULTRA Cartridges <sup>†</sup>	
Phases	50 x 2.1	150 x 2.1	100 x 4.6	150 x 4.6	250 x 4.6	3/pk	3/pk
<b>XB-C18</b>	—	—	<a href="#">00D-4744-E0</a>	<a href="#">00F-4744-E0</a>	—	—	<a href="#">AJ0-8768</a>
<b>PAH</b>	<a href="#">00B-4764-AN</a>	<a href="#">00F-4764-AN</a>	<a href="#">00D-4764-E0</a>	<a href="#">00F-4764-E0</a>	<a href="#">00G-4764-E0</a>	<a href="#">AJ0-9535</a>	<a href="#">AJ0-9533</a>

for 2.1 mm ID      for 4.6 mm ID

5 µm Minibore Columns (mm)					SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
<b>EVO C18</b>	<a href="#">00A-4633-AN</a>	<a href="#">00B-4633-AN</a>	<a href="#">00D-4633-AN</a>	<a href="#">00F-4633-AN</a>	<a href="#">AJ0-9298</a>
<b>F5</b>	—	<a href="#">00B-4724-AN</a>	<a href="#">00D-4724-AN</a>	<a href="#">00F-4724-AN</a>	<a href="#">AJ0-9322</a>
<b>Biphenyl</b>	<a href="#">00A-4627-AN</a>	<a href="#">00B-4627-AN</a>	<a href="#">00D-4627-AN</a>	—	<a href="#">AJ0-9209</a>
<b>XB-C18</b>	<a href="#">00A-4605-AN</a>	<a href="#">00B-4605-AN</a>	<a href="#">00D-4605-AN</a>	—	<a href="#">AJ0-8782</a>
<b>C18</b>	<a href="#">00A-4601-AN</a>	<a href="#">00B-4601-AN</a>	<a href="#">00D-4601-AN</a>	<a href="#">00F-4601-AN</a>	<a href="#">AJ0-8782</a>
<b>C8</b>	—	<a href="#">00B-4608-AN</a>	<a href="#">00D-4608-AN</a>	—	<a href="#">AJ0-8784</a>
<b>Phenyl-Hexyl</b>	—	<a href="#">00B-4603-AN</a>	—	—	<a href="#">AJ0-8788</a>
<b>HILIC</b>	—	<a href="#">00B-4606-AN</a>	—	—	<a href="#">AJ0-8786</a>

for 2.1 mm ID



5 µm MidBore™ Columns (mm)					SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 3.0	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
<b>EVO C18</b>	<a href="#">00A-4633-Y0</a>	<a href="#">00B-4633-Y0</a>	<a href="#">00D-4633-Y0</a>	<a href="#">00F-4633-Y0</a>	<a href="#">AJ0-9297</a>
<b>F5</b>	—	—	<a href="#">00D-4724-Y0</a>	<a href="#">00F-4724-Y0</a>	<a href="#">AJ0-9321</a>
<b>Biphenyl</b>	—	<a href="#">00B-4627-Y0</a>	<a href="#">00D-4627-Y0</a>	<a href="#">00F-4627-Y0</a>	<a href="#">AJ0-9208</a>
<b>XB-C18</b>	—	<a href="#">00B-4605-Y0</a>	<a href="#">00D-4605-Y0</a>	<a href="#">00F-4605-Y0</a>	<a href="#">AJ0-8775</a>
<b>C18</b>	<a href="#">00A-4601-Y0</a>	<a href="#">00B-4601-Y0</a>	<a href="#">00D-4601-Y0</a>	<a href="#">00F-4601-Y0</a>	<a href="#">AJ0-8775</a>
<b>C8</b>	—	<a href="#">00B-4608-Y0</a>	<a href="#">00D-4608-Y0</a>	—	<a href="#">AJ0-8777</a>
<b>Phenyl-Hexyl</b>	—	<a href="#">00B-4603-Y0</a>	<a href="#">00D-4603-Y0</a>	—	<a href="#">AJ0-8781</a>

for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
<b>EVO C18</b>	<a href="#">00B-4633-E0</a>	<a href="#">00D-4633-E0</a>	<a href="#">00F-4633-E0</a>	<a href="#">00G-4633-E0</a>	<a href="#">AJ0-9296</a>
<b>F5</b>	<a href="#">00B-4724-E0</a>	<a href="#">00D-4724-E0</a>	<a href="#">00F-4724-E0</a>	<a href="#">00G-4724-E0</a>	<a href="#">AJ0-9320</a>
<b>Biphenyl</b>	<a href="#">00B-4627-E0</a>	<a href="#">00D-4627-E0</a>	<a href="#">00F-4627-E0</a>	<a href="#">00G-4627-E0</a>	<a href="#">AJ0-9207</a>
<b>XB-C18</b>	<a href="#">00B-4605-E0</a>	<a href="#">00D-4605-E0</a>	<a href="#">00F-4605-E0</a>	<a href="#">00G-4605-E0</a>	<a href="#">AJ0-8768</a>
<b>C18</b>	<a href="#">00B-4601-E0</a>	<a href="#">00D-4601-E0</a>	<a href="#">00F-4601-E0</a>	<a href="#">00G-4601-E0</a>	<a href="#">AJ0-8768</a>
<b>C8</b>	<a href="#">00B-4608-E0</a>	<a href="#">00D-4608-E0</a>	<a href="#">00F-4608-E0</a>	<a href="#">00G-4608-E0</a>	<a href="#">AJ0-8770</a>
<b>Phenyl-Hexyl</b>	<a href="#">00B-4603-E0</a>	<a href="#">00D-4603-E0</a>	<a href="#">00F-4603-E0</a>	<a href="#">00G-4603-E0</a>	<a href="#">AJ0-8774</a>
<b>HILIC</b>	—	—	<a href="#">00F-4606-E0</a>	<a href="#">00G-4606-E0</a>	<a href="#">AJ0-8772</a>

for 4.6 mm ID

5 µm Semi-Preparative Columns (mm)					SecurityGuard SemiPrep Cartridges <sup>**</sup>
Phases	100 x 10	150 x 10	250 x 10	10 x 10 /3pk	
<b>EVO C18</b>	—	<a href="#">00F-4633-N0</a>	<a href="#">00G-4633-N0</a>	<a href="#">AJ0-9306</a>	
<b>F5</b>	—	—	<a href="#">00G-4724-N0</a>	<a href="#">AJ0-9323</a>	
<b>C18</b>	<a href="#">00D-4601-N0</a>	<a href="#">00F-4601-N0</a>	<a href="#">00G-4601-N0</a>	<a href="#">AJ0-9278</a>	
<b>Biphenyl</b>	—	<a href="#">00F-4627-N0</a>	<a href="#">00G-4627-N0</a>	<a href="#">AJ0-9280</a>	
<b>XB-C18</b>	—	<a href="#">00F-4605-N0</a>	<a href="#">00G-4605-N0</a>	<a href="#">AJ0-9278</a>	

for 9-16 mm ID

5 µm Axia™ Packed Preparative Columns (mm)						SecurityGuard PREP Cartridges <sup>**</sup>
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	15 x 21.2 /ea	
<b>EVO C18</b>	<a href="#">00B-4633-PO-AX</a>	<a href="#">00D-4633-PO-AX</a>	<a href="#">00F-4633-PO-AX</a>	<a href="#">00G-4633-PO-AX</a>	<a href="#">AJ0-9304</a>	
<b>F5</b>	—	—	<a href="#">00F-4724-PO-AX</a>	<a href="#">00G-4724-PO-AX</a>	<a href="#">AJ0-9324</a>	
<b>Biphenyl</b>	<a href="#">00B-4627-PO-AX</a>	<a href="#">00D-4627-PO-AX</a>	<a href="#">00F-4627-PO-AX</a>	<a href="#">00G-4627-PO-AX</a>	<a href="#">AJ0-9272</a>	
<b>XB-C18</b>	<a href="#">00B-4605-PO-AX</a>	<a href="#">00D-4605-PO-AX</a>	<a href="#">00F-4605-PO-AX</a>	<a href="#">00G-4605-PO-AX</a>	<a href="#">AJ0-9145</a>	
<b>C18</b>	<a href="#">00B-4601-PO-AX</a>	<a href="#">00D-4601-PO-AX</a>	<a href="#">00F-4601-PO-AX</a>	<a href="#">00G-4601-PO-AX</a>	<a href="#">AJ0-9145</a>	
<b>C8</b>	<a href="#">00B-4608-PO-AX</a>	<a href="#">00D-4608-PO-AX</a>	<a href="#">00F-4608-PO-AX</a>	<a href="#">00G-4608-PO-AX</a>	<a href="#">AJ0-9205</a>	
<b>Phenyl-Hexyl</b>	<a href="#">00B-4603-PO-AX</a>	<a href="#">00D-4603-PO-AX</a>	<a href="#">00F-4603-PO-AX</a>	<a href="#">00G-4603-PO-AX</a>	<a href="#">AJ0-9147</a>	
<b>HILIC</b>	—	<a href="#">00D-4606-PO-AX</a>	<a href="#">00F-4606-PO-AX</a>	<a href="#">00G-4606-PO-AX</a>	<a href="#">AJ0-9277</a>	

for 18-29 mm ID

5 µm Axia Packed Preparative Columns (mm)						SecurityGuard PREP Cartridges <sup>**</sup>
Phases	50 x 30	100 x 30	150 x 30	250 x 30	15 x 30 /ea	
<b>EVO C18</b>	<a href="#">00B-4633-U0-AX</a>	<a href="#">00D-4633-U0-AX</a>	<a href="#">00F-4633-U0-AX</a>	<a href="#">00G-4633-U0-AX</a>	<a href="#">AJ0-9305</a>	
<b>F5</b>	<a href="#">00B-4724-U0-AX</a>	<a href="#">00D-4724-U0-AX</a>	<a href="#">00F-4724-U0-AX</a>	—	<a href="#">AJ0-9325</a>	
<b>Biphenyl</b>	—	—	<a href="#">00F-4627-U0-AX</a>	<a href="#">00G-4627-U0-AX</a>	<a href="#">AJ0-9273</a>	
<b>XB-C18</b>	<a href="#">00B-4605-U0-AX</a>	<a href="#">00D-4605-U0-AX</a>	<a href="#">00F-4605-U0-AX</a>	<a href="#">00G-4605-U0-AX</a>	<a href="#">AJ0-9204</a>	
<b>C18</b>	<a href="#">00B-4601-U0-AX</a>	<a href="#">00D-4601-U0-AX</a>	<a href="#">00F-4601-U0-AX</a>	<a href="#">00G-4601-U0-AX</a>	<a href="#">AJ0-9204</a>	
<b>C8</b>	<a href="#">00B-4608-U0-AX</a>	<a href="#">00D-4608-U0-AX</a>	<a href="#">00F-4608-U0-AX</a>	<a href="#">00G-4608-U0-AX</a>	<a href="#">AJ0-9217</a>	
<b>Phenyl-Hexyl</b>	—	—	<a href="#">00F-4603-U0-AX</a>	<a href="#">00G-4603-U0-AX</a>	<a href="#">AJ0-9216</a>	
<b>HILIC</b>	—	—	<a href="#">00D-4606-U0-AX</a>	—	—	

for 30-49 mm ID

<sup>†</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)  
<sup>\*</sup> PREP SecurityGuard Cartridges require holder, Part No.: [AJ0-8223](#)  
<sup>\*\*</sup> PREP SecurityGuard Cartridges require holder, Part No.: [AJ0-8277](#)  
<sup>\*\*\*</sup> SemiPrep SecurityGuard Cartridges require holder, Part No.: [AJ0-9281](#)

# Ordering Information (cont'd)

2.6 µm Micro LC Columns (mm)						
Phases	30 x 0.3	50 x 0.3	100 x 0.3	150 x 0.3	50 x 0.5	150 x 0.5
Biphenyl	—	00B-4622-AC	—	00F-4622-AC	00B-4622-AF	—
C18	<a href="#">00A-4462-AC</a>	<a href="#">00B-4462-AC</a>	—	<a href="#">00F-4462-AC</a>	<a href="#">00B-4462-AF</a>	—
EVO C18	—	<a href="#">00B-4725-AC</a>	—	<a href="#">00F-4725-AC</a>	<a href="#">00B-4725-AF</a>	—
F5	—	<a href="#">00B-4723-AC</a>	<a href="#">00D-4723-AC</a>	<a href="#">00F-4723-AC</a>	<a href="#">00B-4723-AF</a>	—
XB-C18	<a href="#">00A-4496-AC</a>	<a href="#">00B-4496-AC</a>	<a href="#">00D-4496-AC</a>	<a href="#">00F-4496-AC</a>	<a href="#">00B-4496-AF</a>	<a href="#">00F-4496-AF</a>

2.6µm Minibore Columns (mm)						SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 2.1	50 x 2.1	75 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	<a href="#">00A-4725-AN</a>	<a href="#">00B-4725-AN</a>	—	<a href="#">00D-4725-AN</a>	<a href="#">00F-4725-AN</a>	<a href="#">AJO-9298</a>
PS C18	<a href="#">00A-4780-AN</a>	<a href="#">00B-4780-AN</a>	—	<a href="#">00D-4780-AN</a>	<a href="#">00F-4780-AN</a>	<a href="#">AJO-8951</a>
Polar C18	<a href="#">00A-4759-AN</a>	<a href="#">00B-4759-AN</a>	—	<a href="#">00D-4759-AN</a>	<a href="#">00F-4759-AN</a>	<a href="#">AJO-9532</a>
F5	<a href="#">00A-4723-AN</a>	<a href="#">00B-4723-AN</a>	—	<a href="#">00D-4723-AN</a>	<a href="#">00F-4723-AN</a>	<a href="#">AJO-9322</a>
Biphenyl	<a href="#">00A-4622-AN</a>	<a href="#">00B-4622-AN</a>	—	<a href="#">00D-4622-AN</a>	<a href="#">00F-4622-AN</a>	<a href="#">AJO-9209</a>
XB-C18	<a href="#">00A-4496-AN</a>	<a href="#">00B-4496-AN</a>	<a href="#">00C-4496-AN</a>	<a href="#">00D-4496-AN</a>	<a href="#">00F-4496-AN</a>	<a href="#">AJO-8782</a>
C18	<a href="#">00A-4462-AN</a>	<a href="#">00B-4462-AN</a>	<a href="#">00C-4462-AN</a>	<a href="#">00D-4462-AN</a>	<a href="#">00F-4462-AN</a>	<a href="#">AJO-8782</a>
C8	<a href="#">00A-4497-AN</a>	<a href="#">00B-4497-AN</a>	<a href="#">00C-4497-AN</a>	<a href="#">00D-4497-AN</a>	<a href="#">00F-4497-AN</a>	<a href="#">AJO-8784</a>
HILIC	<a href="#">00A-4461-AN</a>	<a href="#">00B-4461-AN</a>	<a href="#">00C-4461-AN</a>	<a href="#">00D-4461-AN</a>	<a href="#">00F-4461-AN</a>	<a href="#">AJO-8786</a>
Phenyl-Hexyl	<a href="#">00A-4495-AN</a>	<a href="#">00B-4495-AN</a>	<a href="#">00C-4495-AN</a>	<a href="#">00D-4495-AN</a>	<a href="#">00F-4495-AN</a>	<a href="#">AJO-8788</a>

for 2.1 mm ID

2.6µm MidBore™ Columns (mm)						SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 3.0	50 x 3.0	75 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	<a href="#">00A-4725-YO</a>	<a href="#">00B-4725-YO</a>	—	<a href="#">00D-4725-YO</a>	<a href="#">00F-4725-YO</a>	<a href="#">AJO-9297</a>
PS C18	<a href="#">00B-4780-YO</a>	<a href="#">00D-4780-YO</a>	—	<a href="#">00D-4780-YO</a>	<a href="#">00F-4780-YO</a>	<a href="#">AJO-8950</a>
Polar C18	—	<a href="#">00B-4759-YO</a>	—	<a href="#">00D-4759-YO</a>	<a href="#">00F-4759-YO</a>	<a href="#">AJO-9531</a>
F5	—	<a href="#">00B-4723-YO</a>	—	<a href="#">00D-4723-YO</a>	<a href="#">00F-4723-YO</a>	<a href="#">AJO-9321</a>
Biphenyl	—	<a href="#">00B-4622-YO</a>	—	<a href="#">00D-4622-YO</a>	<a href="#">00F-4622-YO</a>	<a href="#">AJO-9208</a>
XB-C18	<a href="#">00A-4496-YO</a>	<a href="#">00B-4496-YO</a>	<a href="#">00C-4496-YO</a>	<a href="#">00D-4496-YO</a>	<a href="#">00F-4496-YO</a>	<a href="#">AJO-8775</a>
C18	<a href="#">00A-4462-YO</a>	<a href="#">00B-4462-YO</a>	<a href="#">00C-4462-YO</a>	<a href="#">00D-4462-YO</a>	<a href="#">00F-4462-YO</a>	<a href="#">AJO-8775</a>
C8	<a href="#">00A-4497-YO</a>	<a href="#">00B-4497-YO</a>	<a href="#">00C-4497-YO</a>	<a href="#">00D-4497-YO</a>	<a href="#">00F-4497-YO</a>	<a href="#">AJO-8777</a>
HILIC	<a href="#">00A-4461-YO</a>	—	—	<a href="#">00D-4461-YO</a>	<a href="#">00F-4461-YO</a>	<a href="#">AJO-8779</a>
Phenyl-Hexyl	—	<a href="#">00B-4495-YO</a>	—	<a href="#">00D-4495-YO</a>	<a href="#">00F-4495-YO</a>	<a href="#">AJO-8781</a>

for 3.0 mm ID

2.6µm Analytical Columns (mm)							SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 4.6	50 x 4.6	75 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	<a href="#">00A-4725-EQ</a>	<a href="#">00B-4725-EQ</a>	—	<a href="#">00D-4725-EQ</a>	<a href="#">00F-4725-EQ</a>	<a href="#">00G-4725-EQ</a>	<a href="#">AJO-9296</a>
PS C18	<a href="#">00A-4780-EQ</a>	<a href="#">00B-4780-EQ</a>	—	<a href="#">00D-4780-EQ</a>	<a href="#">00F-4780-EQ</a>	<a href="#">00G-4780-EQ</a>	<a href="#">AJO-8949</a>
Polar C18	<a href="#">00A-4759-EQ</a>	<a href="#">00B-4759-EQ</a>	—	<a href="#">00D-4759-EQ</a>	<a href="#">00F-4759-EQ</a>	—	<a href="#">AJO-9532</a>
F5	<a href="#">00A-4723-EQ</a>	<a href="#">00B-4723-EQ</a>	—	<a href="#">00D-4723-EQ</a>	<a href="#">00F-4723-EQ</a>	—	<a href="#">AJO-9320</a>
Biphenyl	—	<a href="#">00B-4622-EQ</a>	—	<a href="#">00D-4622-EQ</a>	<a href="#">00F-4622-EQ</a>	—	<a href="#">AJO-9207</a>
XB-C18	—	<a href="#">00B-4496-EQ</a>	<a href="#">00C-4496-EQ</a>	<a href="#">00D-4496-EQ</a>	<a href="#">00F-4496-EQ</a>	—	<a href="#">AJO-8768</a>
C18	<a href="#">00A-4462-EQ</a>	<a href="#">00B-4462-EQ</a>	<a href="#">00C-4462-EQ</a>	<a href="#">00D-4462-EQ</a>	<a href="#">00F-4462-EQ</a>	—	<a href="#">AJO-8768</a>
C8	—	<a href="#">00B-4497-EQ</a>	<a href="#">00C-4497-EQ</a>	<a href="#">00D-4497-EQ</a>	<a href="#">00F-4497-EQ</a>	—	<a href="#">AJO-8770</a>
HILIC	—	<a href="#">00B-4461-EQ</a>	<a href="#">00C-4461-EQ</a>	<a href="#">00D-4461-EQ</a>	<a href="#">00F-4461-EQ</a>	—	<a href="#">AJO-8772</a>
Phenyl-Hexyl	—	<a href="#">00B-4495-EQ</a>	<a href="#">00C-4495-EQ</a>	<a href="#">00D-4495-EQ</a>	<a href="#">00F-4495-EQ</a>	—	<a href="#">AJO-8774</a>

for 4.6 mm ID

1.7µm Minibore Columns (mm)						SecurityGuard™ ULTRA Cartridges <sup>†</sup>
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk	
EVO C18	—	<a href="#">00B-4726-AN</a>	<a href="#">00D-4726-AN</a>	<a href="#">00F-4726-AN</a>	<a href="#">AJO-9298</a>	
F5	—	<a href="#">00B-4722-AN</a>	<a href="#">00D-4722-AN</a>	<a href="#">00F-4722-AN</a>	<a href="#">AJO-9322</a>	
Biphenyl	<a href="#">00A-4628-AN</a>	<a href="#">00B-4628-AN</a>	<a href="#">00D-4628-AN</a>	<a href="#">00F-4628-AN</a>	<a href="#">AJO-9209</a>	
XB-C18	<a href="#">00A-4498-AN</a>	<a href="#">00B-4498-AN</a>	<a href="#">00D-4498-AN</a>	<a href="#">00F-4498-AN</a>	<a href="#">AJO-8782</a>	
C18	<a href="#">00A-4475-AN</a>	<a href="#">00B-4475-AN</a>	<a href="#">00D-4475-AN</a>	<a href="#">00F-4475-AN</a>	<a href="#">AJO-8782</a>	
C8	<a href="#">00A-4499-AN</a>	<a href="#">00B-4499-AN</a>	<a href="#">00D-4499-AN</a>	<a href="#">00F-4499-AN</a>	<a href="#">AJO-8784</a>	
HILIC	<a href="#">00A-4474-AN</a>	<a href="#">00B-4474-AN</a>	<a href="#">00D-4474-AN</a>	—	<a href="#">AJO-8786</a>	
Phenyl-Hexyl	—	<a href="#">00B-4500-AN</a>	<a href="#">00D-4500-AN</a>	<a href="#">00F-4500-AN</a>	<a href="#">AJO-8788</a>	

for 2.1 mm ID

2.6 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
C18	<a href="#">00B-4462-AO</a>	—	—
XB-C18	<a href="#">00B-4496-AO</a>	<a href="#">00D-4496-AO</a>	<a href="#">00F-4496-AO</a>

1.7µm MidBore Columns (mm)					SecurityGuard ULTRA Cartridges <sup>†</sup>
Phases	30 x 3.0	50 x 3.0	100 x 3.0	3/pk	
XB-C18	<a href="#">00A-4498-YO</a>	<a href="#">00B-4498-YO</a>	<a href="#">00D-4498-YO</a>	<a href="#">AJO-8775</a>	
C18	—	<a href="#">00B-4475-YO</a>	<a href="#">00D-4475-YO</a>	<a href="#">AJO-8775</a>	
C8	<a href="#">00A-4499-YO</a>	<a href="#">00B-4499-YO</a>	<a href="#">00D-4499-YO</a>	<a href="#">AJO-8777</a>	
Phenyl	—	—	<a href="#">00D-4500-YO</a>	<a href="#">AJO-8781</a>	
HILIC	—	<a href="#">00B-4474-YO</a>	—	<a href="#">AJO-8779</a>	

for 3.0 mm ID

1.3µm Minibore Columns (mm)		
Phases	30 x 2.1	50 x 2.1
C18	<a href="#">00A-4515-AN</a>	<a href="#">00B-4515-AN</a>

<sup>†</sup> SecurityGuard ULTRA Cartridges require holder, Part No.: [AJO-9000](#)


1.7 µm Microbore Columns (mm)			
Phases	50 x 1.0	100 x 1.0	150 x 1.0
EVO C18	<a href="#">00B-4726-AO</a>	<a href="#">00D-4726-AO</a>	<a href="#">00F-4726-AO</a>
Biphenyl	<a href="#">00B-4628-AO</a>	<a href="#">00D-4628-AO</a>	—

## Setting the Standard for pH Method Development

Rugged reversed phase HPLC columns that offer extended lifetime at extreme pH conditions and excellent stability for reproducible, high efficiency separations.

- Take full advantage of high and low pH conditions (pH 1-12) to manipulate selectivity
- Expect longer column lifetime with patented TWIN-NX™ technology
- For analytical and preparative separations of basic and acidic compounds

Phase	Description	USP Classification
<b>NX-C18</b>	The most rugged Gemini column, offering 5 times the durability of previous generation hybrid columns	L1
<b>C6-Phenyl</b>	A low bleed phenyl phase. For UV and MS detection, which offers an aromatic selectivity complementary to C18 phases	L11
<b>C18</b>	Selectivity, high structural integrity and increased loadability for preparative and purification applications in pre-packed columns and bulk media	L1



# Gemini™ Ordering Information

3 µm Microbore, Minibore and MidBore™ Columns (mm)								SecurityGuard™ Cartridges (mm)			
Phases	50 x 1.0	20 x 2.0	30 x 2.0	50 x 2.0	100 x 2.0	150 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*	
C18	<a href="#">00B-4439-A0</a>	<a href="#">00M-4439-B0</a>	<a href="#">00A-4439-B0</a>	<a href="#">00B-4439-B0</a>	<a href="#">00D-4439-B0</a>	<a href="#">00F-4439-B0</a>	<a href="#">00B-4439-Y0</a>	<a href="#">00D-4439-Y0</a>	<a href="#">00F-4439-Y0</a>	<a href="#">AJ0-7596</a>	
C6-Phenyl	—	—	—	<a href="#">00B-4443-B0</a>	<a href="#">00D-4443-B0</a>	<a href="#">00F-4443-B0</a>	<a href="#">00B-4443-Y0</a>	<a href="#">00D-4443-Y0</a>	<a href="#">00F-4443-Y0</a>	<a href="#">AJ0-7914</a>	
NX-C18	<a href="#">00B-4453-A0</a>	<a href="#">00M-4453-B0</a>	<a href="#">00A-4453-B0</a>	<a href="#">00B-4453-B0</a>	<a href="#">00D-4453-B0</a>	<a href="#">00F-4453-B0</a>	<a href="#">00B-4453-Y0</a>	<a href="#">00D-4453-Y0</a>	<a href="#">00F-4453-Y0</a>	<a href="#">AJ0-8367</a>	

for 2.0-3.0 mm ID

3 µm Analytical Columns (mm)						SecurityGuard Cartridges (mm)	
Phases	20 x 4.0	30 x 4.6	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
C18	<a href="#">00M-4439-D0</a>	<a href="#">00A-4439-E0</a>	<a href="#">00B-4439-E0</a>	<a href="#">00D-4439-E0</a>	<a href="#">00F-4439-E0</a>	<a href="#">00G-4439-E0</a>	<a href="#">AJ0-7597</a>
C6-Phenyl	<a href="#">00A-4443-E0</a>	<a href="#">00A-4443-E0</a>	<a href="#">00B-4443-E0</a>	<a href="#">00D-4443-E0</a>	<a href="#">00F-4443-E0</a>	<a href="#">00G-4443-E0</a>	<a href="#">AJ0-7915</a>
NX-C18	—	<a href="#">00A-4453-E0</a>	<a href="#">00B-4453-E0</a>	<a href="#">00D-4453-E0</a>	<a href="#">00F-4453-E0</a>	<a href="#">00G-4453-E0</a>	<a href="#">AJ0-8368</a>

for 3.2-8.0 mm ID

5 µm Minibore and MidBore Columns (mm)								SecurityGuard Cartridges (mm)	
Phases	30 x 2.0	50 x 2.0	150 x 2.0	250 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	250 x 3.0	4 x 2.0*
C18	<a href="#">00A-4435-B0</a>	<a href="#">00B-4435-B0</a>	<a href="#">00F-4435-B0</a>	<a href="#">00G-4435-B0</a>	<a href="#">00B-4435-Y0</a>	<a href="#">00D-4435-Y0</a>	<a href="#">00F-4435-Y0</a>	<a href="#">00G-4435-Y0</a>	<a href="#">AJ0-7596</a>
C6-Phenyl	—	<a href="#">00B-4444-B0</a>	<a href="#">00F-4444-B0</a>	—	<a href="#">00B-4444-Y0</a>	—	<a href="#">00F-4444-Y0</a>	<a href="#">00G-4444-Y0</a>	<a href="#">AJ0-7914</a>
NX-C18	<a href="#">00A-4454-B0</a>	<a href="#">00B-4454-B0</a>	<a href="#">00F-4454-B0</a>	—	<a href="#">00B-4454-Y0</a>	<a href="#">00D-4454-Y0</a>	<a href="#">00F-4454-Y0</a>	<a href="#">00G-4454-Y0</a>	<a href="#">AJ0-8367</a>

for 2.0-3.0 mm ID

5 µm Analytical Columns (mm)						SecurityGuard Cartridges (mm)
Phases	30 x 4.6	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
C18	<a href="#">00A-4435-E0</a>	<a href="#">00B-4435-E0</a>	<a href="#">00D-4435-E0</a>	<a href="#">00F-4435-E0</a>	<a href="#">00G-4435-E0</a>	<a href="#">AJ0-7597</a>
C6-Phenyl	—	<a href="#">00B-4444-E0</a>	<a href="#">00D-4444-E0</a>	<a href="#">00F-4444-E0</a>	<a href="#">00G-4444-E0</a>	<a href="#">AJ0-7915</a>
NX-C18	—	<a href="#">00B-4454-E0</a>	<a href="#">00D-4454-E0</a>	<a href="#">00F-4454-E0</a>	<a href="#">00G-4454-E0</a>	<a href="#">AJ0-8368</a>

for 3.2-8.0 mm ID

\*SecurityGuard Analytical Cartridges require holder, Part No.: [KJ0-4282](#)

# Vitamin Testing



[www.phenomenex.com](http://www.phenomenex.com)

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**Australia**  
t: +61 (0)2-9428-6444  
[auinfo@phenomenex.com](mailto:auinfo@phenomenex.com)

**Austria**  
t: +43 (0)1-319-1301  
[anfrage@phenomenex.com](mailto:anfrage@phenomenex.com)

**Belgium**  
t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
[beinfo@phenomenex.com](mailto:beinfo@phenomenex.com)

**Canada**  
t: +1 (800) 543-3681  
[info@phenomenex.com](mailto:info@phenomenex.com)

**China**  
t: +86 400-606-8099  
[cninfo@phenomenex.com](mailto:cninfo@phenomenex.com)

**Czech Republic**  
t: +420 272 017 077  
[cz-info@phenomenex.com](mailto:cz-info@phenomenex.com)

**Denmark**  
t: +45 4824 8048  
[nordicinfo@phenomenex.com](mailto:nordicinfo@phenomenex.com)

**Finland**  
t: +358 (0)9 4789 0063  
[nordicinfo@phenomenex.com](mailto:nordicinfo@phenomenex.com)

**France**  
t: +33 (0)1 30 09 21 10  
[franceinfo@phenomenex.com](mailto:franceinfo@phenomenex.com)

**Germany**  
t: +49 (0)6021-58830-0  
[anfrage@phenomenex.com](mailto:anfrage@phenomenex.com)

**Hong Kong**  
t: +852 6012 8162  
[hkinfo@phenomenex.com](mailto:hkinfo@phenomenex.com)

**India**  
t: +91 (0)40-3012 2400  
[indiainfo@phenomenex.com](mailto:indiainfo@phenomenex.com)

**Indonesia**  
t: +62 21 5010 9707  
[indoinfo@phenomenex.com](mailto:indoinfo@phenomenex.com)

**Ireland**  
t: +353 (0)1 247 5405  
[eireinfo@phenomenex.com](mailto:eireinfo@phenomenex.com)

**Italy**  
t: +39 051 6327511  
[italiainfo@phenomenex.com](mailto:italiainfo@phenomenex.com)

**Japan**  
t: +81 (0) 120-149-262  
[jpinfo@phenomenex.com](mailto:jpinfo@phenomenex.com)

**Luxembourg**  
t: +31 (0)30-2418700  
[nlinfo@phenomenex.com](mailto:nlinfo@phenomenex.com)

**Mexico**  
t: 01-800-844-5226  
[tecnicomx@phenomenex.com](mailto:tecnicomx@phenomenex.com)

**The Netherlands**  
t: +31 (0)30-2418700  
[nlinfo@phenomenex.com](mailto:nlinfo@phenomenex.com)

**New Zealand**  
t: +64 (0)9-4780951  
[nzinfo@phenomenex.com](mailto:nzinfo@phenomenex.com)

**Norway**  
t: +47 810 02 005  
[nordicinfo@phenomenex.com](mailto:nordicinfo@phenomenex.com)

**Poland**  
t: +48 22 104 21 72  
[pl-info@phenomenex.com](mailto:pl-info@phenomenex.com)

**Portugal**  
t: +351 221 450 488  
[ptinfo@phenomenex.com](mailto:ptinfo@phenomenex.com)

**Singapore**  
t: +65 800-852-3944  
[sginfo@phenomenex.com](mailto:sginfo@phenomenex.com)

**Slovakia**  
t: +420 272 017 077  
[sk-info@phenomenex.com](mailto:sk-info@phenomenex.com)

**Spain**  
t: +34 91-413-8613  
[espinfo@phenomenex.com](mailto:espinfo@phenomenex.com)

**Sweden**  
t: +46 (0)8 611 6950  
[nordicinfo@phenomenex.com](mailto:nordicinfo@phenomenex.com)

**Switzerland**  
t: +41 (0)61 692 20 20  
[swissinfo@phenomenex.com](mailto:swissinfo@phenomenex.com)

**BE-HAPPY™**  
GUARANTEE

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.  
[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

Terms and Conditions  
Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**  
Axia, Gemini, TWIN-NX, Luna, Kinetex, Strata, Impact, Verex, Novum, β-Gone, BE-HAPPY, Presston, Phree, MidBore, and SecurityGuard are trademarks of Phenomenex. Agilent is a registered trademark of Agilent Technologies, Inc. SCIEX and QTRAP are registered trademarks and Triple Quad, API 4000 and API 4500 are trademarks of AB SCIEX Pte. Ltd. AB SCIEX is being used under license. IMCStzyme is a registered trademark of Integrated Micro-chromatography Systems. Santoprene is a trademark of Exxon Mobil Corporation. Cohesive is a registered trademark of Thermo Fisher Scientific Inc.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145  
Gemini is patented by Phenomenex. U.S. Patent Nos. 7, 563, 367 and 8, 658, 038 and foreign counterparts.  
*FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.*

© 2022 Phenomenex, Inc. All rights reserved.