

Arsenic Species

Column: Hamilton® 10 µm PRP-X100, LC Column 250 x 4.1 mm, Ea

Dimensions: 250 x 4.1 mm ID

Order No: CH0-0687

Elution Type: Gradient

Eluent A: 20 mM Phosphate buffer pH 6.0

Eluent B: 100 mM Phosphate buffer pH 6.0

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	100	0
	2	2	100	0
	3	2.1	50	50
	4	5.1	50	50
	5	5.2	100	0
	6	12.2	100	0

Flow Rate: 0 mL/min

Col. Temp.: ambient

Detection: Unknown @ nm (ambient)

Analyst Note: SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, RP-1 4 x 3.0mm, 10/Pk Part No.: AJ0-5809

- Holder Part No.: KJ0-4282



Products used in this application:

Hamilton®

No Chromatogram Available

ANALYTES:

- 1 As (III)
- 2 DMA
- 3 MMA
- 4 As (V)

