


Inorganic Anions

Column: N/A
Dimensions: N/A
Order No: N/A
Elution Type: N/A
Eluent A: 0.03 M Sodium salicylate, pH 4.0
Gradient

Step No.	Time (min)	Pct A
1	0	100

Profile:
Flow Rate: 0.5 mL/min
Col. Temp.: ambient
Detection: Refractive Index (RI) @ nm (ambient)
Analyst Note: SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, SAX 4 x 3.0mm, 10/Pk Part No.: AJ0-4311
- Holder Part No.: KJ0-4282



ANALYTES:

- 1 Phosphate
- 2 Chloride
- 3 Nitrate
- 4 Sulfate

